



LAPAROSCOPY NEW MODULAR SYSTEM



1. Modular System | 5mm | 10mm

2. Modular System | 3mm

1.1 | MODULAR HANDLES

STANDARD HANDLES
OPTIONAL HANDLES

1.2 | SHEATHS

STANDARD SHEATHS
OPTIONAL SHEATHS

1.3 | INSTRUMENTS

SCISSORS
GRASPING FORCEPS
SPOON FORCEPS
PUNCHES
NEEDLE HOLDERS

2.1 | MODULAR HANDLES

STANDARD HANDLES
OPTIONAL HANDLES

2.2 | INSTRUMENTS

SCISSORS
GRASPING FORCEPS

1. MODULAR SYSTEM | 5MM | 10MM

3-PART-SYSTEM

Ø 5mm

Ø 10mm

1 Handle


2 Sheath

3 Insert




1. MODULAR SYSTEM | 5MM | 10MM

ASSEMBLING



1. MODULAR SYSTEM | 5MM | 10MM

DISASSEMBLING



1.1 MODULAR HANDLES

Ø 5mm

Ø 10mm

HF CONNECTION

ASSEMBLING SHEATH

RATCHET ON/OFF

ROTATION MODULE

RELEASE KEY

Ø 5mm

Ø 10mm

1.1 | MODULAR SYSTEM

STANDARD MODULAR HANDLES

Ordering no.	Description
EC-001-001	Laparoscopic modular handle with HF-connection



Ordering no.	Description
ED-001-001	Laparoscopic modular handle with disengageable ratchet



Ordering no.	Description
ED-001-009	Laparoscopic modular handle with french ratchet



Ø 5mm

Ø 10mm

1.1 | MODULAR SYSTEM

STANDARD MODULAR HANDLES

Ordering no.	Description
ED-001-015	Laparoscopic modular handle with disengageable and stepless ratchet, ratchet release key



Ordering no.	Description
ED-001-016	Laparoscopic modular handle with stepless ratchet ratchet release key



Ø 5mm

Ø 10mm

1.2 | MODULAR SYSTEM

OPTIONAL MODULAR HANDLES


Ordering no.	Description
ED-001-011	Ring handle used for modular laparoscopic system without ratchet




Ordering no.	Description
ED-001-007	Ring handle used for modular laparoscopic system with ratchet



Ordering no.	Description
ED-001-010	Metal Axial handle used for modular system with spring and ratchet



Ordering no.	Description
ED-001-012	Metal Axial handle used for modular system with spring and without ratchet



The handles are **not available** as a complete instrument

Ø 5mm

Ø 10mm

1.2 | MODULAR SYSTEM

OPTIONAL MODULAR HANDLES

Ordering no.	Description
ED-001-013	PISTOL Metal Axial handle used for modular system with spring and ratchet



Ordering no.	Description
ED-001-014	PISTOL Metal Axial handle used for modular system with spring and without ratchet



The handles are **not available** as a complete instrument

1.2 SHEATHS

Ø 5mm

Ø 10mm



Standard sheath
without luer lock



Optional sheath
with luer lock

1.2 | MODULAR SYSTEM

STANDARD SHEATHS FOR LAPAROSCOPIC HANDLES



Ordering no.	Description	Diameter	Length
EC-001-003-25	Insulated sheaths for modular system, without Luer Lock	Ø 5mm	250mm
EC-001-003			360mm
EC-001-005			450mm



Ordering no.	Description	Diameter	Length
EC-001-004	Insulated sheaths for modular system, without Luer Lock	Ø 10mm	360mm
EC-001-006			450mm

Ø 5mm

1.2 | MODULAR SYSTEM

OPTIONAL SHEATHS FOR LAPAROSCOPIC HANDLES



Ordering no.	Description	Diameter	Length
EC-001-003-25-LL	Tube sheaths for modular system, with Luer Lock	Ø 5mm	250mm
EC-001-003-LL			360mm
EC-001-005-LL			450mm

Ø 10mm



Ordering no.	Description	Diameter	Length
EC-001-004-LL	Tube sheaths for modular system, with Luer Lock	Ø 10mm	360mm
EC-001-006-LL			450mm

1.2 | MODULAR SYSTEM

OPTIONAL SHEATHS FOR LAPAROSCOPIC HANDLES



Ordering no.	Description	Diameter	Length
EC-001-012	Tube sheaths for modular system, Ø 10mm reduced to Ø 5mm, without Luer Lock	Ø 10mm	360mm
EC-001-014			450mm



Ordering no.	Description	Diameter	Length
EC-001-012-LL	Tube sheaths for modular system, Ø 10mm reduced to Ø 5mm, with Luer Lock	Ø 10mm	360mm
EC-001-014-LL			450mm

1.3 INSTRUMENTS

Ø 5mm

Ø 10mm

COMPLETE INSTRUMENT
HANDLE **WITHOUT RATCHET**
AND **STANDARD SHEATH**



COMPLETE INSTRUMENT
HANDLE **WITH RATCHET**
AND **STANDARD SHEATH**



1.3 | MODULAR SYSTEM


INSTRUMENTS | SCISSORS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	EC-005-021-25	 10.5mm	Peritoneal scissors, serrated , single action	EC-005-023-25
360mm	EC-005-021			EC-005-023
450mm	EC-005-022			EC-005-024
250mm	EC-005-081-25	 11mm	Peritoneal scissors, serrated , slim , single action	EC-005-083-25
360mm	EC-005-081			EC-005-083
450mm	EC-005-082			EC-005-084
250mm	EC-005-051-25	 5.5mm	Micro scissors, straight , single action	EC-005-053-25
360mm	EC-005-051			EC-005-053
450mm	EC-005-052			EC-005-054
250mm	EC-005-351-25	 5.3mm	Micro scissors, curved , very delicate , single action	EC-005-353-25
360mm	EC-005-351			EC-005-353
450mm	EC-005-352			EC-005-354
250mm	EC-005-031-25	 9mm	Hook scissors, single action	EC-005-033-25
360mm	EC-005-031			EC-005-033
450mm	EC-005-032			EC-005-034
250mm	EC-005-041-25	 18mm	Scissors, straight , double action	EC-005-043-25
360mm	EC-005-041			EC-005-043
450mm	EC-005-042			EC-005-044



Scissors sets can only be ordered with handles
EC-001-001, **ED-001-011** and **ED-001-012**

1.3 | MODULAR SYSTEM

INSTRUMENTS | SCISSORS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	EC-005-091-25	 18mm	Supercut scissors, straight , double action	EC-005-093-25
360mm	EC-005-091			EC-005-093
450mm	EC-005-092			EC-005-094
250mm	EC-005-101-25	 20mm	Supercut METZENBAUM scissors, curved, long , double action	EC-005-103-25
360mm	EC-005-101			EC-005-103
450mm	EC-005-102			EC-005-104
250mm	EC-005-001-25	 15mm	METZENBAUM scissors, curved, short , double action	EC-005-003-25
360mm	EC-005-001			EC-005-003
450mm	EC-005-002			EC-005-004
250mm	EC-005-011-25	 17mm	METZENBAUM scissors, curved, long , double action	EC-005-013-25
360mm	EC-005-011			EC-005-013
450mm	EC-005-012			EC-005-014
250mm	EC-005-061-25	 17mm	METZENBAUM scissors, curved, serrated , double action	EC-005-063-25
360mm	EC-005-061			EC-005-063
450mm	EC-005-062			EC-005-064
250mm	EC-005-311-25	 12mm	Mini METZENBAUM scissors, curved, double action	EC-005-313-25
360mm	EC-005-311			EC-005-313
450mm	EC-005-312			EC-005-314



Scissors sets can only be ordered with handles
EC-001-001, ED-001-011 and ED-001-012

1.3 | MODULAR SYSTEM



INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	EC-005-151-25	 9.5mm	Dissecting- and grasping forceps, double action	EC-005-153-25
360mm	EC-005-151			EC-005-153
450mm	EC-005-152			EC-005-154
250mm	EC-005-111-25	 15mm	Dissecting- and grasping forceps, fenestrated , double action	EC-005-113-25
360mm	EC-005-111			EC-005-113
450mm	EC-005-112			EC-005-114
250mm	EC-005-161-25	 16mm	Dissecting- and grasping forceps, blunt , double action	EC-005-163-25
360mm	EC-005-161			EC-005-163
450mm	EC-005-162			EC-005-164
250mm	EC-005-481-25	 14mm	Dissecting forceps, flat, tapered , double action	EC-005-483-25
360mm	EC-005-481			EC-005-483
450mm	EC-005-482			EC-005-484
250mm	EC-005-211-25	 13.5mm	Dissecting- and grasping forceps, blunt with spoon , double action	EC-005-213-25
360mm	EC-005-211			EC-005-213
450mm	EC-005-212			EC-005-214
250mm	EC-005-501-25	 13.5mm	Dissecting- and grasping forceps, with big reservoir, blunt , double action	EC-005-503-25
360mm	EC-005-501			EC-005-503
450mm	EC-005-502			EC-005-504

1.3 | MODULAR SYSTEM


INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	EC-005-471-25	 17mm	Strong TC , grasping forceps, rough serrated , double action	EC-005-473-25
360mm	EC-005-471			EC-005-473
450mm	EC-005-472			EC-005-474
250mm	EC-005-241-25	 13.5mm	Dissecting forceps, pointed, with spoon , double action	EC-005-243-25
360mm	EC-005-241			EC-005-243
450mm	EC-005-242			EC-005-244
250mm	EC-005-121-25	 17.5mm	MARYLAND dissecting forceps, double action	EC-005-123-25
360mm	EC-005-121			EC-005-123
450mm	EC-005-122			EC-005-124
250mm	EC-005-441-25	 16.5mm	MARYLAND dissecting forceps, conical , double action	EC-005-443-25
360mm	EC-005-441			EC-005-443
450mm	EC-005-442			EC-005-444
250mm	EC-005-451-25	 29mm	MARYLAND dissecting forceps, long , double action	EC-005-453-25
360mm	EC-005-451			EC-005-453
450mm	EC-005-452			EC-005-454
250mm	EC-005-271-25	 20mm	MARYLAND dissecting forceps, cross serrated , double action	EC-005-273-25
360mm	EC-005-271			EC-005-273
450mm	EC-005-272			EC-005-274

1.3 | MODULAR SYSTEM


INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	EC-005-181-25	 11mm	KELLY dissecting- and grasping forceps, double action	EC-005-183-25
360mm	EC-005-181			EC-005-183
450mm	EC-005-182			EC-005-184
250mm	EC-005-391-25	 21.5mm	Dissecting- and grasping forceps, fenestrated , double action	EC-005-393-25
360mm	EC-005-391			EC-005-393
450mm	EC-005-392			EC-005-394
250mm	EC-005-521-25	 17.5mm	DOLPHIN dissecting forceps, double action	EC-005-523-25
360mm	EC-005-521			EC-005-523
450mm	EC-005-522			EC-005-524
250mm	EC-005-131-25	 16mm	DOLPHIN dissecting forceps, with reservoir , double action	EC-005-133-25
360mm	EC-005-131			EC-005-133
450mm	EC-005-132			EC-005-134
250mm	EC-005-511-25	 13.5mm	DOLPHIN dissecting forceps, short , double action	EC-005-513-25
360mm	EC-005-511			EC-005-513
450mm	EC-005-512			EC-005-514
250mm	EC-005-141-25	 26mm	DOLPHIN dissecting forceps, long , double action	EC-005-143-25
360mm	EC-005-141			EC-005-143
450mm	EC-005-142			EC-005-144

1.3 | MODULAR SYSTEM


INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	EC-005-301-25	 13.5mm	REDDICK-OLSEN dissecting- and grasping forceps, double action	EC-005-303-25
360mm	EC-005-301			EC-005-303
450mm	EC-005-302			EC-005-304
250mm	EC-005-231-25	 18mm	Slim dissecting forceps, cross serrated , single action	EC-005-233-25
360mm	EC-005-231			EC-005-233
450mm	EC-005-232			EC-005-234
250mm	EC-005-331-25	 16.5mm	Dissecting forceps, fine serrated, straight , single action	EC-005-333-25
360mm	EC-005-331			EC-005-333
450mm	EC-005-332			EC-005-334
250mm	EC-005-431-25	 13.5mm	Atraumatic dissecting- and grasping forceps, fenestrated, curved , single action	EC-005-433-25
360mm	EC-005-431			EC-005-433
450mm	EC-005-432			EC-005-434
250mm	EC-005-421-25	 20mm	90° MIXTER dissecting- and grasping forceps, double action	EC-005-423-25
360mm	EC-005-421			EC-005-423
450mm	EC-005-422			EC-005-424
250mm	EC-005-171-25	 17.5mm	45° MIXTER dissecting- and grasping forceps, double action	EC-005-173-25
360mm	EC-005-171			EC-005-173
450mm	EC-005-172			EC-005-174

1.3 | MODULAR SYSTEM



INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-011-25	 17.5mm	FUNDUS grasping forceps, double action	ED-005-013-25
360mm	ED-005-011			ED-005-013
450mm	ED-005-012			ED-005-014
250mm	ED-005-471-25	 21.5mm	Maxi FUNDUS grasping forceps, double action	ED-005-473-25
360mm	ED-005-471			ED-005-473
450mm	ED-005-472			ED-005-474
250mm	ED-005-491-25	 19.5mm	Dissecting- and grasping forceps, curved, fenestrated, double action	ED-005-493-25
360mm	ED-005-491			ED-005-493
450mm	ED-005-492			ED-005-494
250mm	ED-005-231-25	 12mm	Strong grasping forceps, fenestrated, double action	ED-005-233-25
360mm	ED-005-231			ED-005-233
450mm	ED-005-232			ED-005-234
250mm	ED-005-041-25	 20mm	ENDO CLINCH grasping forceps, rounded, fenestrated, double action	ED-005-043-25
360mm	ED-005-041			ED-005-043
450mm	ED-005-042			ED-005-044
250mm	ED-005-561-25	 20mm	ENDO CLINCH grasping forceps, fenestrated, double action	ED-005-563-25
360mm	ED-005-561			ED-005-563
450mm	ED-005-562			ED-005-564

1.3 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-731-25	 21mm	BABCOCK grasping forceps, double action	ED-005-733-25
360mm	ED-005-731			ED-005-733
450mm	ED-005-732			ED-005-734
250mm	ED-005-371-25	 21mm	BABCOCK grasping forceps, short , double action	ED-005-373-25
360mm	ED-005-371			ED-005-373
450mm	ED-005-372			ED-005-374
250mm	ED-005-051-25	 31mm	BABCOCK grasping forceps, double action	ED-005-053-25
360mm	ED-005-051			ED-005-053
450mm	ED-005-052			ED-005-054
250mm	ED-005-631-25	 33mm	BABCOCK paddle grasping forceps, double action	ED-005-633-25
360mm	ED-005-631			ED-005-633
450mm	ED-005-632			ED-005-634
250mm	ED-005-081-25	 55mm	DEBAKEY grasping forceps, double action	ED-005-083-25
360mm	ED-005-081			ED-005-083
450mm	ED-005-082			ED-005-084
250mm	ED-005-591-25	 40mm	DEBAKEY grasping forceps, fenestrated, long version , double action	ED-005-593-25
360mm	ED-005-591			ED-005-593
450mm	ED-005-592			ED-005-594

1.3 | MODULAR SYSTEM




INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-601-25	 45mm	DEBAKEY grasping forceps, long, double action	ED-005-603-25
360mm	ED-005-601			ED-005-603
450mm	ED-005-602			ED-005-604
250mm	ED-005-611-25	 28mm	DEBAKEY grasping forceps, short, double action	ED-005-613-25
360mm	ED-005-611			ED-005-613
450mm	ED-005-612			ED-005-614
250mm	ED-005-331-25	 33.5mm	DEBAKEY/MARYLAND grasping forceps, atraumatic, double action	ED-005-333-25
360mm	ED-005-331			ED-005-333
450mm	ED-005-332			ED-005-334
250mm	ED-005-021-25	 13mm	Cobra grasping forceps, 2x4 teeth, double action	ED-005-023-25
360mm	ED-005-021			ED-005-023
450mm	ED-005-022			ED-005-024
250mm	ED-005-341-25	 11mm	Cobra grasping forceps, 2x4 teeth, single action	ED-005-343-25
360mm	ED-005-341			ED-005-343
450mm	ED-005-342			ED-005-344
250mm	ED-005-211-25	 13.5mm	Grasping forceps, 1x2 teeth, double action	ED-005-213-25
360mm	ED-005-211			ED-005-213
450mm	ED-005-212			ED-005-214

1.3 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-521-25	 13.5mm	Grasping forceps, 1x2 teeth , single action	ED-005-523-25
360mm	ED-005-521			ED-005-523
450mm	ED-005-522			ED-005-524
250mm	ED-005-511-25	 25mm	Grasping forceps, double teeth , double action	ED-005-513-25
360mm	ED-005-511			ED-005-513
450mm	ED-005-512			ED-005-514
250mm	ED-005-101-25	 18.5mm	ALLIS grasping forceps, double action	ED-005-103-25
360mm	ED-005-101			ED-005-103
450mm	ED-005-102			ED-005-104
250mm	ED-005-361-25	 12mm	ALLIS atraumatic forceps, 3x4 teeth , double action	ED-005-363-25
360mm	ED-005-361			ED-005-363
450mm	ED-005-362			ED-005-364
250mm	ED-005-621-25	 10mm	ALLIS atraumatic forceps, short , double action	ED-005-623-25
360mm	ED-005-621			ED-005-623
450mm	ED-005-622			ED-005-624
250mm	ED-005-191-25	 11mm	ALLIS grasping forceps, single action	ED-005-193-25
360mm	ED-005-191			ED-005-193
450mm	ED-005-192			ED-005-194

1.3 | MODULAR SYSTEM



INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-381-25	 11mm	Alligator grasping forceps, single action	ED-005-383-25
360mm	ED-005-381			ED-005-383
450mm	ED-005-382			ED-005-384
250mm	ED-005-031-25	 11mm	Alligator grasping- and biopsy forceps, double action	ED-005-033-25
360mm	ED-005-031			ED-005-033
450mm	ED-005-032			ED-005-034
250mm	ED-005-201-25	 13.5mm	Alligator grasping forceps, double action	ED-005-203-25
360mm	ED-005-201			ED-005-203
450mm	ED-005-202			ED-005-204
250mm	ED-005-311-25	 15mm	Duck Bill forceps, short, fenestrated, single action	ED-005-313-25
360mm	ED-005-311			ED-005-313
450mm	ED-005-312			ED-005-314
250mm	ED-005-171-25	 20mm	Duck Bill forceps, fenestrated, single action	ED-005-173-25
360mm	ED-005-171			ED-005-173
450mm	ED-005-172			ED-005-174
250mm	ED-005-321-25	 40mm	Duck Bill forceps, fenestrated, long, single action	ED-005-323-25
360mm	ED-005-321			ED-005-323
450mm	ED-005-322			ED-005-324

1.3 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-181-25	 21.5mm	Wave style grasping forceps, atraumatic , single action	ED-005-183-25
360mm	ED-005-181			ED-005-183
450mm	ED-005-182			ED-005-184
250mm	ED-005-061-25	 26.5mm	Wave style grasping forceps, atraumatic , double action	ED-005-063-25
360mm	ED-005-061			ED-005-063
450mm	ED-005-062			ED-005-064
250mm	ED-005-281-25	 22mm	Waved grasping forceps, with teeth , single action	ED-005-283-25
360mm	ED-005-281			ED-005-283
450mm	ED-005-282			ED-005-284
250mm	ED-005-071-25	 15mm	Slim style forceps, atraumatic, cross serrated , double action	ED-005-073-25
360mm	ED-005-071			ED-005-073
450mm	ED-005-072			ED-005-074
250mm	ED-005-461-25	 14mm	Slim style forceps, atraumatic, short, cross serrated , double action	ED-005-463-25
360mm	ED-005-461			ED-005-463
450mm	ED-005-462			ED-005-464
250mm	ED-005-581-25	 23mm	Grasping forceps, rough tooth, fenestrated , double action	ED-005-583-25
360mm	ED-005-581			ED-005-583
450mm	ED-005-582			ED-005-584

1.3 | MODULAR SYSTEM






INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-421-25		Tenaculum forceps, double action	ED-005-423-25
360mm	ED-005-421			ED-005-423
450mm	ED-005-422			ED-005-424
250mm	ED-005-651-25		Tenaculum forceps, curved tooth, double action	ED-005-653-25
360mm	ED-005-651			ED-005-653
450mm	ED-005-652			ED-005-654
250mm	ED-005-431-25		HUNTER forceps, fenestrated, double action	ED-005-433-25
360mm	ED-005-431			ED-005-433
450mm	ED-005-432			ED-005-434
250mm	ED-005-711-25		HUNTER forceps, fenestrated, fine grooved, double action	ED-005-713-25
360mm	ED-005-711			ED-005-713
450mm	ED-005-712			ED-005-714
250mm	ED-005-161-25		Claw forceps, 2x3 teeth, single action	ED-005-163-25
360mm	ED-005-161			ED-005-163
450mm	ED-005-162			ED-005-164
250mm	ED-005-391-25		Claw forceps, aggressive, double action	ED-005-393-25
360mm	ED-005-391			ED-005-393
450mm	ED-005-392			ED-005-394

1.3 | MODULAR SYSTEM






INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-091-25		BOWEL intestinal grasping forceps, double action	ED-005-093-25
360mm	ED-005-091			ED-005-093
450mm	ED-005-092			ED-005-094
250mm	ED-005-111-25		BOWEL grasping forceps, double action	ED-005-113-25
360mm	ED-005-111			ED-005-113
450mm	ED-005-112			ED-005-114
250mm	ED-005-151-25		DORSEY intestinal grasping forceps, double action	ED-005-153-25
360mm	ED-005-151			ED-005-153
450mm	ED-005-152			ED-005-154
250mm	ED-005-501-25		DORSEY intestinal grasping forceps, short , double action	ED-005-503-25
360mm	ED-005-501			ED-005-503
450mm	ED-005-502			ED-005-504
250mm	ED-005-531-25		DORSEY intestinal grasping forceps, long , double action	ED-005-533-25
360mm	ED-005-531			ED-005-533
450mm	ED-005-532			ED-005-534

1.3 | MODULAR SYSTEM


INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-241-25	 17mm	Strong TC grasping forceps, fine serrated, double action	ED-005-243-25
360mm	ED-005-241			ED-005-243
450mm	ED-005-242			ED-005-244
250mm	ED-005-571-25	 22mm	Retraction forceps, 2x3 teeth, double action	ED-005-573-25
360mm	ED-005-571			ED-005-573
450mm	ED-005-572			ED-005-574
250mm	ED-005-641-25	 13.5mm	REDDICK-OLSEN dissecting forceps, without teeth, double action	ED-005-643-25
360mm	ED-005-641			ED-005-643
450mm	ED-005-642			ED-005-644
250mm	ED-005-671-25	 14mm	Grasping forceps, with fine crosswise teeth, double action	ED-005-673-25
360mm	ED-005-671			ED-005-673
450mm	ED-005-672			ED-005-674
250mm	ED-005-681-25	 23mm	TC grasping forceps, fenestrated, double action	ED-005-683-25
360mm	ED-005-681			ED-005-683
450mm	ED-005-682			ED-005-684
250mm	ED-005-691-25	 15mm	Grasping forceps, 1x2 teeth, double action	ED-005-693-25
360mm	ED-005-691			ED-005-693
450mm	ED-005-692			ED-005-694

1.3 | MODULAR SYSTEM



INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-721-25	 11.5mm	Grasping forceps, distal with row teeth , double action	ED-005-723-25
360mm	ED-005-721			ED-005-723
450mm	ED-005-722			ED-005-724
250mm	ED-005-741-25	 29mm	Grasping forceps, with fine crosswise teeth , double action	ED-005-743-25
360mm	ED-005-741			ED-005-743
450mm	ED-005-742			ED-005-744
250mm	ED-005-751-25	 24mm	Grasping forceps, DeBackey teeth , double action	ED-005-753-25
360mm	ED-005-751			ED-005-753
450mm	ED-005-752			ED-005-754
250mm	ED-005-771-25	 18mm	SWANSTROM grasping forceps, double action	ED-005-773-25
360mm	ED-005-771			ED-005-773
450mm	ED-005-772			ED-005-774
250mm	ED-005-481-25	 12mm	Micro grasping forceps, fenestrated , double action	ED-005-483-25
360mm	ED-005-481			ED-005-483
450mm	ED-005-482			ED-005-484
250mm	ED-005-661-25	 13mm	Spoon biopsy forceps, teethed , double action	ED-005-663-25
360mm	ED-005-661			ED-005-663
450mm	ED-005-662			ED-005-664

1.3 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-541-25	 12mm	Micro grasping forceps, fenestrated , double action	ED-005-543-25
360mm	ED-005-541			ED-005-543
450mm	ED-005-542			ED-005-544
250mm	ED-005-551-25	 10mm	Micro grasping forceps, slim version, fenestrated , double action	ED-005-553-25
360mm	ED-005-551			ED-005-553
450mm	ED-005-552			ED-005-554
250mm	ED-005-761-25	 17mm	Micro Claw grasping forceps, double action	ED-005-763-25
360mm	ED-005-761			ED-005-763
450mm	ED-005-762			ED-005-764

1.3 | MODULAR SYSTEM


INSTRUMENTS | SPOON FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	EC-005-561-25	 15mm	Tumor spoon forceps, single action	EC-005-563-25
360mm	EC-005-561			EC-005-563
450mm	EC-005-562			EC-005-564
250mm	EC-005-251-25	 14.5mm	Spoon biopsy forceps, single action	EC-005-253-25
360mm	EC-005-251			EC-005-253
450mm	EC-005-252			EC-005-254
250mm	EC-005-261-25	 16.5mm	Spoon biopsy forceps, with spike , single action	EC-005-263-25
360mm	EC-005-261			EC-005-263
450mm	EC-005-262			EC-005-264
250mm	EC-005-221-25	 10mm	BLAKESLY biopsy forceps, single action	EC-005-223-25
360mm	EC-005-221			EC-005-223
450mm	EC-005-222			EC-005-224
250mm	EC-005-191-25	 8mm	Spoon biopsy forceps, double action	EC-005-193-25
360mm	EC-005-191			EC-005-193
450mm	EC-005-192			EC-005-194
250mm	EC-005-201-25	 8mm	Spoon biopsy forceps, with spike , double action	EC-005-203-25
360mm	EC-005-201			EC-005-203
450mm	EC-005-202			EC-005-204

1.3 | MODULAR SYSTEM

INSTRUMENTS | PUNCHES

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	ED-005-261-25	 8.5mm	Punch, single action	ED-005-263-25
360mm	ED-005-261			ED-005-263
450mm	ED-005-262			ED-005-264
250mm	ED-005-701-25	 8.5mm	Punch, with spike, single action	ED-005-703-25
360mm	ED-005-701			ED-005-703
450mm	ED-005-702			ED-005-704
250mm	ED-005-451-25	 17mm	Punch, hook shaped, single action	ED-005-453-25
360mm	ED-005-451			ED-005-453
450mm	ED-005-452			ED-005-454

1.3 | MODULAR SYSTEM

INSTRUMENTS | NEEDLE HOLDERS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
250mm	EG-003-062-25		Needleholder, straight, TC-insert, single action	EG-003-066-25
360mm	EG-003-062			EG-003-066
450mm	EG-003-062-45			EG-003-066-45
250mm	EG-003-063-25		Needleholder, curved left, TC-insert, single action	EG-003-067-25
360mm	EG-003-063			EG-003-067
450mm	EG-003-063-45			EG-003-067-45
250mm	EG-003-064-25		Needleholder, curved right, TC-insert, single action	EG-003-068-25
360mm	EG-003-064			EG-003-068
450mm	EG-003-064-45			EG-003-068-45
250mm	EG-003-065-25		Needleholder, straight, selfadjusting, TC-insert, single action	EG-003-069-25
360mm	EG-003-065			EG-003-069
450mm	EG-003-065-45			EG-003-069-45

1.3 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
360mm	EC-010-021		METZENBAUM scissors , straight, long, single action	EC-010-023
450mm	EC-010-022			EC-010-024
360mm	EC-010-121		MARYLAND dissceting forceps, double action	EC-010-123
450mm	EC-010-122			EC-010-124
360mm	EC-010-221		OVERHOLT dissecting forceps, fine teeth, curved, double action	EC-010-223
450mm	EC-010-222			EC-010-224
360mm	EC-010-171		MIXTER dissecting forceps, double action	EC-010-173
450mm	EC-010-172			EC-010-174
360mm	EC-010-191		Spoon biopsy forceps, single action	EC-010-193
450mm	EC-010-192			EC-010-194

1.3 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
360mm	ED-010-261		Claw forceps, 2x3 teeth , single action	ED-010-263
450mm	ED-010-262			ED-010-264
360mm	ED-010-161		Claw forceps, 2x3 teeth , double action	ED-010-163
450mm	ED-010-162			ED-010-164
360mm	ED-010-101		ALLIS grasping forceps, double action	ED-010-103
450mm	ED-010-102			ED-010-104
360mm	ED-010-221		Tenaculum forceps, double action	ED-010-223
450mm	ED-010-222			ED-010-224
360mm	ED-010-111		DORSEY grasping forceps, fenestrated , double action	ED-010-113
450mm	ED-010-112			ED-010-114

1.3 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
360mm	ED-010-051		BABCOCK forceps, double action	ED-010-053
450mm	ED-010-052			ED-010-054
360mm	ED-010-121		BABCOCK grasping forceps, round fenestrated, double action	ED-010-123
450mm	ED-010-122			ED-010-124
360mm	ED-010-131		BABCOCK grasping forceps, cross serrated teeth, double action	ED-010-133
450mm	ED-010-132			ED-010-134
360mm	ED-010-251		DUVALL/COLLIN forceps, double action	ED-010-253
450mm	ED-010-252			ED-010-254
360mm	ED-010-241		PENNINGTON forceps, double action	ED-010-243
450mm	ED-010-242			ED-010-244

1.3 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT ONLY



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
360mm	ED-010-141		KOCHER forceps, jaw length 30mm, double action	ED-010-143
450mm	ED-010-142			ED-010-144
360mm	ED-010-311		KOCHER forceps, jaw length 40mm, double action	ED-010-313
450mm	ED-010-312			ED-010-314
360mm	ED-010-321		KOCHER grasping forceps, jaw length 50mm, double action	ED-010-323
450mm	ED-010-322			ED-010-324
360mm	ED-010-091		DEBAKEY grasping forceps, jaw length 30mm, double action	ED-010-093
450mm	ED-010-092			ED-010-094
360mm	ED-010-081		DEBAKEY grasping forceps, jaw length 40mm, double action	ED-010-083
450mm	ED-010-082			ED-010-084
360mm	ED-010-341		DEBAKEY grasping forceps, jaw length 50mm, double action	ED-010-343
450mm	ED-010-342			ED-010-344

2. MODULAR SYSTEM | 3MM




Ø 3mm

2-PART-SYSTEM





2. MODULAR SYSTEM | 3MM

ASSEMBLING



2. MODULAR SYSTEM | 3MM

DISASSEMBLING



2.1 MODULAR HANDLES

Ø 3mm

ASSEMBLING SHEATH

HF CONNECTION

RATCHET ON/OFF

RELEASE KEY

ROTATION MODULE

2.1 | MODULAR SYSTEM

STANDARD MODULAR HANDLES

Ordering no.	Description
EC-001-001	Laparoscopic modular handle with HF-connection



Ordering no.	Description
ED-001-001	Laparoscopic modular handle with disengageable ratchet



Ordering no.	Description
ED-001-009	Laparoscopic modular handle with french ratchet



2.1 | MODULAR SYSTEM

STANDARD MODULAR HANDLES

Ordering no.	Description
ED-001-015	Laparoscopic modular handle with disengageable and stepless ratchet, ratchet release key



Ordering no.	Description
ED-001-016	Laparoscopic modular handle with stepless ratchet ratchet release key



2.2 | MODULAR SYSTEM

OPTIONAL MODULAR HANDLES


Ordering no.	Description
ED-001-011	Ring handle used for modular laparoscopic system without ratchet




Ordering no.	Description
ED-001-007	Ring handle used for modular laparoscopic system with ratchet



Ordering no.	Description
ED-001-010	Metal Axial handle used for modular system with spring and ratchet



Ordering no.	Description
ED-001-012	Metal Axial handle used for modular system with spring and without ratchet



The handles are **not available** as a complete instrument

2.2 INSTRUMENTS

Ø 3mm

COMPLETE INSTRUMENT
WITH HANDLE **WITHOUT RATCHET**
AND **INSERT WITH SHEATH**



COMPLETE INSTRUMENT
WITH HANDLE **WITH RATCHET**
AND **INSERT WITH SHEATH**



2.2 | MODULAR SYSTEM

INSTRUMENTS | SCISSORS

INSERT WITH SHEATH



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
210mm	EC-003-041-20		Scissors, straight, double action	EC-003-043-20
330mm	EC-003-041			EC-003-043
210mm	EC-003-061-20		Scissors, straight, serrated, single action	EC-003-063-20
330mm	EC-003-061			EC-003-063
210mm	EC-003-001-20		METZENBAUM scissors, short, curved, double action	EC-003-003-20
330mm	EC-003-001			EC-003-003
210mm	EC-003-031-20		Hook scissors, single action	EC-003-033-20
330mm	EC-003-031			EC-003-033
210mm	EC-003-051-20		Micro scissors, straight, single action	EC-003-053-20
330mm	EC-003-051			EC-003-053



Scissors sets can only be ordered with handles
EC-001-001, ED-001-011 and ED-001-012

2.2 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT WITH SHEATH



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
210mm	EC-003-121-20		MARYLAND dissecting forceps, double action	EC-003-123-20
330mm	EC-003-121			EC-003-123
210mm	EC-003-101-20		MARYLAND dissecting forceps, cross serrated teeth, long, double action	EC-003-103-20
330mm	EC-003-101			EC-003-103
210mm	EC-003-111-20		Dissecting- and grasping forceps, fenestrated, double action	EC-003-113-20
330mm	EC-003-111			EC-003-113
210mm	EC-003-081-20		Dissecting- and grasping forceps, fenestrated, curved, double action	EC-003-083-20
330mm	EC-003-081			EC-003-083
210mm	EC-003-151-20		Dissecting- and grasping forceps, with reservoir, blunt, double action	EC-003-153-20
330mm	EC-003-151			EC-003-153
210mm	EC-003-091-20		Dissecting- and grasping forceps, with reservoir, sharpened, double action	EC-003-093-20
330mm	EC-003-091			EC-003-093

2.2 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT WITH SHEATH



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
210mm	EC-003-141-20		DOLPHIN dissecting forceps, long, double action	EC-003-143-20
330mm	EC-003-141			EC-003-143
210mm	EC-003-171-20		MIXTER dissecting- and grasping forceps, 90° , double action	EC-003-173-20
330mm	EC-003-171			EC-003-173
210mm	EC-003-181-20		KELLY dissecting- and grasping forceps, double action	EC-003-183-20
330mm	EC-003-181			EC-003-183
210mm	EC-003-191-20		Spoon biopsy forceps, double action	EC-003-193-20
330mm	EC-003-191			EC-003-193
210mm	EC-003-131-20		Spoon biopsy forceps, with thorn , double action	EC-003-133-20
330mm	EC-003-131			EC-003-133

2.2 | MODULAR SYSTEM

INSTRUMENTS | GRASPING FORCEPS

INSERT WITH SHEATH



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
210mm	ED-003-181-20		Universal grasping forceps, double action	ED-003-183-20
330mm	ED-003-181			ED-003-183
210mm	ED-003-011-20		FUNDUS grasping forceps, double action	ED-003-013-20
330mm	ED-003-011			ED-003-013
210mm	ED-003-071-20		Grasping forceps, cross serrated teeth, sharp , double action	ED-003-073-20
330mm	ED-003-071			ED-003-073
210mm	ED-003-091-20		JOHANN/BOWEL intestinal grasping forceps, double action	ED-003-093-20
330mm	ED-003-091			ED-003-093
210mm	ED-003-161-20		Claw forceps, 2x3 teeth , single action	ED-003-163-20
330mm	ED-003-161			ED-003-163
210mm	ED-003-171-20		JOHANN dissecting- and grasping forceps, single action	ED-003-173-20
330mm	ED-003-171			ED-003-173

2.2 | MODULAR SYSTEM


INSTRUMENTS | GRASPING FORCEPS

INSERT WITH SHEATH



INSTRUMENT COMPLETE



Length	Ordering No.	Tip	Description	Instrument complete
210mm	ED-003-051-20		BABCOCK grasping forceps, double action	ED-003-053-20
330mm	ED-003-051			ED-003-053
210mm	ED-003-421-20		Tenaculum grasping forceps, double action	ED-003-423-20
330mm	ED-003-421			ED-003-423
210mm	ED-003-081-20		DEBAKEY grasping forceps, atraumatic , double action	ED-003-083-20
330mm	ED-003-081			ED-003-083
210mm	ED-003-101-20		ALLIS grasping forceps, double action	ED-003-103-20
330mm	ED-003-101			ED-003-103
210mm	ED-003-031-20		Alligator grasping forceps, double action	ED-003-033-20
330mm	ED-003-031			ED-003-033
210mm	ED-003-111-20		Wave style grasping forceps, atraumatic , single action	ED-003-113-20
330mm	ED-003-111			ED-003-113

ARE YOU INTERESTED
IN OUR PRODUCTS?
WE WILL GLADLY
SUBMIT YOU
AN OFFER.

(if available)

Quantity

WWW.STEMA-MEDIZINTECHNIK.DE

INDEX

EC-001-001	20 , 57	EC-003-131	63	EC-005-051-25	29
EC-001-003	25	EC-003-131-20	63	EC-005-052	29
EC-001-003-25	25	EC-003-133	63	EC-005-053	29
EC-001-003-25-LL	26	EC-003-133-20	63	EC-005-053-25	29
EC-001-003-LL	26	EC-003-141	63	EC-005-054	29
EC-001-004	25	EC-003-141-20	63	EC-005-061	30
EC-001-004-LL	26	EC-003-143	63	EC-005-061-25	30
EC-001-005	25	EC-003-143-20	63	EC-005-062	30
EC-001-005-LL	26	EC-003-151	62	EC-005-063	30
EC-001-006	25	EC-003-151-20	62	EC-005-063-25	30
EC-001-006-LL	26	EC-003-153	62	EC-005-064	30
EC-001-012	27	EC-003-153-20	62	EC-005-081	29
EC-001-012-LL	27	EC-003-171	63	EC-005-081-25	29
EC-001-014	27	EC-003-171-20	63	EC-005-082	29
EC-001-014-LL	27	EC-003-173	63	EC-005-083	29
EC-003-001	61	EC-003-173-20	63	EC-005-083-25	29
EC-003-001-20	61	EC-003-181	63	EC-005-084	29
EC-003-003	61	EC-003-181-20	63	EC-005-091	30
EC-003-003-20	61	EC-003-183	63	EC-005-091-25	30
EC-003-031	61	EC-003-183-20	63	EC-005-092	30
EC-003-031-20	61	EC-003-191	63	EC-005-093	30
EC-003-033	61	EC-003-191-20	63	EC-005-093-25	30
EC-003-033-20	61	EC-003-193	63	EC-005-094	30
EC-003-041	61	EC-003-193-20	63	EC-005-101	30
EC-003-041-20	61	EC-005-001	30	EC-005-101-25	30
EC-003-043	61	EC-005-001-25	30	EC-005-102	30
EC-003-043-20	61	EC-005-002	30	EC-005-103	30
EC-003-051	61	EC-005-003	30	EC-005-103-25	30
EC-003-051-20	61	EC-005-003-25	30	EC-005-104	30
EC-003-053	61	EC-005-004	30	EC-005-111	31
EC-003-053-20	61	EC-005-011	30	EC-005-111-25	31
EC-003-061	61	EC-005-011-25	30	EC-005-112	31
EC-003-061-20	61	EC-005-012	30	EC-005-113	31
EC-003-063	61	EC-005-013	30	EC-005-113-25	31
EC-003-063-20	61	EC-005-013-25	30	EC-005-114	31
EC-003-081	62	EC-005-014	30	EC-005-121	32
EC-003-081-20	62	EC-005-021	29	EC-005-121-25	32
EC-003-083	62	EC-005-021-25	29	EC-005-122	32
EC-003-083-20	62	EC-005-022	29	EC-005-123	32
EC-003-091	62	EC-005-023	29	EC-005-123-25	32
EC-003-091-20	62	EC-005-023-25	29	EC-005-124	32
EC-003-093	62	EC-005-024	29	EC-005-131	33
EC-003-093-20	62	EC-005-031	29	EC-005-131-25	33
EC-003-101	62	EC-005-031-25	29	EC-005-132	33
EC-003-101-20	62	EC-005-032	29	EC-005-133	33
EC-003-103	62	EC-005-033	29	EC-005-133-25	33
EC-003-103-20	62	EC-005-033-25	29	EC-005-134	33
EC-003-111	62	EC-005-034	29	EC-005-141	33
EC-003-111-20	62	EC-005-041	29	EC-005-141-25	33
EC-003-113	62	EC-005-041-25	29	EC-005-142	33
EC-003-113-20	62	EC-005-042	29	EC-005-143	33
EC-003-121	62	EC-005-043	29	EC-005-143-25	33
EC-003-121-20	62	EC-005-043-25	29	EC-005-144	33
EC-003-123	62	EC-005-044	29	EC-005-151	31
EC-003-123-20	62	EC-005-051	29	EC-005-151-25	31

INDEX

EC-005-152	31	EC-005-243	32	EC-005-423-25	34
EC-005-153	31	EC-005-243-25	32	EC-005-424	34
EC-005-153-25	31	EC-005-244	32	EC-005-431	34
EC-005-154	31	EC-005-251	35	EC-005-431-25	34
EC-005-161	31	EC-005-251-25	35	EC-005-432	34
EC-005-161-25	31	EC-005-252	35	EC-005-433	34
EC-005-162	31	EC-005-253	35	EC-005-433-25	34
EC-005-163	31	EC-005-253-25	35	EC-005-434	34
EC-005-163-25	31	EC-005-254	35	EC-005-441	32
EC-005-164	31	EC-005-261	35	EC-005-441-25	32
EC-005-171	34	EC-005-261-25	35	EC-005-442	32
EC-005-171-25	34	EC-005-262	35	EC-005-443	32
EC-005-172	34	EC-005-263	35	EC-005-443-25	32
EC-005-173	34	EC-005-263-25	35	EC-005-444	32
EC-005-173-25	34	EC-005-264	35	EC-005-451	32
EC-005-174	34	EC-005-271	32	EC-005-451-25	32
EC-005-181	33	EC-005-271-25	32	EC-005-452	32
EC-005-181-25	33	EC-005-272	32	EC-005-453	32
EC-005-182	33	EC-005-273	32	EC-005-453-25	32
EC-005-183	33	EC-005-273-25	32	EC-005-454	32
EC-005-183-25	33	EC-005-274	32	EC-005-471	32
EC-005-184	33	EC-005-301	34	EC-005-471-25	32
EC-005-191	35	EC-005-301-25	34	EC-005-472	32
EC-005-191-25	35	EC-005-302	34	EC-005-473	32
EC-005-192	35	EC-005-303	34	EC-005-473-25	32
EC-005-193	35	EC-005-303-25	34	EC-005-474	32
EC-005-193-25	35	EC-005-304	34	EC-005-481	31
EC-005-194	35	EC-005-311	30	EC-005-481-25	31
EC-005-201	35	EC-005-311-25	30	EC-005-482	31
EC-005-201-25	35	EC-005-312	30	EC-005-483	31
EC-005-202	35	EC-005-313	30	EC-005-483-25	31
EC-005-203	35	EC-005-313-25	30	EC-005-484	31
EC-005-203-25	35	EC-005-314	30	EC-005-501	31
EC-005-204	35	EC-005-331	34	EC-005-501-25	31
EC-005-211	31	EC-005-331-25	34	EC-005-502	31
EC-005-211-25	31	EC-005-332	34	EC-005-503	31
EC-005-212	31	EC-005-333	34	EC-005-503-25	31
EC-005-213	31	EC-005-333-25	34	EC-005-504	31
EC-005-213-25	31	EC-005-334	34	EC-005-511	33
EC-005-214	31	EC-005-351	29	EC-005-511-25	33
EC-005-221	35	EC-005-351-25	29	EC-005-512	33
EC-005-221-25	35	EC-005-352	29	EC-005-513	33
EC-005-222	35	EC-005-353	29	EC-005-513-25	33
EC-005-223	35	EC-005-353-25	29	EC-005-514	33
EC-005-223-25	35	EC-005-354	29	EC-005-521	33
EC-005-224	35	EC-005-391	33	EC-005-521-25	33
EC-005-231	34	EC-005-391-25	33	EC-005-522	33
EC-005-231-25	34	EC-005-392	33	EC-005-523	33
EC-005-232	34	EC-005-393	33	EC-005-523-25	33
EC-005-233	34	EC-005-393-25	33	EC-005-524	33
EC-005-233-25	34	EC-005-394	33	EC-005-561	35
EC-005-234	34	EC-005-421	34	EC-005-561-25	35
EC-005-241	32	EC-005-421-25	34	EC-005-562	35
EC-005-241-25	32	EC-005-422	34	EC-005-563	35
EC-005-242	32	EC-005-423	34	EC-005-563-25	35

INDEX

EC-005-564	35
EC-010-021	49
EC-010-022	49
EC-010-023	49
EC-010-024	49
EC-010-121	49
EC-010-122	49
EC-010-123	49
EC-010-124	49
EC-010-171	49
EC-010-172	49
EC-010-173	49
EC-010-174	49
EC-010-191	49
EC-010-192	49
EC-010-193	49
EC-010-194	49
EC-010-221	49
EC-010-222	49
EC-010-223	49
EC-010-224	49
ED-001-001	20, 57
ED-001-007	22, 59
ED-001-009	20, 57
ED-001-010	22, 59
ED-001-011	22, 59
ED-001-012	22, 59
ED-001-013	23
ED-001-014	23
ED-001-015	21, 58
ED-001-016	21, 58
ED-003-011	64
ED-003-011-20	64
ED-003-013	64
ED-003-013-20	64
ED-003-031	65
ED-003-031-20	65
ED-003-033	65
ED-003-033-20	65
ED-003-051	65
ED-003-051-20	65
ED-003-053	65
ED-003-053-20	65
ED-003-071	64
ED-003-071-20	64
ED-003-073	64
ED-003-073-20	64
ED-003-081	65
ED-003-081-20	65
ED-003-083	65
ED-003-083-20	65
ED-003-091	64
ED-003-091-20	64
ED-003-093	64
ED-003-093-20	64
ED-003-101	65
ED-003-101-20	65
ED-003-103	65
ED-003-103-20	65
ED-003-111	65
ED-003-111-20	65
ED-003-113	65
ED-003-113-20	65
ED-003-161	64
ED-003-161-20	64
ED-003-163	64
ED-003-163-20	64
ED-003-171	64
ED-003-171-20	64
ED-003-173	64
ED-003-173-20	64
ED-003-181	64
ED-003-181-20	64
ED-003-183	64
ED-003-183-20	64
ED-003-421	65
ED-003-421-20	65
ED-003-423	65
ED-003-423-20	65
ED-005-011	36
ED-005-011-25	36
ED-005-012	36
ED-005-021	38
ED-005-021-25	38
ED-005-022	38
ED-005-023	38
ED-005-023-25	38
ED-005-024	38
ED-005-031	40
ED-005-031-25	40
ED-005-032	40
ED-005-033	40
ED-005-033-25	40
ED-005-034	40
ED-005-043	36
ED-005-043-25	36
ED-005-044	36
ED-005-051	37
ED-005-051-25	37
ED-005-052	37
ED-005-053	37
ED-005-053-25	37
ED-005-054	37
ED-005-061	41
ED-005-061-25	41
ED-005-062	41
ED-005-063	41
ED-005-063-25	41
ED-005-064	41
ED-005-071	41
ED-005-071-25	41
ED-005-072	41
ED-005-073	41
ED-005-073-25	41
ED-005-074	41
ED-005-081	37
ED-005-081-25	37
ED-005-082	37
ED-005-083	37
ED-005-083-25	37
ED-005-084	37
ED-005-091	43
ED-005-091-25	43
ED-005-092	43
ED-005-093	43
ED-005-093-25	43
ED-005-094	43
ED-005-101	39
ED-005-101-25	39
ED-005-102	39
ED-005-103	39
ED-005-103-25	39
ED-005-104	39
ED-005-111	43
ED-005-111-25	43
ED-005-112	43
ED-005-113	43
ED-005-113-25	43
ED-005-114	43
ED-005-151	43
ED-005-151-25	43
ED-005-152	43
ED-005-153	43
ED-005-153-25	43
ED-005-154	43
ED-005-161	42
ED-005-161-25	42
ED-005-162	42
ED-005-163	42
ED-005-163-25	42
ED-005-164	42
ED-005-171	40
ED-005-171-25	40
ED-005-172	40
ED-005-173	40
ED-005-173-25	40
ED-005-174	40
ED-005-181	41
ED-005-181-25	41
ED-005-182	41
ED-005-183	41
ED-005-183-25	41
ED-005-184	41
ED-005-191	39
ED-005-191-25	39

INDEX

ED-005-192	39	ED-005-341	38	ED-005-471-25	36
ED-005-193	39	ED-005-341-25	38	ED-005-472	36
ED-005-193-25	39	ED-005-342	38	ED-005-481	45
ED-005-194	39	ED-005-343	38	ED-005-481-25	45
ED-005-201	40	ED-005-343-25	38	ED-005-482	45
ED-005-201-25	40	ED-005-344	38	ED-005-483	45
ED-005-202	40	ED-005-361	39	ED-005-483-25	45
ED-005-203	40	ED-005-361-25	39	ED-005-484	45
ED-005-203-25	40	ED-005-362	39	ED-005-491	36
ED-005-204	40	ED-005-363	39	ED-005-491-25	36
ED-005-211	38	ED-005-363-25	39	ED-005-492	36
ED-005-211-25	38	ED-005-364	39	ED-005-493	36
ED-005-212	38	ED-005-371	37	ED-005-493-25	36
ED-005-213	38	ED-005-371-25	37	ED-005-494	36
ED-005-213-25	38	ED-005-372	37	ED-005-501	43
ED-005-214	38	ED-005-373	37	ED-005-501-25	43
ED-005-233	36	ED-005-373-25	37	ED-005-502	43
ED-005-233-25	36	ED-005-374	37	ED-005-503	43
ED-005-234	36	ED-005-381	40	ED-005-503-25	43
ED-005-241	44	ED-005-381-25	40	ED-005-504	43
ED-005-241-25	44	ED-005-382	40	ED-005-511	39
ED-005-242	44	ED-005-383	40	ED-005-511-25	39
ED-005-243	44	ED-005-383-25	40	ED-005-512	39
ED-005-243-25	44	ED-005-384	40	ED-005-513	39
ED-005-244	44	ED-005-391	42	ED-005-513-25	39
ED-005-261	47	ED-005-391-25	42	ED-005-514	39
ED-005-261-25	47	ED-005-392	42	ED-005-521	39
ED-005-262	47	ED-005-393	42	ED-005-521-25	39
ED-005-263	47	ED-005-393-25	42	ED-005-522	39
ED-005-263-25	47	ED-005-394	42	ED-005-523	39
ED-005-264	47	ED-005-421	42	ED-005-523-25	39
ED-005-281	41	ED-005-421-25	42	ED-005-524	39
ED-005-281-25	41	ED-005-422	42	ED-005-531	43
ED-005-282	41	ED-005-423	42	ED-005-531-25	43
ED-005-283	41	ED-005-423-25	42	ED-005-532	43
ED-005-283-25	41	ED-005-424	42	ED-005-533	43
ED-005-284	41	ED-005-431	42	ED-005-533-25	43
ED-005-311	40	ED-005-431-25	42	ED-005-534	43
ED-005-311-25	40	ED-005-432	42	ED-005-541	46
ED-005-312	40	ED-005-433	42	ED-005-541-25	46
ED-005-313	40	ED-005-433-25	42	ED-005-542	46
ED-005-313-25	40	ED-005-434	42	ED-005-543	46
ED-005-314	40	ED-005-451	47	ED-005-543-25	46
ED-005-321	40	ED-005-451-25	47	ED-005-544	46
ED-005-321-25	40	ED-005-452	47	ED-005-551	46
ED-005-322	40	ED-005-453	47	ED-005-551-25	46
ED-005-323	40	ED-005-453-25	47	ED-005-552	46
ED-005-323-25	40	ED-005-454	47	ED-005-553	46
ED-005-324	40	ED-005-461	41	ED-005-553-25	46
ED-005-331	38	ED-005-461-25	41	ED-005-554	46
ED-005-331-25	38	ED-005-462	41	ED-005-563	36
ED-005-332	38	ED-005-463	41	ED-005-563-25	36
ED-005-333	38	ED-005-463-25	41	ED-005-564	36
ED-005-333-25	38	ED-005-464	41	ED-005-571	44
ED-005-334	38	ED-005-471	36	ED-005-571-25	44

INDEX

ED-005-572	44	ED-005-663	45	ED-005-753-25	45
ED-005-573	44	ED-005-663-25	45	ED-005-754	45
ED-005-573-25	44	ED-005-664	45	ED-005-761	46
ED-005-574	44	ED-005-671	44	ED-005-761-25	46
ED-005-581	41	ED-005-671-25	44	ED-005-762	46
ED-005-581-25	41	ED-005-672	44	ED-005-763	46
ED-005-582	41	ED-005-673	44	ED-005-763-25	46
ED-005-583	41	ED-005-673-25	44	ED-005-764	46
ED-005-583-25	41	ED-005-674	44	ED-005-771	45
ED-005-584	41	ED-005-681	44	ED-005-771-25	45
ED-005-591	37	ED-005-681-25	44	ED-005-772	45
ED-005-591-25	37	ED-005-682	44	ED-005-773	45
ED-005-592	37	ED-005-683	44	ED-005-773-25	45
ED-005-593	37	ED-005-683-25	44	ED-005-774	45
ED-005-593-25	37	ED-005-684	44	ED-010-051	51
ED-005-594	37	ED-005-691	44	ED-010-052	51
ED-005-601	38	ED-005-691-25	44	ED-010-053	51
ED-005-601-25	38	ED-005-692	44	ED-010-054	51
ED-005-602	38	ED-005-693	44	ED-010-081	52
ED-005-603	38	ED-005-693-25	44	ED-010-082	52
ED-005-603-25	38	ED-005-694	44	ED-010-083	52
ED-005-604	38	ED-005-701	47	ED-010-084	52
ED-005-611	38	ED-005-701-25	47	ED-010-091	52
ED-005-611-25	38	ED-005-702	47	ED-010-092	52
ED-005-612	38	ED-005-703	47	ED-010-093	52
ED-005-613	38	ED-005-703-25	47	ED-010-094	52
ED-005-613-25	38	ED-005-704	47	ED-010-101	50
ED-005-614	38	ED-005-711	42	ED-010-102	50
ED-005-621	39	ED-005-711-25	42	ED-010-103	50
ED-005-621-25	39	ED-005-712	42	ED-010-104	50
ED-005-622	39	ED-005-713	42	ED-010-111	50
ED-005-623	39	ED-005-713-25	42	ED-010-112	50
ED-005-623-25	39	ED-005-714	42	ED-010-113	50
ED-005-624	39	ED-005-721	45	ED-010-114	50
ED-005-631	37	ED-005-721-25	45	ED-010-121	51
ED-005-631-25	37	ED-005-722	45	ED-010-122	51
ED-005-632	37	ED-005-723	45	ED-010-123	51
ED-005-633	37	ED-005-723-25	45	ED-010-124	51
ED-005-633-25	37	ED-005-724	45	ED-010-131	51
ED-005-634	37	ED-005-731	37	ED-010-132	51
ED-005-641	44	ED-005-731-25	37	ED-010-133	51
ED-005-641-25	44	ED-005-732	37	ED-010-134	51
ED-005-642	44	ED-005-733	37	ED-010-141	52
ED-005-643	44	ED-005-733-25	37	ED-010-142	52
ED-005-643-25	44	ED-005-734	37	ED-010-143	52
ED-005-644	44	ED-005-741	45	ED-010-144	52
ED-005-651	42	ED-005-741-25	45	ED-010-161	50
ED-005-651-25	42	ED-005-742	45	ED-010-162	50
ED-005-652	42	ED-005-743	45	ED-010-163	50
ED-005-653	42	ED-005-743-25	45	ED-010-164	50
ED-005-653-25	42	ED-005-744	45	ED-010-221	50
ED-005-654	42	ED-005-751	45	ED-010-222	50
ED-005-661	45	ED-005-751-25	45	ED-010-223	50
ED-005-661-25	45	ED-005-752	45	ED-010-224	50
ED-005-662	45	ED-005-753	45	ED-010-241	51

INDEX

ED-010-242	51
ED-010-243	51
ED-010-244	51
ED-010-251	51
ED-010-252	51
ED-010-253	51
ED-010-254	51
ED-010-261	50
ED-010-262	50
ED-010-263	50
ED-010-311	52
ED-010-312	52
ED-010-313	52
ED-010-314	52
ED-010-321	52
ED-010-322	52
ED-010-323	52
ED-010-324	52
ED-010-341	52
ED-010-342	52
ED-010-343	52
EG-003-062	48
EG-003-062-25	48
EG-003-062-45	48
EG-003-063	48
EG-003-063-25	48
EG-003-063-45	48
EG-003-064	48
EG-003-064-25	48
EG-003-064-45	48
EG-003-065	48
EG-003-065-25	48
EG-003-065-45	48
EG-003-066	48
EG-003-066-25	48
EG-003-066-45	48
EG-003-067	48
EG-003-067-25	48
EG-003-067-45	48
EG-003-068	48
EG-003-068-25	48
EG-003-068-45	48
EG-003-069	48
EG-003-069-25	48
EG-003-069-45	48

WE COVER THE FOLLOWING AREAS IN SURGERY WITH HIGH QUALITY PRODUCTS

ARTHROSCOPY

CARDIAC MIS

CONTAINER

ENT

GASTROSCOPY

GYNAECOLOGY

HF-SURGERY

LAMINECTOMY

LAPAROSCOPY

PROCTOSCOPY

SPINE SURGERY

TITANIUM PLATING

UNITS

UROLOGY



STEMA MEDIZINTECHNIK GMBH | TAKE-OFF GEWERBEPARK 74 | 78579 NEUHAUSEN OB ECK
INFO@STEMA-MEDIZINTECHNIK.DE | WWW.STEMA-MEDIZINTECHNIK.DE | +49 7467 79 630-0

11/25 | PRINTED IN GERMANY | COPYRIGHT BY STEMA MEDIZINTECHNIK GMBH |
WE RESERVE THE RIGHT TO MAKE ALTERATIONS