




UROLOGY



UROLOGY



ABOUT US



SINCE 20 YEARS



EXPORT IN
+70 COUNTRIES



24H/SERVICE



+ 10.000
PRODUCTS



LOCATED IN THE
WORLD CENTER
FOR MEDICAL
TECHNOLOGY

CERTIFICATED

STEMA Medizintechnik GmbH is certified in accordance with the highest standards:

- EN ISO 13485:2016 + AC:2018 + A11:2021- ISO 13485:2016
- EU 2017/745 (**MDR**)
- Council Directive 93/42/EEC Annex II excluding Section 4
- Council Directive 93/42/EEC Annex V excluding Section 3
- FDA

These standards create an essential base for our company and our partners. They also symbolize the commitment and strive to keep up the quality of our products.



PRODUCT RANGE

The current product range of STEMA Medizintechnik GmbH includes more than 10.000 different items.

We cover the following areas in surgery with high quality products:

ARTHROSCOPY

CARDIAC MIS

CONTAINER-SYSTEMS

GASTROENTEROLOGY

GYNAECOLOGY

HF-SURGERY

LAMINECTOMY

LAPAROSCOPY

MICRO MOTOR SYSTEMS

PROCTOSCOPY

SPINE-SURGERY

TITANIUM PLATING SYSTEM

UNITS

UROLOGY



CONTENT

UROLOGY

PRODUCT

PAGE

1. CYSTO-URETHROSCOPY

1.1 Cysto-Urethroscopy	8
1.2 Slim Cysto-Urethroscopy	25
1.3 Laser Cysto-Urethroscopy	28
1.4 Mini HD Compact Cystoscope	30
1.5 Flexible Cystoscope	31
1.6 Cystoscopy Set	33

2. RESECTOSCOPY

2.1 Monopolar Resectoscopy	37
2.2 Saline Resectoscopy	59
2.3 Laser Resectoscopy	76
2.4 Accessories	87
2.5 Resectoscopy Sets	91

3. URETHROTOMY

3.1 Urethrotome	97
3.2 Optical Urethrotomy	98
3.3 Urethrotomy Set	106

4. NEPHROSCOPY

4.1 HD Nephroscopy	109
4.2 Mini Nephroscopy	111
4.3 Standard Nephroscopy	114
4.4 Accessories	118

5. URETERO-RENOSCOPY

5.1 Uretero-Renoscopes semi-rigid	123
5.2 Flexible Uretero-Renoscope	142

CONTENT

UROLOGY

PRODUCT

PAGE

6. LITHOTRIPSY

6.1 Lithotripsy	146
-------------------	-----

7. MORCELLATION

7.1 Morcellation 26 Fr.	151
7.2 Morcellation 22 Fr.	153

8. PEDIATRIC UROLOGY

8.1 Pediatric Cysto-Urethroscopy	156
8.2 Pediatric Resectoscopy	163
8.3 Pediatric Urethrotomy	171

9. ACCESSORIES

9.1 Instruments	179
9.2 Sterilization	184
9.3 Cleaning and Maintenance	186

10. INDEX

10.1 Index	189
--------------	-----

1. CYSTO-URETHROSCOPY

1.1 | CYSTO-URETHROSCOPY

CYSTOSCOPES (4.0MM)

SHEATHS

OBTURATORS

ENDOSCOPE BRIDGES

WORKING ELEMENTS

FORCEPS, SCISSORS & INSTRUMENTS

STONE EXTRACTORS

MONOPOLAR FLEXIBLE ELECTRODES

MONOPOLAR CONNECTING CABLES

1.2 | SLIM CYSTO-URETHROSCOPY

CYSTOSCOPES (2.9MM)

SHEATHS

ENDOSCOPE BRIDGES

WORKING ELEMENTS

1.3 | LASER CYSTO-URETHROSCOPY

CYSTOSCOPES (4.0MM)

SHEATHS

OBTURATORS

1.4 | HD MINI COMPACT CYSTOSCOPES

1.5 | FLEXIBLE CYSTOSCOPE

CYSTOSCOPE FLEX I

ACCESSORIES


1.6 | CYSTOSCOPY SET

1. CYSTO-URETHROSCOPY


4.0mm

2.9mm

1 CYTOSCOPES



2 SHEATHS & OBTURATORS



3 ACCESSORIES



1.1 | CYSTO-URETHROSCOPY CYSTOSCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-001
	Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-001-W
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-011
	HD Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-021

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-002
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-012
	Ultra HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-022

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-003-W
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013
	HD Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-013-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-023

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 70°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-004
	Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-004-W
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-014
	HD Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-014-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-024

4.0mm

1.1 | CYSTO-URETHROSCOPY SHEATHS



Ordering no.	Description	French	Colour code
AB-001-002	Cysto-Urethroscope sheath with 2 fixed stopcock and standard obturator	17 Fr.	Yellow
AB-001-003		19 Fr.	Green
AB-001-004		21 Fr.	Red
AB-001-005		23 Fr.	Blue
AB-001-006		25 Fr.	White



Ordering no.	Description	French	Colour code
AB-002-001	Cysto-Urethroscope sheath with central valve and standard obturator	17 Fr.	Yellow
AB-002-002		19 Fr.	Green
AB-002-003		21 Fr.	Red
AB-002-004		23 Fr.	Blue
AB-002-005		25 Fr.	White

1.1 | CYSTO-URETHROSCOPY

OBTURATORS | ADAPTER



Ordering no.	Description	French	Colour code
AB-001-012	Standard obturator for Cysto-Urethroscope sheaths	17 Fr.	Yellow
AB-001-013		19 Fr.	Green
AB-001-014		21 Fr.	Red
AB-001-015		23 Fr.	Blue
AB-001-016		25 Fr.	White



Ordering no.	Description	French	Colour code
AB-001-022	Viewing obturator for Cysto-Urethroscope sheaths	17 Fr.	Yellow
AB-001-023		19 Fr.	Green
AB-001-024		21 Fr.	Red
AB-001-025		23 Fr.	Blue
AB-001-026		25 Fr.	White

Ordering no.	Description
AB-001-010	Cystoscope adapter for female urethroscopy




1.1 | CYSTO-URETHROSCOPY ENDOSCOPE BRIDGES



Ordering no.	Description
AB-003-002	Endoscope bridge diagnostic



Ordering no.	Description
AB-003-003	Endoscope bridge with 1 instrument channel



Ordering no.	Description
AB-003-004	Endoscope bridge with 2 instrument channels



Ordering no.	Description
AB-003-015	Endoscope bridge with 2 instrument channels, short



Ordering no.	Description
AJ-010-001	Rubber cap for instrument channel

1.1 | CYSTO-URETHROSCOPY WORKING ELEMENTS



Ordering no.	Description
AB-003-006	Working element with ALBARRAN deflector and 1 instrument channel
AB-003-008	Working element with ALBARRAN deflector and 1 instrument channel, without ratchet



Ordering no.	Description
AB-003-007	Working element with ALBARRAN deflector and 2 instrument channels
AB-003-009	Working element with ALBARRAN deflector and 2 instrument channels, without ratchet

1.1 | CYSTO-URETHROSCOPY

OPTICAL FORCEPS AND SCISSORS



Info: Only for sheaths 19-25 Fr.



	Ordering no.	For scopes	Description
	AB-004-005	12°/30°	Biopsy forceps (down), double action
	AB-004-006	70°	Biopsy forceps (down), double action
	AB-004-001	0°/12°/30°	Biopsy forceps, double action
	AB-004-002	0°/12°/30°	Biopsy forceps, through cutting
	AB-004-004	0°/12°	Hook scissors, double action
	AB-004-010	0°/12°	Forceps, double action, for removal of stents
	AB-004-003	0°/12°	Grasping forceps, double action
	AB-004-008		Adapter for use with resectoscope sheaths

1.1 | CYSTO-URETHROSCOPY

OPTICAL FORCEPS




Info: Only for sheats 25 Fr.



	Ordering no.	For scopes	Description
	AB-004-007	0°/12°/30°	Optical stone punch
	AB-004-009		Adapter for use with resectoscope sheaths 26 - 28.5 Fr.
	AB-004-009-QL		Adapter for use with resectoscope sheaths 26 - 28.5 Fr. - STEMA Click-System

1.1 | CYSTO-URETHROSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps, single action	Ø 2.3mm	7 Fr.	AC-010-004	AC-010-204	
	Grasping forceps, single action	Ø 2.3mm	7 Fr.	AC-011-004	AC-011-204	
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-001	AC-001-201	AC-001-001-R
		Ø 1.2mm	4 Fr.	AC-001-008	AC-001-208	AC-001-008-R
		Ø 1.6mm	5 Fr.	AC-001-002	AC-001-202	AC-001-002-R
		Ø 2.1mm	6 Fr.	AC-001-003	AC-001-203	AC-001-003-R
		Ø 2.3mm	7 Fr.	AC-001-004	AC-001-204	AC-001-004-R
	Biopsy forceps oval, serrated, double action	Ø 3.0mm	9 Fr.	AC-001-005	AC-001-205	AC-001-005-R
		Ø 1.0mm	3 Fr.	AC-002-001	AC-002-201	AC-002-001-R
		Ø 1.6mm	5 Fr.	AC-002-002	AC-002-202	AC-002-002-R
		Ø 2.1mm	6 Fr.	AC-002-003	AC-002-203	AC-002-003-R
		Ø 2.3mm	7 Fr.	AC-002-004	AC-002-204	AC-002-004-R
	Biopsy forceps oval, serrated, teeth interlock, double action	Ø 3.0mm	9 Fr.	AC-002-005	AC-002-205	AC-002-005-R
		Ø 1.6mm	5 Fr.	AC-015-002	AC-015-202	AC-015-002-R
		Ø 2.1mm	6 Fr.		AC-015-203	AC-015-003-R
	Grasping forceps for stone street, double action	Ø 2.3mm	7 Fr.	AC-015-004	AC-015-204	AC-015-004-R
		Ø 1.2mm	4 Fr.	AC-017-001	AC-017-201	AC-017-001-R
	Grasping forceps, alligator jaw, double action	Ø 2.3mm	7 Fr.	AC-017-004	AC-017-204	AC-017-004-R
		Ø 1.0mm	3 Fr.	AC-006-001	AC-006-201	AC-006-001-R
		Ø 1.2mm	4 Fr.	AC-006-003	AC-006-203	AC-006-003-R
		Ø 1.6mm	5 Fr.	AC-006-002	AC-006-202	AC-006-002-R
		Ø 2.3mm	7 Fr.	AC-006-004	AC-006-204	AC-006-004-R
		Ø 3.0mm	9 Fr.	AC-006-005	AC-006-205	AC-006-005-R

1.1 | CYSTO-URETHROSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	Ø 2.3mm	7 Fr.	AC-008-004	AC-008-204	AC-008-004-R
		Ø 3.0mm	9 Fr.	AC-008-005	AC-008-205	AC-008-005-R
	Grasping forceps teeth double action	Ø 1.2mm	4 Fr.	AC-016-001	AC-016-201	AC-016-001-R
		Ø 1.6mm	5 Fr.	AC-016-002	AC-016-202	AC-016-002-R
	Grasping forceps teeth interlock, double action	Ø 1.6mm	5 Fr.	AC-020-002	AC-020-202	AC-020-002-R
		Ø 2.3mm	7 Fr.	AC-020-004	AC-020-204	AC-020-004-R
	Grasping forceps rat tooth, double action	Ø 1.0mm	3 Fr.	AC-019-001	AC-019-201	AC-019-001-R
		Ø 1.6mm	5 Fr.	AC-019-002	AC-019-202	AC-019-002-R
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-001	AC-007-201	AC-007-001-R
		Ø 1.6mm	5 Fr.	AC-007-002	AC-007-202	AC-007-002-R
	Grasping forceps for big stones, double action	Ø 1.2mm	4 Fr.	AC-009-001	AC-009-201	AC-009-001-R
		Ø 1.6mm	5 Fr.	AC-009-002	AC-009-202	AC-009-002-R
		Ø 2.3mm	7 Fr.	AC-009-004	AC-009-204	AC-009-004-R
	3-Prongs retriever	Ø 1.0mm	3 Fr.	AC-012-001		
		Ø 1.6mm	5 Fr.	AC-012-002		
		Ø 2.3mm	7 Fr.	AC-012-004		


1.1 | CYSTO-URETHROSCOPY

SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-004	AC-004-204	AC-004-004-R
		Ø 3.0mm	9 Fr.	AC-004-005	AC-004-205	AC-004-005-R
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-000	AC-003-200	AC-003-000-R
		Ø 1.2mm	4 Fr.	AC-003-001	AC-003-201	AC-003-001-R
		Ø 1.6mm	5 Fr.	AC-003-002	AC-003-202	AC-003-002-R
		Ø 2.1mm	6 Fr.	AC-003-003	AC-003-203	AC-003-003-R
		Ø 2.3mm	7 Fr.	AC-003-004	AC-003-204	AC-003-004-R
		Ø 3.0mm	9 Fr.	AC-003-005	AC-003-205	AC-003-005-R
	Uretero slitting scissors, retrograde cutting single action	Ø 1.6mm	4 Fr.	AC-018-001	AC-018-201	AC-018-001-R
		Ø 2.3mm	5 Fr.	AC-018-002	AC-018-202	AC-018-002-R
	Scissors pointed, single action	Ø 1.6mm	5 Fr.		AC-003-302	
		Ø 2.3mm	7 Fr.		AC-003-304	
	Scissors, serrated flexible, double action	Ø 2.3mm	7 Fr.	AC-005-004	AC-005-204	AC-005-004-R

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL



FEATURES

- 1 QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
- 2 ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
- 3 OVERLOAD PROTECTION AVAILABLE
- 4 SIDE-MOUNTED LUER LOCK

1.1 | CYSTO-URETHROSCOPY

INSTRUMENTS DETACHABLE | SEMI-RIGID



Ordering no.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-200-T	AC-003-200-OT
		Ø 1.2mm	4 Fr.	AC-003-201-T	AC-003-201-OT
		Ø 1.6mm	5 Fr.	AC-003-202-T	AC-003-202-OT
		Ø 2.3mm	7 Fr.	AC-003-204-T	AC-003-204-OT
		Ø 3.0mm	9 Fr.	AC-003-205-T	AC-003-205-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-003-612-T	AC-003-612-OT
		Ø 2.3mm	7 Fr.	AC-003-614-T	AC-003-614-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-204-T	AC-004-204-OT
		Ø 3.0mm	9 Fr.	AC-004-205-T	AC-004-205-OT
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-201-T	AC-001-201-OT
		Ø 1.2mm	4 Fr.	AC-001-206-T	AC-001-206-OT
		Ø 1.6mm	5 Fr.	AC-001-202-T	AC-001-202-OT
		Ø 2.3mm	7 Fr.	AC-001-203-T	AC-001-203-OT
		Ø 3.0mm	9 Fr.	AC-001-204-T	AC-001-204-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-006-201-T	AC-006-201-OT
		Ø 1.2mm	4 Fr.	AC-006-203-T	AC-006-203-OT
		Ø 1.6mm	5 Fr.	AC-006-202-T	AC-006-202-OT
		Ø 2.3mm	7 Fr.	AC-006-204-T	AC-006-204-OT
		Ø 3.0mm	9 Fr.	AC-006-205-T	AC-006-205-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-201-T	AC-007-201-OT
		Ø 1.6mm	5 Fr.	AC-007-202-T	AC-007-202-OT
	Grasping forceps for big stones, double action	Ø 1.2mm	5 Fr.	AC-009-201-T	AC-009-201-OT
		Ø 1.6mm	5 Fr.	AC-009-202-T	AC-009-202-OT

1.1 | CYSTO-URETHROSCOPY

INSTRUMENTS DETACHABLE | RIGID



Ordering no.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-103-200-T	AC-103-200-OT
		Ø 1.2mm	4 Fr.	AC-103-201-T	AC-103-201-OT
		Ø 1.6mm	5 Fr.	AC-103-202-T	AC-103-202-OT
		Ø 2.3mm	7 Fr.	AC-103-204-T	AC-103-204-OT
		Ø 3.0mm	9 Fr.	AC-103-205-T	AC-103-205-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-103-612-T	AC-103-612-OT
		Ø 2.3mm	7 Fr.	AC-103-614-T	AC-103-614-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-104-204-T	AC-104-204-OT
		Ø 3.0mm	9 Fr.	AC-104-205-T	AC-104-205-OT
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-101-201-T	AC-101-201-OT
		Ø 1.2mm	4 Fr.	AC-101-206-T	AC-101-206-OT
		Ø 1.6mm	5 Fr.	AC-101-202-T	AC-101-202-OT
		Ø 2.3mm	7 Fr.	AC-101-203-T	AC-101-203-OT
		Ø 3.0mm	9 Fr.	AC-101-204-T	AC-101-204-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-106-201-T	AC-106-201-OT
		Ø 1.2mm	4 Fr.	AC-106-203-T	AC-106-203-OT
		Ø 1.6mm	5 Fr.	AC-106-202-T	AC-106-202-OT
		Ø 2.3mm	7 Fr.	AC-106-204-T	AC-106-204-OT
		Ø 3.0mm	9 Fr.	AC-106-205-T	AC-106-205-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-107-201-T	AC-107-201-OT
		Ø 1.6mm	5 Fr.	AC-107-202-T	AC-107-202-OT
	Grasping forceps for big stones, double action	Ø 1.2mm	5 Fr.	AC-109-201-T	AC-109-201-OT

1.1 | CYSTO-URETHROSCOPY

STONE EXTRACTORS | 4 WIRES



Ordering no.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering no.	Wires	French	Working length	Description
AF-001-001	4	3 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-002	4	3 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-015	4	3 Fr.	900mm	Dormia stone extractor, basket with tubus only
AF-001-016	4	3 Fr.	900mm	Dormia stone extractor complete with handle
AF-001-023	4	3 Fr.	1200mm	Dormia stone extractor, basket with tubus only
AF-001-022	4	3 Fr.	1200mm	Dormia stone extractor complete with handle
AF-001-003	4	5 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-004	4	5 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-017	4	5 Fr.	900mm	Dormia stone extractor, basket with tubus only
AF-001-019	4	5 Fr.	900mm	Dormia stone extractor complete with handle
AF-001-005	4	7 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-006	4	7 Fr.	600mm	Dormia stone extractor complete with handle

1.1 | CYSTO-URETHROSCOPY

STONE EXTRACTORS | 5-7 WIRES



Ordering no.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering no.	Wires	French	Working length	Description
AF-001-007	5	5 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-008	5	5 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-025	5	5 Fr.	900mm	Dormia stone extractor, basket with tubus only
AF-001-018	5	5 Fr.	900mm	Dormia stone extractor complete with handle
AF-001-009	5	7 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-010	5	7 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-011	6	5 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-012	6	5 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-013	6	7 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-014	7	7 Fr.	600mm	Dormia stone extractor complete with handle

1.1 | CYSTO-URETHROSCOPY

MONOPOLAR FLEXIBLE ELECTRODES




	Description	French	Length	Ordering no.	Length	Ordering no.
	Loop electrode flexible, monopolar	5 Fr.	530mm	AD-003-003	650mm	AD-003-013
		7 Fr.	530mm	AD-003-005	650mm	AD-003-015
		9 Fr.	530mm	AD-003-007	650mm	AD-003-017
		10 Fr.	530mm	AD-003-008	650mm	AD-003-018
	Needle electrode flexible, monopolar	3 Fr.	530mm	AD-002-001	650mm	AD-002-011
		4 Fr.	530mm	AD-002-002	650mm	AD-002-012
		5 Fr.	530mm	AD-002-003	650mm	AD-002-013
		7 Fr.	530mm	AD-002-005	650mm	AD-002-015
	Button electrode flexible, monopolar	3 Fr.	530mm	AD-001-001	650mm	AD-001-011
		4 Fr.	530mm	AD-001-002	650mm	AD-001-012
		5 Fr.	530mm	AD-001-003	650mm	AD-001-013
		6 Fr.	530mm	AD-001-004	650mm	AD-001-014
		7 Fr.	530mm	AD-001-005	650mm	AD-001-015
		8 Fr.	530mm	AD-001-006	650mm	AD-001-016
		9 Fr.	530mm	AD-001-007	650mm	AD-001-017
		10 Fr.	530mm	AD-001-008	650mm	AD-001-018

1.1 | CYSTO-URETHROSCOPY

MONOPOLAR CONNECTING CABLES



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4.0mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4.0mm	
	AD-010-002	3m	Ø 5.0mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5.0mm	
	AD-010-003	3m	Ø 8.0mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8.0mm	
	AD-010-024	4.5m	Ø 2.0mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

1.2 | SLIM CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012
	HD Cystoscope autoclavable wide angle	Ø 2.9mm	302mm	AA-003-012-W



Ordering no.	Description
AB-001-007	Cysto-Urethroscope sheath 15 Fr. with 2 fixed stopcocks, 5 Fr. instrument channel , incl. obturator




Ordering no.	Description
AB-001-000	Cysto-Urethroscope sheath 16 Fr. with 2 fixed stopcocks, Quick-Lock , incl. obturator


1.2 | SLIM CYSTO-URETHROSCOPY

ENDOSCOPE BRIDGES | QUICK-LOCK SYSTEM

Quick-Lock System



Ordering no.	Description
AB-003-010	Endoscope bridge diagnostic used for 16 Fr. sheath



Ordering no.	Description
AB-003-011	Endoscope bridge with 1 instrument channel used for 16 Fr. sheath



Ordering no.	Description
AB-003-012	Endoscope bridge with 2 instrument channels used for 16 Fr. sheath



Ordering no.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0,6-0,8mm

1.2 | SLIM CYSTO-URETHROSCOPY WORKING ELEMENTS



Ordering no.	Description
AB-003-013	Working element with ALBARRAN deflector and 1 instrument channel used for 16 Fr. sheath
AB-003-014	Working element with ALBARRAN deflector and 2 instrument channels used for 16 Fr. sheath



Ordering no.	Description
AB-003-020	Protective rod for working elements with ALBARRAN deflector

1.3 | LASER CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013
	ULTRA HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-023



Ordering no.	Description
AB-001-040	LASER Cysto-Urethroscope outer sheath 23 Fr. with 2 fixed stopcock used for 30° endoscope



Ordering no.	Description
AB-001-041	LASER Cysto-Urethroscope inner sheath with 8 Fr. working channel

1.3 | LASER CYSTO-URETHROSCOPY OBTURATORS



Ordering no.	Description
AB-001-043	Obturator for LASER cysto-urethroscope sheath AB-001-040



Ordering no.	Description
AB-001-042	Viewing obturator for LASER cysto-urethroscope sheath AB-001-040

1.4 | HD MINI COMPACT CYSTOSCOPE



Ordering no.	Product
AA-007-060	HD Mini Compact Cystoscope

Technical data

Eyepiece	30° angled ocular
Sheath	7.5/9/12 Fr.
Direction of view	0° viewing angle
Working length	220mm
Instrument channel	5 Fr.
Luer Lock	3 connections
Optic	40.000 px

Delivery contains

HD Mini compact cystoscope

1.5 | FLEXIBLE CYSTOSCOPE FLEX I

- *High luminance*
- *Exact adjustment of focus*
- *Precise up-/down angulation*



Ordering no.	Product
AA-008-002	Flexible Cystoscope FLEX I

Technical data

Cystoscope

Insertion tube	Ø 5.5mm
Working Channel	Ø 2.3mm
Working length	380mm
Field of view	100°
Pixel	18.000
Angulation up-down	210°- 120°
Accessories	Suction valve, leak tester

Delivery contains

Flexible cystoscope FLEX I, leak tester, suction valve

INCLUDED IN SCOPE OF DELIVERY:

AA-008-103

Leak tester

AA-008-105

Suction valve



AA-008-103



AA-008-105

1.5 | FLEXIBLE CYSTOSCOPE ACCESSORIES



Ordering no.	Description
AA-008-100	Bridge with Luer-Lock and 2 stop cocks for AA-008-002



Ordering no.	Description
AA-008-102	Suitcase grey for STEMA flexible scopes



Ordering no.	Description
AA-008-103	Leak tester



Ordering no.	Description
AA-008-104	Single use Luer-Lock adapter, 25 pcs./box



Ordering no.	Description
AA-008-105	Suction valve

1.6 CYSTO-URETHROSCOPY SET

1 HD HYSTEROSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable, 302mm working length

1 AA-002-013



2 CYSTO-URETHROSCOPE SHEATH

21 Fr. with 2 fixed stopcock and obturator

2 AB-001-004

3 CYSTO-URETHROSCOPE SHEATH

25 Fr. with 2 fixed stopcock and obturator

3 AB-001-006



4 ENDOSCOPE BRIDGE

with 1 instrument channel

4 AB-003-003



5 WORKING ELEMENT

with Albaran deflector and 2 instrument channels

5 AB-003-007



6 BIOPSY FORCEPS

flexible, oval, double action, Ø 1.6mm, WL 40cm

6 AC-001-002



7 SCISSORS

flexible, single action, Ø 1.6mm, WL 40cm

7 AC-003-002



8 GRASPING FORCEPS ALLIGATOR

flexible, double action, Ø 2.3mm, WL 40cm

8 AC-006-004



1.6 CYSTO-URETHROSCOPY SET

9 BUTTON ELECTRODE

flexible, 5 Fr., monopolar, length 530mm



9 AD-001-003

10 NEEDLE ELECTRODE

flexible, 5 Fr., monopolar, length 530mm



10 AD-002-003



11 HF MONOPOLAR CABLE

HF monopolar connecting cable for Berchtold, Martin, Aesculap, Ø 4mm

11 AD-010-003



12 DORMIA STONE EXTRACTOR

5 wires, Ø 1.6mm, complete with handle, 60cm Working length

12 AF-001-008



13 BIOPSY FORCEPS

double action

13 AB-004-001



14 GRASPING FORCEPS

double action

14 AB-004-003



15 OPTICAL STONE PUNCH

15 AB-004-007



2. RESECTOSCOPY

2.1 | MONOPOLAR RESECTOSCOPY

24/26/28.5 FRENCH
19/22 FRENCH (SLIM)

2.2 | SALINE RESECTOSCOPY

24/26 FRENCH
19/22 FRENCH (SLIM)
17.5 FRENCH (MINI)

2.3 | LASER RESECTOSCOPY

24/26 FRENCH
19/22 FRENCH (SLIM)
17.5 FRENCH (MINI)

2.4 | ACCESSORIES

2.5 | RESECTOSCOPY SETS

2. RESECTOSCOPY

24/26 Fr.
28.5 Fr.

19/22 Fr.


17.5 Fr.




1 CYSTOSCOPES




2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS



4 ELECTRODES



2.1 | MONOPOLAR RESECTOSCOPY CYSTOSCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-001
	Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-001-W
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-011
	HD Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-021


Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-002
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-012
	Ultra HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-022

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-003-W
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013
	HD Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-013-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-023

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 70°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-004
	Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-004-W
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-014
	HD Cystoscope autoclavable wide angle	Ø 4.0mm	302mm	AA-002-014-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-024

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | PASSIVE CUTTING



Ordering no.	Description
AG-001-001	Working element used for resectoscopes, passive cutting, handle closed , cutting spring assisted
AG-001-001-HF	Working element used for resectoscopes, passive cutting, handle closed , cutting spring assisted, HF connection on top



Ordering no.	Description
AG-001-002	Working element used for resectoscopes, passive cutting, handle open , cutting spring assisted

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | ACTIVE CUTTING



Ordering no.	Description
AG-001-003	Working element used for resectoscopes, active cutting, handle closed , cutting again spring assistance
AG-001-003-HF	Working element used for resectoscopes, active cutting, handle closed , cutting again spring assistance, HF connection on top



Ordering no.	Description
AG-001-004	Working element used for resectoscopes, active cutting, handle open , cutting spring assisted

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | TITAN



Ordering no.	Description
AG-001-007	Working element used for resectoscopes, passive cutting , handle closed, cutting spring assisted



Ordering no.	Description
AG-001-008	Working element used for resectoscopes, active cutting , handle closed, cutting again spring assistance

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | GYRUS



Ordering no.	Description
AG-001-017	GYRUS Saline Working element used for resectoscopes, active cutting, handle closed TITANIUM VERSION



Ordering no.	Description
AG-001-018	GYRUS Saline Working element used for resectoscopes, passive cutting, handle closed TITANIUM VERSION

2.1 | MONOPOLAR RESECTOSCOPY

CONTINUOUS FLOW SHEATHS FIXED



	Ordering no.	Description	French	Colour code
	AH-001-202	Outer sheath with perforation	26 Fr.	Yellow
	AH-001-212	Outer sheath with perforation	28.5 Fr.	Black
	AH-001-204	Outer sheath with upper/lower perforation	26 Fr.	Yellow
	AH-001-203	Outer sheath with surrounding fenestrated perforation	26 Fr.	Yellow
	AH-001-205	Outer sheath with upper fenestrated perforation	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-201	Inner sheath with ceramic isolation	26 Fr.	Yellow
AH-001-211	Inner sheath with ceramic isolation	28.5 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY

CONTINUOUS FLOW SHEATHS ROTATABLE



	Ordering no.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow
	AH-001-312	Outer sheath with perforation, rotatable	28.5 Fr.	Black
	AH-001-304	Outer sheath with upper/lower perforation, rotatable	26 Fr.	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow
AH-001-311	Inner sheath with ceramic isolation, rotatable	28.5 Fr.	Black



Ordering no.	Description	French	Colour code
AH-001-308	Inner sheath without ceramic isolation, with laser channel up to 7 Fr., for use with AH-001-301	26 Fr.	Yellow

2.1 | MONOPOLAR RESECTOSCOPY

STANDARD SHEATHS WITH INTERMITTENT IRRIGATION



Ordering no.	Description	French	Colour code
AH-001-001	Standard resectoscope sheath with 1 fixed stopcock	24 Fr.	Yellow
AH-001-002	Standard resectoscope sheath with 1 fixed stopcock	27 Fr.	Black



Ordering no.	Description	French	Colour code
AH-001-101	Standard resectoscope sheath with central valve	24 Fr.	Yellow
AH-001-102	Standard resectoscope sheath with central valve	27 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY OBTURATORS



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-012	Standard obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-015	Deflecting obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-016	Deflecting obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-013	Viewing obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-017	Viewing obturator for resectoscope sheaths with 5 Fr. Instrument channel	24 Fr.	26 Fr.	Yellow
AH-001-018	Viewing obturator for resectoscope sheaths with 5 Fr. Instrument channel	27 Fr.	28.5 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY

ELECTRODES FOR MONOPOLAR RESECTOSCOPY




	Ordering no.	French	Description
	AE-001-001	24 Fr. (yellow)	Cutting loop electrode 30° angled
	AE-002-001	27 Fr. (brown)	
	AE-001-002	24 Fr. (yellow)	Cutting loop electrode angled
	AE-002-002	27 Fr. (brown)	
	AE-001-003	24 Fr. (yellow)	Cutting loop electrode straight
	AE-002-003	27 Fr. (brown)	
	AE-001-004	24 Fr. (yellow)	Rollerball electrode Ø 3mm
	AE-002-004	27 Fr. (brown)	
	AE-001-005	24 Fr. (yellow)	Rollerball electrode Ø 5mm
	AE-002-005	27 Fr. (brown)	
	AE-001-010	24 Fr. (yellow)	Roller electrode Ø 3mm
	AE-002-010	27 Fr. (brown)	
	AE-001-011	24 Fr. (yellow)	Roller electrode Ø 5mm
	AE-002-011	27 Fr. (brown)	
	AE-001-012	24 Fr. (yellow)	Spike electrode Ø 3mm
	AE-002-012	27 Fr. (brown)	
	AE-001-013	24 Fr. (yellow)	Spike electrode Ø 5mm
	AE-002-013	27 Fr. (brown)	

2.1 | MONOPOLAR RESECTOSCOPY

ELECTRODES FOR MONOPOLAR RESECTOSCOPY



	Ordering no.	French	Description
	AE-001-014	24 Fr. (yellow)	Vaporisation electrode Ø 3mm
	AE-002-014	27 Fr. (brown)	
	AE-001-015	24 Fr. (yellow)	Vaporisation electrode Ø 5mm
	AE-002-015	27 Fr. (brown)	
	AE-001-009	24 Fr. (yellow)	Conical electrode
	AE-001-016	24 Fr. (yellow)	Knife electrode
	AE-002-016	27 Fr. (brown)	
	AE-001-006	24 Fr. (yellow)	Band electrode angled
	AE-002-006	27 Fr. (brown)	
	AE-001-007	24 Fr. (yellow)	Band electrode 30° angled
	AE-002-007	27 Fr. (brown)	
	AE-001-008	24 Fr. (yellow)	Band vaporization electrode angled
	AE-001-017	24 Fr. (yellow)	Curette sharp, round (not for use with HF)



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.1 | MONOPOLAR SLIM RESECTOSCOPY

CYSTOSCOPES AUTOCLAVABLE Ø 2.9MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS | PASSIVE CUTTING



Ordering no.	Description
AG-002-001	Working element used for resectoscopes, passive cutting, handle closed , cutting against spring assistance



Ordering no.	Description
AG-002-002	Working element used for resectoscopes, passive cutting, handle open , cutting against spring assistance

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS | ACTIVE CUTTING




Ordering no.	Description
AG-002-003	Working element used for slim resectoscopes, active cutting, handle closed , cutting against spring assistance



Ordering no.	Description
AG-002-004	Working element used for slim resectoscopes, active cutting, handle open , cutting against spring assistance

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS | TITAN



Ordering no.	Description
AG-002-006	Working element used for slim resectoscopes, active cutting , handle closed, TITANIUM VERSION

2.1 | MONOPOLAR SLIM RESECTOSCOPY

CONTINUOUS FLOW SHEATHS | FIXED



	Ordering no.	Description	French	Colour code
	AH-002-212	Outer sheath with perforation, fixed	22 Fr.	White



Ordering no.	Description	French	Colour code
AH-002-201	Inner sheath with ceramic isolation, fixed	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY

CONTINOUS FLOW SHEATHS | ROTATABLE



	Ordering no.	Description	French	Colour code
	AH-002-302	Outer sheath with perforation, rotatable	22 Fr.	White



Ordering no.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, rotatable	22 Fr.	White



19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY STANDARD SHEATHS



	Ordering no.	Description	French	Colour code
	AH-002-001	Standard resectoscope sheath, 1 fixed stopcock , without obturator	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY

ELECTRODES FOR MONOPOLAR RESECTOSCOPY



	Ordering no.	French	Description
	AE-003-001	19 Fr.	Cutting loop electrode 30° angled
	AE-003-003	19 Fr.	Cutting loop electrode straight
	AE-003-004	19 Fr.	Rollerball electrode Ø 3mm
	AE-003-010	19 Fr.	Roller electrode Ø 3mm
	AE-003-014	19 Fr.	Vaporization electrode Ø 3mm
	AE-003-016	19 Fr.	Knife electrode



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.1 | MONOPOLAR SLIM RESECTOSCOPY

OBTURATORS



Ordering no.	Description	French	Colour code
AH-002-011	Standard obturator	19 Fr.	White




Ordering no.	Description	French	Colour code
AH-002-013	Viewing obturator	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY

MONOPOLAR CONNECTING CABLES



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4.0mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4.0mm	
	AD-010-002	3m	Ø 5.0mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5.0mm	
	AD-010-003	3m	Ø 8.0mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8.0mm	
	AD-010-024	4.5m	Ø 2.0mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

2.1 | MONOPOLAR SLIM RESECTOSCOPY

CONNECTING CABLES MONOPOLAR



CONNECTION UNIT SIDE



Instrument Side	Ordering no.	Length	Description
	AD-010-025	4.5m	HF Monopolar connecting cable for Storz resectoscopes with 90° angled connector
	AD-010-026	4.5m	HF Monopolar connecting cable for Olympus resectoscopes, diameter Ø 3mm



RECOMMENDED PRODUCTS:

CI-001-030 STEMA HF-U PRO TOUCH

- For all monopolar & bipolar techniques
- with 7" touch screen display
- provided to perform the most advanced electrosurgery
- ADC / APC System

CI-001-040 DOUBLE PEDAL FOOTSWITCH

2.2 | SALINE RESECTOSCOPY CYSTOSCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-012

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013
	ULTRA HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-023

2.2 | SALINE RESECTOSCOPY WORKING ELEMENTS



Ordering no.	Description
AG-001-011	SALINE working element used for resectoscopes, 24 Fr., passive cutting , handle closed, cutting spring assistance




Ordering no.	Description
AG-001-012	SALINE working element used for resectoscopes, 24 Fr. , active cutting , handle closed, cutting spring assistance

2.2 | SALINE RESECTOSCOPY

CONTINUOUS FLOW SHEATHS | ROTATABLE



	Ordering no.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow
	AH-001-312	Outer sheath with perforation, rotatable	28.5 Fr.	Black
	AH-001-304	Outer sheath with upper/lower perforation, rotatable	26 Fr.	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow
AH-001-311	Inner sheath with ceramic isolation, rotatable	28.5 Fr.	Black

2.2 | SALINE RESECTOSCOPY OBTURATORS



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-012	Standard obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-015	Deflecting obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-016	Deflecting obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-013	Viewing obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-017	Viewing obturator for resectoscope sheaths with 5 Fr. instrument channel	24 Fr.	26 Fr.	Yellow
AH-001-018	Viewing obturator for resectoscope sheaths with 5 Fr. instrument channel	27 Fr.	28.5 Fr.	Black

2.2 | SALINE RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY



	Ordering no.	French	Description
	AE-001-030	24 Fr.	Cutting loop electrode 30° angled, strong cut
	AE-001-031	24 Fr.	Cutting loop electrode 30° angled
	AE-001-033	24 Fr.	Cutting loop electrode 90 ° angled
	AE-001-048	24 Fr.	Cutting loop electrode 90 ° angled, for 0° Scope
	AE-001-036	24 Fr.	Cutting loop electrode straight
	AE-001-032	24 Fr.	Knife electrode
	AE-001-041	24 Fr.	Band electrode, 30° angled
	AE-001-038	24 Fr.	Vaporisation spike rollerball electrode Ø 3mm
	AE-001-037	24 Fr.	Vaporisation rollerball electrode Ø 3mm
	AE-001-043	24 Fr.	Vaporisation rollerball electrode Ø 5mm
	AE-001-040	24 Fr.	Roller electrode smooth Ø 3mm
	AE-001-042	24 Fr.	Roller electrode smooth Ø 5mm

2.2 | SALINE RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY



	Ordering no.	French	Description
	AE-001-034	24 Fr.	Rollerball electrode Ø 3mm
	AE-001-035	24 Fr.	Rollerball electrode Ø 5mm
	AE-001-044	24 Fr.	Conical electrode „Mushroom“
	AE-001-045	24 Fr.	Conical electrode „Mushroom“, convex
	AE-001-046	24 Fr.	Conical electrode „Mushroom“, concav
	AE-001-047	24 Fr.	Conical electrode „Mushroom“, grooved
	AE-001-039	24 Fr.	Conical electrode

	Ordering no.	French	Description
	AE-002-031	27 Fr.	Cutting loop electrode, 30° angled



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of eletcrodes and cold knives

2.2 | SALINE SLIM RESECTOSCOPY

CYSTOSCOPES AUTOCLAVABLE Ø 2.9MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-004
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-014

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.2 | SALINE SLIM RESECTOSCOPY

WORKING ELEMENTS



Ordering no.	Description
AG-002-011	SALINE working element used for resectoscopes, 19 Fr., passive cutting , handle closed, cutting spring assistance



Ordering no.	Description
AG-002-012	SALINE working element used for resectoscopes, 19 Fr., active cutting , handle closed, cutting spring assistance

2.2 | SALINE SLIM RESECTOSCOPY SHEATHS



Ordering no.	Description
AH-002-302	Outer sheath 22 Fr. with perforation, rotatable



Ordering no.	Description
AH-002-301	Inner sheath 22 Fr. with ceramic isolation, rotatable

2.2 | SALINE SLIM RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-011	Standard obturator for resectoscope sheaths 19 Fr.



Ordering no.	Description
AH-002-013	Viewing obturator for resectoscope sheaths 19 Fr.

2.2 | SALINE SLIM RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY



	Ordering no.	French	Description
	AE-003-031	19 Fr.	Cutting loop electrode 30° angled
	AE-003-032	19 Fr.	Cutting loop electrode 90° angled
	AE-003-033	19 Fr.	Cutting loop electrode straight
	AE-003-034	19 Fr.	Rollerball electrode Ø 3mm
	AE-003-040	19 Fr.	Roller electrode Ø 3mm
	AE-003-036	19 Fr.	Knife electrode



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives




17.5 Fr.

2.2 | SALINE MINI RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
30°	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.2 | SALINE MINI RESECTOSCOPY WORKING ELEMENTES



Ordering no.	Description
AG-002-015	SALINE working element used for resectoscopes, 18.5 Fr., passive cutting , handle closed, cutting spring assistance



Ordering no.	Description
AG-002-013	SALINE working element used for resectoscopes, 18.5 Fr., active cutting , handle closed, cutting spring assistance

2.2 | SALINE MINI RESECTOSCOPY
SHEATHS



Ordering no.	Description
AH-002-213	Outer sheath 18.5 Fr. with perforation



Ordering no.	Description
AH-002-211	Inner sheath 17.5 Fr. with ceramic isolation, rotatable



Ordering no.	Description
AH-002-216	Inner sheath 17.5 Fr. without ceramic isolation , with laser channel up to 5 Fr., for use with AH-001-213

2.2 | SALINE MINI RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-015	Standard obturator 17.5 Fr.



Ordering no.	Description
AH-002-016	Viewing obturator 17.5 Fr. with 5 fr. instrument channel

2.2 | SALINE MINI RESECTOSCOPY

ELECTRODES FOR BIPOLAR RESECTOSCOPY




	Ordering no.	French	Description
	AE-003-041	17.5 Fr.	Cutting loop electrode 30° angled
	AE-003-043	17.5 Fr.	Cutting loop electrode straight
	AE-003-044	17.5 Fr.	Rollerball electrode Ø 2mm
	AE-003-045	17.5 Fr.	Roller electrode Ø 2mm
	AE-003-046	17.5 Fr.	Knife electrode



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.2 | SALINE RESECTOSCOPY

CONNECTING CABLES | BIPOLAR



CONNECTION INSTRUMENT SIDE



Connection Unit Side	Ordering no.	Length	Description
	AD-010-011	4.5m	Bipolar resectoscope connecting cable for Martin Maxium
	AD-010-012	4.5m	Bipolar resectoscope connecting cable for Erbe Vio
	AD-010-013	4.5m	Bipolar resectoscope connecting cable for Valleylab, Lamidey, EMC
	AD-010-015	4.5m	Bipolar resectoscope connecting cable for Erbe Vio 3

2.3 | LASER RESECTOSCOPY
CYSTOSCOPES AUTOCLAVABLE Ø 4.0MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013

2.3 | LASER RESECTOSCOPY WORKING ELEMENTS



Ordering no.	Description
AG-001-006	Working element used for resectoscopes, titanium, passive cutting, handle closed, 24/26 Fr. with probe channel 0.8mm
AG-001-009	Working element used for resectoscopes, titanium, passive cutting, handle closed, 24/26 Fr. with probe channel 1.2mm

Ordering no.	Description
AG-100-006	Adapter used for laser probes
AG-100-007	TUOHY-Borst adapter

2.3 | LASER RESECTOSCOPY SHEATHS | OBTURATORS



Ordering no.	Description	French	Colour code
AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-308	Inner sheath without ceramic isolation, with laser channel up to 7 Fr. , for use with AH-001-301	26 Fr.	Yellow



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow

2.3 | LASER SLIM RESECTOSCOPY
CYSTOSCOPES AUTOCLAVABLE Ø 2.9MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.3 | LASER SLIM RESECTOSCOPY WORKING ELEMENTS



Ordering no.	Description
AG-002-007	Working element used for resectoscopes, titanium, passive cutting, handle closed, 19/22 Fr. with probe channel 0.8mm
AG-002-008	Working element used for resectoscopes, titanium, passive cutting, handle closed, 19/22 Fr. with probe channel 1.2mm

2.3 | LASER SLIM RESECTOSCOPY SHEATHS



Ordering no.	Description
AH-002-302	Outer sheath 22 Fr. with perforation, rotatable



Ordering no.	Description
AH-002-301	Inner sheath 22 Fr. with ceramic isolation, rotatable

2.3 | LASER SLIM RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-011	Standard obturator for resectoscope sheaths 19 Fr.



Ordering no.	Description
AH-002-013	Viewing obturator for resectoscope sheaths 19 Fr.

2.3 | LASER MINI RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.3 | LASER MINI RESECTOSCOPY WORKING ELEMENTS



Ordering no.	Description
AG-002-016	Working element used for resectoscopes, titanium, passive cutting, handle closed, 17.5 Fr. with probe channel 0.8mm

2.3 | LASER MINI RESECTOSCOPY SHEATHS



Ordering no.	Description
AH-002-213	Outer sheath 18.5 Fr. with perforation



Ordering no.	Description
AH-002-211	Inner sheath 17.5 Fr. with ceramic isolation, rotatable



Ordering no.	Description
AH-002-216	Inner sheath 17.5 Fr. without ceramic isolation , with laser channel up to 5 Fr., for use with AH-001-213

2.3 | LASER MINI RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-015	Standard obturator 17.5 Fr.



Ordering no.	Description
AH-002-016	Viewing obturator 17.5 Fr. with 5 Fr. instrument channel

2.4 | ACCESSORIES

SUPRAPUBIC ASPIRATION



Ordering no.	Description
AB-007-001	Suprapubic aspiration set, complete set



Ordering no.	Description
AB-007-002	Trocar sleeve 16 Fr. with pressure measuring channel



Ordering no.	Description
AB-007-003	Adjustable stabilisation disc, lateral opening



Ordering no.	Description
AB-007-008	Guide sleeve, lateral opening , for inserting a suprapubic catheter

2.4 | ACCESSORIES

SUPRAPUBIC ASPIRATION



Ordering no.	Description
AB-007-004	Trocar with cannula tip



Ordering no.	Description
AB-007-005	Trocar with pyramidal tip



Ordering no.	Description
AB-007-006	Aspiration cannula



Ordering no.	Description
AB-007-007	Overflow cannula for pressure stabilisation

2.4 | ACCESSORIES


EVACUATOR | BLADDER SYRINGES



Ordering no.	Description
AJ-002-001	Ellik evacuator complete with sheath connection and locking device
AJ-002-002	Bowl only
AJ-002-003	Suction Ball only
AJ-002-004	Cone with tube only
AJ-002-005	Sheath connection for Ellik only




Not available in the EU




Ordering no.	Description
AJ-001-000	RECORD syringe with adapter, 50ml
AJ-001-010	Replacement cylinder for RECORD syringe with adapter, 50ml
AJ-001-001	RECORD syringe with adapter, 75ml
AJ-001-011	Replacement cylinder for RECORD syringe with adapter, 75ml
AJ-001-002	RECORD syringe with adapter, 100ml
AJ-001-012	Replacement cylinder for RECORD syringe with adapter, 100ml
AJ-001-003	RECORD syringe with adapter, 150ml
AJ-001-013	Replacement cylinder for RECORD syringe with adapter, 150ml

2.4 | ACCESSORIES

TRAY INSERT



Ordering no.	Description
AJ-010-020	Tray insert for resectoscopy (without instruments, without sterilization box), dimension 370x255x45mm (LxWxH)



RECOMMENDED PRODUCTS:

- AA-002-013** HD autoclavable Scope 30° | 4mm
- AG-001-011** SALINE working element, passive cutting, handle closed
- AH-001-302** Standard resectoscope sheath 24 Fr., with 1 fixed stopcock
- AH-001-301** Inner sheath 26 Fr. with ceramic isolation, rotatable
- AH-001-011** Standard obturator for resectoscope sheaths 24 Fr.
- AJ-010-006** Protection sleeve made of steel for sterilization and storage of electrodes and cold knives
- AE-001-031** Saline Cutting loop electrode 30° angled, 24 Fr. wire thickness 0,35mm
- AE-001-034** Saline Rollerball electrode 3mm, 24 Fr.
+ Connecting cable (depends on unit)

2.5 MONOPOLAR RESECTOSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable,
WL 302mm

1 AA-002-013



2 WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-001



3 OUTER SHEATH

26 Fr. with perforation, rotatable

3 AH-001-302



4 INNER SHEATH

26 Fr. with ceramic isolation, rotatable

4 AH-001-301



5 STANDARD OBTURATOR

24 Fr.

5 AH-001-011



6 CUTTING LOOP ELECTRODE

30° angled, 24 Fr.

6 AE-001-001



7 ROLLERBALL ELECTRODE 3 MM

24 Fr.

7 AE-001-004



8 KNIFE ELECTRODE

24 Fr.

8 AE-001-016



2.5 MONOPOLAR RESECTOSCOPY SET

9 HF MONOPOLAR CONNECTING CABLE
for Berchtold, Martin, Aesculap, Ø 4mm

9 AD-010-001



10 RECORD SYRINGE
100ml

10 AJ-001-002



11 ELLIK EVACUATOR
complete

11 AJ-002-001



12 BIOPSY FORCEPS
double action

12 AB-004-001



13 GRASPING FORCEPS
double action

13 AB-004-003



14 Adapter
for use with resectoscope sheaths

14 AB-004-008



15 PROTECTION SLEEVE
For sterilization and storage of electrodes

15 AJ-010-006



2.5 SALINE RESECTOSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable,
WL 302mm

1 AA-002-013



2 SALINE WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-011



3 OUTER SHEATH

26 Fr. with perforation, rotatable

3 AH-001-302



4 INNER SHEATH

26 Fr. with ceramic isolation, rotatable

4 AH-001-301



5 STANDARD OBTURATOR

24 Fr.

5 AH-001-011



6 CUTTING LOOP ELECTRODE

30° angled, 24 Fr.

6 AE-001-031



7 ROLLERBALL ELECTRODE 3 MM

24 Fr.

7 AE-001-034



8 KNIFE ELECTRODE

24 Fr.

8 AE-001-032



2.5 SALINE RESECTOSCOPY SET

- 9 HF BIPOLAR CONNECTING CABLE**
for Valleylab, Lamidey, EMC, 4.5mm

9 AD-010-013



- 10 RECORD SYRINGE**
100 ml

10 AJ-001-002



- 11 ELLIK EVACUATOR**
complete

11 AJ-002-001



- 12 BIOPSY FORCEPS**
double action

12 AB-004-001



- 13 GRASPING FORCEPS**
double action

13 AB-004-003



- 14 ADAPTER**
for use with resectoscope sheaths

14 AB-004-008



- 15 PROTECTION SLEEVE**
for sterilization and storage of electrodes

15 AJ-010-006



3. URETHROTOMY

3.1 | URETHROTOME

3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES

SHEATHS

OBTURATORS

WORKING ELEMENTS

KNIVES

3.3 | URETHROTOMY SET

3. URETHROTOMY

21 Fr.

15.5 Fr.


1 CYSTOSCOPES



2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS




4 KNIVES



3.1 | URETHROTOME



Ordering no.	Description
 AI-002-001	OTIS urethrotome with 3 exchangeable tips inclusive 2 knives

SPARE KNIVES



Ordering no.	Description
AI-002-002	Spare knives for OTIS urethrotome

3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES | SHEATHS



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-001
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-011
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-021



Ordering no.	Description
AI-001-001	Urethrotome sheath 21 Fr. with instrument channel 5 Fr.
AI-001-003	Endoscope bridge with 1 instrument channel , fit with urethrotome sheath AI-001-001

3.2 | OPTICAL URETHROTOMY SUPPLEMENTARY SHEATHS



Ordering no.	Description
AI-001-002	Supplementary sheath for AI-001-001, open to side



Ordering no.	Description
AI-001-006	Supplementary sheath 25 Fr. for continuous flow

3.2 | OPTICAL URETHROTOMY OBTURATORS



Ordering no.	Description
AI-001-004	Standard obturator for urethrotome sheath 21 Fr.



Ordering no.	Description
AI-001-005	Obturator for urethrotome sheath 21 Fr. with 5 Fr. instrument channel




Ordering no.	Description
AI-001-012	Viewing obturator for urethrotome sheath 21 Fr.

3.2 | OPTICAL URETHROTOMY

WORKING ELEMENTS | PASSIVE CUTTING



Ordering no.	Description
AG-001-001	Working element used for resectoscopes, passive cutting, handle closed , cutting spring assisted
AG-001-001-HF	Working element used for resectoscopes, passive cutting, handle closed , cutting spring assisted, HF connection on top



Ordering no.	Description
AG-001-002	Working element used for resectoscopes, passive cutting, handle open , cutting spring assisted

3.2 | OPTICAL URETHROTOMY

COLD KNIVES



	Ordering no.	Description
	AI-001-007	Cold knife straight
	AI-001-008	Cold knife straight, wave shape
	AI-001-009	Cold knife round
	AI-001-011	Cold knife hook shaped




Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES | SHEATHS | OBTURATORS



Viewing angle	Description	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-001
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-011



Ordering no.	Description
AI-005-001	Urethrotome sheath 15.5 Fr. with instrument channel 3 Fr
AI-005-003	Endoscope bridge with 1 instrument channel , fit with urethrotome sheath AI-005-001



Ordering no.	Description
AI-005-004	Standard obturator for urethrotome sheath 15.5 Fr.

3.2 | OPTICAL URETHROTOMY

WORKING ELEMENT | PASSIVE CUTTING



Ordering no.	Description
AG-001-010	Working element used for resectoscopes, passive cutting, handle closed , without HF-connection

3.2 | OPTICAL URETHROTOMY

COLD KNIVES



	Ordering no.	Description
	AI-005-007	Cold knife straight
	AI-005-008	Cold knife straight, wave shape
	AI-005-011	Cold knife hook shaped



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

3.4 OPTICAL UREHTROTOMY SET

1 HD CYSTOSCOPE

Ø 4.0mm, 0° viewing angle, autoclavable,
302mm working length


1 AA-002-011



2 WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-001



3 URETHROTOME SHEATH

21 Fr. with instrument channel 5 Fr.

3 AI-001-001



4 SUPPLEMENTARY SHEATH

25 Fr. for continous flow

4 AI-001-002



5 STANDARD OBTURATOR

24 Fr.

5 AI-001-004



6 ENDOSCOPE BRIDGE

with 1 instrument channel

6 AI-001-003



7 COLD KNIFE

straight

7 AI-001-007



8 PROTECTION SLEEVE

for sterilization and storage of electrodes

8 AJ-010-006



4. NEPHROSCOPY

4.1 | HD NEPHROSCOPY

HD NEPHROSCOPES
SHEATHS
OBTURATORS
HD MINI NEPHROSCOPE

4.2 | MINI NEPHROSCOPY

MINI NEPHROSCOPES
SHEATHS
DILATATORS
INSTRUMENTS

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 200MM
SHEATHS
OBTURATORS

NEPHROSCOPES 250MM
SHEATHS
OBTURATORS

4.4 | ACCESSORIES

TELESCOPE BOUGIES
FORCEPS
CANNULAS

4. NEPHROSCOPY

1 HD NEPHROSCOPES



2 SHEATHS & OBTURATORS



3 FORCEPS, INSTRUMENTS, KNIVES



4.1 | HD NEPHROSCOPY

HD NEPHROSCOPES | SHEATHS | OBTURATORS | SEALING CAPS



Ordering no.	Working length	Viewing angle	Description
AA-006-069	220mm	0°	HD Nephroscope - 45° angled ocular - Ø 18.8 Fr. - working channel Ø 12.9 Fr.
AA-006-070		15°	



Ordering no.	French	Description
AB-006-170	23 Fr.	Nephroscope sheath for continous flow, rotatable stopcock (fit with nephrosopes WL 220mm / scope 15°)



Ordering no.	French	Description
AB-006-171	23 Fr.	Obturator hollow (fit with nephrosocopy sheat AB-006-170 / scope 15°)


Ordering no.	Description
AA-006-081	Sealing cap, Ø 15mm/1.0mm
AA-006-083	Sealing cap, Ø 20mm/2.0mm

4.1 | HD MINI COMPACT NEPHROSCOPE



Ordering no.	Working length	Viewing angle	Description
AA-006-100	220mm	0°	HD Mini compact nephroscope with 30° angled ocular, sheath 12 Fr. , instrument channel 6.6 Fr. , 3 luerlock connections, 40.000 px

4.2 | MINI NEPHROSCOPY
MINI NEPHROSCOPES




Ordering no.	French	Working length	Viewing angle	Description
AA-006-050	Ø 19 Fr.	220mm	8°	Mini nephroscope with 45° angled ocular, instrument channel 10 Fr.
AA-006-051	Ø 13.5 Fr.	220mm	8°	Mini nephroscope with 45° angled ocular, instrument channel 6 Fr.
AA-006-052	Ø 13.5 Fr.	320mm	8°	Mini nephroscope with 45° angled ocular, instrument channel 6 Fr.

	Ordering no.	Description
	AA-006-053	Instrument bridge used for Mini-Nephroscopes (not included in Mini Nephroscope)
	AJ-010-001 AJ-010-003	Rubber cap for instrument channel Rubber cap for instrument channel and URS, Ø 0,6-0,8mm

4.2 | MINI NEPHROSCOPY

SHEATHS | DILATATORS | APLLICATOR



Ordering no.	French	Description
AB-006-100	15/18 Fr.	Mini nephroscope sheath without perforation, continous flow, stopcock fixed
AB-006-101	16.5/19.5 Fr.	
AB-006-102	21/24 Fr.	




Ordering no.	French	Description
AB-006-110	15/18 Fr.	Dilatator with central cannulation Ø 3 Fr.
AB-006-111	16.5/19.5 Fr.	Dilatator with 2 channels Ø 2x3 Fr.
AB-006-112	21/24 Fr.	Dilatator with 2 channels Ø 2x5 Fr.



Ordering no.	Description
AB-006-115	Applicator used for the mini nephroscope sheaths, incl. plunger

4.2 | MINI NEPHROSCOPY


INSTRUMENTS | FORCEPS FOR MINI NEPHROSCOPES




	Ordering no.	WL	Diameter	Description
	AB-006-242	400mm	Ø 1.8mm	Biopsy forceps serrated with ring handle
	AB-006-342	500mm	Ø 1.8mm	
	AB-006-442	400mm	Ø 3.0mm	
	AB-006-243	400mm	Ø 1.8mm	Biopsy forceps, single action with ring handle
	AB-006-343	500mm	Ø 1.8mm	
	AB-006-443	400mm	Ø 3.0mm	
	AB-006-244	400mm	Ø 1.8mm	Alligator grasping forceps with ring handle
	AB-006-344	500mm	Ø 1.8mm	
	AB-006-444	400mm	Ø 3.0mm	
	AB-006-245	400mm	Ø 1.8mm	Scissors, single action with ring handle
	AB-006-345	500mm	Ø 1.8mm	
	AB-006-445	400mm	Ø 3.0mm	
	AB-006-246	400mm	Ø 1.8mm	3-Prongs retriever for larger stones with coil handle
	AB-006-346	500mm	Ø 1.8mm	
	AB-006-446	400mm	Ø 3.0mm	

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 200MM



Ordering no.	Working length	Viewing angle	Description
AA-006-001	200mm	0°	Nephroscope - 45° angled ocular - Ø 7.45mm - suction channel Ø 1.5mm - working channel Ø 5.0mm
AA-006-021		6°	
AA-006-004		15°	



Ordering no.	Working length	Viewing angle	Description
AA-006-002	200mm	0°	Nephroscope - 90° angled ocular - Ø 7.45mm - suction channel Ø 1.5mm - working channel Ø 5.0mm
AA-006-022		6°	
AA-006-005		15°	



Ordering no.

Description

AA-006-080

Flow valve, Ø 6.1mm, used for STEMA nephroscopes

AA-006-081

Sealing cap, Ø 15mm/1.0mm, **for AA-006-080**

AA-006-083

Sealing cap, Ø 20mm/2.0mm, **for AA-006-001 - AA-006-034**


AJ-010-015

Stopcock left, large (passage 4mm), with LL-connection and tube attachment M6 x 0.35mm

4.3 | STANDARD NEPHROSCOPY SHEATHS | OBTURATORS



Ordering no.	French	Description
AB-006-001	24 Fr.	Nephroscope sheath for continous flow, rotatable stopcock
AB-006-007	26 Fr.	Nephroscope sheath for continous flow, rotatable stopcock
AB-006-005	24 Fr.	Nephroscope sheath for continous flow, fixed stopcock
AB-006-006	26 Fr.	Nephroscope sheath for continous flow, fixed stopcock




Ordering no.	French	Description
AB-006-003	26 Fr.	Nephroscope sheath for continous flow, perforated, rotatable stopcock
AB-006-008	26 Fr.	Nephroscope sheath for continous flow, perforated, fixed stopcock




Ordering no.	French	Description
AB-006-002	24 Fr.	Obturator hollow used for nephroscopy sheaths (200mm)
AB-006-004	26 Fr.	Obturator hollow used for nephroscopy sheaths (200mm)

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 250MM



Ordering no.	Working length	Viewing angle	Description
AA-006-011	250mm	0°	Nephroscope - 45° angled ocular - Ø 7.45mm - suction channel Ø 1.5mm - working channel Ø 5.0mm
AA-006-031		6°	
AA-006-034		15°	



Ordering no.	Working length	Viewing angle	Description
AA-006-012	250mm	0°	Nephroscope - 90° angled ocular - Ø 7.45mm - suction channel Ø 1.5mm - working channel Ø 5.0mm
AA-006-032		6°	
AA-006-006		15°	


	Ordering no.	Description
	AA-006-080	Flow valve, Ø 6.1mm, used for STEMA nephroscopes
	AA-006-081	Sealing cap, Ø 15mm/1.0mm, for AA-006-080
	AA-006-083	Sealing cap, Ø 20mm/2.0mm, for AA-006-001 - AA-006-034
	AJ-010-015	Stopcock left, large (passage 4mm), with LL-connection and tube attachment M6 x 0.35mm

4.3 | STANDARD NEPHROSCOPY

SHEATHS | OBTURATORS



Ordering no.	French	Description
AB-006-061	24 Fr.	Nephroscope sheath for continous flow, rotatable stopcock
AB-006-067	26 Fr.	
AB-006-065	24 Fr.	Nephroscope sheath for continous flow, fixed stopcock
AB-006-066	26 Fr.	



Ordering no.	French	Description
AB-006-063	26 Fr.	Nephroscope sheath for continous flow, perforated, rotatable stopcock
AB-006-068	26 Fr.	Nephroscope sheath for continous flow, perforated, fixed stopcock



Ordering no.	French	Description
AB-006-062	24 Fr.	Obturator hollow used for nephroscopy sheaths (250mm)
AB-006-064	26 Fr.	

4.4 | ACCESSORIES

TELESCOPE BOUGIES




Ordering no.	Description
AB-006-011	Set of telescope bougies consisting of: 6 dilating bougies 9, 12, 15, 18 ,21 and 24 Fr., incl. 1 rigid and 1 flexible guide probe
AB-006-012	Guide probe, rigid
AB-006-013	Guide probe, flexible

Ordering no.	French	Description
AB-006-020	9 Fr.	Dilatating sheath
AB-006-021	12 Fr.	
AB-006-022	15 Fr.	
AB-006-023	18 Fr.	
AB-006-024	21 Fr.	
AB-006-025	24 Fr.	
AB-006-026	27 Fr.	
AB-006-027	30 Fr.	

4.4 | ACCESSORIES

FORCEPS



	Ordering no.	WL	Diameter	Description
	AB-006-046	380mm	Ø 3.5mm	Fenestrated stone grasping forceps for smaller stones with coil handle
	AB-006-040	430mm		
	AB-006-048	380mm	Ø 3.5mm	Fenestrated stone grasping forceps for smaller stones with ring handle
	AB-006-042	430mm		
	AB-006-047	380mm	Ø 3.5mm	3-Prongs retriever for larger stones with coil handle
	AB-006-041	430mm		
	AB-006-049	380mm	Ø 3.5mm	Biopsy forceps, single action with ring handle
	AB-006-043	430mm		
	AB-006-050	380mm	Ø 3.5mm	Alligator grasping forceps with ring handle
	AB-006-044	430mm		
	AB-006-051	380mm	Ø 3.5mm	Scissors, single action with ring handle
	AB-006-045	430mm		



4.4 | ACCESSORIES CANNULAS

Ordering no.	Diameter	Description
AB-006-052		Puncture cannula, inner and outer cannula, 5 Pack
AB-006-053	Ø 3mm	Dilatation cannula, inner and outer cannula

5. URETERO-RENOSCOPY

5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES
URETERO-RENOSCOPES
BICHLER URETERO-RENOSCOPES
MARBERGER URETERO-RENOSCOPES
FORCEPS, SCISSORS, INSTRUMENTS
STONE EXTRACTORS
MONOPOLAR FLEXIBLE ELECTRODES
MONOPOLAR CONNECTING CABLES

5.2 | FLEXIBLE URETERO-RENSCOPE

FLEX II
ACCESSORIES

5. HD URETERO-RENOSCOPES

40.000 px

17.000 px



1 URETERO-RENOSCOPE



2 FORCEPS



5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-003	425mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40.000 px
AA-007-016		5°	
AA-007-017		12°	

Ordering no.	Working length	Viewing angle	Description
AA-007-004	310mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40.000 px
AA-007-018		5°	
AA-007-019		12°	


Ordering no.	Working length	Viewing angle	Description
AA-007-005	150mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40.000 px
AA-007-036		5°	
AA-007-037		12°	

Delivery contains

HD uretero-renaloscope, brush Ø 3.0mm (length 590mm), brush Ø 3.8mm (length 350mm), sealing cap silicone (grey), Duckbill sealing

5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES




Ordering no.	Working length	Viewing angle	Description
AA-007-009	425mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 to 13.2 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40.000 px
AA-007-028		5°	
AA-007-029		12°	

Delivery contains

HD uretero-renaloscope, brush Ø 3.0mm (length 590mm), brush Ø 3.8mm (length 350mm), sealing cap silicone (grey), Duckbill sealing

5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-015	425mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17.000 px
AA-007-034		5°	
AA-007-035		12°	

Ordering no.	Working length	Viewing angle	Description
AA-007-014	310mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17.000 px
AA-007-032		5°	
AA-007-033		12°	


Ordering no.	Working length	Viewing angle	Description
AA-007-013	140mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17.000 px
AA-007-030		5°	
AA-007-031		12°	

Delivery contains

HD uretero-renaloscope, brush Ø 3.0mm (length 590mm), brush Ø 3.8mm (length 350mm), sealing cap silicone (grey), Duckbill sealing

5.1 | URETERO-RENOSCOPY

HD URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-006	425mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 4.9 Fr. to 6.5 Fr. - Working channel 1 x 3 Fr. - 17.000 px
AA-007-022		5°	
AA-007-023		12°	

Ordering no.	Working length	Viewing angle	Description
AA-007-007	310mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 4.9 Fr. to 6.5 Fr. - Working channel 1 x 3 Fr. - 17.000 px
AA-007-024		5°	
AA-007-025		12°	


Ordering no.	Working length	Viewing angle	Description
AA-007-008	150mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 4.9 Fr. to 6.5 Fr. - Working channel 1 x 3 Fr. - 17.000 px
AA-007-026		5°	
AA-007-027		12°	

Delivery contains

HD uretero-renaloscope, brush Ø 3.0mm (length 590mm), brush Ø 3.8mm (length 350mm), sealing cap silicone (grey), Duckbill sealing

5.1 | URETERO-RENOSCOPES SEMI-RIGID

URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-001	310mm	8°	Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 Fr. to 10 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - Field of view: 85°
AA-007-002	425mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-011	310mm	8°	Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 Fr. to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - Field of view: 100°
AA-007-012	425mm		




Ordering no.	Description
AA-007-021	Instrument bridge used for uretero-renaloscopes (included in the delivery content)
AJ-010-001 AJ-010-003	Rubber cap for instrument channel Rubber cap for instrument channel and URS, Ø 0,6-0,8mm



5.1 | URETERO-RENOSCOPES SEMI-RIGID

BICHLER URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-056	310mm	5°	BICHLER Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 4.5 Fr. to 6.5 Fr. - Working channel 1 x 3.3 Fr.
AA-007-055	430mm		


Ordering no.	Working length	Viewing angle	Description
AA-007-043	310mm	12°	BICHLER Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 Fr. to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr.
AA-007-044	430mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-047	310mm	12°	BICHLER Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 8 Fr. to 9.8 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr.
AA-007-048	430mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-051	310mm	12°	BICHLER Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 8.5 Fr. to 11.5 Fr. - Working channel 1 x 6 Fr. or 2 x 4 Fr.
AA-007-052	430mm		

5.1 | URETERO-RENSCOPES SEMI-RIGID

MARBERGER URETERO-RENSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-054	310mm	5°	MARBERGER Uretero-roscope semi-rigid - 45° angled ocular - Working Ø steps by 4.5 Fr. to 6.5 Fr. - Working channel 1 x 3.3 Fr.
AA-007-053	430mm		


Ordering no.	Working length	Viewing angle	Description
AA-007-041	310mm	12°	MARBERGER Uretero-roscope semi-rigid - 45° angled ocular - Working Ø steps by 6 Fr. to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr.
AA-007-042	430mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-045	310mm	12°	MARBERGER Uretero-roscope semi-rigid - 45° angled ocular - Working Ø steps by 8 Fr. to 9.8 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr.
AA-007-046	430mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-049	310mm	12°	MARBERGER Uretero-roscope semi-rigid - 45° angled ocular - Working Ø steps by 8.5 Fr. to 11.5 Fr. - Working channel 1 x 6 Fr. or 2 x 4 Fr.
AA-007-050	430mm		

5.1 | URETERO-RENOSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps, single action	Ø 2.3mm	7 Fr.	AC-010-014	AC-010-214	
	Grasping forceps, single action	Ø 2.3mm	7 Fr.		AC-011-214	
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-011	AC-001-211	AC-001-011-R
		Ø 1.2mm	4 Fr.	AC-001-010	AC-001-210	AC-001-010-R
		Ø 1.6mm	5 Fr.	AC-001-012	AC-001-212	AC-001-012-R
		Ø 2.1mm	6 Fr.	AC-001-013	AC-001-213	AC-001-013-R
		Ø 2.3mm	7 Fr.	AC-001-014	AC-001-214	AC-001-014-R
		Ø 3.0mm	9 Fr.	AC-001-015	AC-001-215	AC-001-015-R
	Biopsy forceps oval, serrated, double action	Ø 1.0mm	3 Fr.	AC-002-011	AC-002-211	AC-002-011-R
		Ø 1.6mm	5 Fr.	AC-002-012	AC-002-212	AC-002-012-R
		Ø 2.1mm	6 Fr.	AC-002-013	AC-002-213	AC-002-013-R
		Ø 2.3mm	7 Fr.	AC-002-014	AC-002-214	AC-002-014-R
		Ø 3.0mm	9 Fr.	AC-002-015	AC-002-215	AC-002-015-R
	Biopsy forceps oval, serrated, teeth interlock, double action	Ø 1.6mm	5 Fr.	AC-015-012	AC-015-212	AC-015-012-R
		Ø 2.1mm	6 Fr.		AC-015-213	AC-015-013-R
		Ø 2.3mm	7 Fr.	AC-015-014	AC-015-214	AC-015-014-R
	Grasping forceps, double action	Ø 1.2mm	4 Fr.	AC-017-011	AC-017-211	AC-017-011-R
		Ø 2.3mm	7 Fr.	AC-017-014	AC-017-214	AC-017-014-R
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-006-011	AC-006-211	AC-006-011-R
		Ø 1.2mm	4 Fr.	AC-006-013	AC-006-213	AC-006-013-R
		Ø 1.6mm	5 Fr.	AC-006-012	AC-006-212	AC-006-012-R
		Ø 2.3mm	7 Fr.	AC-006-014	AC-006-214	AC-006-014-R
		Ø 3.0mm	9 Fr.	AC-006-015	AC-006-215	AC-006-015-R

5.1 | URETERO-RENSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	Ø 2.3mm	7 Fr.	AC-008-014	AC-008-214	AC-008-014-R
		Ø 3.0mm	9 Fr.	AC-008-015	AC-008-215	AC-008-015-R
	Grasping forceps teeth double action	Ø 1.2mm	4 Fr.	AC-016-011	AC-016-211	AC-016-011-R
		Ø 1.6mm	5 Fr.	AC-016-012	AC-016-212	AC-016-012-R
	Grasping forceps teeth interlock, double action	Ø 1.6mm	5 Fr.	AC-020-012	AC-020-212	AC-020-012-R
		Ø 2.3mm	7 Fr.	AC-020-014	AC-020-214	AC-020-014-R
	Grasping forceps rat tooth, double action	Ø 1.0mm	3 Fr.	AC-019-011	AC-019-211	AC-019-011-R
		Ø 1.6mm	5 Fr.	AC-019-012	AC-019-212	AC-019-012-R
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-011	AC-007-211	AC-007-011-R
		Ø 1.6mm	5 Fr.	AC-007-012	AC-007-212	AC-007-012-R
	Grasping forceps for big stones, double action	Ø 1.2mm	4 Fr.	AC-009-011	AC-009-211	AC-009-011-R
		Ø 1.6mm	5 Fr.	AC-009-012	AC-009-212	AC-009-012-R
		Ø 2.3mm	7 Fr.	AC-009-014	AC-009-214	AC-009-014-R
	3-Prongs retriever	Ø 1.0mm	3 Fr.	AC-012-011		
		Ø 1.6mm	5 Fr.	AC-012-012		
		Ø 2.3mm	7 Fr.	AC-012-014		


5.1 | URETERO-RENOSCOPY

SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-011	AC-003-211	AC-003-011-R
		Ø 1.2mm	4 Fr.	AC-003-010	AC-003-210	AC-003-010-R
		Ø 1.6mm	5 Fr.	AC-003-012	AC-003-212	AC-003-012-R
		Ø 2.1mm	6 Fr.	AC-003-013	AC-003-213	AC-003-013-R
		Ø 2.3mm	7 Fr.	AC-003-014	AC-003-214	AC-003-014-R
		Ø 3.0mm	9 Fr.	AC-003-015	AC-003-215	AC-003-015-R
	Uretero slitting scissors, retrograde cutting, single action	Ø 1.6mm	4 Fr.	AC-018-011	AC-018-211	AC-018-011-R
		Ø 2.3mm	5 Fr.	AC-018-012	AC-018-212	AC-018-012-R
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-014	AC-004-214	AC-004-014-R
		Ø 3.0mm	9 Fr.	AC-004-015	AC-004-215	AC-004-015-R
	Scissors, serrated flexible, double action					
		Ø 2.3mm	7 Fr.	AC-005-014	AC-005-214	AC-005-014-R

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




FEATURES

- 1 QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
- 2 ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
- 3 OVERLOAD PROTECTION AVAILABLE
- 4 SIDE-MOUNTED LUER LOCK

5.1 | URETERO-RENOSCOPY

INSTRUMENTS SEMI-RIGID | DETACHABLE




Ordering no.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-210-T	AC-003-210-OT
		Ø 1.2mm	4 Fr.	AC-003-211-T	AC-003-211-OT
		Ø 1.6mm	5 Fr.	AC-003-212-T	AC-003-212-OT
		Ø 2.3mm	7 Fr.	AC-003-214-T	AC-003-214-OT
		Ø 3.0mm	9 Fr.	AC-003-215-T	AC-003-215-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-004-622-T	AC-004-622-OT
		Ø 2.3mm	7 Fr.	AC-004-624-T	AC-004-624-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-214-T	AC-004-214-OT
		Ø 3.0mm	9 Fr.	AC-004-215-T	AC-004-215-OT

5.1 | URETERO-RENOSCOPY

INSTRUMENTS SEMI-RIGID | DETACHABLE




Ordering no.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-211-T	AC-001-211-OT
		Ø 1.2mm	4 Fr.	AC-001-210-T	AC-001-210-OT
		Ø 1.6mm	5 Fr.	AC-001-212-T	AC-001-212-OT
		Ø 2.3mm	7 Fr.	AC-001-214-T	AC-001-214-OT
		Ø 3.0mm	9 Fr.	AC-001-215-T	AC-001-215-OT
	Biopsy forceps oval, serrated, double action	Ø 1.2mm	4 Fr.	AC-016-211-T	AC-016-211-OT
		Ø 1.6mm	5 Fr.	AC-016-212-T	AC-016-212-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-006-211-T	AC-006-211-OT
		Ø 1.2mm	4 Fr.	AC-006-213-T	AC-006-213-OT
		Ø 1.6mm	5 Fr.	AC-006-212-T	AC-006-212-OT
		Ø 2.3mm	7 Fr.	AC-006-214-T	AC-006-214-OT
		Ø 3.0mm	9 Fr.	AC-006-215-T	AC-006-215-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-211-T	AC-007-211-OT
		Ø 1.6mm	5 Fr.	AC-007-212-T	AC-007-212-OT
	Grasping forceps for big stones, double action	Ø 1.2mm	5 Fr.	AC-009-211-T	AC-009-211-OT
		Ø 1.6mm	5 Fr.	AC-009-212-T	AC-009-212-OT

5.1 | URETERO-RENOSCOPY


INSTRUMENTS RIGID | DETACHABLE



Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-103-210-T	AC-103-210-OT
		Ø 1.2mm	4 Fr.	AC-103-211-T	AC-103-211-OT
		Ø 1.6mm	5 Fr.	AC-103-212-T	AC-103-212-OT
		Ø 2.3mm	7 Fr.	AC-103-214-T	AC-103-214-OT
		Ø 3.0mm	9 Fr.	AC-103-215-T	AC-103-215-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-104-622-T	AC-104-622-OT
		Ø 2.3mm	7 Fr.	AC-104-624-T	AC-104-624-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-104-214-T	AC-104-214-OT
		Ø 3.0mm	9 Fr.	AC-104-215-T	AC-104-215-OT
	Biopsy punch forceps, single action	Ø 1.6mm	5 Fr.	AC-114-212-T	AC-114-212-OT

5.1 | URETERO-RENOSCOPY INSTRUMENTS RIGID | DETACHABLE



Ordering no.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Biopsy forceps oval, double action	3 Fr. Ø 1.0mm	3 Fr.	AC-101-211-T	AC-101-211-OT
		Ø 1.2mm	4 Fr.	AC-101-210-T	AC-101-210-OT
		Ø 1.6mm	5 Fr.	AC-101-212-T	AC-101-212-OT
		Ø 2.3mm	7 Fr.	AC-101-214-T	AC-101-214-OT
		Ø 3.0mm	9 Fr.	AC-101-215-T	AC-101-215-OT
	Biopsy forceps oval, double action	Ø 1.2mm	4 Fr.	AC-116-211-T	AC-116-211-OT
		Ø 1.6mm	5 Fr.	AC-116-212-T	AC-116-212-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-106-211-T	AC-106-211-OT
		Ø 1.2mm	4 Fr.	AC-106-213-T	AC-106-213-OT
		Ø 1.6mm	5 Fr.	AC-106-212-T	AC-106-212-OT
		Ø 2.3mm	7 Fr.	AC-106-214-T	AC-106-214-OT
		Ø 3.0mm	9 Fr.	AC-106-215-T	AC-106-215-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-107-211-T	AC-107-211-OT
		Ø 1.6mm	5 Fr.	AC-107-212-T	AC-107-212-OT
	Grasping forceps for big stones, double action	4 Fr. Ø 1.2mm	5 Fr.	AC-109-211-T	AC-109-211-OT
		Ø 1.6mm	5 Fr.	AC-109-212-T	AC-109-212-OT

5.1 | URETERO-RENSCOPY

STONE EXTRACTORS | 4 WIRES



Ordering no.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering no.	Wires	French	Working length	Description
AF-001-001	4	3 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-002	4	3 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-015	4	3 Fr.	900mm	Dormia stone extractor, basket with tubus only
AF-001-016	4	3 Fr.	900mm	Dormia stone extractor complete with handle
AF-001-023	4	3 Fr.	1200mm	Dormia stone extractor, basket with tubus only
AF-001-022	4	3 Fr.	1200mm	Dormia stone extractor complete with handle
AF-001-003	4	5 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-004	4	5 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-017	4	5 Fr.	900mm	Dormia stone extractor, basket with tubus only
AF-001-019	4	5 Fr.	900mm	Dormia stone extractor complete with handle
AF-001-005	4	7 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-006	4	7 Fr.	600mm	Dormia stone extractor complete with handle

5.1 | URETERO-RENOSCOPY

STONE EXTRACTORS | 5-7 WIRES



Ordering no.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering no.	Wires	French	Working length	Description
AF-001-007	5	5 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-008	5	5 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-025	5	5 Fr.	900mm	Dormia stone extractor, basket with tubus only
AF-001-018	5	5 Fr.	900mm	Dormia stone extractor complete with handle
AF-001-009	5	7 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-010	5	7 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-011	6	5 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-012	6	5 Fr.	600mm	Dormia stone extractor complete with handle
AF-001-013	6	7 Fr.	600mm	Dormia stone extractor, basket with tubus only
AF-001-014	7	7 Fr.	600mm	Dormia stone extractor complete with handle

5.1 | URETERO-RENSCOPY

MONOPOLAR FLEXIBLE ELECTRODES




	Description	French	Length	Ordering no.	Length	Ordering no.
	Loop electrode flexible, monopolar	5 Fr.	530mm	AD-003-003	650mm	AD-003-013
		7 Fr.	530mm	AD-003-005	650mm	AD-003-015
		9 Fr.	530mm	AD-003-007	650mm	AD-003-017
		10 Fr.	530mm	AD-003-008	650mm	AD-003-018
	Needle electrode flexible, monopolar	3 Fr.	530mm	AD-002-001	650mm	AD-002-011
		4 Fr.	530mm	AD-002-002	650mm	AD-002-012
		5 Fr.	530mm	AD-002-003	650mm	AD-002-013
		7 Fr.	530mm	AD-002-005	650mm	AD-002-015
	Button electrode flexible, monopolar	3 Fr.	530mm	AD-001-001	650mm	AD-001-011
		4 Fr.	530mm	AD-001-002	650mm	AD-001-012
		5 Fr.	530mm	AD-001-003	650mm	AD-001-013
		6 Fr.	530mm	AD-001-004	650mm	AD-001-014
		7 Fr.	530mm	AD-001-005	650mm	AD-001-015
		8 Fr.	530mm	AD-001-006	650mm	AD-001-016
		9 Fr.	530mm	AD-001-007	650mm	AD-001-017
		10 Fr.	530mm	AD-001-008	650mm	AD-001-018

5.1 | URETERO-RENOSCOPY

MONOPOLAR CONNECTING CABLES



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4mm	
	AD-010-002	3m	Ø 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5mm	
	AD-010-003	3m	Ø 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8mm	
	AD-010-024	4.5m	Ø 2mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

5.2 | FLEXIBLE URETERO-RENSCOPE FLEX II

- *High luminance*
- *Exact adjustment of focus*
- *Precise up-/down angulation*



Ordering no.	Product
AA-008-003	Flexible uretero-renaloscope FLEX II

Technical data

Cystoscope

Insertion tube	Ø 2.8mm
Working Channel	Ø 1.2mm
Working length	680mm
Field of view	88°
Pixel	18.000
Angulation up-down	275°
Accessories	Suction valve, leak tester

Delivery contains

Flexible uretero-renaloscope FLEX II, leak tester, suction valve

INCLUDED IN SCOPE OF DELIVERY:

AA-008-103

Leak tester

AA-008-105

Suction valve



AA-008-103



AA-008-105

5.2 | FLEXIBLE URETERO-RENSCOPE ACCESSORIES



Ordering no.	Description
AA-008-100	Bridge with Luer-Lock and 2 stopcocks for AA-008-002

Ordering no.	Description
AA-008-101	Single use cleaning brushes, 1.0mm, 120cm length, 1 box = 250 pieces



Ordering no.	Description
AA-008-102	Suitcase grey for STEMA flexible scopes



Ordering no.	Description
AA-008-103	Leak tester



Ordering no.	Description
AA-008-104	Single use Luer-Lock adapter, 1 box = 25 pieces



Ordering no.	Description
AA-008-105	Suction valve

6. LITHOTRIPSY


6.1 | LITHOTRIPSY
CYTOSCOPIES
SHEATHS
OBTURATORS
WORKING ELEMENTS

6. LITHOTRIPSY


1 SCOPES



2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS



6.1 | LITHOTRIPSY CYSTOSCOPES | SHEATHS



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
0°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-001
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-011

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013



Ordering no.	Description
AB-002-010	Lithotripsy sheath with central valve used for AB-004-011



Ordering no.	Description
AB-002-013	Lithotripsy sheath used for AB-004-011

6.1 | LITHOTRIPSY

SHEATHS | OBTURATORS



Ordering no.	Description
AB-002-011	Diagnostic sheath with working channel for instruments 7 Fr.



Ordering no.	Description
AB-002-014	Adapter for use with Ellik



Ordering no.	Description
AB-002-012	Obturator for AB-002-010 + AB-002-013

6.1 | LITHOTRIPSY

WORKING ELEMENT | LITHOTRIPTOR



Ordering no.	Description
AB-004-011	Working element used for punch lithotripsy



Ordering no.	Description
AB-004-012	Mechanical lithotripter with handle and turning screw, Ø 24 Fr., 2 stopcocks with luer lock adapter, for use with scope AA-002-004

7. MORCELLATION

7.1 | MORCELLATION 26 FR.
NEPHROSCOPE
SHEATHS
OBTURATORS

7.2 | MORCELLATION 22 FR.
MORCESCOPE
SHEATHS
OBTURATORS


7. MORCELLATION

26 Fr.

22 Fr.


1 HD NEPHROSCOPE

2 SHEATHS
& OBTURATORS



3 MORCESCOPE

4 SHEATHS
& OBTURATORS





26 Fr.

7.1 | MORCELLATION
NEPHROSCOPE



Ordering no.	Working length	Viewing angle	Description
AA-006-012	250mm	0°	Nephroscope with 90° angled ocular - Ø 7.45mm - Suction channel Ø 1.5mm - Working channel Ø 5.0mm

	Ordering no.	Description
	AB-004-010-QL	Adapter for use with nephroscope and resectoscope sheaths 26 Fr. STEMA Click-System
	AA-006-080	Flow valve, Ø 6.1mm, used for STEMA nephroscopes
	AA-006-081	Sealing cap, Ø 15mm/1.0mm, for AA-006-080
	AA-006-083	Sealing cap, Ø 20mm/2.0mm, for AA-006-001 - AA-006-034

7.1 | MORCELLATION

SHEATHS | OBTURATORS



Ordering no.	Description	French	Colour code
AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow




Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow

7.2 | MORCELLATION MORCESCOPE



Ordering no.	Working length	Direction of view	Description
AA-006-072	180mm	0°	Morcescope with 45° angled ocular - 22 Fr. outer tube - Working channel 15 Fr. - Field of view 85° - 40.000 px

Ordering no.	Description
AA-006-084	Sealing cap white, Ø 13.5mm/10 pieces

7.2 | MORCELLATION SHEATHS | OBTURATORS



Ordering no.	Description	French	Colour code
AH-002-302	Outer sheath with perforation, rotatable	22 Fr.	White



Ordering no.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, rotatable	22 Fr.	White



Ordering no.	Description	French	For outer sheath	Colour code
AH-002-011	Standard obturator for resectoscope sheaths	19 Fr.	22 Fr.	White



Ordering no.	Description	French	For outer sheath	Colour code
AH-002-013	Viewing obturator	19 Fr.	22 Fr.	White

8. PEDIATRIC UROLOGY

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
FORCEPS
MONOPOLAR ELECTRODES
CONNECTING CABLES

8.2 | PEDIATRIC RESECTOSCOPY

CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
ELECTRODES
CONNECTING CABLES

8.3 | PEDIATRIC URETHROTOMY


CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
COLD KNIVES AND ACCESSORIES

8. PEDIATRIC UROLOGY

1 PEDIATRIC CYSTOSCOPY



2 PEDIATRIC RESECTOSCOPY



3 PEDIATRIC URETHROTOMY




8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-005-001	Pediatric Cystoscope autoclavable	0°	Ø 2.7mm	187.5mm
AA-005-003	Pediatric Cystoscope autoclavable	30°	Ø 2.7mm	187.5mm
AA-005-004	Pediatric Cystoscope autoclavable	70°	Ø 2.7mm	187.5mm



Ordering no.	Description	French	Colour code
AB-001-001	Cysto-Urethroscope sheath with 2 fixed stopcocks, 1 x 3 Fr. instrument channel incl. obturator	13 Fr.	Red
AB-001-001-13	Cysto-Urethroscope sheath with 2 fixed stopcocks, 1 x 5 Fr. or 2 x 4 Fr. instrument channel incl. obturator	13 Fr.	Red



Ordering no.	Description	For outer sheath	Colour code
AB-001-011	Standard obturator for Cysto-urethroscope sheaths	13 Fr.	Red



13 Fr.

8.1 | PEDIATRIC CYSTO-URETHROSCOPY
WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering no.	Description	For outer sheath	Colour code
AB-003-005	Working element with ALBARRAN deflector	13 Fr.	Red



Ordering no.	Description	For outer sheath	Colour code
AB-003-001	Endoscope bridge	13 Fr.	Red

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-004-011	Pediatric Cystoscope autoclavable	0°	Ø 2.0mm	187mm
AA-004-012	Pediatric Cystoscope autoclavable	30°	Ø 2.0mm	187mm



Ordering no.	Description	French	Colour code
AB-001-001-9	Cysto-Urethroscope sheath with 2 fixed stopcocks, 1x3 FR. instrument channel incl. obturator	9.5 Fr.	Blue
AB-001-001-11	Cysto-Urethroscope sheath with 2 fixed stopcocks, 1x3 FR. instrument channel incl. obturator	11 Fr.	Green



Ordering no.	Description	For outer sheath	Colour code
AB-001-011-9	Standard obturator for Cysto-urethroscope sheaths	9.5 Fr.	Blue
AB-001-011-11	Standard obturator for Cysto-urethroscope sheaths	11 Fr.	Green

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-004-013	Pediatric Cystoscope autoclavable	0°	Ø 1.9mm	145mm
AA-004-014	Pediatric Cystoscope autoclavable	30°	Ø 1.9mm	145mm




Ordering no.	Description	French	Colour code
AB-001-001-7	Cysto-Urethroscope sheath with 2 fixed stopcocks	7.5 Fr.	Yellow



Ordering no.	Description	For outer sheath	Colour code
AB-001-011-7	Standard obturator for Cysto-urethroscope sheaths	7.5 Fr.	Yellow

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

FORCEPS | ELECTRODES



	Description	WL	Diameter	French	Flexible	Semi-rigid
	Biopsy forceps oval, double action	28cm	Ø1.0mm	3 Fr.	AC-001-021	AC-001-221
	Biopsy forceps oval, serrated, teeth interlock, double action	28cm	Ø1.0mm	3 Fr.	AC-015-021	AC-015-221
	Grasping forceps, alligator jaw, double action	28cm	Ø1.0mm	3 Fr.	AC-006-021	AC-006-221
	Grasping forceps, rat tooth, double action	28cm	Ø1.0mm	3 Fr.	AC-019-021	
	Scissors, Single action	28cm	Ø1.0mm	3 Fr.	AC-003-021	AC-003-221



	Ordering no.	French	Diameter	Description
	AD-002-021	3 Fr.	280mm	Needle electrode flexible, monopolar
	AD-001-021	3 Fr.	280mm	Button electrode flexible, monopolar


8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CONNECTING CABLES | MONOPOLAR



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4mm	
	AD-010-002	3m	Ø 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5mm	
	AD-010-003	3m	Ø 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8mm	
	AD-010-024	4.5m	Ø 2mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

8.2 | PEDIATRIC RESECTOSCOPY

RESECTOSCOPES | SHEATHS | OBTURATORS



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-005-001	Pediatric Cystoscope autoclavable	0°	Ø 2.7mm	187.5mm



Ordering no.	Description	French	Colour code
AH-003-001	Standard resectoscope sheath with 1 rotatable stopcock, without obturator	13 Fr.	Red



Ordering no.	Description	For outer sheath	Colour code
AH-003-011	Standard obturator for Cysto-urethroscope sheaths	13 Fr.	Red

8.2 | PEDIATRIC RESECTOSCOPY

WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering no.	Description	For outer sheath
AG-003-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	13 Fr.



Ordering no.	Description	For outer sheath
AI-003-003	Endoscope bridge with 3 Fr. instrument channel	13 Fr.

8.2 | PEDIATRIC RESECTOSCOPY ELECTRODES



	Ordering no.	French	Color code	Description
	AE-004-002	13 Fr.	Red	Cutting loop electrode angled
	AE-004-018	13 Fr.	Red	Cutting loop electrode small angled
	AE-004-004	13 Fr.	Red	Rollerball electrode
	AE-004-016	13 Fr.	Red	Knife electrode angled
	AE-004-019	13 Fr.	Red	Hook electrode without ball



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

8.2 | PEDIATRIC RESECTOSCOPY

RESECTOSCOPES | SHEATHS | OBTURATORS



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-004-011	Pediatric Cystoscope autoclavable	0°	Ø 2.0mm	187.5mm



Ordering no.	Description	French	Colour code
AH-005-001	Standard resectoscope sheath with 1 rotatable stopcock , without obturator	9 Fr.	Blue
AH-004-001	Standard resectoscope sheath with 1 rotatable stopcock , without obturator	11 Fr.	Green



Ordering no.	Description	For outer sheath	Colour code
AH-005-011	Standard obturator for resectoscope sheaths	9 Fr.	Blue
AH-004-011	Standard obturator for resectoscope sheaths	11 Fr.	Green

8.2 | PEDIATRIC RESECTOSCOPY

WORKING ELEMENTS



Ordering no.	Description	For outer sheath
AG-005-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	9 Fr.
AG-004-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	11 Fr.

8.2 | PEDIATRIC RESECTOSCOPY ELECTRODES




	Ordering no.	French	Color code	Description
	AE-006-002	9 Fr.	Blue	Cutting loop electrode angled
	AE-005-002	11 Fr.	Green	
	AE-006-018	9 Fr.	Blue	Cutting loop electrode small angled
	AE-005-018	11 Fr.	Green	
	AE-006-004	9 Fr.	Blue	Rollerball electrode
	AE-005-004	11 Fr.	Green	
	AE-006-016	9 Fr.	Blue	Knife electrode angled
	AE-005-016	11 Fr.	Green	
	AE-006-019	9 Fr.	Blue	Hook electrode without ball
	AE-005-019	11 Fr.	Green	



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives


8.1 | PEDIATRIC RESECTOSCOPY

CONNECTING CABLES | MONOPOLAR



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4mm	HF Monopolar connecting cable for Berchold, Martin, Aesculap
	AD-010-006	5m	Ø 4mm	
	AD-010-002	3m	Ø 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5mm	
	AD-010-003	3m	Ø 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8mm	
	AD-010-024	4.5m	Ø 2mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

8.1 | PEDIATRIC RESECTOSCOPY

CONNECTING CABLES | MONOPOLAR



CONNECTION UNIT SIDE



Instrument Side	Ordering no.	Length	Description
	AD-010-025	4.5m	HF Monopolar connecting cable for Storz resectoscopes, 90° angled connector
	AD-010-026	4.5m	HF Monopolar connecting cable for Olympus resectoscopes, diameter Ø 3mm



RECOMMENDED PRODUCTS:

CI-001-030 STEMA HF-U PRO TOUCH

- For all monopolar & bipolar techniques
- with 7" touch screen display
- provided to perform the most advanced electrosurgery
- ADC / APC System

CI-001-040 DOUBLE PEDAL FOOTSWITCH

8.3 | PEDIATRIC URETHROTOMY

CYSTOSCOPES | SHEATHS | OBTURATORS

Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-005-001	Pediatric Cystoscope autoclavable	0°	Ø2.7mm	187.5mm



Ordering no.	Description	French	Colour code
AI-003-001	Standard urethrotome sheath with 2 fixed stopcocks incl. obturator	13 Fr.	Red



Ordering no.	Description	For outer sheath	Colour code
AI-003-004	Standard obturator for urethrotome sheath	13 Fr.	Red

8.3 | PEDIATRIC URETHROTOMY

WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering no.	Description	For outer sheath
AG-003-001	Working element used for pediatric urethrosopes, passive cutting , handle open	13 Fr.



Ordering no.	Description	For outer sheath
AI-003-003	Endoscope bridge with 3 Fr. instrument channel	13 Fr.

8.3 | PEDIATRIC URETHROTOMY

COLD KNIVES



	Ordering no.	French	Color code	Description
	AI-003-007	13 Fr.	Red	Cold knife straight
	AI-003-009	13 Fr.	Red	Cold knife round
	AI-003-010	13 Fr.	Red	Cold knife ring shaped
	AI-003-011	13 Fr.	Red	Cold knife hook shaped



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

8.3 | PEDIATRIC URETHROTOMY

CYSTOSCOPES | SHEATHS | OBTURATORS

Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-004-011	Pediatric Cystoscope autoclavable	0°	Ø2.0mm	187.5mm



Ordering no.	Description	French	Colour code
AI-004-001	Standard urethrotome sheath with 2 fixed stopcocks incl. obturator	11 Fr.	Green



Ordering no.	Description	For outer sheath	Colour code
AI-004-004	Standard obturator for urethrotome sheath	11 Fr.	Green

8.3 | PEDIATRIC URETHROTOMY

WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering no.	Description	For outer sheath
AG-004-001	Working element used for pediatric urethrosopes, passive cutting , handle open	11 Fr.

8.3 | PEDIATRIC URETHROTOMY

COLD KNIVES



	Ordering no.	French	Color code	Description
	AI-004-007	11 Fr.	Green	Cold knife straight
	AI-004-009	11 Fr.	Green	Cold knife round
	AI-004-010	11 Fr.	Green	Cold knife ring shaped
	AI-004-011	11 Fr.	Green	Cold knife hook shaped



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

9. ACCESSORIES

9.1 | INSTRUMENTS

DILATATING BOUGIES
METAL CATHETERS
ACCESSORIES CATHETERS

9.2 | STERILIZATION


9.3 | CLEANING AND MAINTENANCE

9. ACCESSORIES & INSTRUMENTS

1 DILATATORS



2 STERILIZATION



3 PROTECTION



9.1 | INSTRUMENTS | DILATATOR BOUGIES



VAN BUREN

DITTEL,
STRAIGHTDITTEL,
CURVED

French	Ordering no. VAN BUREN	Ordering no. DITTEL straight	Ordering no. DITTEL curved
8 Fr.	AJ-003-001	AJ-003-021	AJ-003-050
9 Fr.		AJ-003-022	AJ-003-051
10 Fr.	AJ-003-002	AJ-003-023	AJ-003-052
11 Fr.		AJ-003-024	AJ-003-053
12 Fr.	AJ-003-003	AJ-003-025	AJ-003-054
13 Fr.		AJ-003-026	AJ-003-055
14 Fr.	AJ-003-004	AJ-003-027	AJ-003-056
15 Fr.		AJ-003-028	AJ-003-057
16 Fr.	AJ-003-005	AJ-003-029	AJ-003-058
17 Fr.		AJ-003-030	AJ-003-059
18 Fr.	AJ-003-006	AJ-003-031	AJ-003-060
19 Fr.		AJ-003-032	AJ-003-061
20 Fr.	AJ-003-007	AJ-003-033	AJ-003-062
21 Fr.		AJ-003-034	AJ-003-063
22 Fr.	AJ-003-008	AJ-003-035	AJ-003-064
23 Fr.		AJ-003-036	AJ-003-065
24 Fr.	AJ-003-009	AJ-003-037	AJ-003-066
25 Fr.		AJ-003-038	AJ-003-067
26 Fr.	AJ-003-010	AJ-003-039	AJ-003-068
27 Fr.		AJ-003-040	AJ-003-069
28 Fr.	AJ-003-011	AJ-003-041	AJ-003-070
29 Fr.		AJ-003-042	AJ-003-071
30 Fr.	AJ-003-012	AJ-003-043	AJ-003-072
31 Fr.			
32 Fr.	AJ-003-013	AJ-003-044	AJ-003-074
33 Fr.			
34 Fr.	AJ-003-014	AJ-003-045	
35 Fr.			
36 Fr.	AJ-003-015	AJ-003-046	
37 Fr.			
38 Fr.	AJ-003-016		

9.1 | INSTRUMENTS | DILATATOR BOUGIES | UTERINE DILATATORS



GUYON


French	Ordering no. GUYON	French	Ordering no. GUYON
8 Fr.	AJ-003-170	25 Fr.	AJ-003-187
9 Fr.	AJ-003-171	26 Fr.	AJ-003-188
10 Fr.	AJ-003-172	27 Fr.	AJ-003-189
11 Fr.	AJ-003-173	28 Fr.	AJ-003-190
12 Fr.	AJ-003-174	29 Fr.	AJ-003-191
13 Fr.	AJ-003-175	30 Fr.	AJ-003-192
14 Fr.	AJ-003-176	31 Fr.	AJ-003-193
15 Fr.	AJ-003-177	32 Fr.	AJ-003-194
16 Fr.	AJ-003-178	33 Fr.	AJ-003-195
17 Fr.	AJ-003-179	34 Fr.	AJ-003-196
18 Fr.	AJ-003-180	35 Fr.	AJ-003-197
19 Fr.	AJ-003-181	36 Fr.	AJ-003-198
20 Fr.	AJ-003-182	37 Fr.	AJ-003-199
21 Fr.	AJ-003-183	38 Fr.	AJ-003-200
22 Fr.	AJ-003-184	39 Fr.	AJ-003-201
23 Fr.	AJ-003-185	40 Fr.	AJ-003-202
24 Fr.	AJ-003-186	50 Fr.	AJ-003-203



HEGAR

French	Ordering no. HEGAR
3/4 Fr.	AJ-003-151
5/6 Fr.	AJ-003-152
7/8 Fr.	AJ-003-153
9/10 Fr.	AJ-003-154
11/12 Fr.	AJ-003-155
13/14 Fr.	AJ-003-156
15/16 Fr.	AJ-003-157
17/18 Fr.	AJ-003-158

9.1 | INSTRUMENTS | METAL CATHETERS



FOR MEN

FOR WOMEN

COXETER
PROSTATE

French	Ordering no.	Ordering no.	Ordering no.
6 Fr.		AJ-003-100	
7 Fr.	AJ-003-081		
8 Fr.	AJ-003-082	AJ-003-102	AJ-003-120
9 Fr.			
10 Fr.	AJ-003-084	AJ-003-104	AJ-003-121
11 Fr.			
12 Fr.	AJ-003-086	AJ-003-106	AJ-003-122
13 Fr.			
14 Fr.	AJ-003-088	AJ-003-108	AJ-003-123
15 Fr.			
16 Fr.	AJ-003-090	AJ-003-110	AJ-003-124
17 Fr.			
18 Fr.	AJ-003-092	AJ-003-112	AJ-003-125
19 Fr.			
20 Fr.	AJ-003-094	AJ-003-114	AJ-003-126
21 Fr.			
22 Fr.	AJ-003-096	AJ-003-116	AJ-003-127
23 Fr.			
24 Fr.	AJ-003-098	AJ-003-118	AJ-003-128

9.1 | INSTRUMENTS

CATHETER INTRODUCER



	Ordering no.	Length	Description
	AJ-003-142	39cm	Catheter introducer Guyon straight
	AJ-003-143	39cm	Catheter introducer Guyon curved

MEATUS CLAMPS

	Ordering no.	Length	Description
	AJ-003-140	8cm	Meatus clamp, Stockmann
	AJ-003-141	11cm	Meatus clamp, Strauss

SECRETION SPOON




	Ordering no.	Length	Description
	AJ-006-012	26.5cm	Secretion spoon Guttman


9.1 | INSTRUMENTS | TUBING CONNECTORS

	Ordering no.	Description
	AJ-010-011	Luer-Lock tubing connector, male Luer-Lock inside Ø 3.0mm, tubing connector Ø 6.0mm
	AJ-010-012	Luer-Lock tubing connector with stopcock
	AJ-010-013	Luer-Lock tubing connector, male Luer-Lock inside , tubing connector Ø 9.0mm
	AJ-010-018	Stopcock left with Luer-Lock, male, tubing connector Ø 9.0mm

9.2 | STERILIZATION | WIRE BASKET




Ordering no.	Dimension (l x w x h)	Description
FA-010-001	290x80x52mm	Wire basket for 1 endoscope with lid
FA-010-002	460x80x52mm	
FA-010-003	670x80x52mm	
FA-010-004	290x80x52mm	Wire basket for 2 endoscopes with lid
FA-010-005	460x80x52mm	



Ordering no.	Dimension (l x w x h)	Description
FA-010-006	640x150x77mm	Wire basket for 1 endoscope and 1 cold light cable with lid

9.2 | STERILIZATION | ENDOSCOPE CONTAINER



Ordering no.	Colour	Dimension (l x w x h)	Description
FA-001-001	Silver	450x70x70mm	Endoscope container, lid and bottom perforated (without filter)
FA-001-001-R	Red		
FA-001-001-G	Green		
FA-002-001	Silver	250x60x50mm	Endoscope container, lid and bottom perforated (without filter)
FA-002-001-R	Red		
FA-002-001-G	Green		



Ordering no.	Description
FA-001-000	Protection tube for endoscopes, plastic, 10 pcs, 19cm

9.3 | CLEANING AND MAINTENANCE

CLEANING PISTOL SET




Ordering no.	Description
EJ-030-010	Cleaning pistol complete inclusive 8 different tips

CLEANING BRUSH



Ordering no.	Diameter	Description
EJ-030-001	Ø 3mm	Cleaning brush, length 420mm , 1 box = 10 pieces
EJ-030-002	Ø 5mm	Cleaning brush, length 330mm , 1 box = 10 pieces
EJ-030-003	Ø 10mm	

GREASE & PASTE



Ordering no.	Weight	Description
EJ-030-031	3g	Special grease for stopcocks, valves and instruments
EJ-030-032	3g	Polishing paste for scopes, 1 box = 50 pieces

ARE YOU INTERESTED
IN OUR PRODUCTS?
WE WILL GLADLY
SUBMIT YOU
AN OFFER.

(if available)

Company Name:

Do you have a question about our products?

10. INDEX

AA-002-001	8 , 37 , 98 , 146	AA-006-083	109 , 114 , 116 , 151	AA-008-002	31
AA-002-001-W	8 , 37	AA-006-084	153	AA-008-003	142
AA-002-002	8 , 37	AA-006-100	110	AA-008-100	32 , 143
AA-002-003	8 , 28 , 37 , 59 , 146	AA-007-001	127	AA-008-101	143
AA-002-003-W	8 , 37	AA-007-002	127	AA-008-102	32 , 143
AA-002-004	8 , 37	AA-007-003	123	AA-008-103	32 , 143
AA-002-004-W	8 , 37	AA-007-004	123	AA-008-103	31 , 142
AA-002-011	8 , 37 , 98 , 146	AA-007-005	123	AA-008-104	32 , 143
AA-002-011	106	AA-007-006	126	AA-008-105	32 , 143
AA-002-011-W	8 , 37	AA-007-007	126	AA-008-105	31 , 142
AA-002-012	8 , 37 , 59	AA-007-008	126	AB-001-000	25
AA-002-013	8 , 28 , 37 , 59 , 76 , 146	AA-007-009	124	AB-001-001	157
AA-002-013	33 , 90 , 91 , 93	AA-007-011	127	AB-001-001-7	160
AA-002-013-W	8 , 37	AA-007-012	127	AB-001-001-9	159
AA-002-014	8 , 37	AA-007-013	125	AB-001-001-11	159
AA-002-014-W	8 , 37	AA-007-014	125	AB-001-001-13	157
AA-002-021	8 , 37 , 98	AA-007-015	125	AB-001-002	9
AA-002-022	8 , 37	AA-007-016	123	AB-001-003	9
AA-002-023	8 , 28 , 37 , 59	AA-007-017	123	AB-001-004	9
AA-002-024	8 , 37	AA-007-018	123	AB-001-004	33
AA-003-001	103	AA-007-019	123	AB-001-005	9
AA-003-002	25 , 48 , 65 , 70	AA-007-021	127	AB-001-006	9
AA-003-004	65	AA-007-022	126	AB-001-006	33
AA-003-011	103	AA-007-023	126	AB-001-007	25
AA-003-012	25 , 48 , 65 , 70 , 79 , 83	AA-007-024	126	AB-001-010	10
AA-003-012-W	25	AA-007-025	126	AB-001-011	157
AA-003-014	65	AA-007-026	126	AB-001-011-7	160
AA-004-011	159 , 166 , 174	AA-007-027	126	AB-001-011-9	159
AA-004-012	159	AA-007-028	124	AB-001-011-11	159
AA-004-013	160	AA-007-029	124	AB-001-012	10
AA-004-014	160	AA-007-030	125	AB-001-013	10
AA-005-001	157 , 163 , 171	AA-007-031	125	AB-001-014	10
AA-005-003	157	AA-007-032	125	AB-001-015	10
AA-005-004	157	AA-007-033	125	AB-001-016	10
AA-006-001	114	AA-007-034	125	AB-001-022	10
AA-006-002	114	AA-007-035	125	AB-001-023	10
AA-006-004	114	AA-007-036	123	AB-001-024	10
AA-006-005	114	AA-007-037	123	AB-001-025	10
AA-006-006	116	AA-007-041	129	AB-001-026	10
AA-006-011	116	AA-007-042	129	AB-001-040	28
AA-006-012	116 , 151	AA-007-043	128	AB-001-041	28
AA-006-021	114	AA-007-044	128	AB-001-042	29
AA-006-022	114	AA-007-045	129	AB-001-043	29
AA-006-031	116	AA-007-046	129	AB-002-001	9
AA-006-032	116	AA-007-047	128	AB-002-002	9
AA-006-034	116	AA-007-048	128	AB-002-003	9
AA-006-050	111	AA-007-049	129	AB-002-004	9
AA-006-051	111	AA-007-050	129	AB-002-005	9
AA-006-052	111	AA-007-051	128	AB-002-010	146
AA-006-053	111	AA-007-052	128	AB-002-011	147
AA-006-069	109	AA-007-053	129	AB-002-012	147
AA-006-070	109	AA-007-054	129	AB-002-013	146
AA-006-072	153	AA-007-055	128	AB-002-014	147
AA-006-080	114 , 116 , 151	AA-007-056	128	AB-003-001	158
AA-006-081	109 , 114 , 116 , 151	AA-007-060	30	AB-003-002	11

10. INDEX

AB-003-003	11 , 33	AB-006-048	119	AC-001-005-R	15
AB-003-004	11	AB-006-049	119	AC-001-008	15
AB-003-005	158	AB-006-050	119	AC-001-008-R	15
AB-003-006	12	AB-006-051	119	AC-001-010	130
AB-003-007	12 , 33	AB-006-052	120	AC-001-010-R	130
AB-003-008	12	AB-006-053	120	AC-001-011	130
AB-003-009	12	AB-006-061	117	AC-001-011-R	130
AB-003-010	26	AB-006-062	117	AC-001-012	130
AB-003-011	26	AB-006-063	117	AC-001-012-R	130
AB-003-012	26	AB-006-064	117	AC-001-013	130
AB-003-013	27	AB-006-065	117	AC-001-013-R	130
AB-003-014	27	AB-006-066	117	AC-001-014	130
AB-003-015	11	AB-006-067	117	AC-001-014-R	130
AB-003-020	27	AB-006-068	117	AC-001-015	130
AB-004-001	13 , 34 , 92 , 94	AB-006-100	112	AC-001-015-R	130
AB-004-002	13	AB-006-101	112	AC-001-021	161
AB-004-003	13 , 34 , 92 , 94	AB-006-102	112	AC-001-201	15
AB-004-004	13	AB-006-110	112	AC-001-201-OT	19
AB-004-005	13	AB-006-111	112	AC-001-201-T	19
AB-004-006	13	AB-006-112	112	AC-001-202	15
AB-004-007	14 , 34	AB-006-115	112	AC-001-202-OT	19
AB-004-008	13 , 92 , 94	AB-006-170	109	AC-001-202-T	19
AB-004-009	14	AB-006-171	109	AC-001-203	15
AB-004-009-QL	14	AB-006-242	113	AC-001-203-OT	19
AB-004-010	13	AB-006-243	113	AC-001-203-T	19
AB-004-010-QL	151	AB-006-244	113	AC-001-204	15
AB-004-011	148	AB-006-245	113	AC-001-204-OT	19
AB-004-012	148	AB-006-246	113	AC-001-204-T	19
AB-006-001	115	AB-006-342	113	AC-001-205	15
AB-006-002	115	AB-006-343	113	AC-001-206-OT	19
AB-006-003	115	AB-006-344	113	AC-001-206-T	19
AB-006-004	115	AB-006-345	113	AC-001-208	15
AB-006-005	115	AB-006-346	113	AC-001-210	130
AB-006-006	115	AB-006-442	113	AC-001-210-OT	135
AB-006-007	115	AB-006-443	113	AC-001-210-T	135
AB-006-008	115	AB-006-444	113	AC-001-211	130
AB-006-011	118	AB-006-445	113	AC-001-211-OT	135
AB-006-012	118	AB-006-446	113	AC-001-211-T	135
AB-006-013	118	AB-007-001	87	AC-001-212	130
AB-006-020	118	AB-007-002	87	AC-001-212-OT	135
AB-006-021	118	AB-007-003	87	AC-001-212-T	135
AB-006-022	118	AB-007-004	88	AC-001-213	130
AB-006-023	118	AB-007-005	88	AC-001-214	130
AB-006-024	118	AB-007-006	88	AC-001-214-OT	135
AB-006-025	118	AB-007-007	88	AC-001-214-T	135
AB-006-026	118	AB-007-008	87	AC-001-215	130
AB-006-027	118	AC-001-001	15	AC-001-215-OT	135
AB-006-040	119	AC-001-001-R	15	AC-001-215-T	135
AB-006-041	119	AC-001-002	15 , 33	AC-001-221	161
AB-006-042	119	AC-001-002-R	15	AC-002-001	15
AB-006-043	119	AC-001-003	15	AC-002-001-R	15
AB-006-044	119	AC-001-003-R	15	AC-002-002	15
AB-006-045	119	AC-001-004	15	AC-002-002-R	15
AB-006-046	119	AC-001-004-R	15	AC-002-003	15
AB-006-047	119	AC-001-005	15	AC-002-003-R	15

10. INDEX

AC-002-004	15	AC-003-202	17	AC-004-624-OT	134
AC-002-004-R	15	AC-003-202-OT	19	AC-004-624-T	134
AC-002-005	15	AC-003-202-T	19	AC-005-004	17
AC-002-005-R	15	AC-003-203	17	AC-005-004-R	17
AC-002-011	130	AC-003-204	17	AC-005-014	132
AC-002-011-R	130	AC-003-204-OT	19	AC-005-014-R	132
AC-002-012	130	AC-003-204-T	19	AC-005-204	17
AC-002-012-R	130	AC-003-205	17	AC-005-214	132
AC-002-013	130	AC-003-205-OT	19	AC-006-001	15
AC-002-013-R	130	AC-003-205-T	19	AC-006-001-R	15
AC-002-014	130	AC-003-210	132	AC-006-002	15
AC-002-014-R	130	AC-003-210-OT	134	AC-006-002-R	15
AC-002-015	130	AC-003-210-T	134	AC-006-003	15
AC-002-015-R	130	AC-003-211	132	AC-006-003-R	15
AC-002-201	15	AC-003-211-OT	134	AC-006-004	15 , 33
AC-002-202	15	AC-003-211-T	134	AC-006-004-R	15
AC-002-203	15	AC-003-212	132	AC-006-005	15
AC-002-204	15	AC-003-212-OT	134	AC-006-005-R	15
AC-002-205	15	AC-003-212-T	134	AC-006-011	130
AC-002-211	130	AC-003-213	132	AC-006-011-R	130
AC-002-212	130	AC-003-214	132	AC-006-012	130
AC-002-213	130	AC-003-214-OT	134	AC-006-012-R	130
AC-002-214	130	AC-003-214-T	134	AC-006-013	130
AC-002-215	130	AC-003-215	132	AC-006-013-R	130
AC-003-000	17	AC-003-215-OT	134	AC-006-014	130
AC-003-000-R	17	AC-003-215-T	134	AC-006-014-R	130
AC-003-001	17	AC-003-221	161	AC-006-015	130
AC-003-001-R	17	AC-003-302	17	AC-006-015-R	130
AC-003-002	17 , 33	AC-003-304	17	AC-006-021	161
AC-003-002-R	17	AC-003-612-OT	19	AC-006-201	15
AC-003-003	17	AC-003-612-T	19	AC-006-201-OT	19
AC-003-003-R	17	AC-003-614-OT	19	AC-006-201-T	19
AC-003-004	17	AC-003-614-T	19	AC-006-202	15
AC-003-004-R	17	AC-004-004	17	AC-006-202-OT	19
AC-003-005	17	AC-004-004-R	17	AC-006-202-T	19
AC-003-005-R	17	AC-004-005	17	AC-006-203	15
AC-003-010	132	AC-004-005-R	17	AC-006-203-OT	19
AC-003-010-R	132	AC-004-014	132	AC-006-203-T	19
AC-003-011	132	AC-004-014-R	132	AC-006-204	15
AC-003-011-R	132	AC-004-015	132	AC-006-204-OT	19
AC-003-012	132	AC-004-015-R	132	AC-006-204-T	19
AC-003-012-R	132	AC-004-204	17	AC-006-205	15
AC-003-013	132	AC-004-204-OT	19	AC-006-205-OT	19
AC-003-013-R	132	AC-004-204-T	19	AC-006-205-T	19
AC-003-014	132	AC-004-205	17	AC-006-211	130
AC-003-014-R	132	AC-004-205-OT	19	AC-006-211-OT	135
AC-003-015	132	AC-004-205-T	19	AC-006-211-T	135
AC-003-015-R	132	AC-004-214	132	AC-006-212	130
AC-003-021	161	AC-004-214-OT	134	AC-006-212-OT	135
AC-003-200	17	AC-004-214-T	134	AC-006-212-T	135
AC-003-200-OT	19	AC-004-215	132	AC-006-213	130
AC-003-200-T	19	AC-004-215-OT	134	AC-006-213-OT	135
AC-003-201	17	AC-004-215-T	134	AC-006-213-T	135
AC-003-201-OT	19	AC-004-622-OT	134	AC-006-214	130
AC-003-201-T	19	AC-004-622-T	134	AC-006-214-OT	135

10. INDEX

AC-006-214-T	135	AC-009-204	16	AC-017-001	15
AC-006-215	130	AC-009-211	131	AC-017-001-R	15
AC-006-215-OT	135	AC-009-211-OT	135	AC-017-004	15
AC-006-215-T	135	AC-009-211-T	135	AC-017-004-R	15
AC-006-221	161	AC-009-212	131	AC-017-011	130
AC-007-001	16	AC-009-212-OT	135	AC-017-011-R	130
AC-007-001-R	16	AC-009-212-T	135	AC-017-014	130
AC-007-002	16	AC-009-214	131	AC-017-014-R	130
AC-007-002-R	16	AC-010-004	15	AC-017-201	15
AC-007-011	131	AC-010-014	130	AC-017-204	15
AC-007-011-R	131	AC-010-204	15	AC-017-211	130
AC-007-012	131	AC-010-214	130	AC-017-214	130
AC-007-012-R	131	AC-011-004	15	AC-018-001	17
AC-007-201	16	AC-011-204	15	AC-018-001-R	17
AC-007-201-OT	19	AC-011-214	130	AC-018-002	17
AC-007-201-T	19	AC-012-001	16	AC-018-002-R	17
AC-007-202	16	AC-012-002	16	AC-018-011	132
AC-007-202-OT	19	AC-012-004	16	AC-018-011-R	132
AC-007-202-T	19	AC-012-011	131	AC-018-012	132
AC-007-211	131	AC-012-012	131	AC-018-012-R	132
AC-007-211-OT	135	AC-012-014	131	AC-018-201	17
AC-007-211-T	135	AC-015-002	15	AC-018-202	17
AC-007-212	131	AC-015-002-R	15	AC-018-211	132
AC-007-212-OT	135	AC-015-003-R	15	AC-018-212	132
AC-007-212-T	135	AC-015-004	15	AC-019-001	16
AC-008-004	16	AC-015-004-R	15	AC-019-001-R	16
AC-008-004-R	16	AC-015-012	130	AC-019-002	16
AC-008-005	16	AC-015-012-R	130	AC-019-002-R	16
AC-008-005-R	16	AC-015-013-R	130	AC-019-011	131
AC-008-014	131	AC-015-014	130	AC-019-011-R	131
AC-008-014-R	131	AC-015-014-R	130	AC-019-012	131
AC-008-015	131	AC-015-021	161	AC-019-012-R	131
AC-008-015-R	131	AC-015-202	15	AC-019-021	161
AC-008-204	16	AC-015-203	15	AC-019-201	16
AC-008-205	16	AC-015-204	15	AC-019-202	16
AC-008-214	131	AC-015-212	130	AC-019-211	131
AC-008-215	131	AC-015-213	130	AC-019-212	131
AC-009-001	16	AC-015-214	130	AC-020-002	16
AC-009-001-R	16	AC-015-221	161	AC-020-002-R	16
AC-009-002	16	AC-016-001	16	AC-020-004	16
AC-009-002-R	16	AC-016-001-R	16	AC-020-004-R	16
AC-009-004	16	AC-016-002	16	AC-020-012	131
AC-009-004-R	16	AC-016-002-R	16	AC-020-012-R	131
AC-009-011	131	AC-016-011	131	AC-020-014	131
AC-009-011-R	131	AC-016-011-R	131	AC-020-014-R	131
AC-009-012	131	AC-016-012	131	AC-020-202	16
AC-009-012-R	131	AC-016-012-R	131	AC-020-204	16
AC-009-014	131	AC-016-201	16	AC-020-212	131
AC-009-014-R	131	AC-016-202	16	AC-020-214	131
AC-009-201	16	AC-016-211	131	AC-100-001	19 , 20 , 134 , 135 , 136 , 137
AC-009-201-OT	19	AC-016-211-OT	135	AC-101-201-OT	20
AC-009-201-T	19	AC-016-211-T	135	AC-101-201-T	20
AC-009-202	16	AC-016-212	131	AC-101-202-OT	20
AC-009-202-OT	19	AC-016-212-OT	135	AC-101-202-T	20
AC-009-202-T	19	AC-016-212-T	135		

10. INDEX

AC-101-203-OT	20	AD-002-002	23 , 140
AC-101-203-T	20	AD-002-003	23 , 34 , 140
AC-101-204-OT	20	AD-002-005	23 , 140
AC-101-204-T	20	AD-002-011	23 , 140
AC-101-206-OT	20	AD-002-012	23 , 140
AC-101-206-T	20	AD-002-013	23 , 140
AC-101-210-OT	137	AD-002-015	23 , 140
AC-101-210-T	137	AD-002-021	161
AC-101-211-OT	137	AD-003-003	23 , 140
AC-101-211-T	137	AD-003-005	23 , 140
AC-101-212-OT	137	AD-003-007	23 , 140
AC-101-212-T	137	AD-003-008	23 , 140
AC-101-214-OT	137	AD-003-013	23 , 140
AC-101-214-T	137	AD-003-015	23 , 140
AC-101-215-OT	137	AD-003-017	23 , 140
AC-101-215-T	137	AD-003-018	23 , 140
AC-103-200-OT	20	AD-010-001	24 , 57 , 92 , 141 , 162 , 169
AC-103-200-T	20	AD-010-002	24 , 57 , 141 , 162 , 169
AC-103-201-OT	20	AD-010-003	24 , 34 , 57 , 141 , 162 , 169
AC-103-201-T	20	AD-010-006	24 , 57 , 141 , 162 , 169
AC-103-202-OT	20	AD-010-007	24 , 57 , 141 , 162 , 169
AC-103-202-T	20	AD-010-008	24 , 57 , 141 , 162 , 169
AC-103-204-OT	20	AD-010-011	75
AC-103-204-T	20	AD-010-012	75
AC-103-205-OT	20	AD-010-013	75 , 94
AC-103-205-T	20	AD-010-015	75
AC-103-210-OT	136	AD-010-024	24 , 57 , 141 , 162 , 169
AC-103-210-T	136	AD-010-025	58 , 170
AC-103-211-OT	136	AD-010-026	58 , 170
AC-103-211-T	136	AE-001-001	46 , 91
AC-103-212-OT	136	AE-001-002	46
AC-103-212-T	136	AE-001-003	46
AC-103-214-OT	136	AE-001-004	46 , 91
AC-103-214-T	136	AE-001-005	46
AC-103-215-OT	136	AE-001-006	47
AC-103-215-T	136	AE-001-007	47
AC-103-612-OT	20	AE-001-008	47
AC-103-612-T	20	AE-001-009	47
AC-103-614-OT	20	AE-001-010	46
AC-103-614-T	20	AE-001-011	46
AC-104-204-OT	20	AE-001-012	46
AC-104-204-T	20	AE-001-013	46
AC-104-205-OT	20	AE-001-014	47
AC-104-205-T	20	AE-001-015	47
AC-104-214-OT	136	AE-001-016	47 , 91
AC-104-214-T	136	AE-001-017	47
AC-104-215-OT	136	AE-001-030	63
AC-104-215-T	136	AE-001-031	63 , 90 , 93
AC-104-622-OT	136	AE-001-032	63 , 93
AC-104-622-T	136	AE-001-033	63
AC-104-624-OT	136	AE-001-034	64 , 90 , 93
AC-104-624-T	136	AE-001-035	64
AC-106-201-OT	20	AE-001-036	63
AC-106-201-T	20		
AC-106-202-OT	20		
AC-106-202-T	20		
AC-106-203-OT	20		
AC-106-203-T	20		
AC-106-204-OT	20		
AC-106-204-T	20		
AC-106-205-OT	20		
AC-106-205-T	20		
AC-106-211-OT	137		
AC-106-211-T	137		
AC-106-212-OT	137		
AC-106-212-T	137		
AC-106-213-OT	137		
AC-106-213-T	137		
AC-106-214-OT	137		
AC-106-214-T	137		
AC-106-215-OT	137		
AC-106-215-T	137		
AC-107-201-OT	20		
AC-107-201-T	20		
AC-107-202-OT	20		
AC-107-202-T	20		
AC-107-211-OT	137		
AC-107-211-T	137		
AC-107-212-OT	137		
AC-107-212-T	137		
AC-109-201-OT	20		
AC-109-201-T	20		
AC-109-211-OT	137		
AC-109-211-T	137		
AC-109-212-OT	137		
AC-109-212-T	137		
AC-114-212-OT	136		
AC-114-212-T	136		
AC-116-211-OT	137		
AC-116-211-T	137		
AC-116-212-OT	137		
AC-116-212-T	137		
AD-001-001	23 , 140		
AD-001-002	23 , 140		
AD-001-003	23 , 34 , 140		
AD-001-004	23 , 140		
AD-001-005	23 , 140		
AD-001-006	23 , 140		
AD-001-007	23 , 140		
AD-001-008	23 , 140		
AD-001-011	23 , 140		
AD-001-012	23 , 140		
AD-001-013	23 , 140		
AD-001-014	23 , 140		
AD-001-015	23 , 140		
AD-001-016	23 , 140		
AD-001-017	23 , 140		
AD-001-018	23 , 140		
AD-001-021	161		
AD-002-001	23 , 140		

10. INDEX

AE-001-037	63	AE-006-004	168	AG-002-013	71
AE-001-038	63	AE-006-016	168	AG-002-015	71
AE-001-039	64	AE-006-018	168	AG-002-016	84
AE-001-040	63	AE-006-019	168	AG-003-001	164 , 172
AE-001-041	63	AF-001-001	21 , 138	AG-004-001	167 , 175
AE-001-042	63	AF-001-002	21 , 138	AG-005-001	167
AE-001-043	63	AF-001-003	21 , 138	AG-100-006	77
AE-001-044	64	AF-001-004	21 , 138	AG-100-007	77
AE-001-045	64	AF-001-005	21 , 138	AH-001-001	44
AE-001-046	64	AF-001-006	21 , 138	AH-001-002	44
AE-001-047	64	AF-001-007	22 , 139	AH-001-011	45 , 62 , 78 , 90 , 91 , 93 , 152
AE-001-048	63	AF-001-008	22 , 139	AH-001-012	45 , 62
AE-002-001	46	AF-001-008	34	AH-001-013	45 , 62
AE-002-002	46	AF-001-009	22 , 139	AH-001-015	45 , 62
AE-002-003	46	AF-001-010	22 , 139	AH-001-016	45 , 62
AE-002-004	46	AF-001-011	22 , 139	AH-001-017	45 , 62
AE-002-005	46	AF-001-012	22 , 139	AH-001-018	45 , 62
AE-002-006	47	AF-001-013	22 , 139	AH-001-101	44
AE-002-007	47	AF-001-014	22 , 139	AH-001-102	44
AE-002-010	46	AF-001-015	21 , 138	AH-001-201	42
AE-002-011	46	AF-001-016	21 , 138	AH-001-202	42
AE-002-012	46	AF-001-017	21 , 138	AH-001-203	42
AE-002-013	46	AF-001-018	22 , 139	AH-001-204	42
AE-002-014	47	AF-001-019	21 , 138	AH-001-205	42
AE-002-015	47	AF-001-021	21 , 22 , 138 , 139	AH-001-211	42
AE-002-016	47	AF-001-022	21 , 138	AH-001-212	42
AE-002-031	64	AF-001-023	21 , 138	AH-001-301	43 , 61 , 78 , 90 , 91 , 93 , 152
AE-003-001	55	AF-001-024	21 , 22 , 138 , 139	AH-001-302	43 , 61 , 78 , 152
AE-003-003	55	AF-001-025	22 , 139	AH-001-302	90 , 91 , 93
AE-003-004	55	AG-001-001	38 , 101	AH-001-303	43 , 61
AE-003-010	55	AG-001-001	91 , 106	AH-001-304	43 , 61
AE-003-014	55	AG-001-001-HF	38 , 101	AH-001-305	43 , 61
AE-003-016	55	AG-001-002	38 , 101	AH-001-306	43 , 61
AE-003-031	69	AG-001-003	39	AH-001-308	43 , 78
AE-003-032	69	AG-001-003-HF	39	AH-001-311	43 , 61
AE-003-033	69	AG-001-004	39	AH-001-312	43 , 61
AE-003-034	69	AG-001-006	77	AH-002-001	54
AE-003-036	69	AG-001-007	40	AH-002-011	68 , 82
AE-003-040	69	AG-001-008	40	AH-002-011	56
AE-003-041	74	AG-001-009	77	AH-002-011	154
AE-003-043	74	AG-001-010	104	AH-002-013	68 , 82
AE-003-044	74	AG-001-011	60	AH-002-013	56 , 154
AE-003-045	74	AG-001-011	90 , 93	AH-002-015	73 , 86
AE-003-046	74	AG-001-012	60	AH-002-016	73 , 86
AE-004-002	165	AG-001-017	41	AH-002-201	52
AE-004-004	165	AG-001-018	41	AH-002-211	72 , 85
AE-004-016	165	AG-002-001	49	AH-002-212	52
AE-004-018	165	AG-002-002	49	AH-002-213	72 , 85
AE-004-019	165	AG-002-003	50	AH-002-216	72 , 85
AE-005-002	168	AG-002-004	50	AH-002-301	67 , 81
AE-005-004	168	AG-002-006	51	AH-002-301	53 , 154
AE-005-016	168	AG-002-007	80	AH-002-302	67 , 81
AE-005-018	168	AG-002-008	80	AH-002-302	53 , 154
AE-005-019	168	AG-002-011	66		
AE-006-002	168	AG-002-012	66		

10. INDEX

AH-003-001	163	AJ-003-002	179	AJ-003-064	179
AH-003-011	163	AJ-003-003	179	AJ-003-065	179
AH-004-001	166	AJ-003-004	179	AJ-003-066	179
AH-004-011	166	AJ-003-005	179	AJ-003-067	179
AH-005-001	166	AJ-003-006	179	AJ-003-068	179
AH-005-011	166	AJ-003-007	179	AJ-003-069	179
AI-001-001	98 , 106	AJ-003-008	179	AJ-003-070	179
AI-001-002	99 , 106	AJ-003-009	179	AJ-003-071	179
AI-001-003	98 , 106	AJ-003-010	179	AJ-003-072	179
AI-001-004	100 , 106	AJ-003-011	179	AJ-003-074	179
AI-001-005	100	AJ-003-012	179	AJ-003-081	181
AI-001-006	99	AJ-003-013	179	AJ-003-082	181
AI-001-007	102	AJ-003-014	179	AJ-003-084	181
AI-001-007	106	AJ-003-015	179	AJ-003-086	181
AI-001-008	102	AJ-003-016	179	AJ-003-088	181
AI-001-009	102	AJ-003-021	179	AJ-003-090	181
AI-001-011	102	AJ-003-022	179	AJ-003-092	181
AI-001-012	100	AJ-003-023	179	AJ-003-094	181
AI-002-001	97	AJ-003-024	179	AJ-003-096	181
AI-002-002	97	AJ-003-025	179	AJ-003-098	181
AI-003-001	171	AJ-003-026	179	AJ-003-100	181
AI-003-003	164 , 172	AJ-003-027	179	AJ-003-102	181
AI-003-004	171	AJ-003-028	179	AJ-003-104	181
AI-003-007	173	AJ-003-029	179	AJ-003-106	181
AI-003-009	173	AJ-003-030	179	AJ-003-108	181
AI-003-010	173	AJ-003-031	179	AJ-003-110	181
AI-003-011	173	AJ-003-032	179	AJ-003-112	181
AI-004-001	174	AJ-003-033	179	AJ-003-114	181
AI-004-004	174	AJ-003-034	179	AJ-003-116	181
AI-004-007	176	AJ-003-035	179	AJ-003-118	181
AI-004-009	176	AJ-003-036	179	AJ-003-120	181
AI-004-010	176	AJ-003-037	179	AJ-003-121	181
AI-004-011	176	AJ-003-038	179	AJ-003-122	181
AI-005-001	103	AJ-003-039	179	AJ-003-123	181
AI-005-003	103	AJ-003-040	179	AJ-003-124	181
AI-005-004	103	AJ-003-041	179	AJ-003-125	181
AI-005-007	105	AJ-003-042	179	AJ-003-126	181
AI-005-008	105	AJ-003-043	179	AJ-003-127	181
AI-005-011	105	AJ-003-044	179	AJ-003-128	181
AJ-001-000	89	AJ-003-045	179	AJ-003-140	182
AJ-001-001	89	AJ-003-046	179	AJ-003-141	182
AJ-001-002	89	AJ-003-050	179	AJ-003-142	182
AJ-001-002	92 , 94	AJ-003-051	179	AJ-003-143	181 , 182
AJ-001-003	89	AJ-003-052	179	AJ-003-151	180
AJ-001-010	89	AJ-003-053	179	AJ-003-152	180
AJ-001-011	89	AJ-003-054	179	AJ-003-153	180
AJ-001-012	89	AJ-003-055	179	AJ-003-154	180
AJ-001-013	89	AJ-003-056	179	AJ-003-155	180
AJ-002-001	89	AJ-003-057	179	AJ-003-156	180
AJ-002-001	92 , 94	AJ-003-058	179	AJ-003-157	180
AJ-002-002	89	AJ-003-059	179	AJ-003-158	180
AJ-002-003	89	AJ-003-060	179	AJ-003-170	180
AJ-002-004	89	AJ-003-061	179	AJ-003-171	180
AJ-002-005	89	AJ-003-062	179	AJ-003-172	180
AJ-003-001	179	AJ-003-063	179	AJ-003-173	180

10. INDEX

AJ-003-174	180	AJ-010-003	26 , 111 , 127
AJ-003-175	180	AJ-010-006	47 , 55 , 64 , 69 , 74 , 90 , 92 , 94 , 102 , 105 , 106 , 165 , 168 , 173 , 176
AJ-003-176	180	AJ-010-011	183
AJ-003-177	180	AJ-010-012	183
AJ-003-178	180	AJ-010-013	183
AJ-003-179	180	AJ-010-015	114 , 116
AJ-003-180	180	AJ-010-018	183
AJ-003-181	180	AJ-010-020	90
AJ-003-182	180	CI-001-030	58 , 170
AJ-003-183	180	CI-001-040	58 , 170
AJ-003-184	180	EJ-030-001	186
AJ-003-185	180	EJ-030-002	186
AJ-003-186	180	EJ-030-003	186
AJ-003-187	180	EJ-030-010	186
AJ-003-188	180	EJ-030-031	186
AJ-003-189	180	EJ-030-032	186
AJ-003-190	180	EK-015-060	24 , 57 , 141 , 162 , 169
AJ-003-191	180	FA-001-000	185
AJ-003-192	180	FA-001-001	185
AJ-003-193	180	FA-001-001-G	185
AJ-003-194	180	FA-001-001-R	185
AJ-003-195	180	FA-002-001	185
AJ-003-196	180	FA-002-001-G	185
AJ-003-197	180	FA-002-001-R	185
AJ-003-198	180	FA-010-001	184
AJ-003-199	180	FA-010-002	184
AJ-003-200	180	FA-010-003	184
AJ-003-201	180	FA-010-004	184
AJ-003-202	180	FA-010-005	184
AJ-003-203	180	FA-010-006	184
AJ-006-012	182		
AJ-010-001	11 , 26 , 111 , 127		

WE COVER THE FOLLOWING AREAS IN SURGERY WITH HIGH QUALITY PRODUCTS

ARTHROSCOPY

CARDIAC MIS

CONTAINER

ENT

GASTROSCOPY

GYNAECOLOGY

HF-SURGERY

LAMINECTOMY

LAPAROSCOPY

PROCTOSCOPY

SPINE SURGERY

TITANIUM PLATING

UNITS

UROLOGY



STEMA MEDIZINTECHNIK GMBH | TAKE-OFF GEWERBEPARK 74 | 78579 NEUHAUSEN OB ECK
INFO@STEMA-MEDIZINTECHNIK.DE | WWW.STEMA-MEDIZINTECHNIK.DE | +49 7467 79 630-0

11/25 | PRINTED IN GERMANY | COPYRIGHT BY STEMA MEDIZINTECHNIK GMBH |
WE RESERVE THE RIGHT TO MAKE ALTERATIONS