





GYNAECOLOGY



ABOUT US



SINCE 20 YEARS



EXPORT IN
+70 COUNTRIES



24H/SERVICE



+ 10.000
PRODUCTS



LOCATED IN THE
WORLD CENTER
FOR MEDICAL
TECHNOLOGY

CERTIFICATED

STEMA Medizintechnik GmbH is certified in accordance with the highest standards:

- EN ISO 13485:2016 + AC:2018 + A11:2021- ISO 13485:2016
- EU 2017/745 (**MDR**)
- Council Directive 93/42/EEC Annex II excluding Section 4
- Council Directive 93/42/EEC Annex V excluding Section 3
- FDA

These standards create an essential base for our company and our partners. They also symbolize the commitment and strive to keep up the quality of our products.



PRODUCT RANGE

The current product range of STEMA Medizintechnik GmbH includes more than 10.000 different items.

We cover the following areas in surgery with high quality products:

ARTHROSCOPY

CARDIAC MIS

CONTAINER-SYSTEMS

GASTROENTEROLOGY

GYNAECOLOGY

HF-SURGERY

LAMINECTOMY

LAPAROSCOPY

MICRO MOTOR SYSTEMS

PROCTOSCOPY

SPINE-SURGERY

TITANIUM PLATING SYSTEM

UNITS

UROLOGY



CONTENT

GYNAECOLOGY

PRODUCT	PAGE
1. HYSTEROSCOPY 4.0MM	
1.1 Cystoscopes 4.0mm	8
1.2 Diagnostic Hysteroscopy	9
1.3 Operating Hysteroscopy	11
2. SLIM HYSTEROSCOPY 2.9MM	
2.1 Cystoscopes 2.9mm	17
2.2 Diagnostic Hysteroscopy	18
2.3 Operating Hysteroscopy	20
2.4 High-Flow Operating sheaths	24
3. MINI HYSTEROSCOPY 2.0MM	
3.1 Cystoscopes 2.0mm	27
3.2 Diagnostic Hysteroscopy	28
3.3 Operating Hysteroscopy	30
3.4 High-Flow Operating sheaths	31
4. INSTRUMENTS & ELECTRODES	
4.1 Forceps, Scissors, Instruments	34
4.2 Electrodes & Connecting cables	41
4.3 Accessories	46
5. SALINE RESECTOSCOPY	
5.1 Saline Resectoscopy	60
5.2 Saline Slim Resectoscopy	66
5.3 Saline Mini Resectoscopy	71
5.4 Connecting Cables	76

CONTENT

GYNAECOLOGY

PRODUCT

PAGE

6. HYSTEROSCOPY SETS

6.1 Diagnostic Hysteroscopy set	78
6.2 Operating Hysteroscopy set	79
6.3 Out-Patient Atraumatic Operating Hysteroscopy set	80
6.4 Saline Resectoscopy set	81

7. INDEX

7.1 Index	84
-------------	----

1. HYSTEROSCOPY 4.0MM

1.1 | CYSTOSCOPES 4.0MM


1.2 | DIAGNOSTIC HYSTEROSCOPY

1.3 | OPERATING HYSTEROSCOPY

1. HYSTEROSCOPY

COMPLETE DIAGNOSTIC HYSTEROSCOPE

4.0mm



1 CYTOSCOPES



2 SHEATHS



3 OBTURATORS



1.1 | HYSTEROSCOPY

Ø 4.0MM CYSTOSCOPES AUTOCLAVABLE



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-001
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-011

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-002
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-012

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 70°	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-014

1.2 | DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-003-001	Ø 5.1mm	277mm	Diagnostic hysteroscope sheath for Ø 4.0mm, 0° endoscope
AN-003-002	Ø 5.1mm	277mm	Diagnostic hysteroscope sheath for Ø 4.0mm, 30° endoscope



Ordering no.	Description
AN-003-003	Obturator for diagnostic hysteroscope sheath for Ø 4.0mm endoscope

1.2 | CONTINUOUS-FLOW DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR



Ordering no.	Diameter	WL	Description
AN-003-103	∅ 6.5mm	246mm	Outer-Diagnostic-Continuous-Flow hysteroscope sheath for 4.0mm endoscope

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-003-101	∅ 5.4mm	261mm	Inner-Diagnostic-Continuous-Flow hysteroscope sheath for 4.0mm/0° endoscope
AN-003-102	∅ 5.4mm	261mm	Inner-Diagnostic-Continuous-Flow hysteroscope sheath for 4.0mm/30° endoscope



Ordering no.	Description
AN-003-104	Obturator for Outer-Diagnostic-Continuous-Flow hysteroscope sheath for 4.0mm endoscope

1.3 | CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS

Shape of sheath




Ordering no.	Diameter	WL	Description
AN-003-201	∅ 8.0mm	193mm	Standard OP-Continuous-Flow hysteroscope sheath for 0°/12° endoscope and 7 Fr. flexible and semi-rigid instruments
AN-003-202	∅ 8.0mm	193mm	Standard OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 7 Fr. flexible and semi-rigid instruments

1.3 | CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS & OBTURATOR



Ordering no.	Diameter	WL	Description
AN-003-302	∅ 8.0mm	180mm	Outer-OP-Continuous-Flow hysteroscope sheath for 12°/30 endoscope and 7 Fr. flexible and semi-rigid instruments

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-003-301	∅ 8.0mm	180mm	Inner-OP-Continuous-Flow hysteroscope sheath for 12°/30 endoscope and 7 Fr. flexible and semi-rigid instruments



Ordering no.	Diameter	WL	Description
AN-003-303	∅ 8.0mm	180mm	Obturator for Outer-OP-Continuous-Flow hysteroscope sheath

1.3 | CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS



Ordering no.	Diameter	WL	Description
AN-005-302	∅ 6.5mm	192mm	Outer-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 5 Fr. flexible/semi-rigid instruments
AN-005-304	∅ 7.0mm	192mm	Outer-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 7 Fr. flexible/semi-rigid instruments

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-005-301	∅ 6.5mm	192mm	Inner-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 5 Fr. flexible/semi-rigid instruments
AN-005-303	∅ 7.0mm	192mm	Inner-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 7 Fr. flexible/semi-rigid instruments

1.3 | CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS & OBTURATOR



Ordering no.	Diameter	WL	Description
AN-004-401	∅ 8.0mm	180mm	Outer-OP-Continuous-Flow hysteroscope sheath for 7 Fr. flexible and semi-rigid 12°/30° endoscope

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-004-402	∅ 8.0mm	180mm	Working Element with Albarran Deflector



Ordering no.	Diameter	WL	Description
AN-004-403	∅ 8.0mm	180mm	Obturator for Outer-OP-Continuous-Flow hysteroscope sheath for 4.0mm endoscope

2. SLIM HYSTEROSCOPY 2.9MM

2.1 | CYSTOSCOPES 2.9MM


2.2 | DIAGNOSTIC HYSTEROSCOPY

2.3 | OPERATING HYSTEROSCOPY

2.4 | HIGH-FLOW OPERATING SHEATHS

2. SLIM HYSTEROSCOPY 2.9MM


2.9mm



1 CYTOSCOPIES



2 SHEATHS



3 OBTURATORS



2.1 | SLIM HYSTEROSCOPY

Ø 2.9MM CYSTOSCOPES AUTOCLAVABLE



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-001
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-011


Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-014

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.9mm

2.2 | SLIM DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-002-001	Ø 4.0mm	277mm	Diagnostic hysteroscope sheath for Ø 2.9mm, 0° endoscope
AN-002-002	Ø 4.0mm	277mm	Diagnostic hysteroscope sheath for Ø 2.9mm, 30° endoscope



Ordering no.	Description
AN-002-003	Obturator for diagnostic hysteroscope sheath for 2.9mm endoscope

2.2 | SLIM CONTINUOUS-FLOW DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR



Ordering no.	Diameter	WL	Colour	Description
AN-002-103	∅ 5mm	246mm	White	Outer-Diagnostic-Continuous-Flow hysteroscope sheath for 2.9mm endoscope

Shape of sheath




Ordering no.	Diameter	WL	Colour	Description
AN-002-101	∅ 3.7mm	261mm	White	Inner-Diagnostic-Continuous-Flow hysteroscope sheath for 2.9mm/0° endoscope
AN-002-102	∅ 3.7mm	261mm	White	Inner-Diagnostic-Continuous-Flow hysteroscope sheath for 2.9mm/30° endoscope



Ordering no.	Description
AN-002-104	Obturator for Outer-Diagnostic-Continuous-Flow hysteroscope sheath for 2.9mm endoscope

2.3 | SLIM CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS

Shape of sheath




Ordering no.	Diameter	WL	Description
AN-002-201	∅ 6.0mm	193mm	Standard OP-Continuous-Flow hysteroscope sheath 0°/12° endoscope and 5 Fr. flexible and semi-rigid instruments
AN-002-202	∅ 6.0mm	193mm	Standard OP-Continuous-Flow hysteroscope sheath 30° endoscope and 7 Fr. flexible and semi-rigid instruments

2.3 | SLIM CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS & OBTURATOR



Ordering no.	Diameter	WL	Description
AN-002-302	Ø 5.0mm	192mm	Outer-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 5 Fr. flexible/semi-rigid instruments, white
AN-002-308	Ø 6.3mm	192mm	Outer-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 9 Fr. flexible/semi-rigid instruments, yellow

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-002-301	Ø 5.0mm	192mm	Inner-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 5 Fr. flexible/semi-rigid instruments, white
AN-002-307	Ø 6.3mm	192mm	Inner-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 9 Fr. flexible/semi-rigid instruments, yellow


2.9mm

2.3 | SLIM CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS



Ordering no.	Diameter	WL	Description
AN-002-312	∅ 5.0mm	216mm	Outer-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 5 Fr. flexible/semi-rigid instruments

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-002-313	∅ 5.0mm	216mm	Inner-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 5 Fr. flexible/semi-rigid instruments

2.9mm

2.3 | SLIM CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS



Ordering no.	Diameter	WL	Description
AN-002-315	∅ 5.0mm	216mm	Outer-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 7 Fr. flexible/semi-rigid instruments

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-002-314	∅ 5.0mm	216mm	Inner-OP-Continuous-Flow hysteroscope sheath for 30° endoscope and 7 Fr. flexible/semi-rigid instruments

2.9mm


2.4 | SLIM HIGH-FLOW CONTINUOUS-FLOW OPERATING SHEATHS

Shape of sheath



Ordering no.	Diameter	Description
AN-002-310	Ø 5.5mm	High-Flow Continuous-Flow Operating hysteroscope sheath for 5 Fr. flexible and semi-rigid instruments for use with 2.9mm endoscope (12°/30°) , single sheath system

Shape of sheath



Ordering no.	Diameter	Description
AN-002-311	Ø 5.5mm	High-Flow Continuous-Flow Operating hysteroscope sheath, for 5 Fr. flexible and semi-rigid instruments for use with 2.9mm endoscope (12°/30°) , single sheath system with stopcock

3. MINI HYSTEROSCOPY 2.0MM

3.1 | CYSTOSCOPES 2.0MM


3.2 | DIAGNOSTIC SHEATHS & OBTURATORS

3.3 | OPERATING SHEATHS & OBTURATORS

3.4 | HIGH-FLOW OPERATING SHEATHS

3. MINI HYSTEROSCOPY 2.0MM

2.0mm



1 CYTOSCOPES



2 SHEATHS



3 OBTURATORS



3.1 | MINI HYSTEROSCOPY Ø 2.0MM CYSTOSCOPES AUTOCLAVABLE



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 2.0mm	260mm	AA-004-001

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	Cystoscope autoclavable	Ø 2.0mm	260mm	AA-004-003

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.0mm	260mm	AA-004-002

3.2 | MINI DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR

Shape of sheath



Ordering no.	Diameter	WL	Description
AN-001-001	Ø 2.5mm	277mm	Diagnostic hysteroscope sheath for Ø 2.0mm, 0° endoscope
AN-001-002	Ø 2.5mm	277mm	Diagnostic hysteroscope sheath for Ø 2.0mm, 30° endoscope



Ordering no.	Description
AN-001-003	Obturator for diagnostic hysteroscope sheath for 2.0mm endoscope



3.2 | MINI CONTINUOUS-FLOW DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR



Ordering no.	Diameter	WL	Colour	Description
AN-001-103	∅ 3.4mm	204mm	Red	Outer-Diagnostic-Continuous-Flow hysteroscope sheath for ∅ 2.0mm endoscope

Shape of sheath



Ordering no.	Diameter	WL	Colour	Description
AN-001-101	∅ 2.6mm	219mm	Red	Inner-Diagnostic-Continuous-Flow hysteroscope sheath for ∅ 2.0mm, 0° endoscope
AN-001-102	∅ 2.6mm	219mm	Red	Inner-Diagnostic-Continuous-Flow hysteroscope sheath for ∅ 2.0mm, 30° endoscope




Ordering no.	Description
AN-001-104	Obturator for Outer-Diagnostic-Continuous-Flow hysteroscope sheath for 2.0mm endoscope

3.3 | MINI CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS



Ordering no.	Diameter	WL	Description
AN-001-302	∅ 4.3mm	181mm	Outer-OP-Continuous-Flow hysteroscope sheath for 2.0mm/30° endoscope and 5 Fr. flexible and semi-rigid instruments


Shape of sheath



Ordering no.	Diameter	WL	Description
AN-001-301	∅ 4.3mm	181mm	Inner-OP-Continuous-Flow hysteroscope sheath for 2.0mm/30° endoscope and 5 Fr. flexible and semi-rigid instruments


3.4 | HIGH-FLOW CONTINUOUS-FLOW OPERATING SHEATHS

Shape of sheath



Ordering no.	Diameter	Description
AN-001-310	Ø 4.3mm	High-Flow Continuous-Flow Operating hysteroscope sheath for 5 Fr. flexible and semi-rigid instruments for use with 2.0mm endoscope (12°/30°) , single sheath system

Shape of sheath



Ordering no.	Diameter	Description
AN-001-311	Ø 4.3mm	High-Flow Continuous-Flow Operating hysteroscope sheath, for 5 Fr. flexible and semi-rigid instruments for use with 2.0mm endoscope (12°/30°) , single sheath system with stopcock

4. INSTRUMENTS & ACCESSORIES

4.1 | FORCEPS, SCISSORS, INSTRUMENTS

4.2 | ELECTRODES & CONNECTING CABLES

4.3 | ACCESSORIES & HYSTEROSCOPY INSTRUMENTS

4. INSTRUMENTS & ACCESSORIES


1 INSTRUMENTS



2 ELECTRODES



3 INTRAUTERINE PROBE




4 RETRACTORS



4.1 | FORCEPS, SCISSORS, INSTRUMENTS


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	WL	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps oval, double action	40cm	∅ 1.0mm	3 Fr.	AC-001-001	AC-001-201	AC-001-001-R
		40cm	∅ 1.2mm	4 Fr.	AC-001-008	AC-001-208	AC-001-008-R
		40cm	∅ 1.6mm	5 Fr.	AC-001-002	AC-001-202	AC-001-002-R
		40cm	∅ 2.1mm	6 Fr.	AC-001-003	AC-001-203	AC-001-003-R
		40cm	∅ 2.3mm	7 Fr.	AC-001-004	AC-001-204	AC-001-004-R
		40cm	∅ 3.0mm	9 Fr.	AC-001-005	AC-001-205	AC-001-005-R
	Biopsy forceps oval, serrated, double action	40cm	∅ 1.0mm	3 Fr.	AC-002-001	AC-002-201	AC-002-001-R
		40cm	∅ 1.6mm	5 Fr.	AC-002-002	AC-002-202	AC-002-002-R
		40cm	∅ 2.1mm	6 Fr.	AC-002-003	AC-002-203	AC-002-003-R
		40cm	∅ 2.3mm	7 Fr.	AC-002-004	AC-002-204	AC-002-004-R
		40cm	∅ 3.0mm	9 Fr.	AC-002-005	AC-002-205	AC-002-005-R
	Biopsy forceps oval, serrated, teeth interlock, double action	40cm	∅ 1.6mm	5 Fr.	AC-015-002	AC-015-202	AC-015-002-R
		40cm	∅ 2.1mm	6 Fr.		AC-015-203	AC-015-003-R
		40cm	∅ 2.3mm	7 Fr.	AC-015-004	AC-015-204	AC-015-004-R
	Grasping forceps, double action	40cm	∅ 1.2mm	4 Fr.	AC-017-001	AC-017-201	AC-017-001-R
		40cm	∅ 2.3mm	7 Fr.	AC-017-004	AC-017-204	AC-017-004-R
	Grasping forceps, alligator jaw, double action	40cm	∅ 1.0mm	3 Fr.	AC-006-001	AC-006-201	AC-006-001-R
		40cm	∅ 1.2mm	4 Fr.	AC-006-003	AC-006-203	AC-006-003-R
		40cm	∅ 1.6mm	5 Fr.	AC-006-002	AC-006-202	AC-006-002-R
		40cm	∅ 2.3mm	7 Fr.	AC-006-004	AC-006-204	AC-006-004-R
		40cm	∅ 3.0mm	9 Fr.	AC-006-005	AC-006-205	AC-006-005-R

4.1 | FORCEPS, SCISSORS, INSTRUMENTS


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	WL	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	40cm	∅ 2.3mm	7 Fr.	AC-008-004	AC-008-204	AC-008-004-R
		40cm	∅ 3.0mm	9 Fr.	AC-008-005	AC-008-205	AC-008-005-R
	Grasping forceps teeth double action	40cm	∅ 1.2mm	4 Fr.	AC-016-001	AC-016-201	AC-016-001-R
		40cm	∅ 1.6mm	5 Fr.	AC-016-002	AC-016-202	AC-016-002-R
	Grasping forceps teeth interlock, double action	40cm	∅ 1.6mm	5 Fr.	AC-020-002	AC-020-202	AC-020-002-R
		40cm	∅ 2.3mm	7 Fr.	AC-020-004	AC-020-204	AC-020-004-R
	Grasping forceps rat tooth, double action	40cm	∅ 1.0mm	3 Fr.	AC-019-001	AC-019-201	AC-019-001-R
		40cm	∅ 1.6mm	5 Fr.	AC-019-002	AC-019-202	AC-019-002-R
	Grasping forceps for stone street, double action	40cm	∅ 1.2mm	4 Fr.	AC-007-001	AC-007-201	AC-007-001-R
		40cm	∅ 1.6mm	5 Fr.	AC-007-002	AC-007-202	AC-007-002-R
	Grasping forceps for big stones, double action	40cm	∅ 1.2mm	4 Fr.	AC-009-001	AC-009-201	AC-009-001-R
		40cm	∅ 1.6mm	5 Fr.	AC-009-002	AC-009-202	AC-009-002-R
		40cm	∅ 2.3mm	7 Fr.	AC-009-004	AC-009-204	AC-009-004-R
	3-Prongs retriever	40cm	∅ 1.0mm	3 Fr.	AC-012-001		
		40cm	∅ 1.6mm	5 Fr.	AC-012-002		
		40cm	∅ 2.3mm	7 Fr.	AC-012-004		


4.1 | FORCEPS, SCISSORS, INSTRUMENTS

SCISSORS FLEXIBLE | SEMI-RIGID | RIGID




	Description	WL	Diameter	French	Flexible	Semi-rigid	Rigid
	Scissor, serrated flexible, double action	40cm	Ø 2.3mm	7 Fr.	AC-005-004	AC-005-204	AC-005-004-R
	Scissor, single action	40cm	Ø 1.0mm	3 Fr.	AC-003-000	AC-003-200	AC-003-000-R
		40cm	Ø 1.2mm	4 Fr.	AC-003-001	AC-003-201	AC-003-001-R
		40cm	Ø 1.6mm	5 Fr.	AC-003-002	AC-003-202	AC-003-002-R
		40cm	Ø 2.1mm	6 Fr.	AC-003-003	AC-003-203	AC-003-003-R
		40cm	Ø 2.3mm	7 Fr.	AC-003-004	AC-003-204	AC-003-004-R
		40cm	Ø 3.0mm	9 Fr.	AC-003-005	AC-003-205	AC-003-005-R
	Uretero slitting scissors, single action	40cm	Ø 1.6mm	5 Fr.		AC-018-201	AC-018-001-R
		40cm	Ø 2.3mm	7 Fr.		AC-018-202	AC-018-002-R
	Scissors pointed, single action	40cm	Ø 1.6mm	5 Fr.		AC-003-302	
		40cm	Ø 2.3mm	7 Fr.		AC-003-304	
	Hook scissors, single action	40cm	Ø 2.3mm	7 Fr.	AC-004-004	AC-004-204	AC-004-004-R
		40cm	Ø 3.0mm	9 Fr.	AC-004-005	AC-004-205	AC-004-005-R

4.1 | FORCEPS, SCISSORS, INSTRUMENTS INSTRUMENTS SEMI-RIGID



	Description	WL	Diameter	Order No.
	Biopsy forceps oval, double action	340mm	Ø 1.6mm Ø 2.3mm	AC-001-502 AC-001-504
		410mm	Ø 1.6mm Ø 2.3mm	AC-001-512 AC-001-514
	Biopsy forceps oval, serrated, double action	340mm	Ø 1.6mm	AC-002-502
		410mm	Ø 1.6mm	AC-002-512
	Biopsy forceps oval, serrated, single action	340mm	Ø 1.6mm Ø 2.3mm	AC-010-502 AC-010-504
		410mm	Ø 1.6mm Ø 2.3mm	AC-010-512 AC-010-514
	Biopsy punch forceps, single action	340mm	Ø 1.6mm Ø 2.3mm	AC-014-502 AC-014-504
		410mm	Ø 1.6mm Ø 2.3mm	AC-014-512 AC-014-514
	Grasping forceps alligator jaw, double action	340mm	Ø 1.6mm Ø 2.3mm	AC-006-502 AC-006-504
		410mm	Ø 1.6mm Ø 2.3mm	AC-006-512 AC-006-514
	Scissors, single action	340mm	Ø 1.6mm Ø 2.3mm	AC-003-502 AC-003-504
		410mm	Ø 1.6mm Ø 2.3mm	AC-003-512 AC-003-514
	Scissors pointed, single action	340mm	Ø 1.6mm Ø 2.3mm	AC-003-602 AC-003-604
		410mm	Ø 1.6mm Ø 2.3mm	AC-003-612 AC-003-614

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




FEATURES

- 1 QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
- 2 ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
- 3 OVERLOAD PROTECTION AVAILABLE
- 4 SIDE-MOUNTED LUER LOCK

4.1 | FORCEPS, SCISSORS, INSTRUMENTS

INSTRUMENTS SEMI-RIGID | DETACHABLE



Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-200-T	AC-003-200-OT
		Ø 1.2mm	4 Fr.	AC-003-201-T	AC-003-201-OT
		Ø 1.6mm	5 Fr.	AC-003-202-T	AC-003-202-OT
		Ø 2.3mm	7 Fr.	AC-003-204-T	AC-003-204-OT
		Ø 3.0mm	9 Fr.	AC-003-205-T	AC-003-205-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-003-612-T	AC-003-612-OT
		Ø 2.3mm	7 Fr.	AC-003-614-T	AC-003-614-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-204-T	AC-004-204-OT
		Ø 3.0mm	9 Fr.	AC-004-205-T	AC-004-205-OT
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-201-T	AC-001-201-OT
		Ø 1.2mm	4 Fr.	AC-001-206-T	AC-001-206-OT
		Ø 1.6mm	5 Fr.	AC-001-202-T	AC-001-202-OT
		Ø 2.3mm	7 Fr.	AC-001-203-T	AC-001-203-OT
		Ø 3.0mm	9 Fr.	AC-001-204-T	AC-001-204-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-006-201-T	AC-006-201-OT
		Ø 1.2mm	4 Fr.	AC-006-203-T	AC-006-203-OT
		Ø 1.6mm	5 Fr.	AC-006-202-T	AC-006-202-OT
		Ø 2.3mm	7 Fr.	AC-006-204-T	AC-006-204-OT
		Ø 3.0mm	9 Fr.	AC-006-205-T	AC-006-205-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-201-T	AC-007-201-OT
		Ø 1.6mm	5 Fr.	AC-007-202-T	AC-007-202-OT
	Grasping forceps for big stones, double action	Ø 1.2mm	5 Fr.	AC-009-201-T	AC-009-201-OT
		Ø 1.6mm	5 Fr.		AC-009-202-OT

4.1 | FORCEPS, SCISSORS, INSTRUMENTS

INSTRUMENTS RIGID | DETACHABLE



Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-103-200-T	AC-103-200-OT
		Ø 1.2mm	4 Fr.	AC-103-201-T	AC-103-201-OT
		Ø 1.6mm	5 Fr.	AC-103-202-T	AC-103-202-OT
		Ø 2.3mm	7 Fr.	AC-103-204-T	AC-103-204-OT
		Ø 3.0mm	9 Fr.	AC-103-205-T	AC-103-205-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-103-612-T	AC-103-612-OT
		Ø 2.3mm	7 Fr.	AC-103-614-T	AC-103-614-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-104-204-T	AC-104-204-OT
		Ø 3.0mm	9 Fr.	AC-104-205-T	AC-104-205-OT
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-101-201-T	AC-101-201-OT
		Ø 1.2mm	4 Fr.	AC-101-206-T	AC-101-206-OT
		Ø 1.6mm	5 Fr.	AC-101-202-T	AC-101-202-OT
		Ø 2.3mm	7 Fr.	AC-101-203-T	AC-101-203-OT
		Ø 3.0mm	9 Fr.	AC-101-204-T	AC-101-204-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-106-201-T	AC-106-201-OT
		Ø 1.2mm	4 Fr.	AC-106-203-T	AC-106-203-OT
		Ø 1.6mm	5 Fr.	AC-106-202-T	AC-106-202-OT
		Ø 2.3mm	7 Fr.	AC-106-204-T	AC-106-204-OT
		Ø 3.0mm	9 Fr.	AC-106-205-T	AC-106-205-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-107-201-T	AC-107-201-OT
		Ø 1.6mm	5 Fr.	AC-107-202-T	AC-107-202-OT
	Grasping forceps for big stones, double action	Ø 1.2mm	5 Fr.	AC-109-201-T	AC-109-201-OT

4.2 | ELECTRODES & CONNECTING CABLES

MONOPOLAR FLEXIBLE ELECTRODES



	Description	French	Length	Ordering no.	Length	Ordering no.
	Loop electrode flexible, monopolar	5 Fr.	530mm	AD-003-003	650mm	AD-003-013
		7 Fr.	530mm	AD-003-005	650mm	AD-003-015
		9 Fr.	530mm	AD-003-007	650mm	AD-003-017
		10 Fr.	530mm	AD-003-008	650mm	AD-003-018
	Needle electrode flexible, monopolar	3 Fr.	530mm	AD-002-001	650mm	AD-002-011
		4 Fr.	530mm	AD-002-002	650mm	AD-002-012
		5 Fr.	530mm	AD-002-003	650mm	AD-002-013
		7 Fr.	530mm	AD-002-005	650mm	AD-002-015
	Button electrode flexible, monopolar	3 Fr.	530mm	AD-001-001	650mm	AD-001-011
		4 Fr.	530mm	AD-001-002	650mm	AD-001-012
		5 Fr.	530mm	AD-001-003	650mm	AD-001-013
		6 Fr.	530mm	AD-001-004	650mm	AD-001-014
		7 Fr.	530mm	AD-001-005	650mm	AD-001-015
		8 Fr.	530mm	AD-001-006	650mm	AD-001-016
		9 Fr.	530mm	AD-001-007	650mm	AD-001-017
		10 Fr.	530mm	AD-001-008	650mm	AD-001-018


4.2 | ELECTRODES & CONNECTING CABLES

CONNECTING CABLES | MONOPOLAR



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	∅ 4mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	∅ 4mm	
	AD-010-002	3m	∅ 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	∅ 5mm	
	AD-010-003	3m	∅ 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	∅ 8mm	
	AD-010-024	4.5m	∅ 2mm	HF Monopolar connecting cable male for STEMA MASTER

Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-004	5m	∅ 4mm	HF Monopolar connecting cable for units with 4mm banana plug Not available in the EU
	AD-010-005		∅ 4mm	

4.2 | ELECTRODES & CONNECTING CABLES

CONNECTING CABLES | MONOPOLAR



CONNECTION UNIT SIDE

Instrument Side	Ordering no.	Length	Diameter	Description
	AD-010-025	4.5m	90° angled connector	HF Monopolar connecting cable for STORZ resectoscope
	AD-010-026	4.5m	∅ 3mm	HF Monopolar connecting cable for OLYMPUS resectoscope



RECOMMENDED PRODUCTS:

CI-001-030 STEMA HF-U PRO TOUCH

- *For all monopolar & bipolar techniques*
- *with 7" touch screen display*
- *provided to perform the most advanced electrosurgery*
- *ADC / APC System*

CI-001-040 DOUBLE PEDAL FOOTSWITCH

4.2 | ELECTRODES & CONNECTING CABLES

BIPOLAR FLEXIBLE ELECTRODES



	Description	French	Length	Ordering no.
	Button electrode flexible, bipolar, flat plug	5 Fr.	360mm	AD-004-001
	Needle electrode angled, flexible, bipolar, flat plug	5 Fr.	360mm	AD-004-002
	Needle electrode 90° angled, flexible, bipolar, flat plug	5 Fr.	360mm	AD-004-003
	Needle electrode straight, flexible, bipolar, flat plug	5 Fr.	360mm	AD-004-004

4.2 | ELECTRODES & CONNECTING CABLES

CONNECTING CABLES | BIPOLAR



CONNECTION INSTRUMENT SIDE



Connection Unit Side	Ordering no.	Length	Description
	EK-010-010	3m	Bipolar connecting cable, flat plug, for Berchtold, Martin, Aesculap
	EK-010-011	5m	
	EK-010-012	3m	Bipolar connecting cable for, flat plug, Erbe
	EK-010-013	5m	
	EK-010-014	3m	Bipolar connecting cable, flat plug, for STEMA, Valleylab/Conmed, EMC
	EK-010-015	5m	
	EK-010-016	3m	Bipolar connecting cable, flat plug, for units with 4mm banana plug
	EK-010-017	5m	

4.3 | ACCESSORIES

COHEN INTRA UTERINE PROBE




Ordering no.	Description
AM-001-001	COHEN intrauterine probe, complete with 2 portio adapters
AM-001-002	COHEN intrauterine probe, only
AM-001-003	Portio adapter, small
AM-001-004	Portio adapter, large



Ordering no.	Description
AM-004-001	Ø 5mm myoma drill, working length 330mm
AM-004-002	Ø 5mm myoma drill, working length 420mm
AM-004-004	Ø 10mm myoma drill, working length 330mm
AM-004-006	Ø 10mm myoma drill, working length 420mm

4.3 | ACCESSORIES

UTERUS MANIPULATOR




Ordering no.	Description
AM-002-001	Uterus manipulator complete incl. 6 adaptors
AM-002-002	Uterus manipulator, only

	Ordering no.	Description
	AM-002-003	Adaptor extra long, \varnothing 10 x 100mm
	AM-002-004	Adaptor long, \varnothing 8 x 80mm
	AM-002-005	Adaptor thin, long, \varnothing 3 x 55mm
	AM-002-006	Adaptor thin, short, \varnothing 3 x 45mm
	AM-002-007	Adaptor with conus, long
	AM-002-008	Adaptor with conus, short

4.3 | ACCESSORIES

UTERUS MANIPULATOR



Ordering no.	Description
AM-002-010	Uterus manipulator set consisting of: Basic instrument, 1 valve, 3 sealing foams and 11 adaptors (20-120mm)
AM-002-011	Adapter set for AM-002-010 consisting of 11 adaptors (20-120mm)
AM-002-012	Ceramic cap small, \varnothing 24mm inner, \varnothing 32mm outer
AM-002-013	Ceramic cap medium, \varnothing 28mm inner, \varnothing 36mm outer
AM-002-014	Ceramic cap large, \varnothing 36mm inner, \varnothing 40mm outer
AM-002-014-1	Teflon cap large, \varnothing 42mm inner, \varnothing 48mm outer
AM-002-015	Foam padded seal sterile, 3 pieces/box
AM-002-016	Valve tube with hand lever, \varnothing 16mm
AM-002-018	Basic instrument for uterus manipulator

4.3 | ACCESSORIES

UTERUS MANIPULATOR | PALPATION PROBE



Ordering no.	Description
AM-002-020	Uterus manipulator complete
AM-002-021	Handle rotatable
AM-002-022	Sheath with fixation nut
AM-002-023	Inner sheath
AM-002-024	Sliding and fixation sheath
AM-002-025	Tube 35mm
AM-002-026	Tube 40mm
AM-002-027	Tube 45mm
AM-002-028	Cervical tenaculum forceps
AM-002-029	Intrauterin adapter, 4 x 40mm
AM-002-030	Intrauterin adapter, 5 x 50mm
AM-002-031	Intrauterin adapter, 6 x 60mm
AM-002-032	Intrauterin adapter, 7 x 70mm



Ordering no.	Description
AN-002-020	Palpation probe with scaling, semi-rigid, 5 Fr., WL 34cm

4.3 | ACCESSORIES

TENACULUM FORCEPS



AJ-005-001



AJ-005-002

Ordering no.	Description
AJ-005-001	Barret, tenaculum forceps, 18cm
AJ-005-002	Pozzi, tenaculum forceps, strong pattern, 25.5cm

4.3 | ACCESSORIES DILATATOR BOUGIES



VAN BUREN

DITTEL,
STRAIGHT

DITTEL,
CURVED

French	Ordering no. VAN BUREN	Ordering no. DITTEL straight	Ordering no. DITTEL curved
8 Fr.	AJ-003-001	AJ-003-021	AJ-003-050
9 Fr.		AJ-003-022	AJ-003-051
10 Fr.	AJ-003-002	AJ-003-023	AJ-003-052
11 Fr.		AJ-003-024	AJ-003-053
12 Fr.	AJ-003-003	AJ-003-025	AJ-003-054
13 Fr.		AJ-003-026	AJ-003-055
14 Fr.	AJ-003-004	AJ-003-027	AJ-003-056
15 Fr.		AJ-003-028	AJ-003-057
16 Fr.	AJ-003-005	AJ-003-029	AJ-003-058
17 Fr.		AJ-003-030	AJ-003-059
18 Fr.	AJ-003-006	AJ-003-031	AJ-003-060
19 Fr.		AJ-003-032	AJ-003-061
20 Fr.	AJ-003-007	AJ-003-033	AJ-003-062
21 Fr.		AJ-003-034	AJ-003-063
22 Fr.	AJ-003-008	AJ-003-035	AJ-003-064
23 Fr.		AJ-003-036	AJ-003-065
24 Fr.	AJ-003-009	AJ-003-037	AJ-003-066
25 Fr.		AJ-003-038	AJ-003-067
26 Fr.	AJ-003-010	AJ-003-039	AJ-003-068
27 Fr.		AJ-003-040	AJ-003-069
28 Fr.	AJ-003-011	AJ-003-041	AJ-003-070
29 Fr.		AJ-003-042	AJ-003-071
30 Fr.	AJ-003-012	AJ-003-043	AJ-003-072
31 Fr.			
32 Fr.	AJ-003-013	AJ-003-044	AJ-003-074
33 Fr.			
34 Fr.	AJ-003-014	AJ-003-045	
35 Fr.			
36 Fr.	AJ-003-015	AJ-003-046	
37 Fr.			
38 Fr.	AJ-003-016		

4.3 | ACCESSORIES

METAL CATHETERS



FOR MEN


FOR WOMEN

COXETER
PROSTATA

French	Ordering no.	Ordering no.	Ordering no.
6 Fr.		AJ-003-100	
7 Fr.	AJ-003-081		
8 Fr.	AJ-003-082	AJ-003-102	AJ-003-120
9 Fr.			
10 Fr.	AJ-003-084	AJ-003-104	AJ-003-121
11 Fr.			
12 Fr.	AJ-003-086	AJ-003-106	AJ-003-122
13 Fr.			
14 Fr.	AJ-003-088	AJ-003-108	AJ-003-123
15 Fr.			
16 Fr.	AJ-003-090	AJ-003-110	AJ-003-124
17 Fr.			
18 Fr.	AJ-003-092	AJ-003-112	AJ-003-125
19 Fr.			
20 Fr.	AJ-003-094	AJ-003-114	AJ-003-126
21 Fr.			
22 Fr.	AJ-003-096	AJ-003-116	AJ-003-127
23 Fr.			
24 Fr.	AJ-003-098	AJ-003-118	AJ-003-128

4.3 | ACCESSORIES

UTERINE PROBES



Ordering no.	Description
AJ-006-001	Sims uterine probe rigid , 32cm
AJ-006-002	Sims uterine probe malleable , 32cm
AJ-006-003	Martin uterine probe rigid , 32cm
AJ-006-004	Martin uterine probe malleable , 32cm


BIOPSY CURETTES



Ordering no.	Description
AJ-006-005	Novak endometrium biopsy curette, Ø 1mm , 23cm
AJ-006-006	Novak endometrium biopsy curette, Ø 2mm , 23cm
AJ-006-007	Novak endometrium biopsy curette, Ø 3mm , 23cm
AJ-006-008	Novak endometrium biopsy curette, Ø 4mm , 23cm

4.3 | ACCESSORIES

HEGAR UTERINE DILATATOR




Ordering no.	Description
AJ-003-078	Hegar uterine dilatation set (8 pieces), double ended, \varnothing 3/4mm - 7/18mm
AJ-003-079	Hegar uterine dilatation set (13 pieces), double ended, \varnothing 1/2mm - 25/26mm

SECRETION SPOON




Ordering no.	Description
AJ-006-012	Guttman secretion spoon, 26.5cm

RETRACTOR




Ordering no.	Description
AM-002-040	LOWSLEY retractor, straight , \varnothing 21.5Fr., working length 280mm
AM-002-041	LOWSLEY retractor, curved , \varnothing 21.5Fr., working length 280mm

4.3 | ACCESSORIES | STERILIZATION



Ordering no.	Dimension (l x w x h)	Description
FA-010-001	290 x 80 x 52mm	Wire basket for 1 endoscope with lid
FA-010-002	460 x 80 x 52mm	
FA-010-004	290 x 80 x 52mm	Wire basket for 2 endoscopes with lid
FA-010-005	460 x 80 x 52mm	



Ordering no.	Dimension (l x w x h)	Description
FA-010-006	670 x 80 x 52mm	Wire basket for 1 endoscope and 1 cold light cable with lid

4.3 | STERILIZATION



Ordering no.	Colour	Dimension (l x w x h)	Description
FA-001-001	Silver	450 x 70 x 70mm	Endoscope container, lid and bottom perforated (without filter)
FA-001-001-R	Red		
FA-001-001-G	Green		
FA-002-001	Silver	250 x 60 x 50mm	Endoscope container, lid and bottom perforated (without filter)
FA-002-001-R	Red		
FA-002-001-G	Green		



Ordering no.	Description
FA-001-000	Protection tube for endoscopes, plastic, 10 pcs, 19cm

5. SALINE RESECTOSCOPY

5.1 | 24/26 FRENCH

5.2 | 19/22 FRENCH (SLIM)

5.3 | 17.5 FRENCH (MINI)

5.4 | CONNECTING CABLES

5. SALINE RESECTOSCOPY

24/26 Fr.
28 Fr.

19/22 Fr.

17.5 Fr.


COMPLETE RESECTOSCOPE




1 CYSTOSCOPES



2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS



4 ELECTRODES



5.1 | SALINE RESECTOSCOPY CYSTOSCOPES AUTOCLAVABLE | Ø 4.0MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-012

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013
	ULTRA HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-023

5.1 | SALINE RESECTOSCOPY

CONTINUOUS FLOW SHEATHS | ROTATABLE



	Ordering no.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow
	AH-001-312	Outer sheath with perforation, rotatable	28.5 Fr.	Black
	AH-001-304	Outer sheath with upper/lower perforation, rotatable	26 Fr.	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow
AH-001-311	Inner sheath with ceramic isolation, rotatable	28.5 Fr.	Black

5.1 | SALINE RESECTOSCOPY OBTURATORS



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-012	Standard Obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-015	Deflecting Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-016	Deflecting Obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-013	Viewing Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-017	Viewing Obturator for resectoscope sheaths with 5 Fr. Instrument channel	24 Fr.	26 Fr.	Yellow
AH-001-018	Viewing Obturator for resectoscope sheaths with 5 Fr. Instrument channel	27 Fr.	28.5 Fr.	Black

5.1 | SALINE RESECTOSCOPY WORKING ELEMENTS



Ordering no.

Description

AG-001-011

SALINE working element used for resectoscopes, **24 Fr.**, **passive cutting**, handle closed, cutting spring assistance



Ordering no.

Description


AG-001-012

SALINE working element used for resectoscopes, **24 Fr.**, **active cutting**, handle closed, cutting spring assistance

5.1 | SALINE RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY



	Ordering no.	French	Description
	AE-001-030	24 Fr.	Cutting loop electrode 30° angled, strong cut
	AE-001-031	24 Fr.	Cutting loop electrode 30° angled
	AE-001-033	24 Fr.	Cutting loop electrode 90° angled
	AE-001-048	24 Fr.	Cutting loop electrode 90° angled, for 0° Scope
	AE-001-036	24 Fr.	Cutting loop electrode straight
	AE-001-032	24 Fr.	Knife electrode
	AE-001-041	24 Fr.	Band electrode, 30° angled
	AE-001-038	24 Fr.	Vaporisation spike rollerball electrode Ø 3mm
	AE-001-037	24 Fr.	Vaporisation rollerball electrode Ø 3mm
	AE-001-043	24 Fr.	Vaporisation rollerball electrode Ø 5mm
	AE-001-040	24 Fr.	Roller electrode smooth Ø 3mm
	AE-001-042	24 Fr.	Roller electrode smooth Ø 5mm

5.1 | SALINE RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY



	Ordering no.	French	Description
	AE-001-039	24 Fr.	Conical electrode
	AE-001-044	24 Fr.	Conical electrode „Mushroom“
	AE-001-045	24 Fr.	Conical electrode „Mushroom“, convex
	AE-001-046	24 Fr.	Conical electrode „Mushroom“, concav
	AE-001-047	24 Fr.	Conical electrode „Mushroom“, grooved

	Ordering no.	French	Description
	AE-002-031	27 Fr.	Cutting loop electrode, 30° angled



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of eletcrodes and cold knives

5.2 | SALINE SLIM RESECTOSCOPY

CYSTOSCOPES AUTOCLAVABLE | Ø 2.9MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
12°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-004
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-014

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

5.2 | SALINE SLIM RESECTOSCOPY SHEATHS



Ordering no.	Description
AH-002-302	Outer sheath 22 Fr. with perforation, rotatable



Ordering no.	Description
AH-002-301	Inner sheath 22 Fr. with ceramic isolation, rotatable

5.2 | SALINE SLIM RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-011	Standardobturator for resectoscope sheaths 19 Fr.



Ordering no.	Description
AH-002-013	Viewing Obturator for resectoscope sheaths 19 Fr.

5.2 | SALINE SLIM RESECTOSCOPY WORKING ELEMENTS



Ordering no.

Description

AG-002-011

SALINE working element used for resectoscopes, **19 Fr., passive cutting**, handle closed, cutting spring assistance



Ordering no.

Description

AG-002-012

SALINE working element used for resectoscopes, **19 Fr., active cutting**, handle closed, cutting spring assistance

5.2 | SALINE SLIM RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY



	Ordering no.	French	Description
	AE-003-031	19 Fr.	Cutting loop electrode 30° angled
	AE-003-032	19 Fr.	Cutting loop electrode 90° angled
	AE-003-033	19 Fr.	Cutting loop electrode straight
	AE-003-034	19 Fr.	Rollerball electrode Ø 3mm
	AE-003-040	19 Fr.	Roller electrode Ø 3mm
	AE-003-036	19 Fr.	Knife electrode




Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

5.3 | SALINE MINI RESECTOSCOPY CYSTOSCOPES




Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

5.3 | SALINE MINI RESECTOSCOPY WORKING ELEMENTES



Ordering no.	Description
AG-002-015	SALINE working element used for resectoscopes, 18.5 Fr. , passive cutting , handle closed, cutting spring assistance



Ordering no.	Description
AG-002-013	SALINE working element used for resectoscopes, 18.5 Fr. , active cutting , handle closed, cutting spring assistance

5.3 | SALINE MINI RESECTOSCOPY SHEATHS



Ordering no.	Description
AH-002-213	Outer sheath 18.5 Fr. with perforation



Ordering no.	Description
AH-002-211	Inner sheath 17.5 Fr. with ceramic isolation, rotatable



Ordering no.	Description
AH-002-216	Inner sheath 17.5 Fr. without ceramic isolation , with laser channel up to 5 Fr., for use with AH-001-213

5.3 | SALINE MINI RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-015	Standardobturator 17.5 Fr.



Ordering no.	Description
AH-002-016	Viewing Obturator 17.5 Fr. with 5 fr. instrument channel

5.3 | SALINE MINI RESECTOSCOPY

ELECTRODES FOR BIPOLAR RESECTOSCOPY




	Ordering no.	French	Description
	AE-003-041	17.5 Fr.	Cutting loop electrode 30° angled
	AE-003-043	17.5 Fr.	Cutting loop electrode straight
	AE-003-044	17.5 Fr.	Rollerball electrode Ø 2mm
	AE-003-045	17.5 Fr.	Roller electrode Ø 2mm
	AE-003-046	17.5 Fr.	Knife electrode



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

5.4 | SALINE RESECTOSCOPY CONNECTING CABLES | BIPOLAR



CONNECTION INSTRUMENT SIDE

Connection Unit Side	Ordering no.	Length	Description
	AD-010-011	4.5m	Bipolar resectoscope connecting cable for Martin Maxium
	AD-010-012	4.5m	Bipolar resectoscope connecting cable for Erbe Vio
	AD-010-013	4.5m	Bipolar resectoscope connecting cable for Valleylab, Lamidey, EMC
	AD-010-015	4.5m	Bipolar resectoscope connecting cable for Erbe Vio 3

6. HYSTEROSCOPY SETS

6.1 | DIAGNOSTIC HYSTEROSCOPY SET

6.2 | OPERATING HYSTEROSCOPY SET

6.3 | OUT-PATIENT ATRAUMATIC OP-HYSTEROSCOPY SET

6.4 | SALINE RESECTOSCOPY SET

6.1 CONTINUOUS-FLOW DIAGNOSTIC HYSTEROSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable, 302mm working length

1 AA-002-013



2 HYSTEROSCOPY INNER SHEATH

Inner-Diagnostic-Continuous-Flow sheath for 4mm / 30° scope, Ø 5.4mm

2 AN-003-102



3 HYSTEROSCOPY OUTER SHEATH

Outer-Diagnostic-Continuous-Flow sheath for 4mm scopes, Ø 6.5mm

3 AN-003-103



4 OBTURATOR

Used for Diagnostic-Continuous-Flow hysteroscope sheaths for 4mm scopes

4 AN-003-104



6.2 CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable, 302mm working length

1 AA-002-013



2 HYSTEROSCOPY OUTER SHEATH

Outer-Diagnostic-Continuous-Flow sheath for 4mm / 30° scope, Ø 6.5mm

2 AN-003-302



3 HYSTEROSCOPY INNER SHEATH

Inner-Diagnostic-Continuous-Flow sheath for 4mm scopes, Ø 5.4mm

3 AN-003-301



4 BIOPSY FORCEPS

oval, flexible, double action, Ø 2.3mm, WL 40cm

4 AC-001-004



5 SCISSOR FLEXIBLE

single action, Ø 2.3mm, WL 40cm

5 AC-003-004



6 GRASPING FORCEPS

alligator jaw flexible, double action, Ø 2.3mm, WL 40cm

6 AC-006-004



7 BUTTON ELECTRODE FLEXIBLE

7Fr., monopolar, length 530mm

7 AD-001-005



8 NEEDLE ELECTRODE FLEXIBLE

7Fr., monopolar, length 530mm

8 AD-002-005



6.3 OUT-PATIENT ATRAUMATIC OPERATING HYSTEROSCOPY SET

1 HD HYSTEROSCOPE

Ø 2.9mm, 30° viewing angle, autoclavable, 302mm working length

1 AA-003-012



2 HYSTEROSCOPY OUTER SHEATH

Outer-Diagnostic-Continuous-Flow sheath for 2.9mm / 30° scope, Ø 6.0mm

2 AN-002-302



3 HYSTEROSCOPY INNER SHEATH

Inner-Diagnostic-Continuous-Flow sheath for 2.9mm scopes, Ø 5.3mm

3 AN-002-301



4 BIOPSY FORCEPS

oval, flexible, double action, Ø 1.6mm, WL 40cm

4 AC-001-002



5 SCISSOR FLEXIBLE

single action, Ø 1.6mm, WL 40cm

5 AC-003-002



6 GRASPING FORCEPS

alligator jaw flexible, double action, Ø 1.6mm, WL 40cm

6 AC-006-002



7 BUTTON ELECTRODE FLEXIBLE

5Fr., monopolar, length 530mm

7 AD-001-003



8 NEEDLE ELECTRODE FLEXIBLE

5Fr., monopolar, length 530mm

8 AD-002-003



6.4 SALINE MINI RESECTOSCOPY SET

1 HD CYSTOSCOPE

Ø 2.9mm, 30° viewing angle, autoclavable, WL 302mm

1 AA-003-012



2 WORKING ELEMENT

for resectoscopes, 18.5 FR., passive cutting, handle closed, cutting spring assistance

2 AG-002-015



3 OUTER SHEATH

18.5 Fr. with perforation

3 AH-002-213



4 INNER SHEATH

17.5 Fr. with ceramic isolation, rotatable

4 AH-002-211



5 STANDARD OBTURATOR

17.5 Fr.

5 AH-002-015



6 CUTTING LOOP ELECTRODE

30° angled, 17.5 Fr.

6 AE-003-041



7 ROLLERBALL ELECTRODE 3 MM

17.5 Fr.

7 AE-003-044



8 KNIFE ELECTRODE

17.5 Fr.

8 AE-003-046



ARE YOU INTERESTED
IN OUR PRODUCTS?
WE WILL GLADLY
SUBMIT YOU
AN OFFER.

7. INDEX

AA-002-001	8	AC-002-004	33	AC-004-204	35
AA-002-002	8	AC-002-004-R	33	AC-004-204-OT	38
AA-002-003	8, 59	AC-002-005	33	AC-004-204-T	38
AA-002-011	8	AC-002-005-R	33	AC-004-205	35
AA-002-012	8, 59	AC-002-201	33	AC-004-205-OT	38
AA-002-013	8, 59	AC-002-202	33	AC-004-205-T	38
AA-002-013	77, 78	AC-002-203	33	AC-005-004	35
AA-002-014	8	AC-002-204	33	AC-005-004-R	35
AA-002-023	59	AC-002-205	33	AC-005-204	35
AA-003-001	17	AC-002-502	36	AC-006-001	33
AA-003-002	17, 65, 70	AC-002-512	36	AC-006-001-R	33
AA-003-004	65	AC-003-000	35	AC-006-002	33, 79
AA-003-011	17	AC-003-000-R	35	AC-006-002-R	33
AA-003-012	17, 65, 70	AC-003-001	35	AC-006-003	33
AA-003-012	79, 80	AC-003-001-R	35	AC-006-003-R	33
AA-003-014	17, 65	AC-003-002	35, 79	AC-006-004	33, 78
AA-004-001	26	AC-003-002-R	35	AC-006-004-R	33
AA-004-002	26	AC-003-003	35	AC-006-005	33
AA-004-003	26	AC-003-003-R	35	AC-006-005-R	33
AC-001-001	33	AC-003-004	35, 78	AC-006-201	33
AC-001-001-R	33	AC-003-004-R	35	AC-006-201-OT	38
AC-001-002	33, 79	AC-003-005	35	AC-006-201-T	38
AC-001-002-R	33	AC-003-005-R	35	AC-006-202	33
AC-001-003	33	AC-003-200	35	AC-006-202-OT	38
AC-001-003-R	33	AC-003-200-OT	38	AC-006-202-T	38
AC-001-004	33, 78	AC-003-200-T	38	AC-006-203	33
AC-001-004-R	33	AC-003-201	35	AC-006-203-OT	38
AC-001-005	33	AC-003-201-OT	38	AC-006-203-T	38
AC-001-005-R	33	AC-003-201-T	38	AC-006-204	33
AC-001-008	33	AC-003-202	35	AC-006-204-OT	38
AC-001-008-R	33	AC-003-202-OT	38	AC-006-204-T	38
AC-001-201	33	AC-003-202-T	38	AC-006-205	33
AC-001-201-OT	38	AC-003-203	35	AC-006-205-OT	38
AC-001-201-T	38	AC-003-204	35	AC-006-205-T	38
AC-001-202	33	AC-003-204-OT	38	AC-006-502	36
AC-001-202-OT	38	AC-003-204-T	38	AC-006-504	36
AC-001-202-T	38	AC-003-205	35	AC-006-512	36
AC-001-203	33	AC-003-205-OT	38	AC-006-514	36
AC-001-203-OT	38	AC-003-205-T	38	AC-007-001	34
AC-001-203-T	38	AC-003-302	35	AC-007-001-R	34
AC-001-204	33	AC-003-304	35	AC-007-002	34
AC-001-204-OT	38	AC-003-502	36	AC-007-002-R	34
AC-001-204-T	38	AC-003-504	36	AC-007-201	34
AC-001-205	33	AC-003-512	36	AC-007-201-OT	38
AC-001-206-OT	38	AC-003-514	36	AC-007-201-T	38
AC-001-206-T	38	AC-003-602	36	AC-007-202	34
AC-001-208	33	AC-003-604	36	AC-007-202-OT	38
AC-001-502	36	AC-003-612	36	AC-007-202-T	38
AC-001-504	36	AC-003-612-OT	38	AC-008-004	34
AC-001-512	36	AC-003-612-T	38	AC-008-004-R	34
AC-001-514	36	AC-003-614	36	AC-008-005	34
AC-002-001	33	AC-003-614-OT	38	AC-008-005-R	34
AC-002-001-R	33	AC-003-614-T	38	AC-008-204	34
AC-002-002	33	AC-004-004	35	AC-008-205	34
AC-002-002-R	33	AC-004-004-R	35	AC-009-001	34
AC-002-003	33	AC-004-005	35	AC-009-001-R	34
AC-002-003-R	33	AC-004-005-R	35	AC-009-002	34

7. INDEX

AC-009-002-R	34	AC-101-201-T	39	AD-001-017	40
AC-009-004	34	AC-101-202-OT	39	AD-001-018	40
AC-009-004-R	34	AC-101-202-T	39	AD-002-001	40
AC-009-201	34	AC-101-203-OT	39	AD-002-002	40
AC-009-201-OT	38	AC-101-203-T	39	AD-002-003	40, 79
AC-009-201-T	38	AC-101-204-OT	39	AD-002-005	40, 78
AC-009-202	34	AC-101-204-T	39	AD-002-011	40
AC-009-202-OT	38	AC-101-206-OT	39	AD-002-012	40
AC-009-204	34	AC-101-206-T	39	AD-002-013	40
AC-010-502	36	AC-103-200-OT	39	AD-002-015	40
AC-010-504	36	AC-103-200-T	39	AD-003-003	40
AC-010-512	36	AC-103-201-OT	39	AD-003-005	40
AC-010-514	36	AC-103-201-T	39	AD-003-007	40
AC-012-001	34	AC-103-202-OT	39	AD-003-008	40
AC-012-002	34	AC-103-202-T	39	AD-003-013	40
AC-012-004	34	AC-103-204-OT	39	AD-003-015	40
AC-014-502	36	AC-103-204-T	39	AD-003-017	40
AC-014-504	36	AC-103-205-OT	39	AD-003-018	40
AC-014-512	36	AC-103-205-T	39	AD-004-001	43
AC-014-514	36	AC-103-612-OT	39	AD-004-002	43
AC-015-002	33	AC-103-612-T	39	AD-004-003	43
AC-015-002-R	33	AC-103-614-OT	39	AD-004-004	43
AC-015-003-R	33	AC-103-614-T	39	AD-010-001	41
AC-015-004	33	AC-104-204-OT	39	AD-010-002	41
AC-015-004-R	33	AC-104-204-T	39	AD-010-003	41
AC-015-202	33	AC-104-205-OT	39	AD-010-004	41
AC-015-203	33	AC-104-205-T	39	AD-010-005	41
AC-015-204	33	AC-106-201-OT	39	AD-010-006	41
AC-016-001	34	AC-106-201-T	39	AD-010-007	41
AC-016-001-R	34	AC-106-202-OT	39	AD-010-008	41
AC-016-002	34	AC-106-202-T	39	AD-010-011	75
AC-016-002-R	34	AC-106-203-OT	39	AD-010-012	75
AC-016-201	34	AC-106-203-T	39	AD-010-013	75
AC-016-202	34	AC-106-204-OT	39	AD-010-015	75
AC-017-001	33	AC-106-204-T	39	AD-010-024	41
AC-017-001-R	33	AC-106-205-OT	39	AD-010-025	42
AC-017-004	33	AC-106-205-T	39	AD-010-026	42
AC-017-004-R	33	AC-107-201-OT	39	AE-001-030	63
AC-017-201	33	AC-107-201-T	39	AE-001-031	63
AC-017-204	33	AC-107-202-OT	39	AE-001-032	63
AC-018-001-R	35	AC-107-202-T	39	AE-001-033	63
AC-018-002-R	35	AC-109-201-OT	39	AE-001-036	63
AC-018-201	35	AC-109-201-T	39	AE-001-037	63
AC-018-202	35	AD-001-001	40	AE-001-038	63
AC-019-001	34	AD-001-002	40	AE-001-039	64
AC-019-001-R	34	AD-001-003	40, 79	AE-001-040	63
AC-019-002	34	AD-001-004	40	AE-001-041	63
AC-019-002-R	34	AD-001-005	40, 78	AE-001-042	63
AC-019-201	34	AD-001-006	40	AE-001-043	63
AC-019-202	34	AD-001-007	40	AE-001-044	64
AC-020-002	34	AD-001-008	40	AE-001-045	64
AC-020-002-R	34	AD-001-011	40	AE-001-046	64
AC-020-004	34	AD-001-012	40	AE-001-047	64
AC-020-004-R	34	AD-001-013	40	AE-001-048	63
AC-020-202	34	AD-001-014	40	AE-002-031	64
AC-020-204	34	AD-001-015	40	AE-003-031	69
AC-101-201-OT	39	AD-001-016	40	AE-003-032	69

7. INDEX

AE-003-033	69	AJ-003-013	50	AJ-003-082	52
AE-003-034	69	AJ-003-014	50	AJ-003-084	52
AE-003-036	69	AJ-003-015	50	AJ-003-086	52
AE-003-040	69	AJ-003-016	50	AJ-003-088	52
AE-003-041	74, 80	AJ-003-021	50	AJ-003-090	52
AE-003-043	74	AJ-003-022	50	AJ-003-092	52
AE-003-044	74, 80	AJ-003-023	50	AJ-003-094	52
AE-003-045	74	AJ-003-024	50	AJ-003-096	52
AE-003-046	74, 80	AJ-003-025	50	AJ-003-098	52
AG-001-011	62	AJ-003-026	50	AJ-003-100	52
AG-001-012	62	AJ-003-027	50	AJ-003-102	52
AG-002-011	68	AJ-003-028	50	AJ-003-104	52
AG-002-012	68	AJ-003-029	50	AJ-003-106	52
AG-002-013	71	AJ-003-030	50	AJ-003-108	52
AG-002-015	71	AJ-003-031	50	AJ-003-110	52
AG-002-015	80	AJ-003-032	50	AJ-003-112	52
AH-001-011	61	AJ-003-033	50	AJ-003-114	52
AH-001-012	61	AJ-003-034	50	AJ-003-116	52
AH-001-013	61	AJ-003-035	50	AJ-003-118	52
AH-001-015	61	AJ-003-036	50	AJ-003-120	52
AH-001-016	61	AJ-003-037	50	AJ-003-121	52
AH-001-017	61	AJ-003-038	50	AJ-003-122	52
AH-001-018	61	AJ-003-039	50	AJ-003-123	52
AH-001-301	60	AJ-003-040	50	AJ-003-124	52
AH-001-302	60	AJ-003-041	50	AJ-003-125	52
AH-001-303	60	AJ-003-042	50	AJ-003-126	52
AH-001-304	60	AJ-003-043	50	AJ-003-127	52
AH-001-305	60	AJ-003-044	50	AJ-003-128	52
AH-001-306	60	AJ-003-045	50	AJ-003-151	51
AH-001-311	60	AJ-003-046	50	AJ-003-152	51
AH-001-312	60	AJ-003-050	50	AJ-003-153	51
AH-002-011	67	AJ-003-051	50	AJ-003-154	51
AH-002-013	67	AJ-003-052	50	AJ-003-155	51
AH-002-015	73	AJ-003-053	50	AJ-003-156	51
AH-002-015	80	AJ-003-054	50	AJ-003-157	51
AH-002-016	73	AJ-003-055	50	AJ-003-158	51
AH-002-211	72	AJ-003-056	50	AJ-003-170	51
AH-002-211	80	AJ-003-057	50	AJ-003-171	51
AH-002-213	72	AJ-003-058	50	AJ-003-172	51
AH-002-213	80	AJ-003-059	50	AJ-003-173	51
AH-002-216	72	AJ-003-060	50	AJ-003-174	51
AH-002-301	66	AJ-003-061	50	AJ-003-175	51
AH-002-301	79	AJ-003-062	50	AJ-003-176	51
AH-002-302	66	AJ-003-063	50	AJ-003-177	51
AH-003-301	78	AJ-003-064	50	AJ-003-178	51
AJ-003-001	50	AJ-003-065	50	AJ-003-179	51
AJ-003-002	50	AJ-003-066	50	AJ-003-180	51
AJ-003-003	50	AJ-003-067	50	AJ-003-181	51
AJ-003-004	50	AJ-003-068	50	AJ-003-182	51
AJ-003-005	50	AJ-003-069	50	AJ-003-183	51
AJ-003-006	50	AJ-003-070	50	AJ-003-184	51
AJ-003-007	50	AJ-003-071	50	AJ-003-185	51
AJ-003-008	50	AJ-003-072	50	AJ-003-186	51
AJ-003-009	50	AJ-003-074	50	AJ-003-187	51
AJ-003-010	50	AJ-003-078	54	AJ-003-188	51
AJ-003-011	50	AJ-003-079	54	AJ-003-189	51
AJ-003-012	50	AJ-003-081	52	AJ-003-190	51

7. INDEX

AJ-003-187	51	AM-002-020	48	AN-003-001	9
AJ-003-188	51	AM-002-021	48	AN-003-002	9
AJ-003-189	51	AM-002-022	48	AN-003-003	9
AJ-003-190	51	AM-002-023	48	AN-003-101	10
AJ-003-191	51	AM-002-024	48	AN-003-102	10, 77
AJ-003-192	51	AM-002-025	48	AN-003-103	10, 77
AJ-003-193	51	AM-002-026	48	AN-003-104	10, 77
AJ-003-194	51	AM-002-027	48	AN-003-201	11
AJ-003-195	51	AM-002-028	48	AN-003-202	11
AJ-003-196	51	AM-002-029	48	AN-003-301	12
AJ-003-197	51	AM-002-030	48	AN-003-302	12, 78
AJ-003-198	51	AM-002-031	48	AN-003-303	12
AJ-003-199	51	AM-002-032	48	AN-004-401	14
AJ-003-200	51	AM-002-040	54	AN-004-402	14
AJ-003-201	51	AM-002-041	54	AN-004-403	14
AJ-003-202	51	AM-004-001	45	AN-005-301	13
AJ-003-203	51	AM-004-002	45	AN-005-302	13
AJ-005-001	49	AM-004-004	45	AN-005-303	13
AJ-005-002	49	AM-004-006	45	AN-005-304	13
AJ-006-001	53	AN-001-001	27		
AJ-006-002	53	AN-001-002	27		
AJ-006-003	53	AN-001-003	27	CI-001-030	42
AJ-006-004	53	AN-001-101	28	CI-001-040	42
AJ-006-005	53	AN-001-102	28		
AJ-006-006	53	AN-001-103	28		
AJ-006-007	53	AN-001-104	28		
AJ-006-008	53	AN-001-301	29	EK-010-010	44
AJ-006-012	54	AN-001-302	29	EK-010-011	44
AJ-010-006	64, 69, 74	AN-001-310	30	EK-010-012	44
AM-001-001	45	AN-001-311	30	EK-010-013	44
AM-001-002	45	AN-002-001	18	EK-010-014	44
AM-001-003	45	AN-002-002	18	EK-010-015	44
AM-001-004	45	AN-002-003	18	EK-010-016	44
AM-002-001	46	AN-002-020	48	EK-010-017	44
AM-002-002	46	AN-002-101	19		
AM-002-003	46	AN-002-102	19		
AM-002-004	46	AN-002-103	19		
AM-002-005	46	AN-002-104	19	FA-001-000	56
AM-002-006	46	AN-002-201	20	FA-001-001	56
AM-002-007	46	AN-002-202	20	FA-001-001-G	56
AM-002-008	46	AN-002-301	21	FA-001-001-R	56
AM-002-010	47	AN-002-302	21, 79	FA-002-001	56
AM-002-011	47	AN-002-307	21	FA-002-001-G	56
AM-002-012	47	AN-002-308	21	FA-002-001-R	56
AM-002-013	47	AN-002-310	23	FA-010-001	55
AM-002-014	47	AN-002-311	23	FA-010-002	55
AM-002-014-1	47	AN-002-312	22	FA-010-004	55
AM-002-015	47	AN-002-313	22	FA-010-005	55
AM-002-016	47	AN-002-312	23	FA-010-006	55
AM-002-018	47	AN-002-313	23		

WE COVER THE FOLLOWING AREAS IN SURGERY WITH HIGH QUALITY PRODUCTS

ARTHROSCOPY

CARDIAC MIS

CONTAINER

ENT

HF-SURGERY

GASTROSCOPY

GYNAECOLOGY

LAMINECTOMY

LAPAROSCOPY

PROCTOSCOPY

SPINE SURGERY

UNITS

UROLOGY



STEMA MEDIZINTECHNIK GMBH | TAKE-OFF GEWERBEPARK 74 | 78579 NEUHAUSEN OB ECK
INFO@STEMA-MEDIZINTECHNIK.DE | WWW.STEMA-MEDIZINTECHNIK.DE | +49 7467 79 630-0

01/2026 | PRINTED IN GERMANY | COPYRIGHT BY STEMA MEDIZINTECHNIK GMBH |
WE RESERVE THE RIGHT TO MAKE ALTERATIONS