




UROLOGY



ABOUT US



SINCE 20 YEARS



EXPORT IN
+70 COUNTRIES



24H/SERVICE



+ 10.000
PRODUCTS



LOCATED IN THE
WORLD CENTER
FOR MEDICAL
TECHNOLOGY

CERTIFICATED

STEMA Medizintechnik GmbH is certified in accordance with the highest standards:

- EN ISO 13485:2016 + AC:2018 + A11:2021- ISO 13485:2016
- EU 2017/745 (**MDR**)
- Council Directive 93/42/EEC Annex II excluding Section 4
- Council Directive 93/42/EEC Annex V excluding Section 3
- FDA

These standards create an essential base for our company and our partners. They also symbolize the commitment and strive to keep up the quality of our products.



PRODUCT RANGE

The current product range of STEMA Medizintechnik GmbH includes more than 10.000 different items.

We cover the following areas in surgery with high quality products:

ARTHROSCOPY

CARDIAC MIS

CONTAINER-SYSTEMS

GASTROENTEROLOGY

GYNAECOLOGY

HF-SURGERY

LAMINECTOMY

LAPAROSCOPY

MICRO MOTOR SYSTEMS

PROCTOSCOPY

SPINE-SURGERY

TITANIUM PLATING SYSTEM

UNITS

UROLOGY



CONTENT

UROLOGY

PRODUCT

PAGE

1. CYSTO-URETHROSCOPY

1.1 Cysto-Urethroscopy	8
1.2 Slim Cysto-Urethroscopy	25
1.3 Laser Cysto-Urethroscopy	28
1.4 HD Mini Compact Cystoscope	30
1.5 Flexible Cystoscope	31
1.6 Cystoscopy Set	33

2. RESECTOSCOPY

2.1 Monopolar Resectoscopy	37
2.2 Saline Resectoscopy	59
2.3 Laser Resectoscopy	76
2.4 Accessories	87
2.5 Resectoscopy Sets	91

3. URETHROTOMY

3.1 Urethrotome	97
3.2 Optical Urethrotomy	98
3.3 Urethrotomy Set	106

4. NEPHROSCOPY

4.1 HD Nephroscopy	109
4.2 Mini Nephroscopy	111
4.3 Standard Nephroscopy	114
4.4 Accessories	118

5. URETERO-RENOSCOPY

5.1 Uretero-Renoscopes semi-rigid	122
5.2 Flexible Uretero-Renoscope	142

CONTENT

UROLOGY

PRODUCT

PAGE

6. LITHOTRIPSY

6.1 Lithotripsy	147
-------------------	-----

7. MORCELLATION

7.1 Morcellation 26 Fr.	153
7.2 Morcellation 22 Fr.	155

8. PEDIATRIC UROLOGY

8.1 Pediatric Cysto-Urethroscopy	159
8.2 Pediatric Resectoscopy	166
8.3 Pediatric Urethrotomy	174

9. ACCESSORIES

9.1 Instruments	181
9.2 Sterilization	186
9.3 Cleaning and Maintenance	188

10. INDEX

10.1 Index	191
--------------	-----

1. CYSTO-URETHROSCOPY

1.1 | CYSTO-URETHROSCOPY

CYSTOSCOPES (4.0 MM)

SHEATHS

OBTURATORS

ENDOSCOPE BRIDGES

WORKING ELEMENTS

FORCEPS, SCISSORS & INSTRUMENTS

STONE EXTRACTORS

MONOPOLAR FLEXIBLE ELECTRODES

MONOPOLAR CONNECTING CABLES

1.2 | SLIM CYSTO-URETHROSCOPY

CYSTOSCOPES (2.9 MM)

SHEATHS

ENDOSCOPE BRIDGES

WORKING ELEMENTS

1.3 | LASER CYSTO-URETHROSCOPY

CYSTOSCOPES (4.0 MM)

SHEATHS

OBTURATORS

1.4 | HD MINI COMPACT CYSTOSCOPES

1.5 | FLEXIBLE CYSTOSCOPE

CYSTOSCOPE FLEX I

ACCESSORIES


1.6 | CYSTO-URETHROSCOPY SET

1. CYSTO-URETHROSCOPY


4.0 mm

2.9 mm

1 CYTOSCOPIES




2 SHEATHS & OBTURATORS



3 ACCESSORIES



1.1 | CYSTO-URETHROSCOPY CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 0°	Cystoscope autoclavable Cystoscope autoclavable wide angle	Ø 4.0 mm	302 mm	AA-002-001 AA-002-001-W
	HD Cystoscope autoclavable HD Cystoscope autoclavable wide angle			AA-002-011 AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-021

Viewing angle	Description	Diameter	Working length	Ordering No.
 12°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-002
	HD Cystoscope autoclavable			AA-002-012
	Ultra HD (4K) Cystoscope autoclavable			AA-002-022

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable Cystoscope autoclavable wide angle	Ø 4.0 mm	302 mm	AA-002-003 AA-002-003-W
	HD Cystoscope autoclavable HD Cystoscope autoclavable wide angle			AA-002-013 AA-002-013-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-023

Viewing angle	Description	Diameter	Working length	Ordering No.
 70°	Cystoscope autoclavable Cystoscope autoclavable wide angle	Ø 4.0 mm	302 mm	AA-002-004 AA-002-004-W
	HD Cystoscope autoclavable HD Cystoscope autoclavable wide angle			AA-002-014 AA-002-014-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-024

1.1 | CYSTO-URETHROSCOPY SHEATHS



Ordering No.	Description	French	Colour code
AB-001-002	Cysto-Urethroscope sheath with 2 fixed stopcocks and standard obturator	17 Fr.	Yellow
AB-001-003		19 Fr.	Green
AB-001-004		21 Fr.	Red
AB-001-005		23 Fr.	Blue
AB-001-006		25 Fr.	White



Ordering No.	Description	French	Colour code
AB-002-001	Cysto-Urethroscope sheath with central valve and standard obturator	17 Fr.	Yellow
AB-002-002		19 Fr.	Green
AB-002-003		21 Fr.	Red
AB-002-004		23 Fr.	Blue
AB-002-005		25 Fr.	White

1.1 | CYSTO-URETHROSCOPY

OBTURATORS | ADAPTER



Ordering No.	Description	French	Colour code
AB-001-012	Standard obturator for cysto-urethroscope sheaths	17 Fr.	Yellow
AB-001-013		19 Fr.	Green
AB-001-014		21 Fr.	Red
AB-001-015		23 Fr.	Blue
AB-001-016		25 Fr.	White



Ordering No.	Description	French	Colour code
AB-001-022	Viewing obturator for Cysto-Urethroscope sheaths	17 Fr.	Yellow
AB-001-023		19 Fr.	Green
AB-001-024		21 Fr.	Red
AB-001-025		23 Fr.	Blue
AB-001-026		25 Fr.	White

Ordering No.	Description
AB-001-010	Cystoscope adapter for female urethroscopy




1.1 | CYSTO-URETHROSCOPY ENDOSCOPE BRIDGES



Ordering No.	Description
AB-003-002	Endoscope bridge diagnostic



Ordering No.	Description
AB-003-003	Endoscope bridge with 1 instrument channel



Ordering No.	Description
AB-003-004	Endoscope bridge with 2 instrument channels



Ordering No.	Description
AB-003-015	Endoscope bridge with 2 instrument channels, short



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel

1.1 | CYSTO-URETHROSCOPY WORKING ELEMENTS



Ordering No.	Description
AB-003-006	Working element with ALBARRAN deflector and 1 instrument channel
AB-003-008	Working element with ALBARRAN deflector and 1 instrument channel, without ratchet



Ordering No.	Description
AB-003-007	Working element with ALBARRAN deflector and 2 instrument channels
AB-003-009	Working element with ALBARRAN deflector and 2 instrument channels, without ratchet

SHEATH/INSTRUMENT COMBINATIONS

Sheaths	Colour	Albarran		Bridge	
		1 working channel	2 working channels	1 working channel	2 working channels
25 FR.		1 × 10 Fr.	2 × 8 Fr.	1 × 12 Fr.	2 × 8 Fr.
23 FR.		1 × 9 Fr.	2 × 6 Fr.	1 × 10 Fr.	2 × 7 Fr.
21 FR.		1 × 6 Fr.	2 × 5 Fr.	1 × 7 Fr.	2 × 6 Fr.
19 FR.		1 × 5 Fr.	2 × 4 Fr.	1 × 6 Fr.	2 × 5 Fr.
17 FR.		-	-	1 × 5 Fr.	-


1.1 | CYSTO-URETHROSCOPY

OPTICAL FORCEPS AND SCISSORS



Info: Only for sheaths 19-25 Fr.



	Ordering No.	For scopes	Description
	AB-004-005	12°/ 30°	Biopsy forceps (down), double action
	AB-004-006	70°	Biopsy forceps (down), double action
	AB-004-001	0°/ 12°/ 30°	Biopsy forceps, double action
	AB-004-002	0°/ 12°/ 30°	Biopsy forceps, through cutting
	AB-004-004	0°/ 12°	Hook scissors, double action
	AB-004-010	0°/ 12°	Forceps, double action, for removal of stents
	AB-004-003	0°/ 12°	Grasping forceps, double action
	AB-004-008		Adapter for use with resectoscope sheaths

1.1 | CYSTO-URETHROSCOPY

OPTICAL FORCEPS



Info: Only for sheaths 25 Fr.



Ordering No.

For scopes

Description

AB-004-007

0°/ 12°/ 30°

Optical stone punch



AB-004-009

Adapter for use with resectoscope sheaths 26–28.5 Fr.




AB-004-009-QL

Adapter with **STEMA click-system** for use with resectoscope sheaths 26–28.5 Fr.

1.1 | CYSTO-URETHROSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps oval, double action	Ø 1.0 mm	3 Fr.	AC-001-001	AC-001-201	AC-001-001-R
		Ø 1.2 mm	4 Fr.	AC-001-008	AC-001-208	AC-001-008-R
		Ø 1.6 mm	5 Fr.	AC-001-002	AC-001-202	AC-001-002-R
		Ø 2.1 mm	6 Fr.	AC-001-003	AC-001-203	AC-001-003-R
		Ø 2.3 mm	7 Fr.	AC-001-004	AC-001-204	AC-001-004-R
		Ø 3.0 mm	9 Fr.	AC-001-005	AC-001-205	AC-001-005-R
	Biopsy forceps oval, serrated, double action	Ø 1.0 mm	3 Fr.	AC-002-001	AC-002-201	AC-002-001-R
		Ø 1.6 mm	5 Fr.	AC-002-002	AC-002-202	AC-002-002-R
		Ø 2.1 mm	6 Fr.	AC-002-003	AC-002-203	AC-002-003-R
		Ø 2.3 mm	7 Fr.	AC-002-004	AC-002-204	AC-002-004-R
		Ø 3.0 mm	9 Fr.	AC-002-005	AC-002-205	AC-002-005-R
	Biopsy forceps oval, serrated, teeth interlock, double action	Ø 1.6 mm	5 Fr.	AC-015-002	AC-015-202	AC-015-002-R
		Ø 2.1 mm	6 Fr.	-	AC-015-203	AC-015-003-R
		Ø 2.3 mm	7 Fr.	AC-015-004	AC-015-204	AC-015-004-R
	Grasping forceps, stone street double action	Ø 1.2 mm	4 Fr.	AC-017-001	AC-017-201	AC-017-001-R
		Ø 2.3 mm	7 Fr.	AC-017-004	AC-017-204	AC-017-004-R
	Grasping forceps, alligator jaw, double action	Ø 1.0 mm	3 Fr.	AC-006-001	AC-006-201	AC-006-001-R
		Ø 1.2 mm	4 Fr.	AC-006-003	AC-006-203	AC-006-003-R
		Ø 1.6 mm	5 Fr.	AC-006-002	AC-006-202	AC-006-002-R
		Ø 2.3 mm	7 Fr.	AC-006-004	AC-006-204	AC-006-004-R
		Ø 3.0 mm	9 Fr.	AC-006-005	AC-006-205	AC-006-005-R

1.1 | CYSTO-URETHROSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID




	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	Ø 2.3 mm	7 Fr.	AC-008-004	AC-008-204	AC-008-004-R
		Ø 3.0 mm	9 Fr.	AC-008-005	AC-008-205	AC-008-005-R
	Grasping forceps teeth double action	Ø 1.2 mm	4 Fr.	AC-016-001	AC-016-201	AC-016-001-R
		Ø 1.6 mm	5 Fr.	AC-016-002	AC-016-202	AC-016-002-R
	Grasping forceps teeth interlock, double action	Ø 1.6 mm	5 Fr.	AC-020-002	AC-020-202	AC-020-002-R
		Ø 2.3 mm	7 Fr.	AC-020-004	AC-020-204	AC-020-004-R
	Grasping forceps rat tooth, double action	Ø 1.0 mm	3 Fr.	AC-019-001	AC-019-201	AC-019-001-R
		Ø 1.6 mm	5 Fr.	AC-019-002	AC-019-202	AC-019-002-R
	Grasping forceps for stone street, double action	Ø 1.2 mm	4 Fr.	AC-007-001	AC-007-201	AC-007-001-R
		Ø 1.6 mm	5 Fr.	AC-007-002	AC-007-202	AC-007-002-R
	Grasping forceps for big stones, double action	Ø 1.2 mm	4 Fr.	AC-009-001	AC-009-201	AC-009-001-R
		Ø 1.6 mm	5 Fr.	AC-009-002	AC-009-202	AC-009-002-R
		Ø 2.3 mm	7 Fr.	AC-009-004	AC-009-204	AC-009-004-R
	3-Prongs retriever	Ø 1.0 mm	3 Fr.	AC-012-001		
		Ø 1.6 mm	5 Fr.	AC-012-002		
		Ø 2.3 mm	7 Fr.	AC-012-004		


1.1 | CYSTO-URETHROSCOPY

SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Hook scissors, single action	Ø 2.3 mm	7 Fr.	AC-004-004	AC-004-204	AC-004-004-R
		Ø 3.0 mm	9 Fr.	AC-004-005	AC-004-205	AC-004-005-R
	Scissors, single action	Ø 1.0 mm	3 Fr.	AC-003-000	AC-003-200	AC-003-000-R
		Ø 1.2 mm	4 Fr.	AC-003-001	AC-003-201	AC-003-001-R
		Ø 1.6 mm	5 Fr.	AC-003-002	AC-003-202	AC-003-002-R
		Ø 2.1 mm	6 Fr.	AC-003-003	AC-003-203	AC-003-003-R
		Ø 2.3 mm	7 Fr.	AC-003-004	AC-003-204	AC-003-004-R
		Ø 3.0 mm	9 Fr.	AC-003-005	AC-003-205	AC-003-005-R
	Uretero-slitting scissors, retrograde cutting single action	Ø 1.6 mm	5 Fr.	AC-018-001	AC-018-201	AC-018-001-R
		Ø 2.3 mm	7 Fr.	AC-018-002	AC-018-202	AC-018-002-R
	Scissors pointed, single action	Ø 1.6 mm	5 Fr.	-	AC-003-302	-
		Ø 2.3 mm	7 Fr.	-	AC-003-304	-
	Scissors, serrated flexible, dou- ble action	Ø 2.3 mm	7 Fr.	AC-005-004	AC-005-204	AC-005-004-R

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




FEATURES

- 1 QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
- 2 ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
- 3 OVERLOAD PROTECTION AVAILABLE
- 4 SIDE-MOUNTED LUER LOCK

1.1 | CYSTO-URETHROSCOPY


INSTRUMENTS DETACHABLE | SEMI-RIGID



Ordering No.		Description			
AC-100-001		Handle for detachable instruments			
Description		Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0 mm	3 Fr.	AC-003-200-T	AC-003-200-OT
		Ø 1.2 mm	4 Fr.	AC-003-201-T	AC-003-201-OT
		Ø 1.6 mm	5 Fr.	AC-003-202-T	AC-003-202-OT
		Ø 2.3 mm	7 Fr.	AC-003-204-T	AC-003-204-OT
		Ø 3.0 mm	9 Fr.	AC-003-205-T	AC-003-205-OT
	Scissors pointed, single action	Ø 1.6 mm	5 Fr.	AC-003-612-T	AC-003-612-OT
		Ø 2.3 mm	7 Fr.	AC-003-614-T	AC-003-614-OT
	Hook scissors, single action	Ø 2.3 mm	7 Fr.	AC-004-204-T	AC-004-204-OT
		Ø 3.0 mm	9 Fr.	AC-004-205-T	AC-004-205-OT
	Biopsy forceps oval, double action	Ø 1.0 mm	3 Fr.	AC-001-201-T	AC-001-201-OT
		Ø 1.2 mm	4 Fr.	AC-001-206-T	AC-001-206-OT
		Ø 1.6 mm	5 Fr.	AC-001-202-T	AC-001-202-OT
		Ø 2.3 mm	7 Fr.	AC-001-203-T	AC-001-203-OT
		Ø 3.0 mm	9 Fr.	AC-001-204-T	AC-001-204-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0 mm	3 Fr.	AC-006-201-T	AC-006-201-OT
		Ø 1.2 mm	4 Fr.	AC-006-203-T	AC-006-203-OT
		Ø 1.6 mm	5 Fr.	AC-006-202-T	AC-006-202-OT
		Ø 2.3 mm	7 Fr.	AC-006-204-T	AC-006-204-OT
		Ø 3.0 mm	9 Fr.	AC-006-205-T	AC-006-205-OT
	Grasping forceps for stone street, double action	Ø 1.2 mm	4 Fr.	AC-007-201-T	AC-007-201-OT
		Ø 1.6 mm	5 Fr.	AC-007-202-T	AC-007-202-OT
	Grasping forceps for big stones, double action	Ø 1.2 mm	4 Fr.	AC-009-201-T	AC-009-201-OT
		Ø 1.6 mm	5 Fr.	AC-009-202-T	AC-009-202-OT

1.1 | CYSTO-URETHROSCOPY

INSTRUMENTS DETACHABLE | RIGID



Ordering No.		Description			
AC-100-001		Handle for detachable instruments			
	Description	Diameter	French	Flexible	Semi-rigid
	Scissors, single action	Ø 1.0 mm	3 Fr.	AC-103-200-T	AC-103-200-OT
		Ø 1.2 mm	4 Fr.	AC-103-201-T	AC-103-201-OT
		Ø 1.6 mm	5 Fr.	AC-103-202-T	AC-103-202-OT
		Ø 2.3 mm	7 Fr.	AC-103-204-T	AC-103-204-OT
		Ø 3.0 mm	9 Fr.	AC-103-205-T	AC-103-205-OT
	Scissors pointed, single action	Ø 1.6 mm	5 Fr.	AC-103-612-T	AC-103-612-OT
		Ø 2.3 mm	7 Fr.	AC-103-614-T	AC-103-614-OT
	Hook scissors, single action	Ø 2.3 mm	7 Fr.	AC-104-204-T	AC-104-204-OT
		Ø 3.0 mm	9 Fr.	AC-104-205-T	AC-104-205-OT
	Biopsy forceps oval, double action	Ø 1.0 mm	3 Fr.	AC-101-201-T	AC-101-201-OT
		Ø 1.2 mm	4 Fr.	AC-101-206-T	AC-101-206-OT
		Ø 1.6 mm	5 Fr.	AC-101-202-T	AC-101-202-OT
		Ø 2.3 mm	7 Fr.	AC-101-203-T	AC-101-203-OT
		Ø 3.0 mm	9 Fr.	AC-101-204-T	AC-101-204-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0 mm	3 Fr.	AC-106-201-T	AC-106-201-OT
		Ø 1.2 mm	4 Fr.	AC-106-203-T	AC-106-203-OT
		Ø 1.6 mm	5 Fr.	AC-106-202-T	AC-106-202-OT
		Ø 2.3 mm	7 Fr.	AC-106-204-T	AC-106-204-OT
		Ø 3.0 mm	9 Fr.	AC-106-205-T	AC-106-205-OT
	Grasping forceps for stone street, double action	Ø 1.2 mm	4 Fr.	AC-107-201-T	AC-107-201-OT
		Ø 1.6 mm	5 Fr.	AC-107-202-T	AC-107-202-OT
	Grasping forceps for big stones, double action	Ø 1.2 mm	5 Fr.	AC-109-201-T	AC-109-201-OT

1.1 | CYSTO-URETHROSCOPY

STONE EXTRACTORS | 4 WIRES



Ordering No.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering No.	Wires	French	Working length	Description
AF-001-001	4	3 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-002	4	3 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-015	4	3 Fr.	90 cm	Dormia stone extractor, basket with tubus only
AF-001-016	4	3 Fr.	90 cm	Dormia stone extractor complete with handle
AF-001-023	4	3 Fr.	120 cm	Dormia stone extractor, basket with tubus only
AF-001-022	4	3 Fr.	120 cm	Dormia stone extractor complete with handle
AF-001-003	4	5 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-004	4	5 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-017	4	5 Fr.	90 cm	Dormia stone extractor, basket with tubus only
AF-001-019	4	5 Fr.	90 cm	Dormia stone extractor complete with handle
AF-001-005	4	7 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-006	4	7 Fr.	60 cm	Dormia stone extractor complete with handle

1.1 | CYSTO-URETHROSCOPY

STONE EXTRACTORS | 5-7 WIRES



Ordering No.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering No.	Wires	French	Working length	Description
AF-001-007	5	5 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-008	5	5 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-025	5	5 Fr.	90 cm	Dormia stone extractor, basket with tubus only
AF-001-018	5	5 Fr.	90 cm	Dormia stone extractor complete with handle
AF-001-009	5	7 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-010	5	7 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-011	6	5 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-012	6	5 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-013	7	7 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-014	7	7 Fr.	60 cm	Dormia stone extractor complete with handle

1.1 | CYSTO-URETHROSCOPY

MONOPOLAR FLEXIBLE ELECTRODES




Description		French	Length	Ordering No.	French	Length	Ordering No.
	Loop electrode flexible, monopolar	5 Fr.	53 cm	AD-003-003	5 Fr.	65 cm	AD-003-013
		7 Fr.	53 cm	AD-003-005	5 Fr.	65 cm	AD-003-015
		9 Fr.	53 cm	AD-003-007	5 Fr.	65 cm	AD-003-017
		10 Fr.	53 cm	AD-003-008	5 Fr.	65 cm	AD-003-018
	Needle electrode flexible, monopolar	3 Fr.	53 cm	AD-002-001	3 Fr.	65 cm	AD-002-011
		4 Fr.	53 cm	AD-002-002	4 Fr.	65 cm	AD-002-012
		5 Fr.	53 cm	AD-002-003	5 Fr.	65 cm	AD-002-013
		7 Fr.	53 cm	AD-002-005	7 Fr.	65 cm	AD-002-015
	Button electrode flexible, monopolar	3 Fr.	53 cm	AD-001-001	3 Fr.	65 cm	AD-001-011
		4 Fr.	53 cm	AD-001-002	4 Fr.	65 cm	AD-001-012
		5 Fr.	53 cm	AD-001-003	5 Fr.	65 cm	AD-001-013
		6 Fr.	53 cm	AD-001-004	6 Fr.	65 cm	AD-001-014
		7 Fr.	53 cm	AD-001-005	7 Fr.	65 cm	AD-001-015
		8 Fr.	53 cm	AD-001-006	8 Fr.	65 cm	AD-001-016
		9 Fr.	53 cm	AD-001-007	9 Fr.	65 cm	AD-001-017
		10 Fr.	53 cm	AD-001-008	10 Fr.	65 cm	AD-001-018

1.1 | CYSTO-URETHROSCOPY

MONOPOLAR CONNECTING CABLES



Connection instrument side



Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001	3 m	Ø 4.0 mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5 m		
	AD-010-002	3 m	Ø 5.0 mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5 m		
	AD-010-003	3 m	Ø 8.0 mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5 m		
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable male for STEMA MASTER



Ordering No. Description

EK-015-060

Universal adapter used for monopolar HF cables
(For cables: AD-010-001 to AD-010-008)

1.2 | SLIM CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS



	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-002
	HD Cystoscope autoclavable			AA-003-012
	HD Cystoscope autoclavable wide angle			AA-003-012-W



Ordering No.	Description
AB-001-007	Cysto-Urethroscope sheath 15 Fr. with 2 fixed stopcocks, 5 Fr. instrument channel , incl. obturator




Ordering No.	Description
AB-001-000	Cysto-Urethroscope sheath 16 Fr. with 2 fixed stopcocks, Quick-Lock , incl. obturator

1.2 | SLIM CYSTO-URETHROSCOPY

ENDOSCOPE BRIDGES | QUICK-LOCK SYSTEM


Quick-Lock System



Ordering No.	Description
AB-003-010	Endoscope bridge diagnostic used for 16 Fr. sheath



Ordering No.	Description
AB-003-011	Endoscope bridge with 1 instrument channel used for 16 Fr. sheath



Ordering No.	Description
AB-003-004	Endoscope bridge with 2 instrument channels



Ordering No.	Description
AB-003-012	Endoscope bridge with 2 instrument channels used for 16 Fr. sheath



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6–0.8 mm

1.2 | SLIM CYSTO-URETHROSCOPY WORKING ELEMENTS



Ordering No.	Description
AB-003-013	Working element with ALBARRAN deflector and 1 instrument channel used for 16 Fr. sheath
AB-003-014	Working element with ALBARRAN deflector and 2 instrument channels used for 16 Fr. sheath



Ordering No.	Description
AB-003-020	Protective rod for working elements with ALBARRAN deflector

1.3 | LASER CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-003
	HD Cystoscope autoclavable			AA-002-013
	ULTRA HD (4K) Cystoscope autoclavable			AA-002-023



Ordering No.	Description
AB-001-040	LASER Cysto-Urethroscope outer sheath 23 Fr. with 2 fixed stopcock used for 30° endoscope



Ordering No.	Description
AB-001-041	LASER Cysto-Urethroscope inner sheath with 8 Fr. working channel

1.3 | LASER CYSTO-URETHROSCOPY OBTURATORS



Ordering No.	Description
AB-001-043	Obturator for LASER cysto-urethroscope sheath AB-001-040



Ordering No.	Description
AB-001-042	Viewing obturator for LASER cysto-urethroscope sheath AB-001-040

1.4 | HD MINI COMPACT CYSTOSCOPE



Ordering No.	Description
AA-007-060	HD Mini Compact Cystoscope

Technical data

Eyepiece	30° angled ocular
Sheath	7.5/9/12 Fr.
Direction of view	0° viewing angle
Working length	220 mm
Instrument channel.....	6.8 Fr. use with 1 × 5 or 2 × 3 Fr. instruments
Luer Lock	3 connections
Optic	40,000 px

Delivery contains

HD Mini compact cystoscope

1.5 | FLEXIBLE CYSTOSCOPE FLEX I

- High luminance
- Exact adjustment of focus
- Precise up-/down angulation



Ordering No.	Description
AA-008-002	Flexible Cystoscope FLEX I

Technical data

Cystoscope

Insertion tube	Ø 5.5 mm
Working Channel	Ø 2.3 mm
Working length	380 mm
Field of view	100°
Pixel	18,000
Angulation up-down	210/120°
Accessories	Suction valve, leak tester

Delivery contains

Flexible cystoscope FLEX I, leak tester, suction valve

INCLUDED IN SCOPE OF DELIVERY:

AA-008-103

Leak tester

AA-008-105

Suction valve



AA-008-103



AA-008-105

1.5 | FLEXIBLE CYSTOSCOPE ACCESSORIES



Ordering No.	Description
AA-008-100	Bridge with Luer-Lock and 2 stopcocks for AA-008-002



Ordering No.	Description
AA-008-102	Suitcase grey for STEMA flexible scopes



Ordering No.	Description
AA-008-103	Leak tester



Ordering No.	Description
AA-008-104	Single use Luer-Lock adapter, 1 box = 25 pieces



Ordering No.	Description
AA-008-105	Suction valve

1.6 CYSTO-URETHROSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0 mm, 30° viewing angle, autoclavable,
302 mm working length

1 AA-002-013



2 CYSTO-URETHROSCOPE SHEATH

21 Fr. with 2 fixed stopcock and obturator

2 AB-001-004

3 CYSTO-URETHROSCOPE SHEATH

25 Fr. with 2 fixed stopcock and obturator

3 AB-001-006



4 ENDOSCOPE BRIDGE

with 1 instrument channel

4 AB-003-003



5 WORKING ELEMENT

with Albaran deflector and 2 instrument channels

5 AB-003-007



6 BIOPSY FORCEPS

flexible, oval, double action, Ø 1.6 mm,
WL 40 cm

6 AC-001-002



7 SCISSORS

flexible, single action, Ø 1.6 mm, WL 40 cm

7 AC-003-002



8 GRASPING FORCEPS ALLIGATOR

flexible, double action, Ø 2.3 mm, WL 40 cm

8 AC-006-004



1.6 CYSTO-URETHROSCOPY SET

9 BUTTON ELECTRODE

flexible, 5 Fr., monopolar, length 530 mm



9 AD-001-003

10 NEEDLE ELECTRODE

flexible, 5 Fr., monopolar, length 530 mm



10 AD-002-003

11 HF MONOPOLAR CABLE

HF monopolar connecting cable for Berchtold, Martin, Aesculap, Ø 4 mm



11 AD-010-003

12 DORMIA STONE EXTRACTOR

5 wires, Ø 1.6 mm, complete with handle, 60 cm Working length



12 AF-001-008

13 BIOPSY FORCEPS

double action



13 AB-004-001

14 GRASPING FORCEPS

double action



14 AB-004-003

15 OPTICAL STONE PUNCH



15 AB-004-007

2. RESECTOSCOPY

2.1 | MONOPOLAR RESECTOSCOPY

24/26/28.5 FRENCH
19/22 FRENCH (SLIM)

2.2 | SALINE RESECTOSCOPY

24/26 FRENCH
19/22 FRENCH (SLIM)
17.5 FRENCH (MINI)

2.3 | LASER RESECTOSCOPY

24/26 FRENCH
19/22 FRENCH (SLIM)
17.5 FRENCH (MINI)

2.4 | ACCESSORIES

2.5 | RESECTOSCOPY SETS

2. RESECTOSCOPY

24/26 Fr.
28.5 Fr.

19/22 Fr.


17.5 Fr.




1 CYSTOSCOPES



2 WORKING ELEMENTS




3 SHEATHS & OBTURATORS



4 ELECTRODES



2.1 | MONOPOLAR RESECTOSCOPY CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 0°	Cystoscope autoclavable Cystoscope autoclavable wide angle	Ø 4.0 mm	302 mm	AA-002-001 AA-002-001-W
	HD Cystoscope autoclavable HD Cystoscope autoclavable wide angle			AA-002-011 AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-021

Viewing angle	Description	Diameter	Working length	Ordering No.
 12°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-002
	HD Cystoscope autoclavable			AA-002-012
	Ultra HD (4K) Cystoscope autoclavable			AA-002-022

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable Cystoscope autoclavable wide angle	Ø 4.0 mm	302 mm	AA-002-003 AA-002-003-W
	HD Cystoscope autoclavable HD Cystoscope autoclavable wide angle			AA-002-013 AA-002-013-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-023


Viewing angle	Description	Diameter	Working length	Ordering No.
 70°	Cystoscope autoclavable Cystoscope autoclavable wide angle	Ø 4.0 mm	302 mm	AA-002-004 AA-002-004-W
	HD Cystoscope autoclavable HD Cystoscope autoclavable wide angle			AA-002-014 AA-002-014-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-024

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS



Ordering No.	Description
AG-001-001	Working element used for resectoscopes, passive cutting , handle closed
AG-001-001-HF	Working element used for resectoscopes, passive cutting , handle closed, HF connection on top




Ordering No.	Description
AG-001-002	Working element used for resectoscopes, passive cutting , handle open

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS



Ordering No.	Description
AG-001-003	Working element used for resectoscopes, active cutting , handle closed
AG-001-003-HF	Working element used for resectoscopes, active cutting , handle closed, HF connection on top



Ordering No.	Description
AG-001-004	Working element used for resectoscopes, active cutting , handle open, cutting spring assisted

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | TITAN



Ordering No.	Description
AG-001-007	Working element used for resectoscopes, passive cutting , handle closed, titanium version



Ordering No.	Description
AG-001-008	Working element used for resectoscopes, active cutting , handle closed, titanium version

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | GYRUS



Ordering No.	Description
AG-001-017	GYRUS Saline working element used for resectoscopes, active cutting , handle closed, titanium version



Ordering No.	Description
AG-001-018	GYRUS Saline working element used for resectoscopes, passive cutting , handle closed, titanium version

2.1 | MONOPOLAR RESECTOSCOPY SHEATHS




	Ordering No.	Description	French	Colour code
	AH-001-202	Outer sheath with perforation	26 Fr.	Yellow
	AH-001-212	Outer sheath with perforation	28.5 Fr.	Black
	AH-001-204	Outer sheath with upper/lower perforation	26 Fr.	Yellow
	AH-001-203	Outer sheath with surrounding fenestrated perforation	26 Fr.	Yellow
	AH-001-205	Outer sheath with upper fenestrated perforation	26 Fr.	Yellow



Ordering No.	Description	French	Colour code
AH-001-201	Inner sheath with ceramic isolation	26 Fr.	Yellow
AH-001-211	Inner sheath with ceramic isolation	28.5 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY SHEATHS



	Ordering No.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow
	AH-001-312	Outer sheath with perforation, rotatable	28.5 Fr.	Black
	AH-001-304	Outer sheath with upper/lower perforation, rotatable	26 Fr.	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow



Ordering No.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow
AH-001-311	Inner sheath with ceramic isolation, rotatable	28.5 Fr.	Black



Ordering No.	Description	French	Colour code
AH-001-308	Inner sheath without ceramic isolation, with laser channel up to 7 Fr. , for use with AH-001-301	26 Fr.	Yellow

2.1 | MONOPOLAR RESECTOSCOPY SHEATHS



Ordering No.	Description	French	Colour code
AH-001-001	Standard resectoscope sheath with 1 fixed stopcock	24 Fr.	Yellow
AH-001-002	Standard resectoscope sheath with 1 fixed stopcock	27 Fr.	Black



Ordering No.	Description	French	Colour code
AH-001-101	Standard resectoscope sheath with central valve	24 Fr.	Yellow
AH-001-102	Standard resectoscope sheath with central valve	27 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY OBTURATORS



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-011	Standard obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-012	Standard obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-015	Deflecting obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-016	Deflecting obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-013	Viewing obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-017	Viewing obturator for resectoscope sheaths with 5 Fr. Instrument channel	24 Fr.	26 Fr.	Yellow
AH-001-018	Viewing obturator for resectoscope sheaths with 5 Fr. Instrument channel	27 Fr.	28.5 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-001-001	Cutting loop electrode 30° angled	24 Fr.	yellow
	AE-002-001		27 Fr.	brown
	AE-001-002	Cutting loop electrode angled	24 Fr.	yellow
	AE-002-002		27 Fr.	brown
	AE-001-003	Cutting loop electrode straight	24 Fr.	yellow
	AE-002-003		27 Fr.	brown
	AE-001-004	Rollerball electrode Ø 3 mm	24 Fr.	yellow
	AE-002-004		27 Fr.	brown
	AE-001-005	Rollerball electrode Ø 5 mm	24 Fr.	yellow
	AE-002-005		27 Fr.	brown
	AE-001-010	Roller electrode Ø 3 mm	24 Fr.	yellow
	AE-002-010		27 Fr.	brown
	AE-001-011	Roller electrode Ø 5 mm	24 Fr.	yellow
	AE-002-011		27 Fr.	brown
	AE-001-012	Spike electrode Ø 3 mm	24 Fr.	yellow
	AE-002-012		27 Fr.	brown
	AE-001-013	Spike electrode Ø 5 mm	24 Fr.	yellow
	AE-002-013		27 Fr.	brown

2.1 | MONOPOLAR RESECTOSCOPY ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-001-014	Vaporisation electrode Ø 3 mm	24 Fr.	yellow
	AE-002-014		27 Fr.	brown
	AE-001-015	Vaporisation electrode Ø 5 mm	24 Fr.	yellow
	AE-002-015		27 Fr.	brown
	AE-001-009	Conical electrode	24 Fr.	yellow
	AE-001-016	Knife electrode	24 Fr.	yellow
	AE-002-016		27 Fr.	brown
	AE-001-006	Band electrode angled	24 Fr.	yellow
	AE-002-006		27 Fr.	brown
	AE-001-007	Band electrode 30° angled	24 Fr.	yellow
	AE-002-007		27 Fr.	brown
	AE-001-008	Band vaporization electrode angled	24 Fr.	yellow
	AE-001-017	Curette sharp, round (not for use with HF)	24 Fr.	yellow



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.1 | MONOPOLAR SLIM RESECTOSCOPY CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-002
	HD Cystoscope autoclavable			AA-003-012

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS



Ordering No.	Description
AG-002-001	Working element used for slim resectoscopes, passive cutting, handle closed




Ordering No.	Description
AG-002-002	Working element used for slim resectoscopes, passive cutting, handle open

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS




Ordering No.	Description
AG-002-003	Working element used for slim resectoscopes, passive cutting, handle closed



Ordering No.	Description
AG-002-004	Working element used for slim resectoscopes, active cutting, handle open

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS | TITAN



Ordering No.	Description
AG-002-006	Working element used for slim resectoscopes, active cutting , handle closed, titanium version

2.1 |

MONOPOLAR SLIM RESECTOSCOPY

SHEATHS




Ordering No.	Description	French	Colour code
AH-002-212	Outer sheath with perforation, fixed	22 Fr.	White



Ordering No.	Description	French	Colour code
AH-002-201	Inner sheath with ceramic isolation, fixed	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-302	Outer sheath with perforation, rotatable	22 Fr.	White



Ordering No.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, rotatable	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY SHEATHS



Ordering No.	Description	French	Colour code
AH-002-001	Standard resectoscope sheath, 1 fixed stopcock , without obturator	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-003-001	Cutting loop electrode 30° angled	19 Fr.	White
	AE-003-003	Cutting loop electrode straight	19 Fr.	White
	AE-003-004	Rollerball electrode Ø 3 mm	19 Fr.	White
	AE-003-010	Roller electrode Ø 3 mm	19 Fr.	White
	AE-003-014	Vaporisation electrode Ø 3 mm	19 Fr.	White
	AE-003-016	Knife electrode	19 Fr.	White



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.1 | MONOPOLAR SLIM RESECTOSCOPY OBTURATORS



Ordering No.	Description	French	Colour code
AH-002-011	Standard obturator	19 Fr.	White



Ordering No.	Description	French	Colour code
AH-002-013	Viewing obturator	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY

MONOPOLAR CONNECTING CABLES



Connection instrument side



Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001	3 m	Ø 4.0 mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5 m		
	AD-010-002	3 m	Ø 5.0 mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5 m		
	AD-010-003	3 m	Ø 8.0 mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5 m		
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable male for STEMA MASTER



Ordering No.	Description
EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

2.1 | MONOPOLAR SLIM RESECTOSCOPY

MONOPOLAR CONNECTING CABLES



Connection unit side



Unit side	Ordering No.	Length	Description
	AD-010-025	4.5 m	HF Monopolar connecting cable for Storz resectoscopes with 90° angled connector
	AD-010-026	4.5 m	HF Monopolar connecting cable for Olympus resectoscopes, diameter Ø 3 mm



RECOMMENDED PRODUCTS:

- CI-001-030** **STEMA HF-U PRO TOUCH**
- *For all monopolar & bipolar techniques*
 - *with 7" touch screen display*
 - *provided to perform the most advanced electrosurgery*
 - *ADC / APC System*
- CI-001-040** **DOUBLE PEDAL FOOTSWITCH**

2.2 | SALINE RESECTOSCOPY CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 12°	HD Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-012

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-003
	HD Cystoscope autoclavable			AA-002-013
	ULTRA HD Cystoscope autoclavable			AA-002-023

2.2 | SALINE RESECTOSCOPY WORKING ELEMENTS





Ordering No.	Description
AG-001-011	Saline working element used for resectoscopes, 24 Fr., passive cutting , handle closed



Ordering No.	Description
AG-001-012	Saline working element used for resectoscopes, 24 Fr., active cutting , handle closed

2.2 | SALINE RESECTOSCOPY SHEATHS



	Ordering No.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow
	AH-001-312	Outer sheath with perforation, rotatable	28.5 Fr.	Black
	AH-001-304	Outer sheath with upper/lower perforation, rotatable	26 Fr.	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow



Ordering No.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow
AH-001-311	Inner sheath with ceramic isolation, rotatable	28.5 Fr.	Black

2.2 | SALINE RESECTOSCOPY OBTURATORS



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-011	Standard obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-012	Standard obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-015	Deflecting obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-016	Deflecting obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-013	Viewing obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-017	Viewing obturator for resectoscope sheaths with 5 Fr. Instrument channel	24 Fr.	26 Fr.	Yellow
AH-001-018	Viewing obturator for resectoscope sheaths with 5 Fr. Instrument channel	27 Fr.	28.5 Fr.	Black


2.2 | SALINE RESECTOSCOPY ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-001-030	Cutting loop electrode 30° angled, strong cut	24 Fr.	yellow/blue
	AE-001-031	Cutting loop electrode 30° angled	24 Fr.	yellow/blue
	AE-001-033	Cutting loop electrode 90° angled	24 Fr.	yellow/blue
	AE-001-048	Cutting loop electrode 90° angled, for 0° Scope	24 Fr.	yellow/blue
	AE-001-036	Cutting loop electrode straight	24 Fr.	yellow/blue
	AE-001-032	Knife electrode	24 Fr.	yellow/blue
	AE-001-041	Band electrode, 90° angled	24 Fr.	yellow/blue
	AE-001-038	Vaporisation spike rollerball electrode Ø 3 mm	24 Fr.	yellow/blue
	AE-001-037	Vaporisation rollerball electrode Ø 3 mm	24 Fr.	yellow/blue
	AE-001-043	Vaporisation rollerball electrode Ø 5 mm		
	AE-001-040	Roller electrode smooth Ø 3 mm	24 Fr.	yellow/blue
	AE-001-042	Roller electrode smooth Ø 5 mm		

2.2 | SALINE RESECTOSCOPY ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-001-034	Rollerball electrode Ø 3 mm	24 Fr.	yellow/blue
	AE-001-035	Rollerball electrode Ø 5 mm		
	AE-001-044	Conical electrode „Mushroom“	24 Fr.	yellow/blue
	AE-001-045	Conical electrode „Mushroom“, convex	24 Fr.	yellow/blue
	AE-001-046	Conical electrode „Mushroom“, concave	24 Fr.	yellow/blue
	AE-001-047	Conical electrode „Mushroom“, grooved	24 Fr.	yellow/blue
	AE-001-039	Conical electrode	24 Fr.	yellow/blue

	Ordering No.	Description	French	Colour code
	AE-002-031	Cutting loop electrode, 30° angled	27 Fr.	yellow/blue



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.2 | SALINE SLIM RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 12°	HD Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-004
				AA-003-014

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-002
	HD Cystoscope autoclavable			AA-003-012

2.2 | SALINE SLIM RESECTOSCOPY WORKING ELEMENTS




Ordering No.	Description
AG-002-011	Saline working element used for resectoscopes, 19 Fr., passive cutting , handle closed



Ordering No.	Description
AG-002-012	Saline working element used for resectoscopes, 19 Fr., active cutting , handle closed

2.2 | SALINE SLIM RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-302	Outer sheath with perforation, rotatable	22 Fr.	White



Ordering No.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, rotatable	19 Fr.	White

2.2 | SALINE SLIM RESECTOSCOPY OBTURATORS



Ordering No.	Description
AH-002-011	Standard obturator for resectoscope sheaths 19 Fr.



Ordering No.	Description
AH-002-013	Viewing obturator for resectoscope sheaths 19 Fr.

2.2 | SALINE SLIM RESECTOSCOPY ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-003-031	Cutting loop electrode 30° angled	19 Fr.	white/blue
	AE-003-032	Cutting loop electrode 90° angled	19 Fr.	white/blue
	AE-003-033	Cutting loop electrode straight	19 Fr.	white/blue
	AE-003-034	Rollerball electrode Ø 3 mm	19 Fr.	white/blue
	AE-003-040	Roller electrode Ø 3 mm	19 Fr.	white/blue
	AE-003-036	Knife electrode	19 Fr.	white/blue



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives



17.5 Fr.

2.2 | SALINE MINI RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-002
	HD Cystoscope autoclavable			AA-003-012

2.2 | SALINE MINI RESECTOSCOPY

WORKING ELEMENTS



Ordering No.	Description
AG-002-015	Saline working element used for resectoscopes, 18.5 Fr., passive cutting , handle closed



Ordering No.	Description
AG-002-013	Saline working element used for resectoscopes, 18.5 Fr., active cutting , handle closed



17.5 Fr.

2.2 | SALINE MINI RESECTOSCOPY
SHEATHS



Ordering No.	Description	French	Colour code
AH-002-213	Outer sheath 18.5 Fr. with perforation	18.5 Fr.	green



Ordering No.	Description
AH-002-211	Inner sheath 17.5 Fr. with ceramic isolation, rotatable



Ordering No.	Description
AH-002-216	Inner sheath 17.5 Fr. without ceramic isolation , with laser channel up to 5 Fr., for use with AH-001-213

2.2 | SALINE MINI RESECTOSCOPY OBTURATORS



Ordering No.	Description
AH-002-015	Standard obturator 17.5 Fr.



Ordering No.	Description
AH-002-016	Viewing obturator 17.5 Fr. with 5 Fr. instrument channel

2.2 | SALINE MINI RESECTOSCOPY ELECTRODES




	Ordering No.	Description	French	Colour code
	AE-003-041	Cutting loop electrode 30° angled	17.5 Fr.	green
	AE-003-043	Cutting loop electrode straight	17.5 Fr.	green
	AE-003-044	Rollerball electrode Ø 2 mm	17.5 Fr.	green
	AE-003-045	Roller electrode Ø 2 mm	17.5 Fr.	green
	AE-003-046	Knife electrode	17.5 Fr.	green



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.2 | SALINE RESECTOSCOPY

BIPOLAR CONNECTING CABLES



Connection instrument side



Unit side	Ordering No.	Length	Description
	AD-010-011	4.5 m	Bipolar resectoscope connecting cable for Martin Maxium
	AD-010-012	4.5 m	Bipolar resectoscope connecting cable for Erbe Vio
	AD-010-013	4.5 m	Bipolar resectoscope connecting cable for Valleylab, Lamidey, EMC
	AD-010-015	4.5 m	Bipolar resectoscope connecting cable for Erbe Vio 3

2.3 | LASER RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	HD Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-013

2.3 | LASER RESECTOSCOPY WORKING ELEMENTS



Ordering No.	Description
AG-001-006	Working element used for resectoscopes, titanium, passive cutting, handle closed, 24/26 Fr. with probe channel 0.8 mm
AG-001-009	Working element used for resectoscopes, titanium, passive cutting, handle closed, 24/26 Fr. with probe channel 1.2 mm

Ordering No.	Description
AG-100-006	Adapter used for laser probes
AG-100-007	TUOHY-Borst adapter

2.3 | LASER RESECTOSCOPY SHEATHS | OBTURATORS



Ordering No.	Description	French	Colour code
AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow



Ordering No.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-308	Inner sheath without ceramic isolation, with laser channel up to 7 Fr., for use with AH-001-301	24 Fr.	26 Fr.	Yellow



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-011	Standard obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow

2.3 | LASER SLIM RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	HD Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-012


2.3 | LASER SLIM RESECTOSCOPY WORKING ELEMENTS



Ordering No.	Description
AG-002-007	Working element used for resectoscopes, titanium, passive cutting, handle closed, 19/22 Fr. with probe channel 0.8 mm
AG-002-008	Working element used for resectoscopes, titanium, passive cutting, handle closed, 19/22 Fr. with probe channel 1.2 mm

2.3 | LASER SLIM RESECTOSCOPY SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-302	Outer sheath with perforation, rotatable	22 Fr.	White



Ordering No.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, rotatable	19 Fr.	White

2.3 | LASER SLIM RESECTOSCOPY OBTURATORS



Ordering No.	Description
AH-002-011	Standard obturator for resectoscope sheaths 19 Fr.



Ordering No.	Description
AH-002-013	Viewing obturator for resectoscope sheaths 19 Fr.

2.3 | LASER MINI RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	HD Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-012

2.3 | LASER MINI RESECTOSCOPY WORKING ELEMENTS



Ordering No.	Description
AG-002-016	Working element used for resectoscopes, titanium, passive cutting, handle closed, 17.5 Fr. with probe channel 0.8 mm

2.3 | LASER MINI RESECTOSCOPY SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-213	Outer sheath 18.5 Fr. with perforation	18.5 Fr.	green



Ordering No.	Description
AH-002-211	Inner sheath 17.5 Fr. with ceramic isolation, rotatable



Ordering No.	Description
AH-002-216	Inner sheath 17.5 Fr. without ceramic isolation , with laser channel up to 5 Fr., for use with AH-001-213

2.3 | LASER MINI RESECTOSCOPY
OBTURATORS



Ordering No.	Description
AH-002-015	Standard obturator 17.5 Fr.



Ordering No.	Description
AH-002-016	Viewing obturator 17.5 Fr. with 5 Fr. instrument channel

2.4 | ACCESSORIES


SUPRAPUBIC ASPIRATION



Ordering No.	Description
AB-007-001	Suprapubic aspiration set, complete set



Ordering No.	Description
AB-007-002	Trocarn sleeve 16 Fr. with pressure measuring channel



Ordering No.	Description
AB-007-003	Adjustable stabilisation disc, lateral opening



Ordering No.	Description
AB-007-008	Guide sleeve, lateral opening , for inserting a suprapubic catheter

2.4 | ACCESSORIES

SUPRAPUBIC ASPIRATION



Ordering No.	Description
AB-007-004	Trocar with cannula tip



Ordering No.	Description
AB-007-005	Trocar with pyramidal tip



Ordering No.	Description
AB-007-006	Aspiration cannula



Ordering No.	Description
AB-007-007	Overflow cannula for pressure stabilisation

2.4 | ACCESSORIES


EVACUATOR | BLADDER SYRINGES



Ordering No.	Description
AJ-002-001	Ellik evacuator complete with sheath connection and locking device
AJ-002-002	Bowl only
AJ-002-003	Suction Ball only
AJ-002-004	Cone with tube only
AJ-002-005	Sheath connection for Ellik only




Not available in the EU



Ordering No.	Description
AJ-001-000	RECORD syringe with adapter, 50 ml
AJ-001-010	Replacement cylinder for RECORD syringe with adapter, 50 ml
AJ-001-001	RECORD syringe with adapter, 75 ml
AJ-001-011	Replacement cylinder for RECORD syringe with adapter, 75 ml
AJ-001-002	RECORD syringe with adapter, 100 ml
AJ-001-012	Replacement cylinder for RECORD syringe with adapter, 100 ml
AJ-001-003	RECORD syringe with adapter, 150 ml
AJ-001-013	Replacement cylinder for RECORD syringe with adapter, 150 ml

2.4 | ACCESSORIES

TRAY INSERT



Ordering No.	Description
AJ-010-020	Tray insert for resectoscopy (without instruments, without sterilization box) , Dimension: 370 × 255 × 45 mm (L × W × H)



RECOMMENDED PRODUCTS:

- AA-002-013** HD autoclavable scope
30° | 4 mm
- AG-001-011** SALINE working element,
passive cutting, handle closed
- AH-001-302** Outer sheath 26 Fr. with perforation,
rotatable
- AH-001-301** Inner sheath 26 Fr. with ceramic isolation,
rotatable
- AH-001-011** Standard obturator for resectoscope
sheaths 24 Fr.
- AJ-010-006** Protection sleeve made of steel for sterilization
and storage of electrodes and cold knives
- AE-001-031** Saline Cutting loop electrode 30° angled, 24 Fr.
wire thickness 0.35 mm
- AE-001-034** Saline rollerball electrode 3 mm, 24 Fr.
+ connecting cable (depends on unit)

2.5 MONOPOLAR RESECTOSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0 mm, 30° viewing angle, autoclavable,
WL 302 mm

1 AA-002-013



2 WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-001



3 OUTER SHEATH

26 Fr. with perforation, rotatable

3 AH-001-302



4 INNER SHEATH

26 Fr. with ceramic isolation, rotatable

4 AH-001-301



5 STANDARD OBTURATOR

24 Fr.

5 AH-001-011



6 CUTTING LOOP ELECTRODE

30° angled, 24 Fr.

6 AE-001-001



7 ROLLERBALL ELECTRODE 3 MM

24 Fr.

7 AE-001-004



8 KNIFE ELECTRODE

24 Fr.

8 AE-001-016



2.5 MONOPOLAR RESECTOSCOPY SET

9 HF MONOPOLAR CONNECTING CABLE
for Berchtold, Martin, Aesculap, Ø 4 mm

9 AD-010-001



10 RECORD SYRINGE
100 ml

10 AJ-001-002



11 ELLIK EVACUATOR
complete

11 AJ-002-001



12 BIOPSY FORCEPS
double action

12 AB-004-001



13 GRASPING FORCEPS
double action

13 AB-004-003



14 Adapter
for use with resectoscope sheaths

14 AB-004-008



15 PROTECTION SLEEVE
For sterilization and storage of electrodes

15 AJ-010-006



2.5 SALINE RESECTOSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0 mm, 30° viewing angle, autoclavable,
WL 302 mm

1 AA-002-013



2 SALINE WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-011



3 OUTER SHEATH

26 Fr. with perforation, rotatable

3 AH-001-302



4 INNER SHEATH

26 Fr. with ceramic isolation, rotatable

4 AH-001-301



5 STANDARD OBTURATOR

24 Fr.

5 AH-001-011



6 CUTTING LOOP ELECTRODE

30° angled, 24 Fr.

6 AE-001-031



7 ROLLERBALL ELECTRODE 3 MM

24 Fr.

7 AE-001-034



8 KNIFE ELECTRODE

24 Fr.

8 AE-001-032



2.5 SALINE RESECTOSCOPY SET

- 9 HF BIPOLAR CONNECTING CABLE**
for Valleylab, Lamidey, EMC, 4.5 m

9 AD-010-013



- 10 RECORD SYRINGE**
100 ml

10 AJ-001-002



- 11 ELLIK EVACUATOR**
complete

11 AJ-002-001



- 12 BIOPSY FORCEPS**
double action

12 AB-004-001



- 13 GRASPING FORCEPS**
double action

13 AB-004-003



- 14 ADAPTER**
for use with resectoscope sheaths

14 AB-004-008



- 15 PROTECTION SLEEVE**
for sterilization and storage of electrodes

15 AJ-010-006



3. URETHROTOMY

3.1 | URETHROTOME

3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES

SHEATHS

OBTURATORS

WORKING ELEMENTS

KNIVES

3.3 | URETHROTOMY SET

3. URETHROTOMY

21 Fr.

15.5 Fr.


1 CYSTOSCOPES



2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS



4 KNIVES



3.1 | URETHROTOME



Ordering No.	Description
AI-002-001	OTIS urethrotome with 3 exchangeable tips inclusive 2 knives

SPARE KNIVES




Ordering No.	Description
AI-002-002	Spare knives for OTIS urethrotome

3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES | SHEATHS



Viewing angle	Description	Diameter	Working length	Ordering No.
 0°	Cystoscope autoclavable Cystoscope autoclavable wide angle	Ø 4.0 mm	302 mm	AA-002-001 AA-002-001-W
	HD Cystoscope autoclavable HD Cystoscope autoclavable wide angle			AA-002-011 AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-021



Ordering No.	Description
AI-001-001	Urethrotome sheath 21 Fr. with instrument channel 5 Fr.
AI-001-003	Endoscope bridge with 1 instrument channel , fit with urethrotome sheath AI-001-001

3.2 | OPTICAL URETHROTOMY SUPPLEMENTARY SHEATHS



Ordering No.	Description
AI-001-002	Supplementary sheath for AI-001-001, open to side



Ordering No.	Description
AI-001-006	Supplementary sheath 25 Fr. for continuous flow

3.2 | OPTICAL URETHROTOMY OBTURATORS



Ordering No.	Description
AI-001-004	Standard obturator for urethrotome sheath 21 Fr.



Ordering No.	Description
AI-001-005	Obturator for urethrotome sheath 21 Fr. with 5 Fr. instrument channel




Ordering No.	Description
AI-001-012	Viewing obturator for urethrotome sheath 21 Fr.

3.2 | OPTICAL URETHROTOMY

WORKING ELEMENTS | PASSIVE CUTTING



Ordering No.	Description
AG-001-001	Working element used for resectoscopes, passive cutting , handle closed
AG-001-001-HF	Working element used for resectoscopes, passive cutting , handle closed, HF connection on top



Ordering No.	Description
AG-001-002	Working element used for resectoscopes, passive cutting , handle open

3.2 | OPTICAL URETHROTOMY

COLD KNIVES



	Ordering No.	Description
	AI-001-007	Cold knife straight
	AI-001-008	Cold knife straight, wave shape
	AI-001-009	Cold knife round
	AI-001-011	Cold knife hook shaped




Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES | SHEATHS | OBTURATORS



Viewing angle	Description	Diameter	Working length	Ordering No.
 0°	Cystoscope autoclavable	ø 2.9 mm	302 mm	AA-003-001
	HD Cystoscope autoclavable			AA-003-011



AI-005-001

AI-005-003

Ordering No.	Description
AI-005-001	Urethrotome sheath 15.5 Fr. with instrument channel 3 Fr.
AI-005-003	Endoscope bridge with 1 instrument channel , fit with urethrotome sheath AI-005-001



Ordering No.	Description
AI-005-004	Standard obturator for urethrotome sheath 15.5 Fr.

3.2 | OPTICAL URETHROTOMY

WORKING ELEMENT



Ordering No.	Description
AG-001-010	Working element used for resectoscopes, passive cutting, handle closed , without HF-connection

3.2 | OPTICAL URETHROTOMY

COLD KNIVES



	Ordering No.	Description
	AI-005-007	Cold knife straight
	AI-005-008	Cold knife straight, wave shape
	AI-005-011	Cold knife hook shaped



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

3.3 OPTICAL URETHROTOMY SET

1 HD CYSTOSCOPE

Ø 4.0 mm, 0° viewing angle, autoclavable,
302 mm working length

1 AA-002-011



2 WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-001



3 URETHROTOME SHEATH

21 Fr. with instrument channel 5 Fr.

3 AI-001-001



4 SUPPLEMENTARY SHEATH

Supplementary sheath for AI-001-001,
open to side

4 AI-001-002



5 STANDARD OBTURATOR

21 Fr.

5 AI-001-004



6 ENDOSCOPE BRIDGE

with 1 instrument channel

6 AI-001-003



7 COLD KNIFE

straight

7 AI-001-007



8 PROTECTION SLEEVE

for sterilization and storage of electrodes

8 AJ-010-006



4. NEPHROSCOPY

4.1 | HD NEPHROSCOPY

HD NEPHROSCOPES
SHEATHS
OBTURATORS
HD MINI NEPHROSCOPE

4.2 | MINI NEPHROSCOPY

MINI NEPHROSCOPES
SHEATHS
DILATATORS
INSTRUMENTS

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 200 MM
SHEATHS
OBTURATORS

NEPHROSCOPES 250 MM
SHEATHS
OBTURATORS

4.4 | ACCESSORIES

TELESCOPE BOUGIES
FORCEPS
CANNULAS

4. NEPHROSCOPY

1 HD NEPHROSCOPES



2 SHEATHS & OBTURATORS



3 FORCEPS, INSTRUMENTS, KNIVES



4.1 | HD NEPHROSCOPY

HD NEPHROSCOPES | SHEATHS | OBTURATORS | SEALING CAPS



Ordering No.	Description	Working length	Viewing angle
AA-006-069	HD Nephroscope - 45° angled ocular	220 mm	0°
AA-006-070	- Ø 18.8 Fr. - Working channel Ø 12.9 Fr.		15°



Ordering No.	Description	French
AB-006-170	Nephroscope sheath for continuous flow, rotatable stopcock (fit with nephrosopes WL 220mm / scope 15°)	23 Fr.



Ordering No.	Description	French
AB-006-171	Obturator hollow (fit with nephroscope sheath AB-006-170 / scope 15°)	23 Fr.


Ordering No.	Description
AA-006-081	Sealing cap, Ø 15 mm/1.0 mm
AA-006-083	Sealing cap, Ø 20 mm/2.0 mm

4.1 | HD MINI NEPHROSCOPE



Ordering No.	Description	Working length	Viewing angle
AA-006-100	HD Mini compact nephroscope - 30° angled ocular - Sheath 12 Fr. , - Instrument channel 6.6 Fr. , - 3 luer lock connections - 40,000 px	220 mm	0°

4.2 | MINI NEPHROSCOPY
MINI NEPHROSCOPES




Ordering No.	Description	Working length	French	Viewing angle
AA-006-050	Mini nephroscope with 45° angled ocular, instrument channel 10 Fr.	220 mm	Ø 19 Fr.	8°
AA-006-051	Mini nephroscope with 45° angled ocular, instrument channel 6 Fr.	220 mm	Ø 13.5 Fr.	8°
AA-006-052	Mini nephroscope with 45° angled ocular, instrument channel 6 Fr.	320 mm	Ø 13.5 Fr.	8°



Ordering No.	Description
AA-006-053	Instrument bridge used for Mini-Nephroscopes (not included in mini Nephroscope)
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

4.2 | MINI NEPHROSCOPY

SHEATHS | DILATATORS | APPLICATOR



Ordering No.	Description	French
AB-006-100	Mini nephroscope sheath without perforation, continuous flow, stopcock fixed	15/18 Fr.
AB-006-101		16.5/19.5 Fr.
AB-006-102		21/24 Fr.




Ordering No.	Description	French
AB-006-110	Dilatator with central cannulation Ø 3 Fr.	15/18 Fr.
AB-006-111	Dilatator with 2 channels Ø 2x3 Fr.	16.5/19.5 Fr.
AB-006-112	Dilatator with 2 channels Ø 2x5 Fr.	21/24 Fr.



Ordering No.	Description
AB-006-115	Applicator used for the mini nephroscope sheaths, incl. plunger



4.2 | MINI NEPHROSCOPY


INSTRUMENTS | FORCEPS FOR MINI NEPHROSCOPES




	Ordering No.	Description	Length	Diameter
	AB-006-242	Biopsy forceps serrated with ring handle	400 mm	Ø 1.8 mm
	AB-006-342		500 mm	Ø 1.8 mm
	AB-006-442		400 mm	Ø 3.0 mm
	AB-006-243	Biopsy forceps, single action with ring handle	400 mm	Ø 1.8 mm
	AB-006-343		500 mm	Ø 1.8 mm
	AB-006-443		400 mm	Ø 3.0 mm
	AB-006-244	Alligator grasping forceps with ring handle	400 mm	Ø 1.8 mm
	AB-006-344		500 mm	Ø 1.8 mm
	AB-006-444		400 mm	Ø 3.0 mm
	AB-006-245	Scissors, single action with ring handle	400 mm	Ø 1.8 mm
	AB-006-345		500 mm	Ø 1.8 mm
	AB-006-445		400 mm	Ø 3.0 mm
	AB-006-246	3-Prongs retriever for larger stones with coil handle	400 mm	Ø 1.8 mm
	AB-006-346		500 mm	Ø 1.8 mm
	AB-006-446		400 mm	Ø 3.0 mm

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 200 MM



Ordering No.	Description	Working length	Viewing angle
AA-006-001	Nephroscope - 45° angled ocular	200 mm	0°
AA-006-021	- Ø 7.45 mm - Suction channel Ø 1.5 mm		6°
AA-006-004	- Working channel Ø 5.0 mm		15°




Ordering No.	Description	Working length	Viewing angle
AA-006-002	Nephroscope - 90° angled ocular	200 mm	0°
AA-006-022	- Ø 7.45 mm - Suction channel Ø 1.5 mm		6°
AA-006-005	- Working channel Ø 5.0 mm		15°




Ordering No.	Description
AA-006-080	Flow valve, Ø 6.1 mm, used for STEMA nephroscopes
AA-006-081	Sealing cap, Ø 15 mm/1.0 mm, for AA-006-080
AA-006-083	Sealing cap, Ø 20 mm/2.0 mm, for AA-006-001–AA-006-034
AJ-010-015	Stopcock left, large (passage 4 mm), with LL-connection and tube attachment M6 × 0.35 mm

4.3 | STANDARD NEPHROSCOPY

SHEATHS | OBTURATORS



Ordering No.	Description	French
AB-006-001	Nephroscope sheath, without perforation, for continuous flow, rotatable stopcock, 200 mm scope	24 Fr.
AB-006-007		26 Fr.
AB-006-005	Nephroscope sheath, without perforation, for continuous flow, fixed stopcock, 200 mm scope	24 Fr.
AB-006-006		26 Fr.




Ordering No.	Description	French
AB-006-003	Nephroscope sheath for continuous flow, perforated, rotatable stopcock, 200 mm scope	26 Fr.
AB-006-008	Nephroscope sheath for continuous flow, perforated, fixed stopcock, 200 mm scope	26 Fr.




Ordering No.	Description	French
AB-006-002	Obturator hollow used for nephroscopy sheaths, 200 mm scope	24 Fr.
AB-006-004		26 Fr.

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 250 MM



Ordering No.	Description	Working length	Viewing angle
AA-006-011	Nephroscope - 45° angled ocular	250 mm	0°
AA-006-031	- Ø 7.45 mm - Suction channel Ø 1.5 mm		6°
AA-006-034	- Working channel Ø 5.0 mm		15°



Ordering No.	Description	Working length	Viewing angle
AA-006-012	Nephroscope - 90° angled ocular	250 mm	0°
AA-006-032	- Ø 7.45 mm - Suction channel Ø 1.5 mm		6°
AA-006-006	- Working channel Ø 5.0 mm		15°




Ordering No.	Description
AA-006-080	Flow valve, Ø 6.1 mm, used for STEMA nephroscopes
AA-006-081	Sealing cap, Ø 15 mm/1.0 mm, for AA-006-080
AA-006-083	Sealing cap, Ø 20 mm/2.0 mm, for AA-006-001–AA-006-034
AJ-010-015	Stopcock left, large (passage 4 mm), with LL-connection and tube attachment M6 × 0.35 mm

4.3 | STANDARD NEPHROSCOPY

SHEATHS | OBTURATORS



Ordering No.	Description	French
AB-006-061	Nephroscope sheath, without perforation, for continuous flow, rotatable stopcock, 250 mm scope	24 Fr.
AB-006-067		26 Fr.
AB-006-065	Nephroscope sheath, without perforation, for continuous flow, fixed stopcock, 250 mm scope	24 Fr.
AB-006-066		26 Fr.



Ordering No.	Description	French
AB-006-063	Nephroscope sheath for continuous flow, perforated, rotatable stopcock , 250 mm scope	26 Fr.
AB-006-068	Nephroscope sheath for continuous flow, perforated, fixed stopcock , 250 mm scope	26 Fr.



Ordering No.	Description	French
AB-006-062	Obturator hollow used for nephroscopy sheaths, 250 mm scope	24 Fr.
AB-006-064		26 Fr.

4.4 | ACCESSORIES

TELESCOPE BOUGIES | CANNULAS




Ordering No.	Description
AB-006-011	Set of telescope bougies consisting of: 6 dilating bougies 9, 12, 15, 18, 21 and 24 Fr., incl. 1 rigid and 1 flexible guide probe
AB-006-012	Guide probe, rigid
AB-006-013	Guide probe, flexible

Ordering no.	Description	French
AB-006-020	Dilating sheath	9 Fr.
AB-006-021		12 Fr.
AB-006-022		15 Fr.
AB-006-023		18 Fr.
AB-006-024		21 Fr.
AB-006-025		24 Fr.
AB-006-026		27 Fr.
AB-006-027		30 Fr.

Ordering No.	Description
AB-006-052	Puncture cannula, inner and outer cannula, 1 box = 5 pieces
AB-006-053	Dilatation cannula, inner and outer cannula, Ø 3.0 mm

4.4 | ACCESSORIES

FORCEPS



	Ordering No.	Description	Length	Diameter
	AB-006-049	Biopsy forceps, single action with ring handle	400 mm	Ø 1.8 mm
	AB-006-043		500 mm	Ø 1.8 mm
	AB-006-050	Alligator grasping forceps with ring handle	500 mm	Ø 1.8 mm
	AB-006-044		400 mm	Ø 3.0 mm
	AB-006-051	Scissors, single action with ring handle	400 mm	Ø 1.8 mm
	AB-006-045		500 mm	Ø 1.8 mm
	AB-006-047	3-Prongs retriever for larger stones with coil handle	400 mm	Ø 1.8 mm
	AB-006-041		500 mm	Ø 1.8 mm
	AB-006-046	Fenestrated stone grasping forceps for smaller stones, with coil handle	380 mm	Ø 3.5 mm
	AB-006-040		430 mm	
	AB-006-048	Fenestrated stone grasping forceps for smaller stones, with ring handle	380 mm	
	AB-006-042		430 mm	

5. URETERO-RENOSCOPY

5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES
URETERO-RENOSCOPES
BICHLER URETERO-RENOSCOPES
MARBERGER URETERO-RENOSCOPES
FORCEPS, SCISSORS, INSTRUMENTS
STONE EXTRACTORS
MONOPOLAR FLEXIBLE ELECTRODES
MONOPOLAR CONNECTING CABLES

5.2 | FLEXIBLE URETERO-RENOSCOPE

FLEX II
ACCESSORIES
INSTRUMENTS

5. HD URETERO-RENOSCOPES

40,000 px

17,000 px



1 URETERO-RENOSCOPE



2 FORCEPS



5.1 | URETERO-RENSCOPES SEMI-RIGID

HD URETERO-RENSCOPES



Ordering No.	Description	Working length	Viewing angle
AA-007-003	HD Uretero-renaloscope semi-rigid - 45° angled ocular	425 mm	0°
AA-007-016	- Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr.		5°
AA-007-017	- 40,000 px		12°

Ordering No.	Description	Working length	Viewing angle
AA-007-004	HD Uretero-renaloscope semi-rigid - 45° angled ocular	310 mm	0°
AA-007-018	- Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr.		5°
AA-007-019	- 40,000 px		12°

Ordering No.	Description	Working length	Viewing angle
AA-007-005	HD Uretero-renaloscope semi-rigid - 45° angled ocular	150 mm	0°
AA-007-036	- Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr.		5°
AA-007-037	- 40,000 px		12°

Delivery contains

HD uretero-renaloscope, brush Ø 3.0 mm (length 590 mm), brush Ø 3.8 mm (length 350 mm), sealing cap silicone (grey), Duckbill sealing



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES



Ordering No.	Description	Working length	Viewing angle
AA-007-009	HD Uretero-renaloscope semi-rigid	425 mm	0°
AA-007-028	- 45° angled ocular - Working Ø steps by 7.5 to 13.2 Fr.		5°
AA-007-029	- Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40,000 px		12°

Delivery contains

HD uretero-renaloscope, brush Ø 3.0 mm (length 590 mm), brush Ø 3.8 mm (length 350 mm), sealing cap silicone (grey), Duckbill sealing



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

5.1 | URETERO-RENSCOPES SEMI-RIGID

HD URETERO-RENSCOPES



Ordering No.	Description	Working length	Viewing angle
AA-007-015	HD Uretero-renaloscope semi-rigid	425 mm	0°
AA-007-034	- 45° angled ocular		5°
AA-007-035	- Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17,000 px		12°

Ordering No.	Description	Working length	Viewing angle
AA-007-014	HD Uretero-renaloscope semi-rigid	310 mm	0°
AA-007-032	- 45° angled ocular		5°
AA-007-033	- Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17,000 px		12°

Ordering No.	Description	Working length	Viewing angle
AA-007-013	HD Uretero-renaloscope semi-rigid	140 mm	0°
AA-007-030	- 45° angled ocular		5°
AA-007-031	- Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17,000 px		12°

Delivery contains

HD uretero-renaloscope, brush Ø 3.0 mm (length 590 mm), brush Ø 3.8 mm (length 350 mm), sealing cap silicone (grey), Duckbill sealing



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

5.1 | 5.1 | URETERO-RENSCOPIES SEMI-RIGID

HD URETERO-RENSCOPIES



Ordering No.	Description	Working length	Viewing angle
AA-007-006	HD Uretero-renaloscope semi-rigid	425 mm	0°
AA-007-022	- 45° angled ocular		5°
AA-007-023	- Working Ø steps by 4.9 to 6.5 Fr. - Working channel 1 × 3 Fr. - 17,000 px		12°

Ordering No.	Description	Working length	Viewing angle
AA-007-007	HD Uretero-renaloscope semi-rigid	310 mm	0°
AA-007-024	- 45° angled ocular		5°
AA-007-025	- Working Ø steps by 4.9 to 6.5 Fr. - Working channel 1 × 3 Fr. - 17,000 px		12°

Ordering No.	Description	Working length	Viewing angle
AA-007-008	HD Uretero-renaloscope semi-rigid	150 mm	0°
AA-007-026	- 45° angled ocular		5°
AA-007-027	- Working Ø steps by 4.9 to 6.5 Fr. - Working channel 1 × 3 Fr. - 17,000 px		12°

Delivery contains


HD uretero-renaloscope, brush Ø 3.0 mm (length 590 mm), brush Ø 3.8 mm (length 350 mm), sealing cap silicone (grey), Duckbill sealing



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

5.1 | URETERO-RENOSCOPES SEMI-RIGID

URETERO-RENOSCOPES




Ordering No.	Description	Working length	Viewing angle
AA-007-001	Uretero-renaloscope semi-rigid - 45° angled ocular - Ø Working steps • distal 7.5 Fr. • first step 9.8 Fr. • tube 13 Fr.	310 mm	8°
AA-007-002		425 mm	
	- Ø Working channel 1 x 5 Fr. or 2 x 3 Fr. - Field of view: 85° - Image bundles: 30 k		

	Ordering No.	Description
	AA-007-021	Instrument bridge used for uretero-renaloscopes (included in the delivery content)
	AJ-010-001	Rubber cap for instrument channel
	AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

5.1 | URETERO-RENOSCOPES SEMI-RIGID

URETERO-RENOSCOPES



Ordering No.	Description	Working length	Viewing angle
AA-007-011	Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps: • distal 6 Fr. • first step 7.5 Fr. • second step 9.9 Fr. • tube 10.6 Fr.	310 mm	8°
AA-007-012		425 mm	
	- Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - Field of view: 100° - Image bundles: 17 k		




Ordering No.	Description
AA-007-021	Instrument bridge used for uretero-renaloscopes (included in the delivery content)
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm




5.1 | URETERO-RENSCOPES SEMI-RIGID


BICHLER URETERO-RENSCOPES




Ordering No.	Description	Working length	Viewing angle
AA-007-056	BICHLER Uretero-roscope semi-rigid - 45° angled ocular	310 mm	5°
AA-007-055	- Working Ø steps by 4.5 Fr. to 6.5 Fr. - Working channel 1 × 3.3 Fr.	430 mm	




Ordering No.	Description	Working length	Viewing angle
AA-007-043	BICHLER Uretero-roscope semi-rigid - 45° angled ocular	310 mm	12°
AA-007-044	- Working Ø steps by 6 Fr. to 7.5 Fr. - Working channel 1 × 4 Fr. or 2 × 2.4 Fr.	430 mm	



Ordering No.	Description	Working length	Viewing angle
AA-007-047	BICHLER Uretero-roscope semi-rigid - 45° angled ocular	310 mm	12°
AA-007-048	- Working Ø steps by 8 Fr. to 9.8 Fr. - Working channel 1 × 5 Fr. or 2 × 3 Fr.	430 mm	




Ordering No.	Description	Working length	Viewing angle
AA-007-051	BICHLER Uretero-roscope semi-rigid - 45° angled ocular	310 mm	12°
AA-007-052	- Working Ø steps by 8.5 Fr. to 11.5 Fr. - Working channel 1 × 6 Fr. or 2 × 4 Fr.	430 mm	




Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm


5.1 | URETERO-RENOSCOPES SEMI-RIGID MARBERGER URETERO-RENOSCOPES




Ordering No.	Description	Working length	Viewing angle
AA-007-054	MARBERGER Uretero-renaloscope semi-rigid - 45° angled ocular	310 mm	5°
AA-007-053	- Working Ø steps by 4.5 Fr. to 6.5 Fr. - Working channel 1 × 3.3 Fr.	430 mm	




Ordering No.	Description	Working length	Viewing angle
AA-007-041	MARBERGER Uretero-renaloscope semi-rigid - 45° angled ocular	310 mm	12°
AA-007-042	- Working Ø steps by 6 Fr. to 7.5 Fr. - Working channel 1 × 4 Fr. or 2 × 2.4 Fr.	430 mm	



Ordering No.	Description	Working length	Viewing angle
AA-007-045	MARBERGER Uretero-renaloscope semi-rigid - 45° angled ocular	310 mm	12°
AA-007-046	- Working Ø steps by 8 Fr. to 9.8 Fr. - Working channel 1 × 5 Fr. or 2 × 3 Fr.	430 mm	




Ordering No.	Description	Working length	Viewing angle
AA-007-049	MARBERGER Uretero-renaloscope semi-rigid - 45° angled ocular	310 mm	12°
AA-007-050	- Working Ø steps by 8.5 Fr. to 11.5 Fr. - Working channel 1 × 6 Fr. or 2 × 4 Fr.	430 mm	



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

5.1 | URETERO-RENSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID




	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps oval, double action	Ø 1.0 mm	3 Fr.	AC-001-011	AC-001-211	AC-001-011-R
		Ø 1.2 mm	4 Fr.	AC-001-010	AC-001-210	AC-001-010-R
		Ø 1.6 mm	5 Fr.	AC-001-012	AC-001-212	AC-001-012-R
		Ø 2.1 mm	6 Fr.	AC-001-013	AC-001-213	AC-001-013-R
		Ø 2.3 mm	7 Fr.	AC-001-014	AC-001-214	AC-001-014-R
		Ø 3.0 mm	9 Fr.	AC-001-015	AC-001-215	AC-001-015-R
	Biopsy forceps oval, serrated, double action	Ø 1.0 mm	3 Fr.	AC-002-011	AC-002-211	AC-002-011-R
		Ø 1.6 mm	5 Fr.	AC-002-012	AC-002-212	AC-002-012-R
		Ø 2.1 mm	6 Fr.	AC-002-013	AC-002-213	AC-002-013-R
		Ø 2.3 mm	7 Fr.	AC-002-014	AC-002-214	AC-002-014-R
		Ø 3.0 mm	9 Fr.	AC-002-015	AC-002-215	AC-002-015-R
	Biopsy forceps oval, serrated, teeth interlock, double action	Ø 1.6 mm	5 Fr.	AC-015-012	AC-015-212	AC-015-012-R
		Ø 2.1 mm	6 Fr.	-	AC-015-213	AC-015-013-R
		Ø 2.3 mm	7 Fr.	AC-015-014	AC-015-214	AC-015-014-R
	Grasping forceps, double action	Ø 1.2 mm	4 Fr.	AC-017-011	AC-017-211	AC-017-011-R
		Ø 2.3 mm	7 Fr.	AC-017-014	AC-017-214	AC-017-014-R
	Grasping forceps, alligator jaw, double action	Ø 1.0 mm	3 Fr.	AC-006-011	AC-006-211	AC-006-011-R
		Ø 1.2 mm	4 Fr.	AC-006-013	AC-006-213	AC-006-013-R
		Ø 1.6 mm	5 Fr.	AC-006-012	AC-006-212	AC-006-012-R
		Ø 2.3 mm	7 Fr.	AC-006-014	AC-006-214	AC-006-014-R
		Ø 3.0 mm	9 Fr.	AC-006-015	AC-006-215	AC-006-015-R

5.1 | URETERO-RENSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	Ø 2.3 mm	7 Fr.	AC-008-014	AC-008-214	AC-008-014-R
		Ø 3.0 mm	9 Fr.	AC-008-015	AC-008-215	AC-008-015-R
	Grasping forceps teeth double action	Ø 1.2 mm	4 Fr.	AC-016-011	AC-016-211	AC-016-011-R
		Ø 1.6 mm	5 Fr.	AC-016-012	AC-016-212	AC-016-012-R
	Grasping forceps teeth interlock, double action	Ø 1.6 mm	5 Fr.	AC-020-012	AC-020-212	AC-020-012-R
		Ø 2.3 mm	7 Fr.	AC-020-014	AC-020-214	AC-020-014-R
	Grasping forceps rat tooth, double action	Ø 1.0 mm	3 Fr.	AC-019-011	AC-019-211	AC-019-011-R
		Ø 1.6 mm	5 Fr.	AC-019-012	AC-019-212	AC-019-012-R
	Grasping forceps for stone street, double action	Ø 1.2 mm	4 Fr.	AC-007-011	AC-007-211	AC-007-011-R
		Ø 1.6 mm	5 Fr.	AC-007-012	AC-007-212	AC-007-012-R
	Grasping forceps for big stones, double action	Ø 1.2 mm	4 Fr.	AC-009-011	AC-009-211	AC-009-011-R
		Ø 1.6 mm	5 Fr.	AC-009-012	AC-009-212	AC-009-012-R
		Ø 2.3 mm	7 Fr.	AC-009-014	AC-009-214	AC-009-014-R
	3-Prongs retriever	Ø 1.0 mm	3 Fr.	AC-012-011		
		Ø 1.6 mm	5 Fr.	AC-012-012		
		Ø 2.3 mm	7 Fr.	AC-012-014		


5.1 | URETERO-RENOSCOPY

SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Scissors, single action	Ø 1.0 mm	3 Fr.	AC-003-011	AC-003-211	AC-003-011-R
		Ø 1.2 mm	4 Fr.	AC-003-010	AC-003-210	AC-003-010-R
		Ø 1.6 mm	5 Fr.	AC-003-012	AC-003-212	AC-003-012-R
		Ø 2.1 mm	6 Fr.	AC-003-013	AC-003-213	AC-003-013-R
		Ø 2.3 mm	7 Fr.	AC-003-014	AC-003-214	AC-003-014-R
		Ø 3.0 mm	9 Fr.	AC-003-015	AC-003-215	AC-003-015-R
	Hook scissors, single action	Ø 2.3 mm	7 Fr.	AC-004-014	AC-004-214	AC-004-014-R
		Ø 3.0 mm	9 Fr.	AC-004-015	AC-004-215	AC-004-015-R
	Uretero-slitting scissors, retrograde cutting single action	Ø 1.2 mm	4 Fr.	AC-018-011	AC-018-211	AC-018-011-R
		Ø 1.6 mm	5 Fr.	AC-018-012	AC-018-212	AC-018-012-R
	Scissors, serrated flexible, dou- ble action	Ø 2.3 mm	7 Fr.	AC-005-014	AC-005-214	AC-005-014-R

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




FEATURES

- 1 QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
- 2 ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
- 3 OVERLOAD PROTECTION AVAILABLE
- 4 SIDE-MOUNTED LUER LOCK

5.1 | URETERO-RENOGRAPHY


INSTRUMENTS SEMI-RIGID | DETACHABLE



Ordering No.		Description			
AC-100-001		Handle for detachable instruments			
Description		Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0 mm	3 Fr.	AC-003-211-T	AC-003-211-OT
		Ø 1.2 mm	4 Fr.	AC-003-210-T	AC-003-210-OT
		Ø 1.6 mm	5 Fr.	AC-003-212-T	AC-003-212-OT
		Ø 2.3 mm	7 Fr.	AC-003-214-T	AC-003-214-OT
		Ø 3.0 mm	9 Fr.	AC-003-215-T	AC-003-215-OT
	Scissors pointed, single action (Minimum order quantity 5pcs.)	Ø 1.6 mm	5 Fr.	AC-004-622-T	AC-004-622-OT
		Ø 2.3 mm	7 Fr.	AC-004-624-T	AC-004-624-OT
	Hook scissors, single action	Ø 2.3 mm	7 Fr.	AC-004-214-T	AC-004-214-OT
		Ø 3.0 mm	9 Fr.	AC-004-215-T	AC-004-215-OT

5.1 | URETERO-RENOSCOPY


INSTRUMENTS SEMI-RIGID | DETACHABLE



Ordering No.		Description			
AC-100-001		Handle for detachable instruments			
Description		Diameter	French	Complete	Insert
	Biopsy forceps oval, double action	Ø 1.0 mm	3 Fr.	AC-001-211-T	AC-001-211-OT
		Ø 1.2 mm	4 Fr.	AC-001-210-T	AC-001-210-OT
		Ø 1.6 mm	5 Fr.	AC-001-212-T	AC-001-212-OT
		Ø 2.3 mm	7 Fr.	AC-001-214-T	AC-001-214-OT
		Ø 3.0 mm	9 Fr.	AC-001-215-T	AC-001-215-OT
	Biopsy forceps oval, serrated, double action	Ø 1.2 mm	4 Fr.	AC-016-211-T	AC-016-211-OT
		Ø 1.6 mm	5 Fr.	AC-016-212-T	AC-016-212-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0 mm	3 Fr.	AC-006-211-T	AC-006-211-OT
		Ø 1.2 mm	4 Fr.	AC-006-213-T	AC-006-213-OT
		Ø 1.6 mm	5 Fr.	AC-006-212-T	AC-006-212-OT
		Ø 2.3 mm	7 Fr.	AC-006-214-T	AC-006-214-OT
		Ø 3.0 mm	9 Fr.	AC-006-215-T	AC-006-215-OT
	Grasping forceps for stone street, double action	Ø 1.2 mm	4 Fr.	AC-007-211-T	AC-007-211-OT
		Ø 1.6 mm	5 Fr.	AC-007-212-T	AC-007-212-OT
	Grasping forceps for big stones, double action	Ø 1.2 mm	4 Fr.	AC-009-211-T	AC-009-211-OT
		Ø 1.6 mm	5 Fr.	AC-009-212-T	AC-009-212-OT

5.1 | URETERO-RENOSCOPY


INSTRUMENTS RIGID | DETACHABLE



Ordering No.		Description			
AC-100-001		Handle for detachable instruments			
Description		Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0 mm	3 Fr.	AC-103-211-T	AC-103-211-OT
		Ø 1.2 mm	4 Fr.	AC-103-210-T	AC-103-210-OT
		Ø 1.6 mm	5 Fr.	AC-103-212-T	AC-103-212-OT
		Ø 2.3 mm	7 Fr.	AC-103-214-T	AC-103-214-OT
		Ø 3.0 mm	9 Fr.	AC-103-215-T	AC-103-215-OT
	Hook scissors, single action	Ø 2.3 mm	7 Fr.	AC-104-214-T	AC-104-214-OT
		Ø 3.0 mm	9 Fr.	AC-104-215-T	AC-104-215-OT
	Biopsy punch forceps, single action	Ø 1.6 mm	5 Fr.	AC-114-212-T	AC-114-212-OT

5.1 | URETERO-RENOSCOPY

INSTRUMENTS RIGID | DETACHABLE



Ordering No.		Description			
AC-100-001		Handle for detachable instruments			
Description		Diameter	French	Complete	Insert
	Biopsy forceps oval, double action	Ø 1.0 mm	3 Fr.	AC-101-211-T	AC-101-211-OT
		Ø 1.2 mm	4 Fr.	AC-101-210-T	AC-101-210-OT
		Ø 1.6 mm	5 Fr.	AC-101-212-T	AC-101-212-OT
		Ø 2.3 mm	7 Fr.	AC-101-214-T	AC-101-214-OT
		Ø 3.0 mm	9 Fr.	AC-101-215-T	AC-101-215-OT
	Biopsy forceps oval, serrated, double action	Ø 1.2 mm	4 Fr.	AC-116-211-T	AC-116-211-OT
		Ø 1.6 mm	5 Fr.	AC-116-212-T	AC-116-212-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0 mm	3 Fr.	AC-106-211-T	AC-106-211-OT
		Ø 1.2 mm	4 Fr.	AC-106-213-T	AC-106-213-OT
		Ø 1.6 mm	5 Fr.	AC-106-212-T	AC-106-212-OT
		Ø 2.3 mm	7 Fr.	AC-106-214-T	AC-106-214-OT
		Ø 3.0 mm	9 Fr.	AC-106-215-T	AC-106-215-OT
	Grasping forceps for stone street, double action	Ø 1.2 mm	4 Fr.	AC-107-211-T	AC-107-211-OT
		Ø 1.6 mm	5 Fr.	AC-107-212-T	AC-107-212-OT
	Grasping forceps for big stones, double action	Ø 1.2 mm	4 Fr.	AC-109-211-T	AC-109-211-OT
		Ø 1.6 mm	5 Fr.	AC-109-212-T	AC-109-212-OT

5.1 | URETERO-RENOSCOPY

STONE EXTRACTORS | 4 WIRES



Ordering No.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering No.	Wires	French	Working length	Description
AF-001-001	4	3 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-002	4	3 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-015	4	3 Fr.	90 cm	Dormia stone extractor, basket with tubus only
AF-001-016	4	3 Fr.	90 cm	Dormia stone extractor complete with handle
AF-001-023	4	3 Fr.	120 cm	Dormia stone extractor, basket with tubus only
AF-001-022	4	3 Fr.	120 cm	Dormia stone extractor complete with handle
AF-001-003	4	5 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-004	4	5 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-017	4	5 Fr.	90 cm	Dormia stone extractor, basket with tubus only
AF-001-019	4	5 Fr.	90 cm	Dormia stone extractor complete with handle
AF-001-005	4	7 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-006	4	7 Fr.	60 cm	Dormia stone extractor complete with handle

5.1 | URETERO-RENOSCOPY

STONE EXTRACTORS | 5-7 WIRES



Ordering No.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering No.	Wires	French	Working length	Description
AF-001-007	5	5 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-008	5	5 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-025	5	5 Fr.	90 cm	Dormia stone extractor, basket with tubus only
AF-001-018	5	5 Fr.	90 cm	Dormia stone extractor complete with handle
AF-001-009	5	7 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-010	5	7 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-011	6	5 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-012	6	5 Fr.	60 cm	Dormia stone extractor complete with handle
AF-001-013	7	7 Fr.	60 cm	Dormia stone extractor, basket with tubus only
AF-001-014	7	7 Fr.	60 cm	Dormia stone extractor complete with handle

5.1 | URETERO-RENOSCOPY

MONOPOLAR FLEXIBLE ELECTRODES




Description		French	Length	Ordering No.	French	Length	Ordering No.
	Loop electrode flexible, monopolar	5 Fr.	53 cm	AD-003-003	5 Fr.	65 cm	AD-003-013
		7 Fr.	53 cm	AD-003-005	5 Fr.	65 cm	AD-003-015
		9 Fr.	53 cm	AD-003-007	5 Fr.	65 cm	AD-003-017
		10 Fr.	53 cm	AD-003-008	5 Fr.	65 cm	AD-003-018
	Needle electrode flexible, monopolar	3 Fr.	53 cm	AD-002-001	3 Fr.	65 cm	AD-002-011
		4 Fr.	53 cm	AD-002-002	4 Fr.	65 cm	AD-002-012
		5 Fr.	53 cm	AD-002-003	5 Fr.	65 cm	AD-002-013
		7 Fr.	53 cm	AD-002-005	7 Fr.	65 cm	AD-002-015
	Button electrode flexible, monopolar	3 Fr.	53 cm	AD-001-001	3 Fr.	65 cm	AD-001-011
		4 Fr.	53 cm	AD-001-002	4 Fr.	65 cm	AD-001-012
		5 Fr.	53 cm	AD-001-003	5 Fr.	65 cm	AD-001-013
		6 Fr.	53 cm	AD-001-004	6 Fr.	65 cm	AD-001-014
		7 Fr.	53 cm	AD-001-005	7 Fr.	65 cm	AD-001-015
		8 Fr.	53 cm	AD-001-006	8 Fr.	65 cm	AD-001-016
		9 Fr.	53 cm	AD-001-007	9 Fr.	65 cm	AD-001-017
		10 Fr.	53 cm	AD-001-008	10 Fr.	65 cm	AD-001-018

5.1 | URETERO-RENSCOPY MONOPOLAR CONNECTING CABLES



Connection instrument side



Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001	3 m	Ø 4.0 mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5 m		
	AD-010-002	3 m	Ø 5.0 mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5 m		
	AD-010-003	3 m	Ø 8.0 mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5 m		
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable male for STEMA MASTER



Ordering No.	Description
EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

5.2 | FLEXIBLE URETERO-RENSCOPE FLEX II

- High luminance
- Exact adjustment of focus
- Precise up-/down angulation



Ordering No.	Description
AA-008-003	Flexible uretero-roscope FLEX II

Technical data

Cystoscope

Insertion tube	Ø 2.8 mm
Working Channel.....	Ø 1.2 mm
Working length	680 mm
Field of view	88°
Pixel	18,000
Angulation up-down	275°
Accessories	Suction valve, leak tester

Delivery contains

Flexible uretero-roscope FLEX II, leak tester, suction valve

INCLUDED IN SCOPE OF DELIVERY:

AA-008-103

Leak tester

AA-008-105

Suction valve



AA-008-103



AA-008-105

5.2 | FLEXIBLE URETERO-RENSCOPE ACCESSORIES



Ordering No.	Description
AA-008-100	Bridge with Luer-Lock and 2 stopcocks for AA-008-002

Ordering No.	Description
AA-008-101	Single use cleaning brushes, 1.0 mm, 120 cm length, 1 box = 250 pieces



Ordering No.	Description
AA-008-102	Suitcase grey for STEMA flexible scopes



Ordering No.	Description
AA-008-103	Leak tester




Ordering No.	Description
AA-008-104	Single use Luer-Lock adapter, 1 box = 25 pieces



Ordering No.	Description
AA-008-105	Suction valve

5.2 | FLEXIBLE URETERO-RENSCOPE

INSTRUMENTS FLEXIBLE



	Description	Diameter	French	Complete	
	Biopsy forceps oval, flexible, double action	Ø 1.0 mm	3 Fr.	AC-001-031	
	Grasping forceps alligator jaw, flexible, double action	Ø 1.0 mm	3 Fr.	AC-006-031	
	3-Prongs retriever, flexible, with handle	Ø 1.0 mm	3 Fr.	AC-012-015	

6. LITHOTRIPSY


6.1 | LITHOTRIPSY
CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS

6. LITHOTRIPSY


1 SCOPES



2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS



6.1 | LITHOTRIPSY CYSTOSCOPES | SHEATHS



Viewing angle	Description	Diameter	Working length	Ordering No.
	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-001
	HD Cystoscope autoclavable			AA-002-011

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-003
	HD Cystoscope autoclavable			AA-002-013

6.1 | LITHOTRIPSY CYSTOSCOPES | SHEATHS



Ordering No.	Description
AB-002-010	Lithotripsy sheath with central valve used for AB-004-011



Ordering No.	Description
AB-002-013	Lithotripsy sheath used for AB-004-011

6.1 | LITHOTRIPSY

SHEATHS | OBTURATORS



Ordering No.	Description
AB-002-011	Diagnostic sheath with working channel for instruments 7 Fr.



Ordering No.	Description
AB-002-014	Adapter for use with Ellik



Ordering No.	Description
AB-002-012	Obturator for AB-002-010 + AB-002-013

6.1 | LITHOTRIPSY

WORKING ELEMENT | LITHOTRIPTOR



Ordering No.	Description
AB-004-011	Working element used for punch lithotripsy



Ordering No.	Description
AB-004-012	Mechanical lithotripter with handle and turning screw, Ø 24 Fr., 2 stopcocks with luer lock adapter, for use with scope AA-002-004

7. MORCELLATION

7.1 | MORCELLATION 26 FR.
NEPHROSCOPE
SHEATHS
OBTURATORS

7.2 | MORCELLATION 22 FR.
MORCESCOPE
SHEATHS
OBTURATORS


7. MORCELLATION

26 Fr.

22 Fr.


1 HD NEPHROSCOPE

2 SHEATHS
& OBTURATORS



3 MORCESCOPE


4 SHEATHS
& OBTURATORS





26 Fr.

7.1 | MORCELLATION
NEPHROSCOPE




Ordering No.	Description	Working length	Viewing angle
AA-006-012	Nephroscope - 90° angled ocular - Ø 7.45 mm - Suction channel Ø 1.5 mm - Working channel Ø 5.0 mm	250 mm	0°



Ordering No.	Description
AA-006-080	Flow valve, Ø 6.1 mm, used for STEMA nephroscopes
AA-006-081	Sealing cap, Ø 15 mm/1.0 mm, for AA-006-080
AA-006-083	Sealing cap, Ø 20 mm/2.0 mm, for AA-006-001–AA-006-034
AB-004-010-QL	Adapter for use with nephroscope and resectoscope sheaths 26 Fr. STEMA Click-System

7.1 | MORCELLATION SHEATHS



	Ordering No.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow



Ordering No.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow

7.1 | MORCELLATION OBTURATORS



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-011	Standard obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow




Ordering No.	Description	French	For outer sheath	Colour code
AH-001-013	Viewing obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow



Ordering No.	Description	French	Colour code
AH-001-308	Inner sheath without ceramic isolation, with laser channel up to 7 Fr. , for use with AH-001-301	26 Fr.	Yellow

22 Fr.

7.2 | MORCELLATION
MORCESCOPE



Ordering No.	Description	Working length	Viewing angle
AA-006-072	Morcescope - 22 Fr. outer tube - Working channel 15 Fr. - Field of view 85° - 40,000 px	180 mm	0°

Ordering No.	Description
AA-006-084	Sealing cap white, Ø 13.5 mm, 1 box = 10 pieces

7.2 | MORCELLATION SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-302	Outer sheath with perforation, rotatable	22 Fr.	White



Ordering No.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, rotatable	19 Fr.	White

7.2 | MORCELLATION OBTURATORS



Ordering No.	Description	French	For outer sheath	Colour code
AH-002-011	Standard obturator for resectoscope sheaths	19 Fr.	22 Fr.	White



Ordering No.	Description	French	For outer sheath	Colour code
AH-002-013	Viewing obturator	19 Fr.	22 Fr.	White

8. PEDIATRIC UROLOGY

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
FORCEPS
MONOPOLAR ELECTRODES
CONNECTING CABLES

8.2 | PEDIATRIC RESECTOSCOPY


CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
ELECTRODES
CONNECTING CABLES

8.3 | PEDIATRIC URETHROTOMY


CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
COLD KNIVES AND ACCESSORIES

8. PEDIATRIC UROLOGY


1 PEDIATRIC CYSTOSCOPY



2 PEDIATRIC RESECTOSCOPY



3 PEDIATRIC URETHROTOMY



8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-005-001	Pediatric cystoscope autoclavable	0°	Ø 2.7 mm	187.5 mm
AA-005-003		30°		
AA-005-004		70°		



Ordering No.	Description	French	Colour code
AB-001-001	Cysto-Urethroscope sheath with 2 fixed stopcocks , 1 × 3 Fr. instrument channel incl. obturator	13 Fr.	Red
AB-001-001-13	Cysto-Urethroscope sheath with 2 fixed stopcocks , 1 × 5 Fr. or 2×4 Fr. instrument channel incl. obturator		



Ordering No.	Description	French	Colour code
AB-001-011	Standard obturator	13 Fr.	Red



13 Fr.

8.1 | PEDIATRIC CYSTO-URETHROSCOPY WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering No.	Description	For sheath	Colour code
AB-003-005	Working element with ALBARRAN deflector	13 Fr.	Red



Ordering No.	Description	For sheath	Colour code
AB-003-001	Endoscope bridge	13 Fr.	Red

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-004-011	Pediatric cystoscope autoclavable	0°	Ø 2.0 mm	187 mm
AA-004-012		30°		



Ordering No.	Description	French	Colour code
AB-001-001-9	Cysto-Urethroscope sheath with 2 fixed stopcocks, 1 x 3 Fr. instrument channel incl. obturator	9.5 Fr.	Blue
AB-001-001-11		11 Fr.	Green



Ordering No.	Description	French	Colour code
AB-001-011-9	Standard obturator	9.5 Fr.	Blue
AB-001-011-11		11 Fr.	Green

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-004-013	Pediatric cystoscope autoclavable	0°	Ø 1.9 mm	145 mm
AA-004-014		30°		




Ordering No.	Description	French	Colour code
AB-001-001-7	Cysto-Urethroscope sheath with 2 fixed stopcocks	7.5 Fr.	Yellow



Ordering No.	Description	French	Colour code
AB-001-011-7	Standard obturator for cysto-urethroscope sheaths	7.5 Fr.	Yellow

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

FORCEPS | ELECTRODES



	Description	WL	Diameter	French	Flexible	Semi-rigid
	Biopsy forceps oval, double action	28 cm	Ø 1.0 mm	3 Fr.	AC-001-021	AC-001-221
	Biopsy forceps oval, serrated, teeth interlock, double action	28 cm	Ø 1.0 mm	3 Fr.	AC-015-021	AC-015-221
	Grasping forceps, alligator jaw, double action	28 cm	Ø 1.0 mm	3 Fr.	AC-006-021	AC-006-221
	Grasping forceps, rat tooth, double action	28 cm	Ø 1.0 mm	3 Fr.	AC-019-021	-
	Scissors, single action	28 cm	Ø 1.0 mm	3 Fr.	AC-003-021	AC-003-221




	Ordering No.	Description	French	Diameter
	AD-002-021	Needle electrode flexible, monopolar	3.0 Fr.	280 mm
	AD-001-021	Button electrode flexible, monopolar		

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CONNECTING CABLES | MONOPOLAR



Connection instrument side



Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001	3 m	Ø 4.0 mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5 m		
	AD-010-002	3 m	Ø 5.0 mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5 m		
	AD-010-003	3 m	Ø 8.0 mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5 m		
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable male for STEMA MASTER



Ordering No.	Description
EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

8.2 | PEDIATRIC RESECTOSCOPY

RESECTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-005-001	Pediatric cystoscope autoclavable	0°	Ø 2.7 mm	187.5 mm



Ordering No.	Description	French	Colour code
AH-003-001	Standard resectoscope sheath with 1 rotatable stopcock, without obturator	13 Fr.	Red



Ordering No.	Description	French	Colour code
AH-003-011	Standard obturator	13 Fr.	Red

8.2 | PEDIATRIC RESECTOSCOPY

WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering No.	Description	For sheath
AG-003-001	Working element used for pediatric resectoscopes, passive cutting , handle open	13 Fr.



Ordering No.	Description	For sheath
AI-003-003	Endoscope bridge with 3 Fr. instrument channel	13 Fr.

8.2 | PEDIATRIC RESECTOSCOPY ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-004-002	Cutting loop electrode angled	13 Fr.	Red
	AE-004-018	Cutting loop electrode small angled	13 Fr.	Red
	AE-004-004	Rollerball electrode	13 Fr.	Red
	AE-004-016	Knife electrode angled	13 Fr.	Red
	AE-004-019	Hook electrode without ball	13 Fr.	Red



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

8.2 | PEDIATRIC RESECTOSCOPY

RESECTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-004-011	Pediatric cystoscope autoclavable	0°	Ø 2.0 mm	187 mm



Ordering No.	Description	French	Colour code
AH-005-001	Standard resectoscope sheath with 1 rotatable stopcock , without obturator	9 Fr.	Blue
AH-004-001		11 Fr.	Green



Ordering No.	Description	French	Colour code
AH-005-011	Standard obturator	9 Fr.	Blue
AH-004-011		11 Fr.	Green

8.2 | PEDIATRIC RESECTOSCOPY

WORKING ELEMENTS



Ordering No.	Description	For sheath
AG-005-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	9 Fr.
AG-004-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	11 Fr.

8.2 | PEDIATRIC RESECTOSCOPY ELECTRODES




	Ordering No.	Description	French	Colour code
	AE-006-002	Cutting loop electrode angled	9 Fr.	Blue
	AE-005-002		11 Fr.	Green
	AE-006-018	Cutting loop electrode small angled	9 Fr.	Blue
	AE-005-018		11 Fr.	Green
	AE-006-004	Rollerball electrode	9 Fr.	Blue
	AE-005-004		11 Fr.	Green
	AE-006-016	Knife electrode angled	9 Fr.	Blue
	AE-005-016		11 Fr.	Green
	AE-006-019	Hook electrode without ball	9 Fr.	Blue
	AE-005-019		11 Fr.	Green



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

8.2 | PEDIATRIC RESECTOSCOPY

CONNECTING CABLES | MONOPOLAR



Connection instrument side



Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001	3 m	Ø 4.0 mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5 m		
	AD-010-002	3 m	Ø 5.0 mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5 m		
	AD-010-003	3 m	Ø 8.0 mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5 m		
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable male for STEMA MASTER



Ordering No.	Description
EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

8.2 | PEDIATRIC RESECTOSCOPY

CONNECTING CABLES | MONOPOLAR



Connection unit side



Unit side	Ordering No.	Length	Description
	AD-010-025	4.5 m	HF Monopolar connecting cable for Storz resectoscopes, 90° angled connector
	AD-010-026	4.5 m	HF Monopolar connecting cable for Olympus resectoscopes, diameter Ø 3 mm



RECOMMENDED PRODUCTS:

CI-001-030 STEMA HF-U PRO TOUCH

- For all monopolar & bipolar techniques
- with 7" touch screen display
- provided to perform the most advanced electrosurgery
- ADC / APC System

CI-001-040 DOUBLE PEDAL FOOTSWITCH

8.3 | PEDIATRIC URETHROTOMY

CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-005-001	Pediatric cystoscope autoclavable	0°	Ø 2.7 mm	187.5 mm



Ordering No.	Description	French	Colour code
AI-003-001	Standard urethrotome sheath with 2 fixed stopcocks incl. obturator	13 Fr.	Red



Ordering No.	Description	French	Colour code
AI-003-004	Standard obturator	13 Fr.	Red

8.3 | PEDIATRIC URETHROTOMY

WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering No.	Description	For sheath
AG-003-001	Working element used for pediatric resectoscopes, passive cutting , handle open	13 Fr.



Ordering No.	Description	For sheath
AI-003-003	Endoscope bridge with 3 Fr. instrument channel	13 Fr.

8.3 | PEDIATRIC URETHROTOMY

COLD KNIVES



	Ordering No.	Description	French	Colour code
	AI-003-007	Cold knife straight	13 Fr.	Red
	AI-003-009	Cold knife round	13 Fr.	Red
	AI-003-010	Cold knife ring shaped	13 Fr.	Red
	AI-003-011	Cold knife hook shaped	13 Fr.	Red



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

8.3 | PEDIATRIC URETHROTOMY

CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-004-011	Pediatric cystoscope autoclavable	0°	Ø 2.0 mm	187 mm



Ordering No.	Description	French	Colour code
AI-004-001	Standard urethrotome sheath with 2 fixed stopcocks incl. obturator	11 Fr.	Green



Ordering No.	Description	French	Colour code
AI-004-004	Standard obturator	11 Fr.	Green

8.3 | PEDIATRIC URETHROTOMY

WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering No.	Description	For sheath
AG-004-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	11 Fr.

8.3 | PEDIATRIC URETHROTOMY

COLD KNIVES



	Ordering No.	Description	French	Colour code
	AI-004-007	Cold knife straight	11 Fr.	Green
	AI-004-009	Cold knife round	11 Fr.	Green
	AI-004-010	Cold knife ring shaped	11 Fr.	Green
	AI-004-011	Cold knife hook shaped	11 Fr.	Green



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

9. ACCESSORIES

9.1 | INSTRUMENTS

DILATING BOUGIES
METAL CATHETERS
ACCESSORIES CATHETERS

9.2 | STERILIZATION


9.3 | CLEANING AND MAINTENANCE

9. ACCESSORIES & INSTRUMENTS

1 DILATATORS



2 STERILIZATION



3 PROTECTION



9.1 | INSTRUMENTS | DILATATOR BOUGIES



VAN BUREN

DITTEL,
STRAIGHTDITTEL,
CURVED

French	Ordering no. VAN BUREN	Ordering no. DITTEL straight	Ordering no. DITTEL curved
8 Fr.	AJ-003-001	AJ-003-021	AJ-003-050
9 Fr.		AJ-003-022	AJ-003-051
10 Fr.	AJ-003-002	AJ-003-023	AJ-003-052
11 Fr.		AJ-003-024	AJ-003-053
12 Fr.	AJ-003-003	AJ-003-025	AJ-003-054
13 Fr.		AJ-003-026	AJ-003-055
14 Fr.	AJ-003-004	AJ-003-027	AJ-003-056
15 Fr.		AJ-003-028	AJ-003-057
16 Fr.	AJ-003-005	AJ-003-029	AJ-003-058
17 Fr.		AJ-003-030	AJ-003-059
18 Fr.	AJ-003-006	AJ-003-031	AJ-003-060
19 Fr.		AJ-003-032	AJ-003-061
20 Fr.	AJ-003-007	AJ-003-033	AJ-003-062
21 Fr.		AJ-003-034	AJ-003-063
22 Fr.	AJ-003-008	AJ-003-035	AJ-003-064
23 Fr.		AJ-003-036	AJ-003-065
24 Fr.	AJ-003-009	AJ-003-037	AJ-003-066
25 Fr.		AJ-003-038	AJ-003-067
26 Fr.	AJ-003-010	AJ-003-039	AJ-003-068
27 Fr.		AJ-003-040	AJ-003-069
28 Fr.	AJ-003-011	AJ-003-041	AJ-003-070
29 Fr.		AJ-003-042	AJ-003-071
30 Fr.	AJ-003-012	AJ-003-043	AJ-003-072
31 Fr.			
32 Fr.	AJ-003-013	AJ-003-044	AJ-003-074
33 Fr.			
34 Fr.	AJ-003-014	AJ-003-045	
35 Fr.			
36 Fr.	AJ-003-015	AJ-003-046	
37 Fr.			
38 Fr.	AJ-003-016		

9.1 | INSTRUMENTS | DILATATOR BOUGIES | UTERINE DILATATORS



GUYON


French	Ordering no. GUYON	French	Ordering no. GUYON
8 Fr.	AJ-003-170	25 Fr.	AJ-003-187
9 Fr.	AJ-003-171	26 Fr.	AJ-003-188
10 Fr.	AJ-003-172	27 Fr.	AJ-003-189
11 Fr.	AJ-003-173	28 Fr.	AJ-003-190
12 Fr.	AJ-003-174	29 Fr.	AJ-003-191
13 Fr.	AJ-003-175	30 Fr.	AJ-003-192
14 Fr.	AJ-003-176	31 Fr.	AJ-003-193
15 Fr.	AJ-003-177	32 Fr.	AJ-003-194
16 Fr.	AJ-003-178	33 Fr.	AJ-003-195
17 Fr.	AJ-003-179	34 Fr.	AJ-003-196
18 Fr.	AJ-003-180	35 Fr.	AJ-003-197
19 Fr.	AJ-003-181	36 Fr.	AJ-003-198
20 Fr.	AJ-003-182	37 Fr.	AJ-003-199
21 Fr.	AJ-003-183	38 Fr.	AJ-003-200
22 Fr.	AJ-003-184	39 Fr.	AJ-003-201
23 Fr.	AJ-003-185	40 Fr.	AJ-003-202
24 Fr.	AJ-003-186	50 Fr.	AJ-003-203



HEGAR

French	Ordering no. HEGAR
3/4 Fr.	AJ-003-151
5/6 Fr.	AJ-003-152
7/8 Fr.	AJ-003-153
9/10 Fr.	AJ-003-154
11/12 Fr.	AJ-003-155
13/14 Fr.	AJ-003-156
15/16 Fr.	AJ-003-157
17/18 Fr.	AJ-003-158

9.1 | INSTRUMENTS | METAL CATHETERS



	FOR MEN	FOR WOMEN	COXETER PROSTATE
French	Ordering no.	Ordering no.	Ordering no.
6 Fr.		AJ-003-100	
7 Fr.	AJ-003-081		
8 Fr.	AJ-003-082	AJ-003-102	AJ-003-120
9 Fr.			
10 Fr.	AJ-003-084	AJ-003-104	AJ-003-121
11 Fr.			
12 Fr.	AJ-003-086	AJ-003-106	AJ-003-122
13 Fr.			
14 Fr.	AJ-003-088	AJ-003-108	AJ-003-123
15 Fr.			
16 Fr.	AJ-003-090	AJ-003-110	AJ-003-124
17 Fr.			
18 Fr.	AJ-003-092	AJ-003-112	AJ-003-125
19 Fr.			
20 Fr.	AJ-003-094	AJ-003-114	AJ-003-126
21 Fr.			
22 Fr.	AJ-003-096	AJ-003-116	AJ-003-127
23 Fr.			
24 Fr.	AJ-003-098	AJ-003-118	AJ-003-128

9.1 | INSTRUMENTS

CATHETER INTRODUCER



	Ordering No.	Description	Length
	AJ-003-142	Catheter introducer Guyon straight	39 cm
	AJ-003-143	Catheter introducer Guyon curved	39 cm

MEATUS CLAMPS

	Ordering No.	Description	Length
	AJ-003-140	Meatus clamp, Stockmann	8 cm
	AJ-003-141	Meatus clamp, Strauss	11 cm

SECRETION SPOON




	Ordering No.	Description	Length
	AJ-006-012	Secretion spoon Guttman	26.5 cm


9.1 | INSTRUMENTS | TUBING CONNECTORS

	Ordering No. Description	
	Ordering No.	Description
	AJ-010-011	Luer-Lock tubing connector, male Luer-Lock inside Ø 3.0 mm, tubing connector Ø 6.0 mm
	AJ-010-012	Luer-Lock tubing connector with stopcock
	AJ-010-013	Luer-Lock tubing connector, male Luer-Lock inside, tubing connector Ø 9.0 mm
	AJ-010-018	Stopcock left with Luer-Lock, male, tubing connector Ø 9.0 mm

9.2 | STERILIZATION | WIRE BASKETS



Ordering No.	Description	Dimension
FA-010-001	Wire basket for 1 endoscope with lid	290 × 80 × 52 mm
FA-010-002		460 × 80 × 52 mm
FA-010-003		670 × 80 × 52 mm
FA-010-004	Wire basket for 2 endoscopes with lid	290 × 80 × 52 mm
FA-010-005		460 × 80 × 52 mm



Ordering No.	Description	Dimension
FA-010-006	Wire basket for 1 endoscope and 1 cold light cable with lid	640 × 150 × 77 mm

9.2 | STERILIZATION

ENDOSCOPE CONTAINER



Ordering No.	Description	Colour code	Dimension
FA-001-001	Endoscope container, lid and bottom perforated (without filter)	Silver	450 × 70 × 72 mm
FA-001-001-R		Red	
FA-001-001-G		Green	
FA-002-001	Endoscope container, lid and bottom perforated (without filter)	Silver	250 × 60 × 50 mm
FA-002-001-R		Red	
FA-002-001-G		Green	

PROTECTION TUBE



Ordering No.	Description
FA-001-000	Protection tube for endoscopes, plastic, 1 box = 10 pieces, 19 cm

CAMERA SLEEVE



Ordering No.	Description
ES-007-001	Camera sleeve, 13 × 250 cm, PE hose film, sterile, 1 box = 100 pieces

9.3 | CLEANING AND MAINTENANCE

CLEANING PISTOL SET



Ordering No.	Description
EJ-030-010	Cleaning pistol complete inclusive 8 different tips

CLEANING BRUSH



Ordering No.	Description	Diameter
EJ-030-001	Cleaning brush, length 420 mm , 1 box = 10 pieces	Ø 3 mm
EJ-030-002	Cleaning brush, length 330 mm , 1 box = 10 pieces	Ø 5 mm
EJ-030-003		Ø 10 mm

GREASE & PASTE



Ordering No.	Description	Weight
EJ-030-031	Special grease for stopcocks, valves and instruments	3 g
EJ-030-032	Polishing paste for scopes, 1 box = 50 pieces	3 g

ARE YOU INTERESTED
IN OUR PRODUCTS?
WE WILL GLADLY
SUBMIT YOU
AN OFFER.

(if available)

Company Name:

Do you have a question about our products?

10. INDEX

AA-002-001	8 , 37 , 98 , 147	AA-007-006	125	AB-001-004	9 , 33
AA-002-001-W	8 , 37 , 98	AA-007-007	125	AB-001-005	9
AA-002-002	8 , 37	AA-007-008	125	AB-001-006	9 , 33
AA-002-003	8 , 28 , 37 , 59 , 147	AA-007-009	123	AB-001-007	25
AA-002-003-W	8 , 37	AA-007-011	127	AB-001-010	10
AA-002-004	8 , 37 , 150	AA-007-012	127	AB-001-011	161
AA-002-004-W	8 , 37	AA-007-013	124	AB-001-011-7	164
AA-002-011	8 , 37 , 98 , 106 , 147	AA-007-014	124	AB-001-011-9	163
AA-002-011-W	8 , 37 , 98	AA-007-015	124	AB-001-011-11	163
AA-002-012	8 , 37 , 59	AA-007-016	122	AB-001-012	10
AA-002-013	8 , 28 , 33 , 37 , 59 , 76 , 90 , 91 , 93 , 147	AA-007-017	122	AB-001-013	10
AA-002-013-W	8 , 37	AA-007-018	122	AB-001-014	10
AA-002-014	8 , 37	AA-007-019	122	AB-001-015	10
AA-002-014-W	8 , 37	AA-007-022	125	AB-001-016	10
AA-002-021	8 , 37 , 98	AA-007-023	125	AB-001-022	10
AA-002-022	8 , 37	AA-007-024	125	AB-001-023	10
AA-002-023	8 , 28 , 37 , 59	AA-007-025	125	AB-001-024	10
AA-002-024	8 , 37	AA-007-026	125	AB-001-025	10
AA-003-001	103	AA-007-027	125	AB-001-026	10
AA-003-002	25 , 48 , 65 , 70	AA-007-028	123	AB-001-040	28
AA-003-004	65	AA-007-029	123	AB-001-041	28
AA-003-011	103	AA-007-030	124	AB-001-042	29
AA-003-012	25 , 48 , 65 , 70 , 79 , 83	AA-007-031	124	AB-001-043	29
AA-003-012-W	25	AA-007-032	124	AB-002-001	9
AA-003-014	65	AA-007-033	124	AB-002-002	9
AA-004-011	163 , 170 , 178	AA-007-034	124	AB-002-003	9
AA-004-012	163	AA-007-035	124	AB-002-004	9
AA-004-013	164	AA-007-036	122	AB-002-005	9
AA-004-014	164	AA-007-037	122	AB-002-010	148
AA-005-001	161 , 167 , 175	AA-007-041	129	AB-002-011	149
AA-005-003	161	AA-007-042	129	AB-002-012	149
AA-005-004	161	AA-007-043	128	AB-002-013	148
AA-006-001	114	AA-007-044	128	AB-002-014	149
AA-006-002	114	AA-007-045	129	AB-003-001	162
AA-006-004	114	AA-007-046	129	AB-003-002	11
AA-006-005	114	AA-007-047	128	AB-003-003	11 , 33
AA-006-006	116	AA-007-048	128	AB-003-004	11 , 26
AA-006-011	116	AA-007-049	129	AB-003-005	162
AA-006-012	116 , 153	AA-007-050	129	AB-003-006	12
AA-006-021	114	AA-007-051	128	AB-003-007	12 , 33
AA-006-022	114	AA-007-052	128	AB-003-008	12
AA-006-031	116	AA-007-053	129	AB-003-009	12
AA-006-032	116	AA-007-054	129	AB-003-010	26
AA-006-034	116	AA-007-055	128	AB-003-011	26
AA-006-050	111	AA-007-056	128	AB-003-012	26
AA-006-051	111	AA-007-060	30	AB-003-013	27
AA-006-052	111	AA-008-002	31	AB-003-014	27
AA-006-053	111	AA-008-003	142	AB-003-015	11
AA-006-069	109	AA-008-100	32 , 143	AB-003-020	27
AA-006-070	109	AA-008-101	143	AB-004-001	13 , 34 , 92 , 94
AA-006-072	156	AA-008-102	32 , 143	AB-004-002	13
AA-006-080	114 , 116 , 153	AA-008-103	31 , 32 , 142	AB-004-003	13 , 34 , 92 , 94
AA-006-081	109 , 114 , 116 , 153	AA-008-104	32 , 143	AB-004-004	13
AA-006-083	109 , 114 , 116 , 153	AA-008-105	31 , 32 , 142	AB-004-005	13
AA-006-084	156	AB-001-000	25	AB-004-006	13
AA-006-100	110	AB-001-001	161	AB-004-007	14 , 34
AA-007-001	126	AB-001-001-7	164	AB-004-008	13 , 92 , 94
AA-007-002	126	AB-001-001-9	163	AB-004-009	14
AA-007-003	122	AB-001-001-11	163	AB-004-009-QL	14
AA-007-004	122	AB-001-001-13	161	AB-004-010	13
AA-007-005	122	AB-001-002	9	AB-004-010-QL	153
		AB-001-003	9	AB-004-011	150

10. INDEX

AB-004-012	150	AB-006-443	113	AC-001-212-T	135
AB-006-001	115	AB-006-444	113	AC-001-213	130
AB-006-002	115	AB-006-445	113	AC-001-214	130
AB-006-003	115	AB-006-446	113	AC-001-214-OT	135
AB-006-004	115	AB-007-001	87	AC-001-214-T	135
AB-006-005	115	AB-007-002	87	AC-001-215	130
AB-006-006	115	AB-007-003	87	AC-001-215-OT	135
AB-006-007	115	AB-007-004	88	AC-001-215-T	135
AB-006-008	115	AB-007-005	88	AC-001-221	165
AB-006-011	118	AB-007-006	88	AC-002-001	15
AB-006-012	118	AB-007-007	88	AC-002-001-R	15
AB-006-013	118	AB-007-008	87	AC-002-002	15
AB-006-020	118	AC-001-001	15	AC-002-002-R	15
AB-006-021	118	AC-001-001-R	15	AC-002-003	15
AB-006-022	118	AC-001-002	15, 33	AC-002-003-R	15
AB-006-023	118	AC-001-002-R	15	AC-002-004	15
AB-006-024	118	AC-001-003	15	AC-002-004-R	15
AB-006-025	118	AC-001-003-R	15	AC-002-005	15
AB-006-026	118	AC-001-004	15	AC-002-005-R	15
AB-006-027	118	AC-001-004-R	15	AC-002-011	130
AB-006-040	119	AC-001-005	15	AC-002-011-R	130
AB-006-041	119	AC-001-005-R	15	AC-002-012	130
AB-006-042	119	AC-001-008	15	AC-002-012-R	130
AB-006-043	119	AC-001-008-R	15	AC-002-013	130
AB-006-044	119	AC-001-010	130	AC-002-013-R	130
AB-006-045	119	AC-001-010-R	130	AC-002-014	130
AB-006-046	119	AC-001-011	130	AC-002-014-R	130
AB-006-047	119	AC-001-011-R	130	AC-002-015	130
AB-006-048	119	AC-001-012	130	AC-002-015-R	130
AB-006-049	119	AC-001-012-R	130	AC-002-201	15
AB-006-050	119	AC-001-013	130	AC-002-202	15
AB-006-051	119	AC-001-013-R	130	AC-002-203	15
AB-006-052	118	AC-001-014	130	AC-002-204	15
AB-006-053	118	AC-001-014-R	130	AC-002-205	15
AB-006-061	117	AC-001-015	130	AC-002-211	130
AB-006-062	117	AC-001-015-R	130	AC-002-212	130
AB-006-063	117	AC-001-021	165	AC-002-213	130
AB-006-064	117	AC-001-031	144	AC-002-214	130
AB-006-065	117	AC-001-201	15, 130	AC-002-215	130
AB-006-066	117	AC-001-201-OT	19, 135, 137	AC-003-000	17
AB-006-067	117	AC-001-201-T	19, 135, 137	AC-003-000-R	17
AB-006-068	117	AC-001-202	15, 130	AC-003-001	17
AB-006-100	112	AC-001-202-OT	19, 135, 137	AC-003-001-R	17
AB-006-101	112	AC-001-202-T	19, 135, 137	AC-003-002	17
AB-006-102	112	AC-001-203	15, 130	AC-003-002-R	17
AB-006-110	112	AC-001-203-OT	19, 135, 137	AC-003-003	17
AB-006-111	112	AC-001-203-T	19, 135, 137	AC-003-003-R	17
AB-006-112	112	AC-001-204	15, 130	AC-003-004	17
AB-006-115	112	AC-001-204-OT	19, 135, 137	AC-003-004-R	17
AB-006-170	109	AC-001-204-T	19, 135, 137	AC-003-005	17
AB-006-171	109	AC-001-205	15, 130	AC-003-005-R	17
AB-006-242	113	AC-001-206-OT	19, 135, 137	AC-003-010	132
AB-006-243	113	AC-001-206-T	19, 135, 137	AC-003-010-R	132
AB-006-244	113	AC-001-208	15, 130	AC-003-011	132
AB-006-245	113	AC-001-210	130	AC-003-011-R	132
AB-006-246	113	AC-001-210-OT	135	AC-003-012	132
AB-006-342	113	AC-001-210-T	135	AC-003-012-R	132
AB-006-343	113	AC-001-211	130	AC-003-013	132
AB-006-344	113	AC-001-211-OT	135	AC-003-013-R	132
AB-006-345	113	AC-001-211-T	135	AC-003-014	132
AB-006-346	113	AC-001-212	130	AC-003-014-R	132
AB-006-442	113	AC-001-212-OT	135	AC-003-015	132

10. INDEX

AC-003-015-R	132	AC-004-624-T	134	AC-007-002	16
AC-003-021	165	AC-005-004	17	AC-007-002-R	16
AC-003-200	17	AC-005-004-R	17	AC-007-011	131
AC-003-200-OT	19	AC-005-014	132	AC-007-011-R	131
AC-003-200-T	19	AC-005-014-R	132	AC-007-012	131
AC-003-201	17	AC-005-204	17	AC-007-012-R	131
AC-003-201-OT	19	AC-005-214	132	AC-007-201	16
AC-003-201-T	19	AC-006-001	15	AC-007-201-OT	19
AC-003-202	17	AC-006-001-R	15	AC-007-201-T	19
AC-003-202-OT	19	AC-006-002	15	AC-007-202	16
AC-003-202-T	19	AC-006-002-R	15	AC-007-202-OT	19
AC-003-203	17	AC-006-003	15	AC-007-202-T	19
AC-003-204	17	AC-006-003-R	15	AC-007-211	131
AC-003-204-OT	19	AC-006-004	15	AC-007-211-OT	135
AC-003-204-T	19	AC-006-004-R	15	AC-007-211-T	135
AC-003-205	17	AC-006-005	15	AC-007-212	131
AC-003-205-OT	19	AC-006-005-R	15	AC-007-212-OT	135
AC-003-205-T	19	AC-006-011	130	AC-007-212-T	135
AC-003-210	132	AC-006-011-R	130	AC-008-004	16
AC-003-210-OT	134	AC-006-012	130	AC-008-004-R	16
AC-003-210-T	134	AC-006-012-R	130	AC-008-005	16
AC-003-211	132	AC-006-013	130	AC-008-005-R	16
AC-003-211-OT	134	AC-006-013-R	130	AC-008-014	131
AC-003-211-T	134	AC-006-014	130	AC-008-014-R	131
AC-003-212	132	AC-006-014-R	130	AC-008-015	131
AC-003-212-OT	134	AC-006-015	130	AC-008-015-R	131
AC-003-212-T	134	AC-006-015-R	130	AC-008-204	16
AC-003-213	132	AC-006-021	165	AC-008-205	16
AC-003-214	132	AC-006-031	144	AC-008-214	131
AC-003-214-OT	134	AC-006-201	15	AC-008-215	131
AC-003-214-T	134	AC-006-201-OT	19	AC-009-001	16
AC-003-215	132	AC-006-201-T	19	AC-009-001-R	16
AC-003-215-OT	134	AC-006-202	15	AC-009-002	16
AC-003-215-T	134	AC-006-202-OT	19	AC-009-002-R	16
AC-003-221	165	AC-006-202-T	19	AC-009-004	16
AC-003-302	17	AC-006-203	15	AC-009-004-R	16
AC-003-304	17	AC-006-203-OT	19	AC-009-011	131
AC-003-612-OT	19	AC-006-203-T	19	AC-009-011-R	131
AC-003-612-T	19	AC-006-204	15	AC-009-012	131
AC-003-614-OT	19	AC-006-204-OT	19	AC-009-012-R	131
AC-003-614-T	19	AC-006-204-T	19	AC-009-014	131
AC-004-004	17	AC-006-205	15	AC-009-014-R	131
AC-004-004-R	17	AC-006-205-OT	19	AC-009-201	16
AC-004-005	17	AC-006-205-T	19	AC-009-201-OT	19
AC-004-005-R	17	AC-006-211	130	AC-009-201-T	19
AC-004-014	132	AC-006-211-OT	135	AC-009-202	16
AC-004-015	132	AC-006-211-T	135	AC-009-202-OT	19
AC-004-204	17	AC-006-212	130	AC-009-202-T	19
AC-004-204-OT	19	AC-006-212-OT	135	AC-009-204	16
AC-004-204-T	19	AC-006-212-T	135	AC-009-211	131
AC-004-205	17	AC-006-213	130	AC-009-211-OT	135
AC-004-205-OT	19	AC-006-213-OT	135	AC-009-211-T	135
AC-004-205-T	19	AC-006-213-T	135	AC-009-212	131
AC-004-214	132	AC-006-214	130	AC-009-212-OT	135
AC-004-214-OT	134	AC-006-214-OT	135	AC-009-212-T	135
AC-004-214-T	134	AC-006-214-T	135	AC-009-214	131
AC-004-215	132	AC-006-215	130	AC-012-001	16
AC-004-215-OT	134	AC-006-215-OT	135	AC-012-002	16
AC-004-215-T	134	AC-006-215-T	135	AC-012-004	16
AC-004-622-OT	134	AC-006-221	165	AC-012-011	131
AC-004-622-T	134	AC-007-001	16	AC-012-012	131
AC-004-624-OT	134	AC-007-001-R	16	AC-012-014	131

10. INDEX

AC-012-015	144	AC-019-011	131	AC-103-215-T	136
AC-015-002	15	AC-019-011-R	131	AC-103-612-OT	20
AC-015-002-R	15	AC-019-012	131	AC-103-612-T	20
AC-015-003-R	15	AC-019-012-R	131	AC-103-614-OT	20
AC-015-004	15	AC-019-021	165	AC-103-614-T	20
AC-015-004-R	15	AC-019-201	16	AC-104-204-OT	20
AC-015-012	130	AC-019-202	16	AC-104-204-T	20
AC-015-012-R	130	AC-019-211	131	AC-104-205-OT	20
AC-015-013-R	130	AC-019-212	131	AC-104-205-T	20
AC-015-014	130	AC-020-002	16	AC-104-214-OT	136
AC-015-014-R	130	AC-020-002-R	16	AC-104-214-T	136
AC-015-021	165	AC-020-004	16	AC-104-215-OT	136
AC-015-202	15	AC-020-004-R	16	AC-104-215-T	136
AC-015-203	15	AC-020-012	131	AC-106-201-OT	20
AC-015-204	15	AC-020-012-R	131	AC-106-201-T	20
AC-015-212	130	AC-020-014	131	AC-106-202-OT	20
AC-015-213	130	AC-020-014-R	131	AC-106-202-T	20
AC-015-214	130	AC-020-202	16	AC-106-203-OT	20
AC-015-221	165	AC-020-204	16	AC-106-203-T	20
AC-016-001	16	AC-020-212	131	AC-106-204-OT	20
AC-016-001-R	16	AC-020-214	131	AC-106-204-T	20
AC-016-002	16	AC-100-001	19, 20, 134, 135, 136, 137	AC-106-205-OT	20
AC-016-002-R	16			AC-106-205-T	20
AC-016-011	131	AC-101-201-OT	20	AC-106-211-OT	137
AC-016-011-R	131	AC-101-201-T	20	AC-106-211-T	137
AC-016-012	131	AC-101-202-OT	20	AC-106-212-OT	137
AC-016-012-R	131	AC-101-202-T	20	AC-106-212-T	137
AC-016-201	16	AC-101-203-OT	20	AC-106-213-OT	137
AC-016-202	16	AC-101-203-T	20	AC-106-213-T	137
AC-016-211	131	AC-101-204-OT	20	AC-106-214-OT	137
AC-016-211-OT	135	AC-101-204-T	20	AC-106-214-T	137
AC-016-211-T	135	AC-101-206-OT	20	AC-106-215-OT	137
AC-016-212	131	AC-101-206-T	20	AC-106-215-T	137
AC-016-212-OT	135	AC-101-210-OT	137	AC-107-201-OT	20
AC-016-212-T	135	AC-101-210-T	137	AC-107-201-T	20
AC-017-001	15	AC-101-211-OT	137	AC-107-202-OT	20
AC-017-001-R	15	AC-101-211-T	137	AC-107-202-T	20
AC-017-004	15	AC-101-212-OT	137	AC-107-211-OT	137
AC-017-004-R	15	AC-101-212-T	137	AC-107-211-T	137
AC-017-011	130	AC-101-214-OT	137	AC-107-212-OT	137
AC-017-014	130	AC-101-214-T	137	AC-107-212-T	137
AC-017-014-R	130	AC-101-215-OT	137	AC-109-201-OT	20
AC-017-201	15	AC-101-215-T	137	AC-109-201-T	20
AC-017-204	15	AC-103-200-OT	20	AC-109-211-OT	137
AC-017-211	130	AC-103-200-T	20	AC-109-211-T	137
AC-017-214	130	AC-103-201-OT	20	AC-109-212-OT	137
AC-018-001	17	AC-103-201-T	20	AC-109-212-T	137
AC-018-001-R	17	AC-103-202-OT	20	AC-114-212-OT	136
AC-018-002	17	AC-103-202-T	20	AC-114-212-T	136
AC-018-002-R	17	AC-103-204-OT	20	AC-116-211-OT	137
AC-018-011	132	AC-103-204-T	20	AC-116-211-T	137
AC-018-011-R	132	AC-103-205-OT	20	AC-116-212-OT	137
AC-018-012	132	AC-103-205-T	20	AC-116-212-T	137
AC-018-012-R	132	AC-103-210-OT	136	AD-001-001	23, 140
AC-018-201	17	AC-103-210-T	136	AD-001-002	23, 140
AC-018-202	17	AC-103-211-OT	136	AD-001-003	23, 34, 140
AC-018-211	132	AC-103-211-T	136	AD-001-004	23, 140
AC-018-212	132	AC-103-212-OT	136	AD-001-005	23, 140
AC-019-001	16	AC-103-212-T	136	AD-001-006	23, 140
AC-019-001-R	16	AC-103-214-OT	136	AD-001-007	23, 140
AC-019-002	16	AC-103-214-T	136	AD-001-008	23, 140
AC-019-002-R	16	AC-103-215-OT	136	AD-001-011	23, 140

10. INDEX

AD-001-012	23, 140	AE-001-037	63	AF-001-004	21, 138
AD-001-013	23, 140	AE-001-038	63	AF-001-005	21, 138
AD-001-014	23, 140	AE-001-039	64	AF-001-006	21, 138
AD-001-015	23, 140	AE-001-040	63	AF-001-007	22, 139
AD-001-016	23, 140	AE-001-041	63	AF-001-008	22, 34, 139
AD-001-017	23, 140	AE-001-042	63	AF-001-009	22, 139
AD-001-018	23, 140	AE-001-043	63	AF-001-010	22, 139
AD-001-021	165	AE-001-044	64	AF-001-011	22, 139
AD-002-001	23, 140	AE-001-045	64	AF-001-012	22, 139
AD-002-002	23, 140	AE-001-046	64	AF-001-013	22, 139
AD-002-003	23, 34, 140	AE-001-047	64	AF-001-014	22, 139
AD-002-005	23, 140	AE-001-048	63	AF-001-015	21, 138
AD-002-011	23, 140	AE-002-001	46	AF-001-016	21, 138
AD-002-012	23, 140	AE-002-002	46	AF-001-017	21, 138
AD-002-013	23, 140	AE-002-003	46	AF-001-018	22, 139
AD-002-015	23, 140	AE-002-004	46	AF-001-019	21, 138
AD-002-021	165	AE-002-005	46	AF-001-021	21, 22, 138, 139
AD-003-003	23, 140	AE-002-006	47	AF-001-022	21, 138
AD-003-005	23, 140	AE-002-007	47	AF-001-023	21, 138
AD-003-007	23, 140	AE-002-010	46	AF-001-024	21, 22, 138, 139
AD-003-008	23, 140	AE-002-011	46	AF-001-025	22, 139
AD-003-013	23, 140	AE-002-012	46	AG-001-001	38, 91, 101, 106
AD-003-015	23, 140	AE-002-013	46	AG-001-001-HF	38, 101
AD-003-017	23, 140	AE-002-014	47	AG-001-002	38, 101
AD-003-018	23, 140	AE-002-015	47	AG-001-003	39
AD-010-001	24, 57, 92, 141, 166, 173	AE-002-016	47	AG-001-003-HF	39
AD-010-002	24, 57, 141, 166, 173	AE-002-031	64	AG-001-004	39
AD-010-003	24, 34, 57, 141, 166, 173	AE-003-001	55	AG-001-006	77
AD-010-006	24, 57, 141, 166, 173	AE-003-003	55	AG-001-007	40
AD-010-007	24, 57, 141, 166, 173	AE-003-004	55	AG-001-008	40
AD-010-008	24, 57, 141, 166, 173	AE-003-010	55	AG-001-009	77
AD-010-011	75	AE-003-014	55	AG-001-010	104
AD-010-012	75	AE-003-016	55	AG-001-011	60, 93
AD-010-013	75, 94	AE-003-031	69	AG-001-012	60
AD-010-015	75	AE-003-032	69	AG-001-017	41
AD-010-024	24, 57, 141, 166, 173	AE-003-033	69	AG-001-018	41
AD-010-025	58, 174	AE-003-034	69	AG-002-001	49
AD-010-026	58, 174	AE-003-036	69	AG-002-002	49
AE-001-001	46, 91	AE-003-040	69	AG-002-003	50
AE-001-002	46	AE-003-041	74	AG-002-004	50
AE-001-003	46	AE-003-043	74	AG-002-006	51
AE-001-004	46, 91	AE-003-044	74	AG-002-007	80
AE-001-005	46	AE-003-045	74	AG-002-008	80
AE-001-006	47	AE-003-046	74	AG-002-011	66
AE-001-007	47	AE-004-002	169	AG-002-012	66
AE-001-008	47	AE-004-004	169	AG-002-013	71
AE-001-009	47	AE-004-016	169	AG-002-015	71
AE-001-010	46	AE-004-018	169	AG-002-016	84
AE-001-011	46	AE-004-019	169	AG-003-001	168, 176
AE-001-012	46	AE-005-002	172	AG-004-001	171, 179
AE-001-013	46	AE-005-004	172	AG-005-001	171
AE-001-014	47	AE-005-016	172	AG-100-006	77
AE-001-015	47	AE-005-018	172	AG-100-007	77
AE-001-016	47, 91	AE-005-019	172	AH-001-001	44
AE-001-017	47	AE-006-002	172	AH-001-002	44
AE-001-030	63	AE-006-004	172	AH-001-011	45, 62, 78, 90, 91, 93,
AE-001-031	63, 93	AE-006-016	172		155
AE-001-032	63, 93	AE-006-018	172	AH-001-012	45, 62
AE-001-033	63	AE-006-019	172	AH-001-013	45, 62, 155
AE-001-034	64, 93	AF-001-001	21, 138	AH-001-015	45, 62
AE-001-035	64	AF-001-002	21, 138	AH-001-016	45, 62
AE-001-036	63	AF-001-003	21, 138	AH-001-017	45, 62

10. INDEX

AH-001-018	45, 62	AI-004-010	180	AJ-003-046	183
AH-001-101	44	AI-004-011	180	AJ-003-050	183
AH-001-102	44	AI-005-001	103	AJ-003-051	183
AH-001-201	42	AI-005-003	103	AJ-003-052	183
AH-001-202	42	AI-005-004	103	AJ-003-053	183
AH-001-203	42	AI-005-007	105	AJ-003-054	183
AH-001-204	42	AI-005-008	105	AJ-003-055	183
AH-001-205	42	AI-005-011	105	AJ-003-056	183
AH-001-211	42	AJ-001-000	89	AJ-003-057	183
AH-001-212	42	AJ-001-001	89	AJ-003-058	183
AH-001-301	43, 61, 78, 90, 91, 93, 154, 155	AJ-001-002	89, 92, 94	AJ-003-059	183
AH-001-302	43, 61, 78, 91, 154	AJ-001-003	89	AJ-003-060	183
AH-001-303	43, 61	AJ-001-010	89	AJ-003-061	183
AH-001-304	43, 61	AJ-001-011	89	AJ-003-062	183
AH-001-305	43, 61	AJ-001-012	89	AJ-003-063	183
AH-001-306	43, 61	AJ-001-013	89	AJ-003-064	183
AH-001-308	43, 78, 155	AJ-002-001	89, 92	AJ-003-065	183
AH-001-311	43, 61	AJ-002-002	89	AJ-003-066	183
AH-001-312	43, 61	AJ-002-003	89	AJ-003-067	183
AH-002-001	54	AJ-002-004	89	AJ-003-068	183
AH-002-011	56, 68, 82, 158	AJ-002-005	89	AJ-003-069	183
AH-002-013	56, 68, 82, 158	AJ-003-001	183	AJ-003-070	183
AH-002-015	73, 86	AJ-003-002	183	AJ-003-071	183
AH-002-016	73, 86	AJ-003-003	183	AJ-003-072	183
AH-002-201	52	AJ-003-004	183	AJ-003-074	183
AH-002-211	72, 85	AJ-003-005	183	AJ-003-081	185
AH-002-212	52	AJ-003-006	183	AJ-003-082	185
AH-002-213	72, 85	AJ-003-007	183	AJ-003-084	185
AH-002-216	72, 85	AJ-003-008	183	AJ-003-086	185
AH-002-301	53, 67, 81, 157	AJ-003-009	183	AJ-003-088	185
AH-002-302	53, 67, 81, 157	AJ-003-010	183	AJ-003-090	185
AH-003-001	167	AJ-003-011	183	AJ-003-092	185
AH-003-011	167	AJ-003-012	183	AJ-003-094	185, 190
AH-004-001	170	AJ-003-013	183	AJ-003-096	185
AH-004-011	170	AJ-003-014	183	AJ-003-098	185
AH-005-001	170	AJ-003-015	183	AJ-003-100	185
AH-005-011	170	AJ-003-016	183	AJ-003-102	185
AI-001-001	98, 106	AJ-003-021	183	AJ-003-104	185
AI-001-002	99, 106	AJ-003-022	183	AJ-003-106	185
AI-001-003	98, 106	AJ-003-023	183	AJ-003-108	185
AI-001-004	100, 106	AJ-003-024	183	AJ-003-110	185
AI-001-005	100	AJ-003-025	183	AJ-003-112	185
AI-001-006	99	AJ-003-026	183	AJ-003-114	185
AI-001-007	102, 106	AJ-003-027	183	AJ-003-116	185
AI-001-008	102	AJ-003-028	183	AJ-003-118	185
AI-001-009	102	AJ-003-029	183	AJ-003-120	185
AI-001-011	102	AJ-003-030	183	AJ-003-121	185
AI-001-012	100	AJ-003-031	183	AJ-003-122	185
AI-002-001	97	AJ-003-032	183	AJ-003-123	185
AI-002-002	97	AJ-003-033	183	AJ-003-124	185
AI-003-001	175	AJ-003-034	183	AJ-003-125	185
AI-003-003	168, 176	AJ-003-035	183	AJ-003-126	185
AI-003-004	175	AJ-003-036	183	AJ-003-127	185
AI-003-007	177	AJ-003-037	183	AJ-003-128	185
AI-003-009	177	AJ-003-038	183	AJ-003-140	186
AI-003-010	177	AJ-003-039	183	AJ-003-141	186
AI-003-011	177	AJ-003-040	183	AJ-003-142	186
AI-004-001	178	AJ-003-041	183	AJ-003-143	186
AI-004-004	178	AJ-003-042	183	AJ-003-151	184
AI-004-007	180	AJ-003-043	183	AJ-003-152	184
AI-004-009	180	AJ-003-044	183	AJ-003-153	184
		AJ-003-045	183	AJ-003-154	184

10. INDEX

AJ-003-155	184	AJ-003-192	184	AJ-010-020	90
AJ-003-156	184	AJ-003-193	184	CI-001-030	58, 174
AJ-003-157	184	AJ-003-194	184	CI-001-040	58, 174
AJ-003-158	184	AJ-003-195	184	EJ-030-001	190
AJ-003-170	184	AJ-003-196	184	EJ-030-002	190
AJ-003-171	184	AJ-003-197	184	EJ-030-003	190
AJ-003-172	184	AJ-003-198	184	EJ-030-010	190
AJ-003-173	184	AJ-003-199	184	EJ-030-031	190
AJ-003-174	184	AJ-003-200	184	EJ-030-032	190
AJ-003-175	184	AJ-003-201	184	EK-015-060	24, 57, 141, 166, 173
AJ-003-176	184	AJ-003-202	184	ES-007-001	189
AJ-003-177	184	AJ-003-203	184	FA-001-000	189
AJ-003-178	184	AJ-006-012	186	FA-001-001	189
AJ-003-179	184	AJ-010-001	11, 26, 111, 122, 123, 124, 125, 126, 127, 128, 129	FA-001-001-G	189
AJ-003-180	184			FA-001-001-R	189
AJ-003-181	184	AJ-010-003	26, 111, 122, 123, 124, 125, 126, 127, 128, 129	FA-002-001	189
AJ-003-182	184			FA-002-001-G	189
AJ-003-183	184	AJ-010-006	47, 55, 64, 69, 74, 90, 92, 94, 102, 105, 106, 169, 172, 177, 180	FA-002-001-R	189
AJ-003-184	184			FA-010-001	188
AJ-003-185	184			FA-010-002	188
AJ-003-186	184			FA-010-003	188
AJ-003-187	184	AJ-010-011	187	FA-010-004	188
AJ-003-188	184	AJ-010-012	187	FA-010-005	188
AJ-003-189	184	AJ-010-013	187	FA-010-006	188
AJ-003-190	184	AJ-010-015	114, 116		
AJ-003-191	184	AJ-010-018	187		

WE COVER THE FOLLOWING AREAS IN SURGERY WITH HIGH QUALITY PRODUCTS

ARTHROSCOPY

CARDIAC MIS

CONTAINER

ENT

HF-SURGERY

GASTROSCOPY

GYNAECOLOGY

LAMINECTOMY

LAPAROSCOPY

PROCTOSCOPY

SPINE SURGERY

UNITS

UROLOGY



STEMA Medizintechnik GmbH | Take-off Gewerbepark 74 | 78579 Neuhausen ob Eck
info@stema-medizintechnik.de | www.stema-medizintechnik.de | +49 7467 79630 0

© 2026 STEMA Medizintechnik GmbH. All rights reserved.
Specifications are subject to change without notice.