
STEMA

UROLOGY





UROLOGY



WWW.STEMA-MEDIZINTECHNIK.DE

ABOUT US

SINCE 20 YEARS

EXPORT IN +70 COUNTRIES

24H/SERVICE

+ 10.000 PRODUCTS

LOCATED IN THE WORLD CENTER FOR MEDICAL TECHNOLOGY

CERTIFICATED

STEMA Medizintechnik GmbH is certified in accordance with the highest standards:
- EN ISO 13485:2016 + AC:2018 + A11:2021- ISO 13485:2016
- EU 2017/745 **(MDR)**
- Council Directive 93/42/EEC Annex II excluding Section 4
- Council Directive 93/42/EEC Annex V excluding Section 3
- FDA

These standards create an essential base for our company and our partners. They also symbolize the commitment and strive to keep up the quality of our products.

STEMA

PRODUCT RANGE

The current product range of STEMA Medizintechnik GmbH includes more than 10.000 different items.
We cover the following areas in surgery with high quality products:

ARTHROSCOPY
CARDIAC MIS
CONTAINER-SYSTEMS
GASTROENTEROLOGY
GYNAECOLOGY
HF-SURGERY
LAMINECTOMY
LAPAROSCOPY
MICRO MOTOR SYSTEMS
PROCTOSCOPY
SPINE-SURGERY
TITANIUM PLATING SYSTEM
UNITS
UROLOGY



STEMA

CONTENT
UROLOGY

PRODUCT	PAGE
1. CYSTO-URETHROSCOPY	
1.1 | Cysto-Urethroscopy	8
1.2 | Slim Cysto-Urethroscopy	25
1.3 | Laser Cysto-Urethroscopy	28
1.4 | HD Mini Compact Cystoscope	30
1.5 | Flexible Cystoscope	31
1.6 | Cystoscopy Set	33
2. RESECTOSCOPY	
2.1 | Monopolar Resectoscopy	37
2.2 | Saline Resectoscopy	59
2.3 | Laser Resectoscopy	76
2.4 | Accessories	87
2.5 | Resectoscopy Sets	91
3. URETHROTOMY	
3.1 | Urethrotome	97
3.2 | Optical Urethrotomy	98
3.3 | Urethrotomy Set	106
4. NEPHROSCOPY	
4.1 | HD Nephroscopy	109
4.2 | Mini Nephroscopy	111
4.3 | Standard Nephroscopy	114
4.4 | Accessories	118
5. URETERO-RENOSCOPY	
5.1 | Uretero-Renoscopes semi-rigid	122
5.2 | Flexible Uretero-Renoscope	142

CONTENT
UROLOGY

PRODUCT	PAGE
6. LITHOTRIPSY	
6.1 | Lithotripsy	147
7. MORCELLATION	
7.1 | Morcellation 26 Fr.	153
7.2 | Morcellation 22 Fr.	155
8. PEDIATRIC UROLOGY	
8.1 | Pediatric Cysto-Urethroscopy	159
8.2 | Pediatric Resectoscopy	166
8.3 | Pediatric Urethrotomy	174
9. ACCESSORIES	
9.1 | Instruments	181
9.2 | Sterilization	186
9.3 | Cleaning and Maintenance	188
10. INDEX	
10.1 | Index	191

1. CYSTO-URETHROSCOPY

1.1 | CYSTO-URETHROSCOPY
CYSTOSCOPES (4.0 MM)
SHEATHS
OBTURATORS
ENDOSCOPE BRIDGES
WORKING ELEMENTS
FORCEPS, SCISSORS & INSTRUMENTS
STONE EXTRACTORS
MONOPOLAR FLEXIBLE ELECTRODES
MONOPOLAR CONNECTING CABLES

1.2 | SLIM CYSTO-URETHROSCOPY
CYSTOSCOPES (2.9 MM)
SHEATHS
ENDOSCOPE BRIDGES
WORKING ELEMENTS

1.3 | LASER CYSTO-URETHROSCOPY
CYSTOSCOPES (4.0 MM)
SHEATHS
OBTURATORS


1.4 | HD MINI COMPACT CYSTOSCOPES

1.5 | FLEXIBLE CYSTOSCOPE
CYSTOSCOPE FLEX I
ACCESSORIES

1.6 | CYSTO-URETHROSCOPY SET


1. CYSTO-URETHROSCOPY

4.0 mm 2.9 mm



1 CYSTOSCOPES

2 SHEATHS & OBTURATORS



3 ACCESSORIES




STEMA

4.0 mm

1.1 | CYSTO-URETHROSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 0°	Cystoscope autoclavable Cystoscope autoclavable **wide angle**	Ø 4.0 mm	302 mm	AA-002-001 AA-002-001-W
	HD Cystoscope autoclavable **HD** Cystoscope autoclavable **wide angle**			AA-002-011 AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-021

Viewing angle	Description	Diameter	Working length	Ordering No.
 12°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-002
	HD Cystoscope autoclavable			AA-002-012
	Ultra HD (4K) Cystoscope autoclavable			AA-002-022

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable Cystoscope autoclavable **wide angle**	Ø 4.0 mm	302 mm	AA-002-003 AA-002-003-W
	HD Cystoscope autoclavable **HD** Cystoscope autoclavable **wide angle**			AA-002-013 AA-002-013-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-023

Viewing angle	Description	Diameter	Working length	Ordering No.
 70°	Cystoscope autoclavable Cystoscope autoclavable **wide angle**	Ø 4.0 mm	302 mm	AA-002-004 AA-002-004-W
	HD Cystoscope autoclavable **HD** Cystoscope autoclavable **wide angle**			AA-002-014 AA-002-014-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-024

4.0 mm

1.1 | CYSTO-URETHROSCOPY
SHEATHS



Ordering No.	Description	French	Colour code
AB-001-002	Cysto-Urethroscope sheath **with 2 fixed stopcocks** and standard obturator	**17 Fr.**	Yellow
AB-001-003		**19 Fr.**	Green
AB-001-004		**21 Fr.**	Red
AB-001-005		**23 Fr.**	Blue
AB-001-006		**25 Fr.**	White



Ordering No.	Description	French	Colour code
AB-002-001	Cysto-Urethroscope sheath **with central valve** and standard obturator	**17 Fr.**	Yellow
AB-002-002		**19 Fr.**	Green
AB-002-003		**21 Fr.**	Red
AB-002-004		**23 Fr.**	Blue
AB-002-005		**25 Fr.**	White

4.0 mm

1.1 | CYSTO-URETHROSCOPY
OBTURATORS | ADAPTER




Ordering No.	Description	French	Colour code
AB-001-012	Standard obturator for cysto-urethroscope sheaths	**17 Fr.**	Yellow
AB-001-013		**19 Fr.**	Green
AB-001-014		**21 Fr.**	Red
AB-001-015		**23 Fr.**	Blue
AB-001-016		**25 Fr.**	White



Ordering No.	Description	French	Colour code
AB-001-022	Viewing obturator for Cysto-Urethroscope sheaths	**17 Fr.**	Yellow
AB-001-023		**19 Fr.**	Green
AB-001-024		**21 Fr.**	Red
AB-001-025		**23 Fr.**	Blue
AB-001-026		**25 Fr.**	White

Ordering No.	Description
AB-001-010	Cystoscope **adapter** for female urethroscopy



4.0 mm

1.1 | CYSTO-URETHROSCOPY
ENDOSCOPE BRIDGES




Ordering No.	Description
AB-003-002	Endoscope bridge **diagnostic**



Ordering No.	Description
AB-003-003	Endoscope bridge **with 1 instrument channel**



Ordering No.	Description
AB-003-004	Endoscope bridge **with 2 instrument channels**



Ordering No.	Description
AB-003-015	Endoscope bridge **with 2 instrument channels, short**



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel

4.0 mm

1.1 | CYSTO-URETHROSCOPY
WORKING ELEMENTS



Ordering No.	Description
AB-003-006	Working element with ALBARRAN deflector and **1 instrument channel**
AB-003-008	Working element with ALBARRAN deflector and **1 instrument channel, without ratchet**



Ordering No.	Description
AB-003-007	Working element with ALBARRAN deflector and **2 instrument channels**
AB-003-009	Working element with ALBARRAN deflector and **2 instrument channels, without ratchet**

SHEATH/INSTRUMENT COMBINATIONS

Sheaths	Colour	Albarran		Bridge	
		1 working channel	2 working channels	1 working channel	2 working channels
25 FR.	○	1 × 10 Fr.	2 × 8 Fr.	1 × 12 Fr.	2 × 8 Fr.
23 FR.		1 × 9 Fr.	2 × 6 Fr.	1 × 10 Fr.	2 × 7 Fr.
21 FR.		1 × 6 Fr.	2 × 5 Fr.	1 × 7 Fr.	2 × 6 Fr.
19 FR.		1 × 5 Fr.	2 × 4 Fr.	1 × 6 Fr.	2 × 5 Fr.
17 FR.		-	-	1 × 5 Fr.	-

4.0 mm

1.1 | CYSTO-URETHROSCOPY
OPTICAL FORCEPS AND SCISSORS



Info: Only for sheaths 19-25 Fr.

	Ordering No.	For scopes	Description
	AB-004-005	**12°/ 30°**	Biopsy forceps (down), double action
	AB-004-006	**70°**	Biopsy forceps (down), double action
	AB-004-001	**0°/ 12°/ 30°**	Biopsy forceps, double action
	AB-004-002	**0°/ 12°/ 30°**	Biopsy forceps, through cutting
	AB-004-004	**0°/ 12°**	Hook scissors, double action
	AB-004-010	**0°/ 12°**	Forceps, double action, for removal of stents
	AB-004-003	**0°/ 12°**	Grasping forceps, double action
	AB-004-008		Adapter for use with resectoscope sheaths

4.0 mm

1.1 | CYSTO-URETHROSCOPY
OPTICAL FORCEPS




Info: Only for sheaths 25 Fr.

	Ordering No.	For scopes	Description
	AB-004-007	**0°/ 12°/ 30°**	Optical stone punch
	AB-004-009		Adapter for use with resectoscope sheaths 26–28.5 Fr.
	AB-004-009-QL		Adapter with **STEMA click-system** for use with resectoscope sheaths 26–28.5 Fr.

4.0 mm


1.1 | CYSTO-URETHROSCOPY
FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps oval, double action	**Ø 1.0 mm**	3 Fr.	AC-001-001	AC-001-201	AC-001-001-R
		Ø 1.2 mm	4 Fr.	AC-001-008	AC-001-208	AC-001-008-R
		Ø 1.6 mm	5 Fr.	AC-001-002	AC-001-202	AC-001-002-R
		Ø 2.1 mm	6 Fr.	AC-001-003	AC-001-203	AC-001-003-R
		Ø 2.3 mm	7 Fr.	AC-001-004	AC-001-204	AC-001-004-R
		Ø 3.0 mm	9 Fr.	AC-001-005	AC-001-205	AC-001-005-R
	Biopsy forceps oval, serrated, double action	**Ø 1.0 mm**	3 Fr.	AC-002-001	AC-002-201	AC-002-001-R
		Ø 1.6 mm	5 Fr.	AC-002-002	AC-002-202	AC-002-002-R
		Ø 2.1 mm	6 Fr.	AC-002-003	AC-002-203	AC-002-003-R
		Ø 2.3 mm	7 Fr.	AC-002-004	AC-002-204	AC-002-004-R
		Ø 3.0 mm	9 Fr.	AC-002-005	AC-002-205	AC-002-005-R
	Biopsy forceps oval, serrated, teeth interlock, double action	**Ø 1.6 mm**	5 Fr.	AC-015-002	AC-015-202	AC-015-002-R
		Ø 2.1 mm	6 Fr.	-	AC-015-203	AC-015-003-R
		Ø 2.3 mm	7 Fr.	AC-015-004	AC-015-204	AC-015-004-R
	Grasping forceps, stone street double action	**Ø 1.2 mm**	4 Fr.	AC-017-001	AC-017-201	AC-017-001-R
		Ø 2.3 mm	7 Fr.	AC-017-004	AC-017-204	AC-017-004-R
	Grasping forceps, alligator jaw, double action	**Ø 1.0 mm**	3 Fr.	AC-006-001	AC-006-201	AC-006-001-R
		Ø 1.2 mm	4 Fr.	AC-006-003	AC-006-203	AC-006-003-R
		Ø 1.6 mm	5 Fr.	AC-006-002	AC-006-202	AC-006-002-R
		Ø 2.3 mm	7 Fr.	AC-006-004	AC-006-204	AC-006-004-R
		Ø 3.0 mm	9 Fr.	AC-006-005	AC-006-205	AC-006-005-R

4.0 mm


1.1 | CYSTO-URETHROSCOPY
FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	**Ø 2.3 mm**	7 Fr.	AC-008-004	AC-008-204	AC-008-004-R
		Ø 3.0 mm	9 Fr.	AC-008-005	AC-008-205	AC-008-005-R
	Grasping forceps teeth double action	**Ø 1.2 mm**	4 Fr.	AC-016-001	AC-016-201	AC-016-001-R
		Ø 1.6 mm	5 Fr.	AC-016-002	AC-016-202	AC-016-002-R
	Grasping forceps teeth interlock, double action	**Ø 1.6 mm**	5 Fr.	AC-020-002	AC-020-202	AC-020-002-R
		Ø 2.3 mm	7 Fr.	AC-020-004	AC-020-204	AC-020-004-R
	Grasping forceps rat tooth, double action	**Ø 1.0 mm**	3 Fr.	AC-019-001	AC-019-201	AC-019-001-R
		Ø 1.6 mm	5 Fr.	AC-019-002	AC-019-202	AC-019-002-R
	Grasping forceps for stone street, double action	**Ø 1.2 mm**	4 Fr.	AC-007-001	AC-007-201	AC-007-001-R
		Ø 1.6 mm	5 Fr.	AC-007-002	AC-007-202	AC-007-002-R
	Grasping forceps for big stones, double action	**Ø 1.2 mm**	4 Fr.	AC-009-001	AC-009-201	AC-009-001-R
		Ø 1.6 mm	5 Fr.	AC-009-002	AC-009-202	AC-009-002-R
		Ø 2.3 mm	7 Fr.	AC-009-004	AC-009-204	AC-009-004-R
	3-Prongs retriever	**Ø 1.0 mm**	3 Fr.	AC-012-001		
		Ø 1.6 mm	5 Fr.	AC-012-002		
		Ø 2.3 mm	7 Fr.	AC-012-004		


4.0 mm

1.1 | CYSTO-URETHROSCOPY
SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Hook scissors, single action	**Ø 2.3 mm**	7 Fr.	AC-004-004	AC-004-204	AC-004-004-R
		Ø 3.0 mm	9 Fr.	AC-004-005	AC-004-205	AC-004-005-R
	Scissors, single action	**Ø 1.0 mm**	3 Fr.	AC-003-000	AC-003-200	AC-003-000-R
		Ø 1.2 mm	4 Fr.	AC-003-001	AC-003-201	AC-003-001-R
		Ø 1.6 mm	5 Fr.	AC-003-002	AC-003-202	AC-003-002-R
		Ø 2.1 mm	6 Fr.	AC-003-003	AC-003-203	AC-003-003-R
		Ø 2.3 mm	7 Fr.	AC-003-004	AC-003-204	AC-003-004-R
		Ø 3.0 mm	9 Fr.	AC-003-005	AC-003-205	AC-003-005-R
	Uretero-slitting scissors, retrograde cutting single action	**Ø 1.2 mm**	4 Fr.	AC-018-001	AC-018-201	AC-018-001-R
		Ø 1.6 mm	5 Fr.	AC-018-002	AC-018-202	AC-018-002-R
	Scissors pointed, single action	**Ø 1.6 mm**	5 Fr.	-	AC-003-302	-
		Ø 2.3 mm	7 Fr.	-	AC-003-304	-
	Scissors, serrated flexible, double action	**Ø 2.3 mm**	7 Fr.	AC-005-004	AC-005-204	AC-005-004-R

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




FEATURES

1. QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
2. ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
3. OVERLOAD PROTECTION AVAILABLE
4. SIDE-MOUNTED LUER LOCK


STEMA

4.0 mm

1.1 | CYSTO-URETHROSCOPY
INSTRUMENTS DETACHABLE | SEMI-RIGID




Ordering No.	Description
AC-100-001	Handle for detachable instruments


	Description	Diameter	French	Complete	Insert
	Scissors, single action	**Ø 1.0 mm**	3 Fr.	AC-003-200-T	AC-003-200-OT
		Ø 1.2 mm	4 Fr.	AC-003-201-T	AC-003-201-OT
		Ø 1.6 mm	5 Fr.	AC-003-202-T	AC-003-202-OT
		Ø 2.3 mm	7 Fr.	AC-003-204-T	AC-003-204-OT
		Ø 3.0 mm	9 Fr.	AC-003-205-T	AC-003-205-OT
	Scissors pointed, single action	**Ø 1.6 mm**	5 Fr.	AC-003-612-T	AC-003-612-OT
		Ø 2.3 mm	7 Fr.	AC-003-614-T	AC-003-614-OT
	Hook scissors, single action	**Ø 2.3 mm**	7 Fr.	AC-004-204-T	AC-004-204-OT
		Ø 3.0 mm	9 Fr.	AC-004-205-T	AC-004-205-OT
	Biopsy forceps oval, double action	**Ø 1.0 mm**	3 Fr.	AC-001-201-T	AC-001-201-OT
		Ø 1.2 mm	4 Fr.	AC-001-206-T	AC-001-206-OT
		Ø 1.6 mm	5 Fr.	AC-001-202-T	AC-001-202-OT
		Ø 2.3 mm	7 Fr.	AC-001-203-T	AC-001-203-OT
		Ø 3.0 mm	9 Fr.	AC-001-204-T	AC-001-204-OT
	Grasping forceps, alligator jaw, double action	**Ø 1.0 mm**	3 Fr.	AC-006-201-T	AC-006-201-OT
		Ø 1.2 mm	4 Fr.	AC-006-203-T	AC-006-203-OT
		Ø 1.6 mm	5 Fr.	AC-006-202-T	AC-006-202-OT
		Ø 2.3 mm	7 Fr.	AC-006-204-T	AC-006-204-OT
		Ø 3.0 mm	9 Fr.	AC-006-205-T	AC-006-205-OT
	Grasping forceps for stone street, double action	**Ø 1.2 mm**	4 Fr.	AC-007-201-T	AC-007-201-OT
		Ø 1.6 mm	5 Fr.	AC-007-202-T	AC-007-202-OT
	Grasping forceps for big stones, double action	**Ø 1.2 mm**	4 Fr.	AC-009-201-T	AC-009-201-OT
		Ø 1.6 mm	5 Fr.	AC-009-202-T	AC-009-202-OT

4.0 mm

1.1 | CYSTO-URETHROSCOPY
INSTRUMENTS DETACHABLE | RIGID



Ordering No.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Flexible	Semi-rigid
	Scissors, single action	**Ø 1.0 mm**	3 Fr.	AC-103-200-T	AC-103-200-OT
		Ø 1.2 mm	4 Fr.	AC-103-201-T	AC-103-201-OT
		Ø 1.6 mm	5 Fr.	AC-103-202-T	AC-103-202-OT
		Ø 2.3 mm	7 Fr.	AC-103-204-T	AC-103-204-OT
		Ø 3.0 mm	9 Fr.	AC-103-205-T	AC-103-205-OT
	Scissors pointed, single action	**Ø 1.6 mm**	5 Fr.	AC-103-612-T	AC-103-612-OT
		Ø 2.3 mm	7 Fr.	AC-103-614-T	AC-103-614-OT
	Hook scissors, single action	**Ø 2.3 mm**	7 Fr.	AC-104-204-T	AC-104-204-OT
		Ø 3.0 mm	9 Fr.	AC-104-205-T	AC-104-205-OT
	Biopsy forceps oval, double action	**Ø 1.0 mm**	3 Fr.	AC-101-201-T	AC-101-201-OT
		Ø 1.2 mm	4 Fr.	AC-101-206-T	AC-101-206-OT
		Ø 1.6 mm	5 Fr.	AC-101-202-T	AC-101-202-OT
		Ø 2.3 mm	7 Fr.	AC-101-203-T	AC-101-203-OT
		Ø 3.0 mm	9 Fr.	AC-101-204-T	AC-101-204-OT
	Grasping forceps, alligator jaw, double action	**Ø 1.0 mm**	3 Fr.	AC-106-201-T	AC-106-201-OT
		Ø 1.2 mm	4 Fr.	AC-106-203-T	AC-106-203-OT
		Ø 1.6 mm	5 Fr.	AC-106-202-T	AC-106-202-OT
		Ø 2.3 mm	7 Fr.	AC-106-204-T	AC-106-204-OT
		Ø 3.0 mm	9 Fr.	AC-106-205-T	AC-106-205-OT
	Grasping forceps for stone street, double action	**Ø 1.2 mm**	4 Fr.	AC-107-201-T	AC-107-201-OT
		Ø 1.6 mm	5 Fr.	AC-107-202-T	AC-107-202-OT
	Grasping forceps for big stones, double action	**Ø 1.2 mm**	5 Fr.	AC-109-201-T	AC-109-201-OT

4.0 mm

1.1 | CYSTO-URETHROSCOPY
STONE EXTRACTORS | 4 WIRES



Ordering No.	Description
AF-001-021	Handle for dormia stone extractor **3 and 5 Fr.**
AF-001-024	Handle for dormia stone extractor **7 Fr.**

Ordering No.	Wires	French	Working length	Description
AF-001-001	4	**3 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-002	4	**3 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-015	4	**3 Fr.**	90 cm	Dormia stone extractor, basket with **tubus only**
AF-001-016	4	**3 Fr.**	90 cm	Dormia stone extractor complete with **handle**
AF-001-023	4	**3 Fr.**	120 cm	Dormia stone extractor, basket with **tubus only**
AF-001-022	4	**3 Fr.**	120 cm	Dormia stone extractor complete with **handle**
AF-001-003	4	**5 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-004	4	**5 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-017	4	**5 Fr.**	90 cm	Dormia stone extractor, basket with **tubus only**
AF-001-019	4	**5 Fr.**	90 cm	Dormia stone extractor complete with **handle**
AF-001-005	4	**7 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-006	4	**7 Fr.**	60 cm	Dormia stone extractor complete with **handle**

4.0 mm

1.1 | CYSTO-URETHROSCOPY
STONE EXTRACTORS | 5-7 WIRES



Ordering No.	Description
AF-001-021	Handle for dormia stone extractor **3 and 5 Fr.**
AF-001-024	Handle for dormia stone extractor **7 Fr.**

Ordering No.	Wires	French	Working length	Description
AF-001-007	5	**5 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-008	5	**5 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-025	5	**5 Fr.**	90 cm	Dormia stone extractor, basket with **tubus only**
AF-001-018	5	**5 Fr.**	90 cm	Dormia stone extractor complete with **handle**
AF-001-009	5	**7 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-010	5	**7 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-011	6	**5 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-012	6	**5 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-013	7	**7 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-014	7	**7 Fr.**	60 cm	Dormia stone extractor complete with **handle**

4.0 mm

1.1 | CYSTO-URETHROSCOPY
MONOPOLAR FLEXIBLE ELECTRODES



	Description	French	Length	Ordering No.	French	Length	Ordering No.
	Loop electrode flexible, monopolar	**5 Fr.**	53 cm	AD-003-003	**5 Fr.**	65 cm	AD-003-013
		7 Fr.	53 cm	AD-003-005	**5 Fr.**	65 cm	AD-003-015
		9 Fr.	53 cm	AD-003-007	**5 Fr.**	65 cm	AD-003-017
		10 Fr.	53 cm	AD-003-008	**5 Fr.**	65 cm	AD-003-018
	Needle electrode flexible, monopolar	**3 Fr.**	53 cm	AD-002-001	**3 Fr.**	65 cm	AD-002-011
		4 Fr.	53 cm	AD-002-002	**4 Fr.**	65 cm	AD-002-012
		5 Fr.	53 cm	AD-002-003	**5 Fr.**	65 cm	AD-002-013
		7 Fr.	53 cm	AD-002-005	**7 Fr.**	65 cm	AD-002-015
	Button electrode flexible, monopolar	**3 Fr.**	53 cm	AD-001-001	**3 Fr.**	65 cm	AD-001-011
		4 Fr.	53 cm	AD-001-002	**4 Fr.**	65 cm	AD-001-012
		5 Fr.	53 cm	AD-001-003	**5 Fr.**	65 cm	AD-001-013
		6 Fr.	53 cm	AD-001-004	**6 Fr.**	65 cm	AD-001-014
		7 Fr.	53 cm	AD-001-005	**7 Fr.**	65 cm	AD-001-015
		8 Fr.	53 cm	AD-001-006	**8 Fr.**	65 cm	AD-001-016
		9 Fr.	53 cm	AD-001-007	**9 Fr.**	65 cm	AD-001-017
		10 Fr.	53 cm	AD-001-008	**10 Fr.**	65 cm	AD-001-018

4.0 mm

1.1 | CYSTO-URETHROSCOPY
MONOPOLAR CONNECTING CABLES



Connection instrument side





2 mm

Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001 AD-010-006	3 m 5 m	Ø 4.0 mm	HF Monopolar connecting cable for **Berchtold, Martin, Aesculap**
	AD-010-002 AD-010-007	3 m 5 m	Ø 5.0 mm	HF Monopolar connecting cable for **Erbe, ACC/ICC, Storz**
	AD-010-003 AD-010-008	3 m 5 m	Ø 8.0 mm	HF Monopolar connecting cable for **STEMA HF-U Series, Valleylab/Conmed**
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable **male** for **STEMA MASTER**






Ordering No.	Description
EK-015-060	**Universal adapter** used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

2.9 mm

1.2 | SLIM CYSTO-URETHROSCOPY
CYSTOSCOPES | SHEATHS



	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-002
	HD Cystoscope autoclavable			AA-003-012
	HD Cystoscope autoclavable wide angle			AA-003-012-W



Ordering No.	Description
AB-001-007	Cysto-Urethroscope sheath **15 Fr. with 2 fixed stopcocks, 5 Fr. instrument channel,** incl. obturator




Ordering No.	Description
AB-001-000	Cysto-Urethroscope sheath **16 Fr. with 2 ixed stopcocks, Quick-Lock,** incl. obturator, 1x5 fr. or 2x3 fr. instruments

STEMA

2.9 mm

1.2 | SLIM CYSTO-URETHROSCOPY
ENDOSCOPE BRIDGES | QUICK-LOCK SYSTEM


Quick-Lock System



Ordering No.	Description
AB-003-010	Endoscope bridge **diagnostic used for 16 Fr. sheath**



Ordering No.	Description
AB-003-011	Endoscope bridge **with 1 instrument channel used for 16 Fr. sheath**



Ordering No.	Description
AB-003-004	Endoscope bridge **with 2 instrument channels**



Ordering No.	Description
AB-003-012	Endoscope bridge **with 2 instrument channels used for 16 Fr. sheath**



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6–0.8 mm

2.9 mm

1.2 | SLIM CYSTO-URETHROSCOPY
WORKING ELEMENTS



Ordering No.	Description
AB-003-013	Working element with ALBARRAN deflector and **1 instrument channel** used for 16 Fr. sheath
AB-003-014	Working element with ALBARRAN deflector and **2 instrument channels** used for 16 Fr. sheath



Ordering No.	Description
AB-003-020	Protective rod for working elements with ALBARRAN deflector

1.3 | LASER CYSTO-URETHROSCOPY
CYSTOSCOPES | SHEATHS



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-003
	HD Cystoscope autoclavable			AA-002-013
	ULTRA HD (4K) Cystoscope autoclavable			AA-002-023



Ordering No.	Description
AB-001-040	**LASER** Cysto-Urethroscope outer sheath 23 Fr. with 2 fixed stopcock used for 30° endoscope



Ordering No.	Description
AB-001-041	**LASER** Cysto-Urethroscope inner sheath with 8 Fr. working channel

1.3 | LASER CYSTO-URETHROSCOPY
OBTURATORS



Ordering No.	Description
AB-001-043	Obturator for LASER cysto-urethroscope sheath **AB-001-040**



Ordering No.	Description
AB-001-042	Viewing obturator for LASER cysto-urethroscope sheath **AB-001-040**

1.4 | HD MINI COMPACT CYSTOSCOPE



Ordering No.	Description
AA-007-060	**HD Mini** Compact Cystoscope

Technical data

Eyepiece	30° angled ocular
Sheath	7.5/9/12 Fr.
Direction of view	0° viewing angle
Working length	220 mm
Instrument channel	6.8 Fr. | use with 1 × 5 or 2 × 3 Fr. instruments
Luer Lock	3 connections
Optic	40,000 px

Delivery contains
HD Mini compact cystoscope

1.5 | FLEXIBLE CYSTOSCOPE
FLEX I

- **High luminance**
- **Exact adjustment of focus**
- **Precise up-/down angulation**



Ordering No.	Description
AA-008-002	**Flexible** Cystoscope FLEX I

Technical data

Cystoscope

Insertion tube	Ø 5.5 mm
Working Channel	Ø 2.3 mm
Working length	380 mm
Field of view	100°
Pixel	18,000
Angulation up-down	210/120°
Accessories	Suction valve, leak tester

Delivery contains
Flexible cystoscope FLEX I, leak tester, suction valve

INCLUDED IN SCOPE OF DELIVERY:

AA-008-103
Leak tester

AA-008-105
Suction valve



AA-008-103



AA-008-105

1.5 | FLEXIBLE CYSTOSCOPE
ACCESSORIES



Ordering No.	Description
AA-008-100	Bridge with Luer-Lock and 2 stopcocks for **AA-008-002**



Ordering No.	Description
AA-008-102	Suitcase grey for STEMA flexible scopes



Ordering No.	Description
AA-008-103	Leak tester



Ordering No.	Description
AA-008-104	Single use Luer-Lock adapter, 1 box = 25 pieces



Ordering No.	Description
AA-008-105	Suction valve

1.6 CYSTO-URETHROSCOPY SET

1 **HD CYSTOSCOPE**
Ø 4.0 mm, 30° viewing angle, autoclavable, 302 mm working length

1 AA-002-013



2 **CYSTO-URETHROSCOPE SHEATH**
21 Fr. with 2 fixed stopcock and obturator

3 **CYSTO-URETHROSCOPE SHEATH**
25 Fr. with 2 fixed stopcock and obturator

2 AB-001-004

3 AB-001-006



4 **ENDOSCOPE BRIDGE**
with 1 instrument channel

4 AB-003-003



5 **WORKING ELEMENT**
with Albaran deflector and 2 instrument channels

5 AB-003-007





6 **BIOPSY FORCEPS**
flexible, oval, double action, Ø 1.6 mm, WL 40 cm

6 AC-001-002

7 **SCISSORS**
flexible, single action, Ø 1.6 mm, WL 40 cm

7 AC-003-002



8 **GRASPING FORCEPS ALLIGATOR**
flexible, double action, Ø 2.3 mm, WL 40 cm

8 AC-006-004



1.6 CYSTO-URETHROSCOPY SET

9 BUTTON ELECTRODE
flexible, 5 Fr., monopolar, length 530 mm

10 NEEDLE ELECTRODE
flexible, 5 Fr., monopolar, length 530 mm



9 AD-001-003



10 AD-002-003



11 HF MONOPOLAR CABLE
HF monopolar connecting cable for Berchtold, Martin, Aesculap, Ø 4 mm

11 AD-010-003



12 DORMIA STONE EXTRACTOR
5 wires, Ø 1.6 mm, complete with handle, 60 cm Working length

12 AF-001-008



13 BIOPSY FORCEPS
double action

14 GRASPING FORCEPS
double action



13 AB-004-001

14 AB-004-003



15 OPTICAL STONE PUNCH

15 AB-004-007



2. RESECTOSCOPY

2.1 | MONOPOLAR RESECTOSCOPY
24/26/28.5 FRENCH
19/22 FRENCH (SLIM)

2.2 | SALINE RESECTOSCOPY
24/26 FRENCH
19/22 FRENCH (SLIM)
17.5 FRENCH (MINI)

2.3 | LASER RESECTOSCOPY
24/26 FRENCH
19/22 FRENCH (SLIM)
17.5 FRENCH (MINI)

2.4 | ACCESSORIES

2.5 | RESECTOSCOPY SETS

2. RESECTOSCOPY

24/26 Fr.
28.5 Fr.

19/22 Fr.


17.5 Fr.



1 CYSTOSCOPES



2 WORKING ELEMENTS



3 SHEATHS &
OBTURATORS



4 ELECTRODES





STEMA

4.0 mm


2.1 | MONOPOLAR RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 0°	Cystoscope autoclavable Cystoscope autoclavable **wide angle**	Ø 4.0 mm	302 mm	AA-002-001 AA-002-001-W
	HD Cystoscope autoclavable **HD** Cystoscope autoclavable **wide angle**			AA-002-011 AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-021

Viewing angle	Description	Diameter	Working length	Ordering No.
 12°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-002
	HD Cystoscope autoclavable			AA-002-012
	Ultra HD (4K) Cystoscope autoclavable			AA-002-022

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable Cystoscope autoclavable **wide angle**	Ø 4.0 mm	302 mm	AA-002-003 AA-002-003-W
	HD Cystoscope autoclavable **HD** Cystoscope autoclavable **wide angle**			AA-002-013 AA-002-013-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-023

Viewing angle	Description	Diameter	Working length	Ordering No.
 70°	Cystoscope autoclavable Cystoscope autoclavable wide angle	Ø 4.0 mm	302 mm	AA-002-004 AA-002-004-W
	HD Cystoscope autoclavable **HD** Cystoscope autoclavable wide angle			AA-002-014 AA-002-014-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-024

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
WORKING ELEMENTS



Ordering No.	Description
AG-001-001	Working element used for resectoscopes, **passive cutting**, handle closed
AG-001-001-HF	Working element used for resectoscopes, **passive cutting**, handle closed**, HF connection on top**




Ordering No.	Description
AG-001-002	Working element used for resectoscopes, **passive cutting**, handle open

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
WORKING ELEMENTS



Ordering No.	Description
AG-001-003	Working element used for resectoscopes, **active cutting**, handle closed
AG-001-003-HF	Working element used for resectoscopes, **active cutting**, handle closed, **HF connection on top**




Ordering No.	Description
AG-001-004	Working element used for resectoscopes, **active cutting**, handle open, cutting spring assisted

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
WORKING ELEMENTS | TITAN



Ordering No.	Description
AG-001-007	Working element used for resectoscopes, **passive cutting,** handle closed, titanium version



Ordering No.	Description
AG-001-008	Working element used for resectoscopes, **active cutting,** handle closed, titanium version

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
WORKING ELEMENTS | GYRUS



Ordering No.	Description
AG-001-017	**GYRUS** Saline working element used for resectoscopes, **active cutting**, handle closed, titanium version



Ordering No.	Description
AG-001-018	**GYRUS** Saline working element used for resectoscopes, **passive cutting**, handle closed, titanium version

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-001-202	Outer sheath with **perforation**	**26 Fr.**	Yellow
	AH-001-212	Outer sheath with **perforation**	**28.5 Fr.**	Black
	AH-001-204	Outer sheath with **upper/lower perforation**	**26 Fr.**	Yellow
	AH-001-203	Outer sheath with **surrounding fenestrated perforation**	**26 Fr.**	Yellow
	AH-001-205	Outer sheath with **upper fenestrated perforation**	**26 Fr.**	Yellow




Ordering No.	Description	French	Colour code
AH-001-201	Inner sheath with **ceramic isolation**	**26 Fr.**	Yellow
AH-001-211	Inner sheath with **ceramic isolation**	**28.5 Fr.**	Black

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, **rotatable**	**26 Fr.**	Yellow
	AH-001-312	Outer sheath with perforation, **rotatable**	**28.5 Fr.**	Black
	AH-001-304	Outer sheath with upper/lower perforation, **rotatable**	**26 Fr.**	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, **rotatable**	**26 Fr.**	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, **rotatable**	**26 Fr.**	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, **rotatable**	**26 Fr.**	Yellow



Ordering No.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, **rotatable**	**26 Fr.**	Yellow
AH-001-311	Inner sheath with ceramic isolation, **rotatable**	**28.5 Fr.**	Black



Ordering No.	Description	French	Colour code
AH-001-308	Inner sheath without ceramic isolation, **with laser channel up** to 7 Fr., for use with AH-001-301	**26 Fr.**	Yellow

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
SHEATHS



Ordering No.	Description	French	Colour code
AH-001-001	Standard resectoscope sheath with **1 fixed stopcock**	**24 Fr.**	Yellow
AH-001-002	Standard resectoscope sheath with **1 fixed stopcock**	**27 Fr.**	Black



Ordering No.	Description	French	Colour code
AH-001-101	Standard resectoscope sheath with **central valve**	**24 Fr.**	Yellow
AH-001-102	Standard resectoscope sheath with **central valve**	**27 Fr.**	Black

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
OBTURATORS



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-011	**Standard** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow
AH-001-012	**Standard** obturator for resectoscope sheaths	**27 Fr.**	**28.5 Fr.**	Black



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-015	**Deflecting** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow
AH-001-016	**Deflecting** obturator for resectoscope sheaths	**27 Fr.**	**28.5 Fr.**	Black



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-013	**Viewing** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-017	**Viewing** obturator for resectoscope sheaths with 5 Fr. **Instrument channel**	**24 Fr.**	**26 Fr.**	Yellow
AH-001-018	**Viewing** obturator for resectoscope sheaths with 5 Fr. **Instrument channel**	**27 Fr.**	**28.5 Fr.**	Black

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-001-001	Cutting loop electrode **30° angled**	**24 Fr.**	yellow
	AE-002-001		**27 Fr.**	brown
	AE-001-002	Cutting loop electrode **angled**	**24 Fr.**	yellow
	AE-002-002		**27 Fr.**	brown
	AE-001-003	Cutting loop electrode **straight**	**24 Fr.**	yellow
	AE-002-003		**27 Fr.**	brown
	AE-001-004	Rollerball electrode **Ø 3 mm**	**24 Fr.**	yellow
	AE-002-004		**27 Fr.**	brown
	AE-001-005	Rollerball electrode **Ø 5 mm**	**24 Fr.**	yellow
	AE-002-005		**27 Fr.**	brown
	AE-001-010	Roller electrode **Ø 3 mm**	**24 Fr.**	yellow
	AE-002-010		**27 Fr.**	brown
	AE-001-011	Roller electrode **Ø 5 mm**	**24 Fr.**	yellow
	AE-002-011		**27 Fr.**	brown
	AE-001-012	Spike electrode **Ø 3 mm**	**24 Fr.**	yellow
	AE-002-012		**27 Fr.**	brown
	AE-001-013	Spike electrode **Ø 5 mm**	**24 Fr.**	yellow
	AE-002-013		**27 Fr.**	brown

24/26 Fr.
28.5 Fr.

2.1 | MONOPOLAR RESECTOSCOPY
ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-001-014	Vaporisation electrode **Ø 3 mm**	**24 Fr.**	yellow
	AE-002-014		**27 Fr.**	brown
	AE-001-015	Vaporisation electrode **Ø 5 mm**	**24 Fr.**	yellow
	AE-002-015		**27 Fr.**	brown
	AE-001-009	Conical electrode	**24 Fr.**	yellow
	AE-001-016	Knife electrode	**24 Fr.**	yellow
	AE-002-016		**27 Fr.**	brown
	AE-001-006	Band electrode **angled**	**24 Fr.**	yellow
	AE-002-006		**27 Fr.**	brown
	AE-001-007	Band electrode **30° angled**	**24 Fr.**	yellow
	AE-002-007		**27 Fr.**	brown
	AE-001-008	Band vaporization electrode **angled**	**24 Fr.**	yellow
	AE-001-017	Curette sharp, **round** (not for use with HF)	**24 Fr.**	yellow



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-002
	HD Cystoscope autoclavable			AA-003-012

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
WORKING ELEMENTS



Ordering No.	Description
AG-002-001	Working element used for slim resectoscopes, **passive cutting**, **handle closed**




Ordering No.	Description
AG-002-002	Working element used for slim resectoscopes, **passive cutting**, **handle open**

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
WORKING ELEMENTS




Ordering No.	Description
AG-002-003	Working element used for slim resectoscopes, **passive cutting**, **handle closed**



Ordering No.	Description
AG-002-004	Working element used for slim resectoscopes, **active cutting**, **handle open**

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
WORKING ELEMENTS | TITAN




Ordering No.	Description
AG-002-006	Working element used for slim resectoscopes, **active cutting**, handle closed, titanium version

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-212	Outer sheath with perforation, **fixed**	**22 Fr.**	White




Ordering No.	Description	French	Colour code
AH-002-201	Inner sheath with ceramic isolation, **fixed**	**19 Fr.**	White

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
SHEATHS





Ordering No.	Description	French	Colour code
AH-002-302	Outer sheath with perforation, **rotatable**	**22 Fr.**	White



Ordering No.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, **rotatable**	**19 Fr.**	White

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-001	Standard resectoscope sheath, **1 fixed stopcock**, without obturator	**19 Fr.**	White

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-003-001	Cutting loop electrode **30° angled**	**19 Fr.**	White
	AE-003-003	Cutting loop electrode **straight**	**19 Fr.**	White
	AE-003-004	Rollerball electrode **Ø 3 mm**	**19 Fr.**	White
	AE-003-010	Roller electrode **Ø 3 mm**	**19 Fr.**	White
	AE-003-014	Vaporisation electrode **Ø 3 mm**	**19 Fr.**	White
	AE-003-016	Knife electrode	**19 Fr.**	White



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

19/22 Fr.

2.1 | MONOPOLAR SLIM RESECTOSCOPY
OBTURATORS




Ordering No.	Description	French	Colour code
AH-002-011	Standard obturator	**19 Fr.**	White




Ordering No.	Description	French	Colour code
AH-002-013	Viewing obturator	**19 Fr.**	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY
MONOPOLAR CONNECTING CABLES



Connection instrument side





Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001 AD-010-006	3 m 5 m	Ø 4.0 mm	HF Monopolar connecting cable for **Berchtold, Martin, Aesculap**
	AD-010-002 AD-010-007	3 m 5 m	Ø 5.0 mm	HF Monopolar connecting cable for **Erbe, ACC/ICC, Storz**
	AD-010-003 AD-010-008	3 m 5 m	Ø 8.0 mm	HF Monopolar connecting cable for **STEMA HF-U Series, Valleylab/Conmed**
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable **male** for **STEMA MASTER**





Ordering No.	Description
EK-015-060	**Universal adapter** used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

2.1 | MONOPOLAR SLIM RESECTOSCOPY
MONOPOLAR CONNECTING CABLES



Connection unit side



Unit side	Ordering No.	Length	Description
	AD-010-025	4.5 m	HF Monopolar connecting cable for **Storz resectoscopes with 90° angled connector**
	AD-010-026	4.5 m	HF Monopolar connecting cable for **Olympus resectoscopes, diameter Ø 3 mm**



RECOMMENDED PRODUCTS:

CI-001-030 STEMA HF-U PRO TOUCH

- *For all monopolar & bipolar techniques*
- *with 7“ touch screen display*
- *provided to perform the most advanced electrosurgery*
- *ADC / APC System*

CI-001-040 DOUBLE PEDAL FOOTSWITCH

24/26 Fr.
28.5 Fr.

2.2 | SALINE RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 12°	**HD** Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-012

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-003
	HD Cystoscope autoclavable			AA-002-013
	ULTRA HD Cystoscope autoclavable			AA-002-023

24/26 Fr.
28.5 Fr.

2.2 | SALINE RESECTOSCOPY
WORKING ELEMENTS





Ordering No.	Description
AG-001-011	Saline working element used for resectoscopes, **24 Fr., passive cutting**, handle closed







Ordering No.	Description
AG-001-012	Saline working element used for resectoscopes, **24 Fr., active cutting**, handle closed

24/26 Fr.
28.5 Fr.

2.2 | SALINE RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, **rotatable**	**26 Fr.**	Yellow
	AH-001-312	Outer sheath with perforation, **rotatable**	**28.5 Fr.**	Black
	AH-001-304	Outer sheath with upper/lower perforation, **rotatable**	**26 Fr.**	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, **rotatable**	**26 Fr.**	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, **rotatable**	**26 Fr.**	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, **rotatable**	**26 Fr.**	Yellow



Ordering No.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, **rotatable**	**26 Fr.**	Yellow
AH-001-311	Inner sheath with ceramic isolation, **rotatable**	**28.5 Fr.**	Black

24/26 Fr.
28.5 Fr.

2.2 | SALINE RESECTOSCOPY
OBTURATORS



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-011	**Standard** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow
AH-001-012	**Standard** obturator for resectoscope sheaths	**27 Fr.**	**28.5 Fr.**	Black



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-015	**Deflecting** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow
AH-001-016	**Deflecting** obturator for resectoscope sheaths	**27 Fr.**	**28.5 Fr.**	Black



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-013	**Viewing** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow




Ordering No.	Description	French	For outer sheath	Colour code
AH-001-017	**Viewing** obturator for resectoscope sheaths with 5 Fr. **Instrument channel**	**24 Fr.**	**26 Fr.**	Yellow
AH-001-018	**Viewing** obturator for resectoscope sheaths with 5 Fr. **Instrument channel**	**27 Fr.**	**28.5 Fr.**	Black

24/26 Fr.
28.5 Fr.

2.2 | SALINE RESECTOSCOPY
ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-001-030	Cutting loop electrode **30° angled, strong cut**	**24 Fr.**	yellow/blue
	AE-001-031	Cutting loop electrode **30° angled**	**24 Fr.**	yellow/blue
	AE-001-033	Cutting loop electrode **90° angled**	**24 Fr.**	yellow/blue
	AE-001-048	Cutting loop electrode **90° angled, for 0° Scope**	**24 Fr.**	yellow/blue
	AE-001-036	Cutting loop electrode **straight**	**24 Fr.**	yellow/blue
	AE-001-032	Knife electrode	**24 Fr.**	yellow/blue
	AE-001-041	Band electrode, **90° angled**	**24 Fr.**	yellow/blue
	AE-001-038	Vaporisation spike rollerball electrode **Ø 3 mm**	**24 Fr.**	yellow/blue
	AE-001-037	Vaporisation rollerball electrode **Ø 3 mm**	**24 Fr.**	yellow/blue
	AE-001-043	Vaporisation rollerball electrode **Ø 5 mm**		
	AE-001-040	Roller electrode smooth **Ø 3 mm**	**24 Fr.**	yellow/blue
	AE-001-042	Roller electrode smooth **Ø 5 mm**		

24/26 Fr.
28.5 Fr.

2.2 | SALINE RESECTOSCOPY
ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-001-034	Rollerball electrode **Ø 3 mm**	**24 Fr.**	yellow/blue
	AE-001-035	Rollerball electrode **Ø 5 mm**		
	AE-001-044	Conical electrode „Mushroom“	**24 Fr.**	yellow/blue
	AE-001-045	Conical electrode „Mushroom“, **convex**	**24 Fr.**	yellow/blue
	AE-001-046	Conical electrode „Mushroom“, **concave**	**24 Fr.**	yellow/blue
	AE-001-047	Conical electrode „Mushroom“, **grooved**	**24 Fr.**	yellow/blue
	AE-001-039	Conical electrode	**24 Fr.**	yellow/blue

	Ordering No.	Description	French	Colour code
	AE-002-031	Cutting loop electrode, **30° angled**	**27 Fr.**	yellow/blue



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

19/22 Fr.

2.2 | SALINE SLIM RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 12°	**HD** Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-004
				AA-003-014

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-002
	HD Cystoscope autoclavable			AA-003-012

19/22 Fr.

2.2 | SALINE SLIM RESECTOSCOPY
WORKING ELEMENTS





Ordering No.	Description
AG-002-011	Saline working element used for resectoscopes, **19 Fr., passive cutting,** handle closed






Ordering No.	Description
AG-002-012	Saline working element used for resectoscopes, **19 Fr., active cutting,** handle closed

19/22 Fr.

2.2 | SALINE SLIM RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-302	Outer sheath with perforation, **rotatable**	**22 Fr.**	White



Ordering No.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, **rotatable**	**19 Fr.**	White

19/22 Fr.

2.2 | SALINE SLIM RESECTOSCOPY
OBTURATORS



Ordering No.	Description
AH-002-011	Standard obturator for resectoscope sheaths **19 Fr.**



Ordering No.	Description
AH-002-013	Viewing obturator for resectoscope sheaths **19 Fr.**

19/22 Fr.

2.2 | SALINE SLIM RESECTOSCOPY
ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-003-031	Cutting loop electrode **30° angled**	**19 Fr.**	white/blue
	AE-003-032	Cutting loop electrode **90° angled**	**19 Fr.**	white/blue
	AE-003-033	Cutting loop electrode **straight**	**19 Fr.**	white/blue
	AE-003-034	Rollerball electrode **Ø 3 mm**	**19 Fr.**	white/blue
	AE-003-040	Roller electrode **Ø 3 mm**	**19 Fr.**	white/blue
	AE-003-036	Knife electrode	**19 Fr.**	white/blue



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

17.5 Fr.

2.2 | SALINE MINI RESECTOSCOPY
CYSTOSCOPES




Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-002
	HD Cystoscope autoclavable			AA-003-012

17.5 Fr.

2.2 | SALINE MINI RESECTOSCOPY
WORKING ELEMENTS



Ordering No.	Description
AG-002-015	Saline working element used for resectoscopes, **18.5 Fr., passive cutting**, handle closed




Ordering No.	Description
AG-002-013	Saline working element used for resectoscopes, **18.5 Fr., active cutting**, handle closed

17.5 Fr.

2.2 | SALINE MINI RESECTOSCOPY
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-002-213	Outer sheath 18.5 Fr. **with perforation**	**18.5 Fr.**	green



Ordering No.	Description
AH-002-211	Inner sheath 17.5 Fr. **with ceramic isolation, rotatable**



Ordering No.	Description
AH-002-216	Inner sheath 17.5 Fr. **without ceramic isolation**, with laser channel up to 5 Fr., for use with AH-001-213

17.5 Fr.

2.2 | SALINE MINI RESECTOSCOPY
OBTURATORS



Ordering No.	Description
AH-002-015	Standard obturator **17.5 Fr.**



Ordering No.	Description
AH-002-016	Viewing obturator **17.5 Fr. with 5 Fr. instrument channel**

17.5 Fr.

2.2 | SALINE MINI RESECTOSCOPY
ELECTRODES




	Ordering No.	Description	French	Colour code
	AE-003-041	Cutting loop electrode **30° angled**	**17.5 Fr.**	green
	AE-003-043	Cutting loop electrode **straight**	**17.5 Fr.**	green
	AE-003-044	Rollerball electrode **Ø 2 mm**	**17.5 Fr.**	green
	AE-003-045	Roller electrode **Ø 2 mm**	**17.5 Fr.**	green
	AE-003-046	Knife electrode	**17.5 Fr.**	green



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.2 | SALINE RESECTOSCOPY
BIPOLAR CONNECTING CABLES



Connection instrument side



Unit side	Ordering No.	Length	Description
	AD-010-011	4.5 m	Bipolar resectoscope connecting cable for **Martin Maxium**
	AD-010-012	4.5 m	Bipolar resectoscope connecting cable for **Erbe Vio**
	AD-010-013	4.5 m	Bipolar resectoscope connecting cable for **Valleylab, Lamidey, EMC**
	AD-010-015	4.5 m	Bipolar resectoscope connecting cable for **Erbe Vio 3**

24/26 Fr.

2.3 | LASER RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	**HD** Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-013

24/26 Fr.

2.3 | LASER RESECTOSCOPY
WORKING ELEMENTS




Ordering No.	Description
AG-001-006	Working element used for resectoscopes, titanium, passive cutting, handle closed, **24/26 Fr.** with probe channel **0.8 mm**
AG-001-009	Working element used for resectoscopes, titanium, passive cutting, handle closed, **24/26 Fr.** with probe channel **1.2 mm**

Ordering No.	Description
AG-100-006	Adapter used for laser probes
AG-100-007	TUOHY-Borst adapter

24/26 Fr.

2.3 | LASER RESECTOSCOPY
SHEATHS | OBTURATORS





Ordering No.	Description	French	Colour code
AH-001-302	Outer sheath with perforation, **rotatable**	**26 Fr.**	Yellow



Ordering No.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, **rotatable**	**26 Fr.**	Yellow



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-308	Inner sheath without ceramic isolation, **with laser channel up** to 7 Fr., for use with AH-001-301	**24 Fr.**	**26 Fr.**	Yellow



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-011	**Standard** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow

19/22 Fr.

2.3 | LASER SLIM RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	**HD** Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-012

19/22 Fr.

2.3 | LASER SLIM RESECTOSCOPY
WORKING ELEMENTS




Ordering No.	Description
AG-002-007	Working element used for resectoscopes, titanium, passive cutting, handle closed, **19/22 Fr.** with probe channel **0.8 mm**
AG-002-008	Working element used for resectoscopes, titanium, passive cutting, handle closed, **19/22 Fr.** with probe channel **1.2 mm**

19/22 Fr.

2.3 | LASER SLIM RESECTOSCOPY
SHEATHS





Ordering No.	Description	French	Colour code
AH-002-302	Outer sheath with perforation, **rotatable**	**22 Fr.**	White



Ordering No.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, **rotatable**	**19 Fr.**	White

19/22 Fr.

2.3 | LASER SLIM RESECTOSCOPY
OBTURATORS



Ordering No.	Description
AH-002-011	Standard obturator for resectoscope sheaths **19 Fr.**



Ordering No.	Description
AH-002-013	Viewing obturator for resectoscope sheaths **19 Fr.**

17.5 Fr.

2.3 | LASER MINI RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	**HD** Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-012

17.5 Fr.

2.3 | LASER MINI RESECTOSCOPY
WORKING ELEMENTS




Ordering No.	Description
AG-002-016	Working element used for resectoscopes, titanium, passive cutting, handle closed, **17.5 Fr.** with probe channel **0.8 mm**

17.5 Fr.

2.3 | LASER MINI RESECTOSCOPY
SHEATHS





Ordering No.	Description	French	Colour code
AH-002-213	Outer sheath 18.5 Fr. **with perforation**	**18.5 Fr.**	green



Ordering No.	Description
AH-002-211	Inner sheath 17.5 Fr. **with ceramic isolation, rotatable**



Ordering No.	Description
AH-002-216	Inner sheath 17.5 Fr. **without ceramic isolation**, with laser channel up to 5 Fr., for use with AH-001-213

17.5 Fr.

2.3 | LASER MINI RESECTOSCOPY
OBTURATORS



Ordering No.	Description
AH-002-015	Standard obturator **17.5 Fr.**



Ordering No.	Description
AH-002-016	Viewing obturator **17.5 Fr. with 5 Fr. instrument channel**

2.4 | ACCESSORIES
SUPRAPUBIC ASPIRAT



Ordering No.	Description
AB-007-001	Suprapubic aspiration set, complete set



Ordering No.	Description
AB-007-002	Trocar sleeve **16 Fr.** with pressure measuring channel



Ordering No.	Description
AB-007-003	Adjustable stabilisation disc, **lateral opening**



Ordering No.	Description
AB-007-008	Guide sleeve, **lateral opening**, for inserting a suprapubic catheter

2.4 | ACCESSORIES
SUPRAPUBIC ASPIRATION



Ordering No.	Description
AB-007-004	Trocar with **cannula tip**



Ordering No.	Description
AB-007-005	Trocar with **pyramidal tip**



Ordering No.	Description
AB-007-006	Aspiration cannula




Ordering No.	Description
AB-007-007	Overflow cannula for pressure stabilisation

STEMA

2.4 | ACCESSORIES
EVACUATOR | BLADDER SYRINGES




Ordering No.	Description
AJ-002-001	Ellik evacuator complete with sheath connection and locking device
AJ-002-002	Bowl only
AJ-002-003	Suction Ball only
AJ-002-004	Cone with tube only
AJ-002-005	Sheath connection for Ellik only



Not available *in the EU*

Ordering No.	Description
AJ-001-000	RECORD syringe with adapter, **50 ml**
AJ-001-010	Replacement cylinder for RECORD syringe with adapter, **50 ml**
AJ-001-001	RECORD syringe with adapter, **75 ml**
AJ-001-011	Replacement cylinder for RECORD syringe with adapter, **75 ml**
AJ-001-002	RECORD syringe with adapter, **100 ml**
AJ-001-012	Replacement cylinder for RECORD syringe with adapter, **100 ml**
AJ-001-003	RECORD syringe with adapter, **150 ml**
AJ-001-013	Replacement cylinder for RECORD syringe with adapter, **150 ml**

2.4 | ACCESSORIES
TRAY INSERT



Ordering No.	Description
AJ-010-020	Tray insert for resectoscopy **(without instruments, without sterilization box)**, Dimension: **370 × 255 × 45 mm (L × W × H)**



RECOMMENDED PRODUCTS:

AA-002-013	HD autoclavable scope 30° | 4 mm
AG-001-011	SALINE working element, passive cutting, handle closed
AH-001-302	Outer sheath 26 Fr. with perforation, rotatable
AH-001-301	Inner sheath 26 Fr. with ceramic isolation, rotatable
AH-001-011	Standard obturator for resectoscope sheaths 24 Fr.
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives
AE-001-031	Saline Cutting loop electrode 30° angled, 24 Fr. wire thickness 0.35 mm
AE-001-034	Saline rollerball electrode 3 mm, 24 Fr.
	+ connecting cable (depends on unit)

2.5 MONOPOLAR RESECTOSCOPY SET

1 **HD CYSTOSCOPE**
Ø 4.0 mm, 30° viewing angle, autoclavable, WL 302 mm

1 AA-002-013



2 **WORKING ELEMENT**
for resectoscopes, passive cutting, handle closed

2 AG-001-001



3 **OUTER SHEATH**
26 Fr. with perforation, rotatable

3 AH-001-302



4 **INNER SHEATH**
26 Fr. with ceramic isolation, rotatable

4 AH-001-301



5 **STANDARD OBTURATOR**
24 Fr.

5 AH-001-011





6 **CUTTING LOOP ELECTRODE**
30° angled, 24 Fr.

6 AE-001-001



7 **ROLLERBALL ELECTRODE 3 MM**
24 Fr.

7 AE-001-004



8 **KNIFE ELECTRODE**
24 Fr.

8 AE-001-016



2.5 MONOPOLAR RESECTOSCOPY SET

9 **HF MONOPOLAR CONNECTING CABLE**
for Berchtold, Martin, Aesculap, Ø 4 mm

9 AD-010-001



10 **RECORD SYRINGE**
100 ml

10 AJ-001-002



11 **ELLIK EVACUATOR**
complete

11 AJ-002-001





12 **BIOPSY FORCEPS**
double action

12 AB-004-001

13 **GRASPING FORCEPS**
double action

13 AB-004-003



14 **Adapter**
for use with resectoscope sheaths

14 AB-004-008



15 **PROTECTION SLEEVE**
For sterilization and storage of electrodes

15 AJ-010-006



2.5 SALINE RESECTOSCOPY SET

1 **HD CYSTOSCOPE**
Ø 4.0 mm, 30° viewing angle, autoclavable, WL 302 mm

2 **SALINE WORKING ELEMENT**
for resectoscopes, passive cutting, handle closed

3 **OUTER SHEATH**
26 Fr. with perforation, rotatable

4 **INNER SHEATH**
26 Fr. with ceramic isolation, rotatable

5 **STANDARD OBTURATOR**
24 Fr.

6 **CUTTING LOOP ELECTRODE**
30° angled, 24 Fr.

7 **ROLLERBALL ELECTRODE 3 MM**
24 Fr.

8 **KNIFE ELECTRODE**
24 Fr.

1 AA-002-013



2 AG-001-011



3 AH-001-302



4 AH-001-301



5 AH-001-011





6 AE-001-031



7 AE-001-034



8 AE-001-032



2.5 SALINE RESECTOSCOPY SET

9 **HF BIPOLAR CONNECTING CABLE**
for Valleylab, Lamidey, EMC, 4.5 m

9 AD-010-013



10 **RECORD SYRINGE**
100 ml

10 AJ-001-002



11 **ELLIK EVACUATOR**
complete

11 AJ-002-001





12 **BIOPSY FORCEPS**
double action

12 AB-004-001

13 **GRASPING FORCEPS**
double action

13 AB-004-003



14 **ADAPTER**
for use with resectoscope sheaths

14 AB-004-008



15 **PROTECTION SLEEVE**
for sterilization and storage of electrodes

15 AJ-010-006



3. URETHROTOMY

3.1 | URETHROTOME

3.2 | OPTICAL URETHROTOMY
CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
KNIVES

3.3 | URETHROTOMY SET

3. URETHROTOMY


21 Fr. 15.5 Fr.




1 CYSTOSCOPES



2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS




4 KNIVES



STEMA

3.1 | URETHROTOME



	Ordering No.	Description
	AI-002-001	OTIS urethrotome with **3 exchangeable tips inclusive 2 knives**


SPARE KNIVES



Ordering No.	Description
AI-002-002	Spare knives for OTIS urethrotome

3.2 | OPTICAL URETHROTOMY
CYSTOSCOPES | SHEATHS



Viewing angle	Description	Diameter	Working length	Ordering No.
 0°	Cystoscope autoclavable Cystoscope autoclavable **wide angle**	Ø 4.0 mm	302 mm	AA-002-001 AA-002-001-W
	HD Cystoscope autoclavable **HD** Cystoscope autoclavable **wide angle**			AA-002-011 AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable			AA-002-021



AI-001-001 AI-001-003

Ordering No.	Description
AI-001-001	Urethrotome sheath **21 Fr. with instrument channel 5 Fr.**
AI-001-003	Endoscope bridge **with 1 instrument channel**, fit with urethrotome sheath AI-001-001

3.2 | OPTICAL URETHROTOMY
SUPPLEMENTARY SHEATHS



Ordering No.	Description
AI-001-002	Supplementary sheath for AI-001-001, open to side



Ordering No.	Description
AI-001-006	Supplementary sheath 25 Fr. **for continuous flow**

21 Fr.

3.2 | OPTICAL URETHROTOMY
OBTURATORS



Ordering No.	Description
AI-001-004	**Standard** obturator for urethrotome sheath 21 Fr.



Ordering No.	Description
AI-001-005	Obturator for urethrotome sheath 21 Fr. **with 5 Fr. instrument channel**



Ordering No.	Description
AI-001-012	**Viewing** obturator for urethrotome sheath 21 Fr.

21 Fr.

3.2 | OPTICAL URETHROTOMY
WORKING ELEMENTS | PASSIVE CUTTING



Ordering No.	Description
AG-001-001	Working element used for resectoscopes, **passive cutting**, handle closed
AG-001-001-HF	Working element used for resectoscopes, **passive cutting**, handle closed**, HF connection on top**



Ordering No.	Description
AG-001-002	Working element used for resectoscopes, **passive cutting**, handle open

21 Fr.

3.2 | OPTICAL URETHROTOMY
COLD KNIVES



	Ordering No.	Description
	AI-001-007	Cold knife **straight**
	AI-001-008	Cold knife straight, **wave shape**
	AI-001-009	Cold knife **round**
	AI-001-011	Cold knife **hook shaped**




Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

15.5 Fr.

3.2 | OPTICAL URETHROTOMY
CYSTOSCOPES | SHEATHS | OBTURATORS



Viewing angle	Description	Diameter	Working length	Ordering No.
 0°	Cystoscope autoclavable	Ø 2.9 mm	302 mm	AA-003-001
	HD Cystoscope autoclavable			AA-003-011



AI-005-001 AI-005-003

Ordering No.	Description
AI-005-001	Urethrotome sheath **15.5 Fr. with instrument channel 3 Fr.**
AI-005-003	Endoscope bridge **with 1 instrument channel**, fit with urethrotome sheath AI-005-001



Ordering No.	Description
AI-005-004	Standard obturator for urethrotome **sheath 15.5 Fr.**

15.5 Fr.

3.2 | OPTICAL URETHROTOMY
WORKING ELEMENT



Ordering No.	Description
AG-001-010	Working element used for resectoscopes, passive cutting, **handle closed**, without HF-connection

15.5 Fr.

3.2 | OPTICAL URETHROTOMY
COLD KNIVES



	Ordering No.	Description
	AI-005-007	Cold knife **straight**
	AI-005-008	Cold knife straight, **wave shape**
	AI-005-011	Cold knife **hook shaped**



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

3.3 OPTICAL URETHROTOMY SET

1 **HD CYSTOSCOPE**
Ø 4.0 mm, 0° viewing angle, autoclavable, 302 mm working length

1 AA-002-011



2 **WORKING ELEMENT**
for resectoscopes, passive cutting, handle closed

2 AG-001-001



3 **URETHROTOME SHEATH**
21 Fr. with instrument channel 5 Fr.

3 AI-001-001



4 **SUPPLEMENTARY SHEATH**
Supplementary sheath for AI-001-001, open to side

4 AI-001-002



5 **STANDARD OBTURATOR**
21 Fr.

5 AI-001-004



6 **ENDOSCOPE BRIDGE**
with 1 instrument channel

6 AI-001-003



7 **COLD KNIFE**
straight

7 AI-001-007

8 **PROTECTION SLEEVE**
for sterilization and storage of electrodes

8 AJ-010-006



4. NEPHROSCOPY

4.1 | HD NEPHROSCOPY
HD NEPHROSCOPES
SHEATHS
OBTURATORS
HD MINI NEPHROSCOPE

4.2 | MINI NEPHROSCOPY
MINI NEPHROSCOPES
SHEATHS
DILATATORS
INSTRUMENTS

4.3 | STANDARD NEPHROSCOPY
NEPHROSCOPES 200 MM
SHEATHS
OBTURATORS

NEPHROSCOPES 250 MM
SHEATHS
OBTURATORS

4.4 | ACCESSORIES
TELESCOPE BOUGIES
FORCEPS
CANNULAS

4. NEPHROSCOPY

1 HD NEPHROSCOPES



2 SHEATHS & OBTURATORS



3 FORCEPS, INSTRUMENTS, KNIVES



STEMA

STEMA

4.1 | HD NEPHROSCOPY

HD NEPHROSCOPES | SHEATHS | OBTURATORS | SEALING CAPS



Ordering No.	Description	Working length	Viewing angle
AA-006-069	**HD Nephroscope** - **45°** angled ocular - **Ø 18.8 Fr.** - Working channel **Ø 12.9 Fr.**	**220 mm**	**0°**
AA-006-070			**15°**




Ordering No.	Description	French
AB-006-170	Nephroscope sheath for continuous flow, rotatable stopcock **(fit with nephroscopes WL 220mm / scope 15°)**	**23 Fr.**



Ordering No.	Description	French
AB-006-171	Obturator hollow **(fit with nephroscopy sheath AB-006-170 / scope 15°)**	**23 Fr.**


Ordering No.	Description
AA-006-081	Sealing cap, **Ø 15 mm/1.0 mm**
AA-006-083	Sealing cap, **Ø 20 mm/2.0 mm**

4.1 | HD MINI NEPHROSCOPE



Ordering No.	Description	Working length	Viewing angle
AA-006-100	**HD Mini compact nephroscope** - 30° angled ocular - Sheath **12 Fr.**, - Instrument channel **6.6 Fr**., - 3 luer lock connections **- 40,000 px**	**220 mm**	**0°**

4.2 | MINI NEPHROSCOPY
MINI NEPHROSCOPES




Ordering No.	Description	Working length	French	Viewing angle
AA-006-050	Mini nephroscope with **45°** angled ocular, instrument channel **10 Fr.**	**220 mm**	**Ø 19 Fr.**	**8°**
AA-006-051	Mini nephroscope with **45°** angled ocular, instrument channel **6 Fr.**	**220 mm**	**Ø 13.5 Fr.**	**8°**
AA-006-052	Mini nephroscope with **45°** angled ocular, instrument channel **6 Fr.**	**320 mm**	**Ø 13.5 Fr.**	**8°**



Ordering No.	Description
AA-006-053	Instrument bridge used for Mini-Nephroscopes (not included in mini Nephroscope)
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

4.2 | MINI NEPHROSCOPY
SHEATHS | DILATATORS | APPLICATOR



Ordering No.	Description	French
AB-006-100	Mini nephroscope sheath without perforation, continuous flow, stopcock fixed	**15/18 Fr.**
AB-006-101		**16.5/19.5 Fr.**
AB-006-102		**21/24 Fr.**




Ordering No.	Description	French
AB-006-110	Dilatator with central cannulation **Ø 3 Fr.**	**15/18 Fr.**
AB-006-111	Dilatator with 2 channels **Ø 2×3 Fr.**	**16.5/19.5 Fr.**
AB-006-112	Dilatator with 2 channels **Ø 2×5 Fr.**	**21/24 Fr.**





Ordering No.	Description
AB-006-115	Applicator used for the mini nephroscope sheaths, incl. plunger



STEMA


4.2 | MINI NEPHROSCOPY
INSTRUMENTS | FORCEPS FOR MINI NEPHROSCOPES




	Ordering No.	Description	Length	Diameter
	AB-006-242	Biopsy forceps serrated **with ring handle**	400 mm	**Ø 1.8 mm**
	AB-006-342		500 mm	**Ø 1.8 mm**
	AB-006-442		400 mm	**Ø 3.0 mm**
	AB-006-243	Biopsy forceps, single action **with ring handle**	400 mm	**Ø 1.8 mm**
	AB-006-343		500 mm	**Ø 1.8 mm**
	AB-006-443		400 mm	**Ø 3.0 mm**
	AB-006-244	Alligator grasping forceps **with ring handle**	400 mm	**Ø 1.8 mm**
	AB-006-344		500 mm	**Ø 1.8 mm**
	AB-006-444		400 mm	**Ø 3.0 mm**
	AB-006-245	Scissors, single action **with ring handle**	400 mm	**Ø 1.8 mm**
	AB-006-345		500 mm	**Ø 1.8 mm**
	AB-006-445		400 mm	**Ø 3.0 mm**
	AB-006-246	3-Prongs retriever for larger stones **with coil handle**	400 mm	**Ø 1.8 mm**
	AB-006-346		500 mm	**Ø 1.8 mm**
	AB-006-446		400 mm	**Ø 3.0 mm**

200 mm

4.3 | STANDARD NEPHROSCOPY
NEPHROSCOPES 200 MM



Ordering No.	Description	Working length	Viewing angle
AA-006-001	**Nephroscope** - **45°** angled ocular - **Ø 7.45 mm** - Suction channel **Ø 1.5 mm** - Working channel **Ø 5.0 mm**	**200 mm**	**0°**
AA-006-021			**6°**
AA-006-004			**15°**




Ordering No.	Description	Working length	Viewing angle
AA-006-002	**Nephroscope** - **90°** angled ocular - **Ø 7.45 mm** - Suction channel **Ø 1.5 mm** - Working channel **Ø 5.0 mm**	**200 mm**	**0°**
AA-006-022			**6°**
AA-006-005			**15°**




Ordering No.	Description
AA-006-080	Flow valve, Ø 6.1 mm, used for STEMA nephroscopes
AA-006-081	Sealing cap, Ø 15 mm/1.0 mm, **for AA-006-080**
AA-006-083	Sealing cap, Ø 20 mm/2.0 mm, **for AA-006-001–AA-006-034**
AJ-010-015	Stopcock left, large (passage 4 mm), with LL-connection and tube attachment M6 × 0.35 mm

200 mm

4.3 | STANDARD NEPHROSCOPY
SHEATHS | OBTURATORS



Ordering No.	Description	French
AB-006-001	Nephroscope sheath, without perforation, for continuous flow, **rotatable** stopcock, 200 mm scope	24 Fr.
AB-006-007		26 Fr.
AB-006-005	Nephroscope sheath, without perforation, for continuous flow, **fixed** stopcock, 200 mm scope	24 Fr.
AB-006-006		26 Fr.




Ordering No.	Description	French
AB-006-003	Nephroscope sheath for continuous flow, perforated, **rotatable** stopcock, 200 mm scope	26 Fr.
AB-006-008	Nephroscope sheath for continuous flow, perforated, **fixed** stopcock, 200 mm scope	26 Fr.



Ordering No.	Description	French
AB-006-002	Obturator hollow used for nephroscopy sheaths, 200 mm scope	24 Fr.
AB-006-004		26 Fr.


250 mm

4.3 | STANDARD NEPHROSCOPY
NEPHROSCOPES 250 MM



AA-006-080

Ordering No.	Description	Working length	Viewing angle
AA-006-011	**Nephroscope** - **45°** angled ocular - **Ø 7.45 mm** - Suction channel **Ø 1.5 mm** - Working channel **Ø 5.0 mm**	**250 mm**	**0°**
AA-006-031			**6°**
AA-006-034			**15°**



AA-006-080


Ordering No.	Description	Working length	Viewing angle
AA-006-012	**Nephroscope** - **90°** angled ocular - **Ø 7.45 mm** - Suction channel **Ø 1.5 mm** - Working channel **Ø 5.0 mm**	**250 mm**	**0°**
AA-006-032			**6°**
AA-006-006			**15°**




Ordering No.	Description
AA-006-080	Flow valve, Ø 6.1 mm, used for STEMA nephroscopes
AA-006-081	Sealing cap, Ø 15 mm/1.0 mm, **for AA-006-080**
AA-006-083	Sealing cap, Ø 20 mm/2.0 mm, **for AA-006-001–AA-006-034**
AJ-010-015	Stopcock left, large (passage 4 mm), with LL-connection and tube attachment M6 × 0.35 mm

250 mm

4.3 | STANDARD NEPHROSCOPY
SHEATHS | OBTURATORS



Ordering No.	Description	French
AB-006-061	Nephroscope sheath, without perforation, for continuous flow, **rotatable** stopcock, 250 mm scope	24 Fr.
AB-006-067		26 Fr.
AB-006-065	Nephroscope sheath, without perforation, for continuous flow, **fixed** stopcock, 250 mm scope	24 Fr.
AB-006-066		26 Fr.



Ordering No.	Description	French
AB-006-063	Nephroscope sheath for continuous flow, perforated, **rotatable stopcock,** 250 mm scope	26 Fr.
AB-006-068	Nephroscope sheath for continuous flow, perforated, **fixed stopcock,** 250 mm scope	26 Fr.



Ordering No.	Description	French
AB-006-062	Obturator hollow used for nephroscopy sheaths, 250 mm scope	24 Fr.
AB-006-064		26 Fr.

4.4 | ACCESSORIES
TELESCOPE BOUGIES | CANNULAS




Ordering No.	Description
AB-006-011	Set of telescope bougies consisting of: **6 dilating bougies** 9, 12, 15, 18, 21 and 24 Fr., incl. **1 rigid and 1 flexible guide probe**
AB-006-012	Guide probe, **rigid**
AB-006-013	Guide probe, **flexible**


Ordering no.	Description	French
AB-006-020	Dilating sheath	**9 Fr.**
AB-006-021		**12 Fr.**
AB-006-022		**15 Fr.**
AB-006-023		**18 Fr.**
AB-006-024		**21 Fr.**
AB-006-025		**24 Fr.**
AB-006-026		**27 Fr.**
AB-006-027		**30 Fr.**


Ordering No.	Description
AB-006-052	Puncture cannula, inner and outer cannula, 1 box = 5 pieces
AB-006-053	Dilatation cannula, inner and outer cannula, **Ø 3.0 mm**

STEMA

4.4 | ACCESSORIES
FORCEPS





	Ordering No.	Description	Length	Diameter
	AB-006-049	Biopsy forceps, single action **with ring handle**	380 mm	**Ø 3.5 mm**
	AB-006-043		430 mm	
	AB-006-050	Alligator grasping forceps **with ring handle**	380 mm	**Ø 3.5mm**
	AB-006-044		430 mm	
	AB-006-051	Scissors, single action **with ring handle**	380 mm	**Ø 3.5 mm**
	AB-006-045		430 mm	
	AB-006-047	3-Prongs retriever for larger stones **with coil handle**	380 mm	**Ø 3.5 mm**
	AB-006-041		430 mm	
	AB-006-046	Fenestrated stone grasping forceps for smaller stones, with **coil** handle	380 mm	**Ø 3.5 mm**
	AB-006-040		430 mm	
	AB-006-048	Fenestrated stone grasping forceps for smaller stones, **with ring handle**	380 mm	
	AB-006-042		430 mm	

5. URETERO-RENOSCOPY

5.1 | URETERO-RENOSCOPES SEMI-RIGID
HD URETERO-RENOSCOPES
URETERO-RENOSCOPES
BICHLER URETERO-RENOSCOPES
MARBERGER URETERO-RENOSCOPES
FORCEPS, SCISSORS, INSTRUMENTS
STONE EXTRACTORS
MONOPOLAR FLEXIBLE ELECTRODES
MONOPOLAR CONNECTING CABLES

5.2 | FLEXIBLE URETERO-RENOSCOPE
FLEX II
ACCESSORIES
INSTRUMENTS

5. HD URETERO-RENOSCOPES

40,000 px

17,000 px



1 URETERO-RENOSCOPE



2 FORCEPS



STEMA

40,000 px

5.1 | URETERO-RENOSCOPES SEMI-RIGID
HD URETERO-RENOSCOPES



Ordering No.	Description	Working length	Viewing angle
AA-007-003	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **7.5 to 12 Fr.** - Working channel **1 × 5 Fr. or 2 × 3 Fr.** - 40,000 px	**425 mm**	0°
AA-007-016			5°
AA-007-017			12°

Ordering No.	Description	Working length	Viewing angle
AA-007-004	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **7.5 to 12 Fr.** - Working channel **1 × 5 Fr. or 2 × 3 Fr.** - 40,000 px	**310 mm**	0°
AA-007-018			5°
AA-007-019			12°

Ordering No.	Description	Working length	Viewing angle
AA-007-005	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **7.5 to 12 Fr.** - Working channel **1 × 5 Fr. or 2 ×3 Fr.** - 40,000 px	**150 mm**	0°
AA-007-036			5°
AA-007-037			12°


Delivery contains
HD uretero-renoscope, brush Ø 3.0 mm (length 590 mm), brush Ø 3.8 mm (length 350 mm), sealing cap silicone (grey), Duckbill sealing



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

40,000 px

5.1 | URETERO-RENOSCOPES SEMI-RIGID
HD URETERO-RENOSCOPES



Ordering No.	Description	Working length	Viewing angle
AA-007-009	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **7.5 to 13.2 Fr.** - Working channel **1 × 5 Fr. or 2 × 3 Fr.** - 40,000 px	**425 mm**	0°
AA-007-028			5°
AA-007-029			12°

Delivery contains
HD uretero-renoscope, brush Ø 3.0 mm (length 590 mm), brush Ø 3.8 mm (length 350 mm), sealing cap silicone (grey), Duckbill sealing



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

17,000 px

5.1 | URETERO-RENOSCOPES SEMI-RIGID
HD URETERO-RENOSCOPES



Ordering No.	Description	Working length	Viewing angle
AA-007-015	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **6 to 7.5 Fr.** - Working channel **1 × 4 Fr. or 2 × 2.4 Fr.** - 17,000 px	**425 mm**	0°
AA-007-034			5°
AA-007-035			12°

Ordering No.	Description	Working length	Viewing angle
AA-007-014	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **6 to 7.5 Fr.** - Working channel **1 × 4 Fr. or 2 × 2.4 Fr.** - 17,000 px	**310 mm**	0°
AA-007-032			5°
AA-007-033			12°

Ordering No.	Description	Working length	Viewing angle
AA-007-013	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **6 to 7.5 Fr.** - Working channel **1 × 4 Fr. or 2 × 2.4 Fr.** - 17,000 px	**140 mm**	0°
AA-007-030			5°
AA-007-031			12°

Delivery contains
HD uretero-renoscope, brush Ø 3.0 mm (length 590 mm), brush Ø 3.8 mm (length 350 mm), sealing cap silicone (grey), Duckbill sealing



Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

17,000 px

5.1 | 5.1 | URETERO-RENOSCOPES SEMI-RIGID
HD URETERO-RENOSCOPES



Ordering No.	Description	Working length	Viewing angle
AA-007-006	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **4.9 to 6.5 Fr.** - Working channel **1 × 3 Fr.** - 17,000 px	**425 mm**	0°
AA-007-022			5°
AA-007-023			12°

Ordering No.	Description	Working length	Viewing angle
AA-007-007	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **4.9 to 6.5 Fr.** - Working channel **1 × 3 Fr.** - 17,000 px	**310 mm**	0°
AA-007-024			5°
AA-007-025			12°


Ordering No.	Description	Working length	Viewing angle
AA-007-008	**HD Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **4.9 to 6.5 Fr.** - Working channel **1 × 3 Fr.** - 17,000 px	**150 mm**	0°
AA-007-026			5°
AA-007-027			12°


Delivery contains
HD uretero-renoscope, brush Ø 3.0 mm (length 590 mm), brush Ø 3.8 mm (length 350 mm), sealing cap silicone (grey), Duckbill sealing




Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

5.1 | URETERO-RENOSCOPES SEMI-RIGID
URETERO-RENOSCOPES







Ordering No.	Description	Working length	Viewing angle
AA-007-001	**Uretero-renoscope semi-rigid** - 45° angled ocular - Ø Working steps • **distal 7.5 Fr.** • **first step 9.8 Fr.** • **tube 13 Fr.** - Ø Working channel **1 × 5 Fr. or 2 × 3 Fr.** - Field of view: **85°** - Image bundles: 30 k	**310 mm**	8°
AA-007-002		**425 mm**	





Ordering No.	Description
AA-007-021	Instrument bridge used for uretero-renoscopes **(included in the delivery content)**




AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

5.1 | URETERO-RENOSCOPES SEMI-RIGID
URETERO-RENOSCOPES







Ordering No.	Description	Working length	Viewing angle
AA-007-011	**Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps: • **distal 6 Fr.** • **first step 7.5 Fr.** • **second step 9.9 Fr.** • **tube 10.6 Fr.** - Working channel **1 × 4 Fr. or 2 × 2.4 Fr.** - Field of view: **100°** - Image bundles: 17 k	**310 mm**	8°
AA-007-012		**425 mm**	




Ordering No.	Description
AA-007-021	Instrument bridge used for uretero-renoscopes **(included in the delivery content)**




AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm

STEMA


5.1 | URETERO-RENOSCOPES SEMI-RIGID
BICHLER URETERO-RENOSCOPES




Ordering No.	Description	Working length	Viewing angle
AA-007-056	**BICHLER Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **4.5 Fr. to 6.5 Fr.** - Working channel **1 × 3.3 Fr.**	**310 mm**	5°
AA-007-055		**430 mm**	




Ordering No.	Description	Working length	Viewing angle
AA-007-043	**BICHLER Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **6 Fr. to 7.5 Fr.** - Working channel **1 × 4 Fr. or 2 × 2.4 Fr.**	**310 mm**	12°
AA-007-044		**430 mm**	



Ordering No.	Description	Working length	Viewing angle
AA-007-047	**BICHLER Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **8 Fr. to 9.8 Fr.** - Working channel **1 × 5 Fr. or 2 × 3 Fr.**	**310 mm**	12°
AA-007-048		**430 mm**	




Ordering No.	Description	Working length	Viewing angle
AA-007-051	**BICHLER Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **8.5 Fr. to 11.5 Fr.** - Working channel **1 × 6 Fr. or 2 × 4 Fr.**	**310 mm**	12°
AA-007-052		**430 mm**	






Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm


5.1 | URETERO-RENOSCOPES SEMI-RIGID
MARBERGER URETERO-RENOSCOPES




Ordering No.	Description	Working length	Viewing angle
AA-007-054	**MARBERGER Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **4.5 Fr. to 6.5 Fr.** - Working channel **1 × 3.3 Fr.**	**310 mm**	5°
AA-007-053		**430 mm**	




Ordering No.	Description	Working length	Viewing angle
AA-007-041	**MARBERGER Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **6 Fr. to 7.5 Fr.** - Working channel **1 × 4 Fr. or 2 × 2.4 Fr.**	**310 mm**	12°
AA-007-042		**430 mm**	



Ordering No.	Description	Working length	Viewing angle
AA-007-045	**MARBERGER Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **8 Fr. to 9.8 Fr.** - Working channel **1 × 5 Fr. or 2 × 3 Fr.**	**310 mm**	12°
AA-007-046		**430 mm**	



Ordering No.	Description	Working length	Viewing angle
AA-007-049	**MARBERGER Uretero-renoscope semi-rigid** - 45° angled ocular - Working Ø steps by **8.5 Fr. to 11.5 Fr.** - Working channel **1 × 6 Fr. or 2 × 4 Fr.**	**310 mm**	12°
AA-007-050		**430 mm**	






Ordering No.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0.6-0.8 mm


STEMA


5.1 | URETERO-RENOSCOPY
FORCEPS FLEXIBLE | SEMI-RIGID | RIGID




	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps oval, double action	**Ø 1.0 mm**	3 Fr.	AC-001-011	AC-001-211	AC-001-011-R
		Ø 1.2 mm	4 Fr.	AC-001-010	AC-001-210	AC-001-010-R
		Ø 1.6 mm	5 Fr.	AC-001-012	AC-001-212	AC-001-012-R
		Ø 2.1 mm	6 Fr.	AC-001-013	AC-001-213	AC-001-013-R
		Ø 2.3 mm	7 Fr.	AC-001-014	AC-001-214	AC-001-014-R
		Ø 3.0 mm	9 Fr.	AC-001-015	AC-001-215	AC-001-015-R
	Biopsy forceps oval, serrated, double action	**Ø 1.0 mm**	3 Fr.	AC-002-011	AC-002-211	AC-002-011-R
		Ø 1.6 mm	5 Fr.	AC-002-012	AC-002-212	AC-002-012-R
		Ø 2.1 mm	6 Fr.	AC-002-013	AC-002-213	AC-002-013-R
		Ø 2.3 mm	7 Fr.	AC-002-014	AC-002-214	AC-002-014-R
		Ø 3.0 mm	9 Fr.	AC-002-015	AC-002-215	AC-002-015-R
	Biopsy forceps oval, serrated, teeth interlock, double action	**Ø 1.6 mm**	5 Fr.	AC-015-012	AC-015-212	AC-015-012-R
		Ø 2.1 mm	6 Fr.	-	AC-015-213	AC-015-013-R
		Ø 2.3 mm	7 Fr.	AC-015-014	AC-015-214	AC-015-014-R
	Grasping forceps, double action	**Ø 1.2 mm**	4 Fr.	AC-017-011	AC-017-211	AC-017-011-R
		Ø 2.3 mm	7 Fr.	AC-017-014	AC-017-214	AC-017-014-R
	Grasping forceps, alligator jaw, double action	**Ø 1.0 mm**	3 Fr.	AC-006-011	AC-006-211	AC-006-011-R
		Ø 1.2 mm	4 Fr.	AC-006-013	AC-006-213	AC-006-013-R
		Ø 1.6 mm	5 Fr.	AC-006-012	AC-006-212	AC-006-012-R
		Ø 2.3 mm	7 Fr.	AC-006-014	AC-006-214	AC-006-014-R
		Ø 3.0 mm	9 Fr.	AC-006-015	AC-006-215	AC-006-015-R

5.1 | URETERO-RENOSCOPY
FORCEPS FLEXIBLE | SEMI-RIGID | RIGID




	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	**Ø 2.3 mm**	7 Fr.	AC-008-014	AC-008-214	AC-008-014-R
		Ø 3.0 mm	9 Fr.	AC-008-015	AC-008-215	AC-008-015-R
	Grasping forceps teeth double action	**Ø 1.2 mm**	4 Fr.	AC-016-011	AC-016-211	AC-016-011-R
		Ø 1.6 mm	5 Fr.	AC-016-012	AC-016-212	AC-016-012-R
	Grasping forceps teeth interlock, double action	**Ø 1.6 mm**	5 Fr.	AC-020-012	AC-020-212	AC-020-012-R
		Ø 2.3 mm	7 Fr.	AC-020-014	AC-020-214	AC-020-014-R
	Grasping forceps rat tooth, double action	**Ø 1.0 mm**	3 Fr.	AC-019-011	AC-019-211	AC-019-011-R
		Ø 1.6 mm	5 Fr.	AC-019-012	AC-019-212	AC-019-012-R
	Grasping forceps for stone street, double action	**Ø 1.2 mm**	4 Fr.	AC-007-011	AC-007-211	AC-007-011-R
		Ø 1.6 mm	5 Fr.	AC-007-012	AC-007-212	AC-007-012-R
	Grasping forceps for big stones, double action	**Ø 1.2 mm**	4 Fr.	AC-009-011	AC-009-211	AC-009-011-R
		Ø 1.6 mm	5 Fr.	AC-009-012	AC-009-212	AC-009-012-R
		Ø 2.3 mm	7 Fr.	AC-009-014	AC-009-214	AC-009-014-R
	3-Prongs retriever	**Ø 1.0 mm**	3 Fr.	AC-012-011		
		Ø 1.6 mm	5 Fr.	AC-012-012		
		Ø 2.3 mm	7 Fr.	AC-012-014		

5.1 | URETERO-RENOSCOPY
SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Scissors, single action	**Ø 1.0 mm**	3 Fr.	AC-003-011	AC-003-211	AC-003-011-R
		Ø 1.2 mm	4 Fr.	AC-003-010	AC-003-210	AC-003-010-R
		Ø 1.6 mm	5 Fr.	AC-003-012	AC-003-212	AC-003-012-R
		Ø 2.1 mm	6 Fr.	AC-003-013	AC-003-213	AC-003-013-R
		Ø 2.3 mm	7 Fr.	AC-003-014	AC-003-214	AC-003-014-R
		Ø 3.0 mm	9 Fr.	AC-003-015	AC-003-215	AC-003-015-R
	Hook scissors, single action	**Ø 2.3 mm**	7 Fr.	AC-004-014	AC-004-214	AC-004-014-R
		Ø 3.0 mm	9 Fr.	AC-004-015	AC-004-215	AC-004-015-R
	Uretero-slitting scissors, retrograde cutting single action	**Ø 1.2 mm**	4 Fr.	AC-018-011	AC-018-211	AC-018-011-R
		Ø 1.6 mm	5 Fr.	AC-018-012	AC-018-212	AC-018-012-R
	Scissors, serrated flexible, double action	**Ø 2.3 mm**	7 Fr.	AC-005-014	AC-005-214	AC-005-014-R

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




FEATURES

1. QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
2. ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
3. OVERLOAD PROTECTION AVAILABLE
4. SIDE-MOUNTED LUER LOCK

STEMA

5.1 | URETERO-RENOSCOPY
INSTRUMENTS SEMI-RIGID | DETACHABLE




Ordering No.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Scissors, single action	**Ø 1.0 mm**	3 Fr.	AC-003-211-T	AC-003-211-OT
		Ø 1.2 mm	4 Fr.	AC-003-210-T	AC-003-210-OT
		Ø 1.6 mm	5 Fr.	AC-003-212-T	AC-003-212-OT
		Ø 2.3 mm	7 Fr.	AC-003-214-T	AC-003-214-OT
		Ø 3.0 mm	9 Fr.	AC-003-215-T	AC-003-215-OT
	Scissors pointed, single action **(Minimum order quantity 5pcs.)**	**Ø 1.6 mm**	5 Fr.	AC-004-622-T	AC-004-622-OT
		Ø 2.3 mm	7 Fr.	AC-004-624-T	AC-004-624-OT
	Hook scissors, single action	**Ø 2.3 mm**	7 Fr.	AC-004-214-T	AC-004-214-OT
		Ø 3.0 mm	9 Fr.	AC-004-215-T	AC-004-215-OT


5.1 | URETERO-RENOSCOPY
INSTRUMENTS SEMI-RIGID | DETACHABLE



Ordering No.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Biopsy forceps oval, double action	**Ø 1.0 mm**	3 Fr.	AC-001-211-T	AC-001-211-OT
		Ø 1.2 mm	4 Fr.	AC-001-210-T	AC-001-210-OT
		Ø 1.6 mm	5 Fr.	AC-001-212-T	AC-001-212-OT
		Ø 2.3 mm	7 Fr.	AC-001-214-T	AC-001-214-OT
		Ø 3.0 mm	9 Fr.	AC-001-215-T	AC-001-215-OT
	Biopsy forceps oval, serrated, double action	**Ø 1.2 mm**	4 Fr.	AC-016-211-T	AC-016-211-OT
		Ø 1.6 mm	5 Fr.	AC-016-212-T	AC-016-212-OT
	Grasping forceps, alligator jaw, double action	**Ø 1.0 mm**	3 Fr.	AC-006-211-T	AC-006-211-OT
		Ø 1.2 mm	4 Fr.	AC-006-213-T	AC-006-213-OT
		Ø 1.6 mm	5 Fr.	AC-006-212-T	AC-006-212-OT
		Ø 2.3 mm	7 Fr.	AC-006-214-T	AC-006-214-OT
		Ø 3.0 mm	9 Fr.	AC-006-215-T	AC-006-215-OT
	Grasping forceps for stone street, double action	**Ø 1.2 mm**	4 Fr.	AC-007-211-T	AC-007-211-OT
		Ø 1.6 mm	5 Fr.	AC-007-212-T	AC-007-212-OT
	Grasping forceps for big stones, double action	**Ø 1.2 mm**	4 Fr.	AC-009-211-T	AC-009-211-OT
		Ø 1.6 mm	5 Fr.	AC-009-212-T	AC-009-212-OT

5.1 | URETERO-RENOSCOPY
INSTRUMENTS RIGID | DETACHABLE




Ordering No.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Scissors, single action	**Ø 1.0 mm**	3 Fr.	AC-103-211-T	AC-103-211-OT
		Ø 1.2 mm	4 Fr.	AC-103-210-T	AC-103-210-OT
		Ø 1.6 mm	5 Fr.	AC-103-212-T	AC-103-212-OT
		Ø 2.3 mm	7 Fr.	AC-103-214-T	AC-103-214-OT
		Ø 3.0 mm	9 Fr.	AC-103-215-T	AC-103-215-OT
	Hook scissors, single action	**Ø 2.3 mm**	7 Fr.	AC-104-214-T	AC-104-214-OT
		Ø 3.0 mm	9 Fr.	AC-104-215-T	AC-104-215-OT
	Biopsy punch forceps, single action	**Ø 1.6 mm**	5 Fr.	AC-114-212-T	AC-114-212-OT

5.1 | URETERO-RENOSCOPY
INSTRUMENTS RIGID | DETACHABLE



Ordering No.	Description
AC-100-001	Handle for detachable instruments

	Description	Diameter	French	Complete	Insert
	Biopsy forceps oval, double action	**Ø 1.0 mm**	3 Fr.	AC-101-211-T	AC-101-211-OT
		Ø 1.2 mm	4 Fr.	AC-101-210-T	AC-101-210-OT
		Ø 1.6 mm	5 Fr.	AC-101-212-T	AC-101-212-OT
		Ø 2.3 mm	7 Fr.	AC-101-214-T	AC-101-214-OT
		Ø 3.0 mm	9 Fr.	AC-101-215-T	AC-101-215-OT
	Biopsy forceps oval, serrated, double action	**Ø 1.2 mm**	4 Fr.	AC-116-211-T	AC-116-211-OT
		Ø 1.6 mm	5 Fr.	AC-116-212-T	AC-116-212-OT
	Grasping forceps, alligator jaw, double action	**Ø 1.0 mm**	3 Fr.	AC-106-211-T	AC-106-211-OT
		Ø 1.2 mm	4 Fr.	AC-106-213-T	AC-106-213-OT
		Ø 1.6 mm	5 Fr.	AC-106-212-T	AC-106-212-OT
		Ø 2.3 mm	7 Fr.	AC-106-214-T	AC-106-214-OT
		Ø 3.0 mm	9 Fr.	AC-106-215-T	AC-106-215-OT
	Grasping forceps for stone street, double action	**Ø 1.2 mm**	4 Fr.	AC-107-211-T	AC-107-211-OT
		Ø 1.6 mm	5 Fr.	AC-107-212-T	AC-107-212-OT
	Grasping forceps for big stones, double action	**Ø 1.2 mm**	4 Fr.	AC-109-211-T	AC-109-211-OT
		Ø 1.6 mm	5 Fr.	AC-109-212-T	AC-109-212-OT

STEMA

5.1 | URETERO-RENOSCOPY
STONE EXTRACTORS | 4 WIRES



Ordering No.	Description
AF-001-021	Handle for dormia stone extractor **3 and 5 Fr.**
AF-001-024	Handle for dormia stone extractor **7 Fr.**

Ordering No.	Wires	French	Working length	Description
AF-001-001	4	**3 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-002	4	**3 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-015	4	**3 Fr.**	90 cm	Dormia stone extractor, basket with **tubus only**
AF-001-016	4	**3 Fr.**	90 cm	Dormia stone extractor complete with **handle**
AF-001-023	4	**3 Fr.**	120 cm	Dormia stone extractor, basket with **tubus only**
AF-001-022	4	**3 Fr.**	120 cm	Dormia stone extractor complete with **handle**
AF-001-003	4	**5 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-004	4	**5 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-017	4	**5 Fr.**	90 cm	Dormia stone extractor, basket with **tubus only**
AF-001-019	4	**5 Fr.**	90 cm	Dormia stone extractor complete with **handle**
AF-001-005	4	**7 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-006	4	**7 Fr.**	60 cm	Dormia stone extractor complete with **handle**

STEMA

5.1 | URETERO-RENOSCOPY
STONE EXTRACTORS | 5-7 WIRES



Ordering No.	Description
AF-001-021	Handle for dormia stone extractor **3 and 5 Fr.**
AF-001-024	Handle for dormia stone extractor **7 Fr.**

Ordering No.	Wires	French	Working length	Description
AF-001-007	5	**5 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-008	5	**5 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-025	5	**5 Fr.**	90 cm	Dormia stone extractor, basket with **tubus only**
AF-001-018	5	**5 Fr.**	90 cm	Dormia stone extractor complete with **handle**
AF-001-009	5	**7 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-010	5	**7 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-011	6	**5 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-012	6	**5 Fr.**	60 cm	Dormia stone extractor complete with **handle**
AF-001-013	7	**7 Fr.**	60 cm	Dormia stone extractor, basket with **tubus only**
AF-001-014	7	**7 Fr.**	60 cm	Dormia stone extractor complete with **handle**

5.1 | URETERO-RENOSCOPY
MONOPOLAR FLEXIBLE ELECTRODES



	Description	French	Length	Ordering No.	French	Length	Ordering No.
	Loop electrode flexible, monopolar	**5 Fr.**	53 cm	AD-003-003	**5 Fr.**	65 cm	AD-003-013
		7 Fr.	53 cm	AD-003-005	**5 Fr.**	65 cm	AD-003-015
		9 Fr.	53 cm	AD-003-007	**5 Fr.**	65 cm	AD-003-017
		10 Fr.	53 cm	AD-003-008	**5 Fr.**	65 cm	AD-003-018
	Needle electrode flexible, monopolar	**3 Fr.**	53 cm	AD-002-001	**3 Fr.**	65 cm	AD-002-011
		4 Fr.	53 cm	AD-002-002	**4 Fr.**	65 cm	AD-002-012
		5 Fr.	53 cm	AD-002-003	**5 Fr.**	65 cm	AD-002-013
		7 Fr.	53 cm	AD-002-005	**7 Fr.**	65 cm	AD-002-015
	Button electrode flexible, monopolar	**3 Fr.**	53 cm	AD-001-001	**3 Fr.**	65 cm	AD-001-011
		4 Fr.	53 cm	AD-001-002	**4 Fr.**	65 cm	AD-001-012
		5 Fr.	53 cm	AD-001-003	**5 Fr.**	65 cm	AD-001-013
		6 Fr.	53 cm	AD-001-004	**6 Fr.**	65 cm	AD-001-014
		7 Fr.	53 cm	AD-001-005	**7 Fr.**	65 cm	AD-001-015
		8 Fr.	53 cm	AD-001-006	**8 Fr.**	65 cm	AD-001-016
		9 Fr.	53 cm	AD-001-007	**9 Fr.**	65 cm	AD-001-017
		10 Fr.	53 cm	AD-001-008	**10 Fr.**	65 cm	AD-001-018

5.1 | URETERO-RENOSCOPY
MONOPOLAR CONNECTING CABLES



Connection instrument side





Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001 AD-010-006	3 m 5 m	Ø 4.0 mm	HF Monopolar connecting cable for **Berchtold, Martin, Aesculap**
	AD-010-002 AD-010-007	3 m 5 m	Ø 5.0 mm	HF Monopolar connecting cable for **Erbe, ACC/ICC, Storz**
	AD-010-003 AD-010-008	3 m 5 m	Ø 8.0 mm	HF Monopolar connecting cable for **STEMA HF-U Series, Valleylab/Conmed**
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable **male** for **STEMA MASTER**





Ordering No.	Description
EK-015-060	**Universal adapter** used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

5.2 | FLEXIBLE URETERO-RENOSCOPE
FLEX II



- **High luminance**
- **Exact adjustment of focus**
- **Precise up-/down angulation**

Ordering No.	Description
AA-008-003	Flexible uretero-renoscope **FLEX II**

Technical data

Cystoscope

Insertion tube	Ø 2.8 mm
Working Channel	Ø 1.2 mm
Working length	680 mm
Field of view	88°
Pixel	18,000
Angulation up-down	275°
Accessories	Suction valve, leak tester

Delivery contains
Flexible uretero-renoscope FLEX II, leak tester, suction valve

INCLUDED IN SCOPE OF DELIVERY:

AA-008-103
Leak tester

AA-008-105
Suction valve



AA-008-103



AA-008-105

5.2 | FLEXIBLE URETERO-RENOSCOPE
ACCESSORIES



Ordering No.	Description
AA-008-100	Bridge with Luer-Lock and 2 stopcocks for **AA-008-002**

Ordering No.	Description
AA-008-101	Single use cleaning brushes, 1.0 mm, 120 cm length, **1 box = 250 pieces**



Ordering No.	Description
AA-008-102	Suitcase grey for STEMA flexible scopes



Ordering No.	Description
AA-008-103	Leak tester





Ordering No.	Description
AA-008-104	Single use Luer-Lock adapter, 1 box = 25 pieces



Ordering No.	Description
AA-008-105	Suction valve

5.2 | FLEXIBLE URETERO-RENOSCOPE
INSTRUMENTS FLEXIBLE



	Description	Diameter	French	Complete
	Biopsy forceps oval, flexible, double action	**Ø 1.0 mm**	3 Fr.	**AC-001-031**
	Grasping forceps alligator jaw, flexible, double action	**Ø 1.0 mm**	3 Fr.	**AC-006-031**
	3-Prongs retriever, flexible, with handle	**Ø 1.0 mm**	3 Fr.	**AC-012-015**




6. LITHOTRIPSY


6.1 | LITHOTRIPSY
CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS

6. LITHOTRIPSY


1 SCOPES



2 WORKING ELEMENTS



3 SHEATHS
& OBTURATORS



STEMA

6.1 | LITHOTRIPSY
CYSTOSCOPES | SHEATHS



Viewing angle	Description	Diameter	Working length	Ordering No.
	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-001
	HD Cystoscope autoclavable			AA-002-011

Viewing angle	Description	Diameter	Working length	Ordering No.
 30°	Cystoscope autoclavable	Ø 4.0 mm	302 mm	AA-002-003
	HD Cystoscope autoclavable			AA-002-013

6.1 | LITHOTRIPSY
CYSTOSCOPES | SHEATHS



Ordering No.	Description
AB-002-010	Lithotripsy sheath with central valve used for **AB-004-011**



Ordering No.	Description
AB-002-013	Lithotripsy sheath used for **AB-004-011**

6.1 | LITHOTRIPSY
SHEATHS | OBTURATORS



Ordering No.	Description
AB-002-011	Diagnostic sheath **with working channel for instruments 7 Fr.**



Ordering No.	Description
AB-002-014	Adapter for use with **Ellik**



Ordering No.	Description
AB-002-012	Obturator for **AB-002-010 + AB-002-013**

6.1 | LITHOTRIPSY
WORKING ELEMENT | LITHOTRIPTOR



Ordering No.	Description
AB-004-011	Working element used for punch lithotripsy



Ordering No.	Description
AB-004-012	Mechanical lithotriptor with handle and turning screw, **Ø 24 Fr., 2 stopcocks** with luer lock adapter, for use with scope **AA-002-004**

7. MORCELLATION


7.1 | MORCELLATION 26 FR.
NEPHROSCOPE
SHEATHS
OBTURATORS

7.2 | MORCELLATION 22 FR.
MORCESCOPE
SHEATHS
OBTURATORS

7. MORCELLATION


26 Fr. 22 Fr.

1 HD NEPHROSCOPE



2 SHEATHS & OBTURATORS

3 MORCESCOPE



4 SHEATHS & OBTURATORS

STEMA

26 Fr.

7.1 | MORCELLATION
NEPHROSCOPE




Ordering No.	Description	Working length	Viewing angle
AA-006-012	**Nephroscope** - **90°** angled ocular - **Ø 7.45 mm** - Suction channel **Ø 1.5 mm** - Working channel **Ø 5.0 mm**	**250 mm**	**0°**

	Ordering No.	Description
	AA-006-080	Flow valve, Ø 6.1 mm, used for STEMA nephroscopes
	AA-006-081	Sealing cap, Ø 15 mm/1.0 mm, **for AA-006-080**
	AA-006-083	Sealing cap, Ø 20 mm/2.0 mm, **for AA-006-001–AA-006-034**
	AB-004-010-QL	Adapter for use with nephroscope and resectoscope sheaths 26 Fr. STEMA Click-System

26 Fr.

7.1 | MORCELLATION
SHEATHS



	Ordering No.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, **rotatable**	**26 Fr.**	Yellow



Ordering No.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, **rotatable**	**26 Fr.**	Yellow

26 Fr.

7.1 | MORCELLATION
OBTURATORS



Ordering No.	Description	French	For outer sheath	Colour code
AH-001-011	**Standard** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow




Ordering No.	Description	French	For outer sheath	Colour code
AH-001-013	**Viewing** obturator for resectoscope sheaths	**24 Fr.**	**26 Fr.**	Yellow



Ordering No.	Description	French	Colour code
AH-001-308	Inner sheath without ceramic isolation, **with laser channel up** to 7 Fr., for use with AH-001-301	**26 Fr.**	Yellow

22 Fr.

7.2 | MORCELLATION
MORCESCOPE



Ordering No.	Description	Working length	Viewing angle
AA-006-072	**Morcescope** - 22 Fr. outer tube - Working channel 15 Fr. - Field of view 85° - 40,000 px	**180 mm**	**0°**

Ordering No.	Description
AA-006-084	Sealing cap white, **Ø 13.5 mm,** 1 box = 10 pieces

22 Fr.

7.2 | MORCELLATION
SHEATHS





Ordering No.	Description	French	Colour code
AH-002-302	Outer sheath with perforation, **rotatable**	**22 Fr.**	White



Ordering No.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, **rotatable**	**19 Fr.**	White

22 Fr.

7.2 | MORCELLATION
OBTURATORS



Ordering No.	Description	French	For outer sheath	Colour code
AH-002-011	**Standard** obturator for resectoscope sheaths	**19 Fr.**	**22 Fr.**	White



Ordering No.	Description	French	For outer sheath	Colour code
AH-002-013	**Viewing** obturator	**19 Fr.**	**22 Fr.**	White

8. PEDIATRIC UROLOGY


8.1 | PEDIATRIC CYSTO-URETHROSCOPY
CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
FORCEPS
MONOPOLAR ELECTRODES
CONNECTING CABLES

8.2 | PEDIATRIC RESECTOSCOPY
CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
ELECTRODES
CONNECTING CABLES


8.3 | PEDIATRIC URETHROTOMY
CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
COLD KNIVES AND ACCESSORIES

8. PEDIATRIC UROLOGY


1 PEDIATRIC CYSTOSCOPY



2 PEDIATRIC RESECTOSCOPY



3 PEDIATRIC URETHROTOMY



STEMA

13 Fr.

8.1 | PEDIATRIC CYSTO-URETHROSCOPY
CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-005-001	Pediatric cystoscope autoclavable	**0°**	Ø 2.7 mm	187.5 mm
AA-005-003		**30°**		
AA-005-004		**70°**		



Ordering No.	Description	French	Colour code
AB-001-001	Cysto-Urethroscope sheath with **2 fixed stopcocks,** **1 × 3 Fr. instrument channel** incl. obturator	**13 Fr.**	Red
AB-001-001-13	Cysto-Urethroscope sheath with **2 fixed stopcocks,** **1 × 5 Fr. or 2×4 Fr. instrument channel** incl. obturator		



Ordering No.	Description	French	Colour code
AB-001-011	Standard obturator	**13 Fr.**	Red

13 Fr.

8.1 | PEDIATRIC CYSTO-URETHROSCOPY
WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering No.	Description	For sheath	Colour code
AB-003-005	Working element with **ALBARRAN deflector**	**13 Fr.**	Red



Ordering No.	Description	For sheath	Colour code
AB-003-001	Endoscope bridge	**13 Fr.**	Red

9.5 Fr.
11 Fr.

8.1 | PEDIATRIC CYSTO-URETHROSCOPY
CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-004-011	Pediatric cystoscope autoclavable	**0°**	Ø 2.0 mm	187 mm
AA-004-012		**30°**		



Ordering No.	Description	French	Colour code
AB-001-001-9	Cysto-Urethroscope sheath with **2 fixed stopcocks,** **1 × 3 Fr. instrument channel** incl. obturator	**9.5 Fr.**	Blue
AB-001-001-11		**11 Fr.**	Green




Ordering No.	Description	French	Colour code
AB-001-011-9	Standard obturator	**9.5 Fr.**	Blue
AB-001-011-11		**11 Fr.**	Green

7.5 Fr.

8.1 | PEDIATRIC CYSTO-URETHROSCOPY
CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-004-013	Pediatric cystoscope autoclavable	**0°**	Ø 1.9 mm	145 mm
AA-004-014		**30°**		




Ordering No.	Description	French	Colour code
AB-001-001-7	Cysto-Urethroscope sheath with **2 fixed stopcocks**	**7.5 Fr.**	Yellow



Ordering No.	Description	French	Colour code
AB-001-011-7	Standard obturator for cysto-urethroscope sheaths	**7.5 Fr.**	Yellow

8.1 | PEDIATRIC CYSTO-URETHROSCOPY
FORCEPS | ELECTRODES



	Description	WL	Diameter	French	Flexible	Semi-rigid
	Biopsy forceps oval, double action	28 cm	Ø 1.0 mm	3 Fr.	AC-001-021	AC-001-221
	Biopsy forceps oval, serrated, teeth interlock, double action	28 cm	Ø 1.0 mm	3 Fr.	AC-015-021	AC-015-221
	Grasping forceps, alligator jaw, double action	28 cm	Ø 1.0 mm	3 Fr.	AC-006-021	AC-006-221
	Grasping forceps, rat tooth, double action	28 cm	Ø 1.0 mm	3 Fr.	AC-019-021	-
	Scissors, single action	28 cm	Ø 1.0 mm	3 Fr.	AC-003-021	AC-003-221



	Ordering No.	Description	French	Diameter
	AD-002-021	Needle electrode flexible, monopolar	3.0 Fr.	280 mm
	AD-001-021	Button electrode flexible, monopolar		

8.1 | PEDIATRIC CYSTO-URETHROSCOPY
CONNECTING CABLES | MONOPOLAR



Connection instrument side





2 mm

Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001	3 m	Ø 4.0 mm	HF Monopolar connecting cable for **Berchtold, Martin, Aesculap**
	AD-010-006	5 m		
	AD-010-002	3 m	Ø 5.0 mm	HF Monopolar connecting cable for **Erbe, ACC/ICC, Storz**
	AD-010-007	5 m		
	AD-010-003	3 m	Ø 8.0 mm	HF Monopolar connecting cable for **STEMA HF-U Series, Valleylab/Conmed**
	AD-010-008	5 m		
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable **male** for **STEMA MASTER**





Ordering No.	Description
EK-015-060	**Universal adapter** used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

13 Fr.

8.2 | PEDIATRIC RESECTOSCOPY
RESECTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-005-001	Pediatric cystoscope autoclavable	**0°**	Ø 2.7 mm	187.5 mm



Ordering No.	Description	French	Colour code
AH-003-001	Standard resectoscope sheath **with 1 rotatable stopcock,** without obturator	**13 Fr.**	Red



Ordering No.	Description	French	Colour code
AH-003-011	Standard obturator	**13 Fr.**	Red

13 Fr.

8.2 | PEDIATRIC RESECTOSCOPY
WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering No.	Description	For sheath
AG-003-001	Working element used for pediatric resectoscopes, **passive cutting**, handle open	**13 Fr.**



Ordering No.	Description	For sheath
AI-003-003	Endoscope bridge **with 3 Fr. instrument channel**	**13 Fr.**

13 Fr.

8.2 | PEDIATRIC RESECTOSCOPY
ELECTRODES



	Ordering No.	Description	French	Colour code
	AE-004-002	Cutting loop electrode angled	**13 Fr.**	Red
	AE-004-018	Cutting loop electrode small angled	**13 Fr.**	Red
	AE-004-004	Rollerball electrode	**13 Fr.**	Red
	AE-004-016	Knife electrode angled	**13 Fr.**	Red
	AE-004-019	Hook electrode without ball	**13 Fr.**	Red



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

9 Fr.
11 Fr.

8.2 | PEDIATRIC RESECTOSCOPY
RESECTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-004-011	Pediatric cystoscope autoclavable	**0°**	Ø 2.0 mm	187 mm



Ordering No.	Description	French	Colour code
AH-005-001	Standard resectoscope sheath **with 1 rotatable stopcock**, without obturator	**9 Fr.**	Blue
AH-004-001		**11 Fr.**	Green



Ordering No.	Description	French	Colour code
AH-005-011	Standard obturator	**9 Fr.**	Blue
AH-004-011		**11 Fr.**	Green

9 Fr.
11 Fr.

8.2 | PEDIATRIC RESECTOSCOPY
WORKING ELEMENTS



Ordering No.	Description	For sheath
AG-005-001	Working element used for pediatric resectoscopes, **passive cutting**, handle closed	**9 Fr.**
AG-004-001	Working element used for pediatric resectoscopes, **passive cutting**, handle closed	**11 Fr.**

9 Fr.
11 Fr.

8.2 | PEDIATRIC RESECTOSCOPY
ELECTRODES




	Ordering No.	Description	French	Colour code
	AE-006-002	Cutting loop electrode angled	9 Fr.	Blue
	AE-005-002		11 Fr.	Green
	AE-006-018	Cutting loop electrode small angled	9 Fr.	Blue
	AE-005-018		11 Fr.	Green
	AE-006-004	Rollerball electrode	9 Fr.	Blue
	AE-005-004		11 Fr.	Green
	AE-006-016	Knife electrode angled	9 Fr.	Blue
	AE-005-016		11 Fr.	Green
	AE-006-019	Hook electrode without ball	9 Fr.	Blue
	AE-005-019		11 Fr.	Green



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

8.2 | PEDIATRIC RESECTOSCOPY
CONNECTING CABLES | MONOPOLAR



Connection instrument side





Unit side	Ordering No.	Length	Diameter	Description
	AD-010-001 AD-010-006	3 m 5 m	Ø 4.0 mm	HF Monopolar connecting cable for **Berchtold, Martin, Aesculap**
	AD-010-002 AD-010-007	3 m 5 m	Ø 5.0 mm	HF Monopolar connecting cable for **Erbe, ACC/ICC, Storz**
	AD-010-003 AD-010-008	3 m 5 m	Ø 8.0 mm	HF Monopolar connecting cable for **STEMA HF-U Series, Valleylab/Conmed**
	AD-010-024	4.5 m	Ø 2.0 mm	HF Monopolar connecting cable **male** for **STEMA MASTER**





Ordering No.	Description
EK-015-060	**Universal adapter** used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

8.2 | PEDIATRIC RESECTOSCOPY
CONNECTING CABLES | MONOPOLAR



Connection unit side



Unit side	Ordering No.	Length	Description
	AD-010-025	4.5 m	HF Monopolar connecting cable for **Storz** resectoscopes, 90° angled connector
	AD-010-026	4.5 m	HF Monopolar connecting cable for **Olympus** resectoscopes, diameter Ø 3 mm



RECOMMENDED PRODUCTS:

CI-001-030 STEMA HF-U PRO TOUCH

- *For all monopolar & bipolar techniques*
- *with 7“ touch screen display*
- *provided to perform the most advanced electrosurgery*
- *ADC / APC System*

CI-001-040 DOUBLE PEDAL FOOTSWITCH

13 Fr.

8.3 | PEDIATRIC URETHROTOMY
CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-005-001	Pediatric cystoscope autoclavable	**0°**	Ø 2.7 mm	187.5 mm



Ordering No.	Description	French	Colour code
AI-003-001	Standard urethrotome sheath **with 2 fixed stopcocks** incl. obturator	**13 Fr.**	Red



Ordering No.	Description	French	Colour code
AI-003-004	Standard obturator	**13 Fr.**	Red

13 Fr.

8.3 | PEDIATRIC URETHROTOMY
WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering No.	Description	For sheath
AG-003-001	Working element used for pediatric resectoscopes, **passive cutting**, handle open	**13 Fr.**



Ordering No.	Description	For sheath
AI-003-003	Endoscope bridge **with 3 Fr. instrument channel**	**13 Fr.**

13 Fr.

8.3 | PEDIATRIC URETHROTOMY
COLD KNIVES



	Ordering No.	Description	French	Colour code
	AI-003-007	Cold knive straight	**13 Fr.**	Red
	AI-003-009	Cold knive round	**13 Fr.**	Red
	AI-003-010	Cold knive ring shaped	**13 Fr.**	Red
	AI-003-011	Cold knive hook shaped	**13 Fr.**	Red



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

11 Fr.

8.3 | PEDIATRIC URETHROTOMY
CYSTOSCOPES | SHEATHS | OBTURATORS



Ordering No.	Description	Viewing angle	Diameter	Working length
AA-004-011	Pediatric cystoscope autoclavable	**0°**	Ø 2.0 mm	187 mm



Ordering No.	Description	French	Colour code
AI-004-001	Standard urethrotome sheath **with 2 fixed stopcocks** incl. obturator	**11 Fr.**	Green



Ordering No.	Description	French	Colour code
AI-004-004	Standard obturator	**11 Fr.**	Green

11 Fr.

8.3 | PEDIATRIC URETHROTOMY
WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering No.	Description	For sheath
AG-004-001	Working element used for pediatric resectoscopes, **passive cutting**, handle closed	**11 Fr.**

11 Fr.

8.3 | PEDIATRIC URETHROTOMY
COLD KNIVES



	Ordering No.	Description	French	Colour code
	AI-004-007	Cold knive straight	**11 Fr.**	Green
	AI-004-009	Cold knive round	**11 Fr.**	Green
	AI-004-010	Cold knive ring shaped	**11 Fr.**	Green
	AI-004-011	Cold knive hook shaped	**11 Fr.**	Green



Ordering No.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

9. ACCESSORIES

9.1 | INSTRUMENTS
DILATING BOUGIES
METAL CATHETERS
ACCESSORIES CATHETERS

9.2 | STERILIZATION


9.3 | CLEANING AND MAINTENANCE

9. ACCESSORIES & INSTRUMENTS

1 DILATATORS



2 STERILIZATION



3 PROTECTION



STEMA

9.1 | INSTRUMENTS | DILATATOR BOUGIES







	VAN BUREN	DITTEL, STRAIGHT	DITTEL, CURVED
French	**Ordering no. VAN BUREN**	**Ordering no. DITTEL straight**	**Ordering no. DITTEL curved**
8 Fr.	AJ-003-001	AJ-003-021	AJ-003-050
9 Fr.		AJ-003-022	AJ-003-051
10 Fr.	AJ-003-002	AJ-003-023	AJ-003-052
11 Fr.		AJ-003-024	AJ-003-053
12 Fr.	AJ-003-003	AJ-003-025	AJ-003-054
13 Fr.		AJ-003-026	AJ-003-055
14 Fr.	AJ-003-004	AJ-003-027	AJ-003-056
15 Fr.		AJ-003-028	AJ-003-057
16 Fr.	AJ-003-005	AJ-003-029	AJ-003-058
17 Fr.		AJ-003-030	AJ-003-059
18 Fr.	AJ-003-006	AJ-003-031	AJ-003-060
19 Fr.		AJ-003-032	AJ-003-061
20 Fr.	AJ-003-007	AJ-003-033	AJ-003-062
21 Fr.		AJ-003-034	AJ-003-063
22 Fr.	AJ-003-008	AJ-003-035	AJ-003-064
23 Fr.		AJ-003-036	AJ-003-065
24 Fr.	AJ-003-009	AJ-003-037	AJ-003-066
25 Fr.		AJ-003-038	AJ-003-067
26 Fr.	AJ-003-010	AJ-003-039	AJ-003-068
27 Fr.		AJ-003-040	AJ-003-069
28 Fr.	AJ-003-011	AJ-003-041	AJ-003-070
29 Fr.		AJ-003-042	AJ-003-071
30 Fr.	AJ-003-012	AJ-003-043	AJ-003-072
31 Fr.			
32 Fr.	AJ-003-013	AJ-003-044	AJ-003-074
33 Fr.			
34 Fr.	AJ-003-014	AJ-003-045	
35 Fr.			
36 Fr.	AJ-003-015	AJ-003-046	
37 Fr.			
38 Fr.	AJ-003-016		

9.1 | INSTRUMENTS | DILATATOR BOUGIES | UTERINE DILATATORS





GUYON

French	Ordering no. GUYON	French	Ordering no. GUYON
8 Fr.	AJ-003-170	25 Fr.	AJ-003-187
9 Fr.	AJ-003-171	26 Fr.	AJ-003-188
10 Fr.	AJ-003-172	27 Fr.	AJ-003-189
11 Fr.	AJ-003-173	28 Fr.	AJ-003-190
12 Fr.	AJ-003-174	29 Fr.	AJ-003-191
13 Fr.	AJ-003-175	30 Fr.	AJ-003-192
14 Fr.	AJ-003-176	31 Fr.	AJ-003-193
15 Fr.	AJ-003-177	32 Fr.	AJ-003-194
16 Fr.	AJ-003-178	33 Fr.	AJ-003-195
17 Fr.	AJ-003-179	34 Fr.	AJ-003-196
18 Fr.	AJ-003-180	35 Fr.	AJ-003-197
19 Fr.	AJ-003-181	36 Fr.	AJ-003-198
20 Fr.	AJ-003-182	37 Fr.	AJ-003-199
21 Fr.	AJ-003-183	38 Fr.	AJ-003-200
22 Fr.	AJ-003-184	39 Fr.	AJ-003-201
23 Fr.	AJ-003-185	40 Fr.	AJ-003-202
24 Fr.	AJ-003-186	50 Fr.	AJ-003-203

HEGAR

French	Ordering no. HEGAR
3/4 Fr.	AJ-003-151
5/6 Fr.	AJ-003-152
7/8 Fr.	AJ-003-153
9/10 Fr.	AJ-003-154
11/12 Fr.	AJ-003-155
13/14 Fr.	AJ-003-156
15/16 Fr.	AJ-003-157
17/18 Fr.	AJ-003-158

9.1 | INSTRUMENTS | METAL CATHETERS







	FOR MEN	FOR WOMEN	COXETER PROSTATE
French	**Ordering no.**	**Ordering no.**	**Ordering no.**
6 Fr.		AJ-003-100	
7 Fr.	AJ-003-081		
8 Fr.	AJ-003-082	AJ-003-102	AJ-003-120
9 Fr.			
10 Fr.	AJ-003-084	AJ-003-104	AJ-003-121
11 Fr.			
12 Fr.	AJ-003-086	AJ-003-106	AJ-003-122
13 Fr.			
14 Fr.	AJ-003-088	AJ-003-108	AJ-003-123
15 Fr.			
16 Fr.	AJ-003-090	AJ-003-110	AJ-003-124
17 Fr.			
18 Fr.	AJ-003-092	AJ-003-112	AJ-003-125
19 Fr.			
20 Fr.	AJ-003-094	AJ-003-114	AJ-003-126
21 Fr.			
22 Fr.	AJ-003-096	AJ-003-116	AJ-003-127
23 Fr.			
24 Fr.	AJ-003-098	AJ-003-118	AJ-003-128

9.1 | INSTRUMENTS

CATHETER INTRODUCER



	Ordering No.	Description	Length
	AJ-003-142	Catheter introducer Guyon **straight**	39 cm
	AJ-003-143	Catheter introducer Guyon **curved**	39 cm

MEATUS CLAMPS

	Ordering No.	Description	Length
	AJ-003-140	Meatus clamp, **Stockmann**	8 cm
	AJ-003-141	Meatus clamp, **Strauss**	11 cm

SECRETION SPOON




	Ordering No.	Description	Length
	AJ-006-012	Secretion spoon **Guttmann**	26.5 cm


9.1 | INSTRUMENTS | TUBING CONNECTORS

	Ordering No.	Description
	AJ-010-011	Luer-Lock tubing connector, male Luer-Lock **inside Ø 3.0 mm**, **tubing connector Ø 6.0 mm**
	AJ-010-012	Luer-Lock tubing connector with **stopcock**
	AJ-010-013	Luer-Lock tubing connector, male Luer-Lock **inside**, **tubing connector Ø 9.0 mm**
	AJ-010-018	Stopcock left with Luer-Lock, male, **tubing connector Ø 9.0 mm**

9.2 | STERILIZATION | WIRE BASKETS



Ordering No.	Description	Dimension
FA-010-001	Wire basket for **1 endoscope** with lid	**290 × 80 × 52 mm**
FA-010-002		**460 × 80 × 52 mm**
FA-010-003		**670 × 80 × 52 mm**
FA-010-004	Wire basket for **2 endoscopes** with lid	**290 × 80 × 52 mm**
FA-010-005		**460 × 80 × 52 mm**



Ordering No.	Description	Dimension
FA-010-006	Wire basket for **1 endoscope** and **1 cold light cable** with lid	**640 × 150 × 77 mm**

9.2 | STERILIZATION

ENDOSCOPE CONTAINER



Ordering No.	Description	Colour code	Dimension
FA-001-001	Endoscope container, lid and bottom perforated (without filter)	Silver	**450 × 70 × 72 mm**
FA-001-001-R		Red	
FA-001-001-G		Green	
FA-002-001	Endoscope container, lid and bottom perforated (without filter	Silver	**250 × 60 × 50 mm**
FA-002-001-R		Red	
FA-002-001-G		Green	

PROTECTION TUBE





Ordering No.	Description
FA-001-000	Protection tube for endoscopes, plastic, 1 box = 10 pieces, **19 cm**

CAMERA SLEEVE



Ordering No.	Description
ES-007-001	Camera sleeve, 13 × 250 cm, PE hose film, sterile, 1 box = 100 pieces

9.3 | CLEANING AND MAINTENANCE

CLEANING PISTOL SET



Ordering No.	Description
EJ-030-010	Cleaning pistol complete inclusive **8 different tips**

CLEANING BRUSH



Ordering No.	Description	Diameter
EJ-030-001	Cleaning brush, **length 420 mm**, 1 box = 10 pieces	**Ø 3 mm**
EJ-030-002	Cleaning brush, **length 330 mm**, 1 box = 10 pieces	**Ø 5 mm**
EJ-030-003		**Ø 10 mm**

GREASE & PASTE





Ordering No.	Description	Weight
EJ-030-031	Special grease for **stopcocks, valves and instruments**	**3 g**
EJ-030-032	Polishing paste for scopes, 1 box = 50 pieces	**3 g**

ARE YOU INTERESTED IN OUR PRODUCTS? WE WILL GLADLY SUBMIT YOU AN OFFER.

OFFER FORM

Client No.:
(if available)

Company Name:

Order No.	Product	Quantity

Do you have a question about our products?

10. INDEX

AA-002-001 8, 37, 98, 147
AA-002-001-W 8, 37, 98
AA-002-002 8, 37
AA-002-003 8, 28, 37, 59, 147
AA-002-003-W 8, 37
AA-002-004 8, 37, 150
AA-002-004-W 8, 37
AA-002-011 8, 37, 98, 106, 147
AA-002-011-W 8, 37, 98
AA-002-012 8, 37, 59
AA-002-013 8, 28, 33, 37, 59, 76, 90, 91, 93, 147
AA-002-013-W 8, 37
AA-002-014 8, 37
AA-002-014-W 8, 37
AA-002-021 8, 37, 98
AA-002-022 8, 37
AA-002-023 8, 28, 37, 59
AA-002-024 8, 37
AA-003-001 103
AA-003-002 25, 48, 65, 70
AA-003-004 65
AA-003-011 103
AA-003-012 25, 48, 65, 70, 79, 83
AA-003-012-W 25
AA-003-014 65
AA-004-011 163, 170, 178
AA-004-012 163
AA-004-013 164
AA-004-014 164
AA-005-001 161, 167, 175
AA-005-003 161
AA-005-004 161
AA-006-001 114
AA-006-002 114
AA-006-004 114
AA-006-005 114
AA-006-006 116
AA-006-011 116
AA-006-012 116, 153
AA-006-021 114
AA-006-022 114
AA-006-031 116
AA-006-032 116
AA-006-034 116
AA-006-050 111
AA-006-051 111
AA-006-052 111
AA-006-053 111
AA-006-069 109
AA-006-070 109
AA-006-072 156
AA-006-080 114, 116, 153
AA-006-081 109, 114, 116, 153
AA-006-083 109, 114, 116, 153
AA-006-084 156
AA-006-100 110
AA-007-001 126
AA-007-002 126
AA-007-003 122
AA-007-004 122
AA-007-005 122
AA-007-006 125
AA-007-007 125
AA-007-008 125
AA-007-009 123
AA-007-011 127
AA-007-012 127
AA-007-013 124
AA-007-014 124
AA-007-015 124
AA-007-016 122
AA-007-017 122
AA-007-018 122
AA-007-019 122
AA-007-022 125
AA-007-023 125
AA-007-024 125
AA-007-025 125
AA-007-026 125
AA-007-027 125
AA-007-028 123
AA-007-029 123
AA-007-030 124
AA-007-031 124
AA-007-032 124
AA-007-033 124
AA-007-034 124
AA-007-035 124
AA-007-036 122
AA-007-037 122
AA-007-041 129
AA-007-042 129
AA-007-043 128
AA-007-044 128
AA-007-045 129
AA-007-046 129
AA-007-047 128
AA-007-048 128
AA-007-049 129
AA-007-050 129
AA-007-051 128
AA-007-052 128
AA-007-053 129
AA-007-054 129
AA-007-055 128
AA-007-056 128
AA-007-060 30
AA-008-002 31
AA-008-003 142
AA-008-100 32, 143
AA-008-101 143
AA-008-102 32, 143
AA-008-103 31, 32, 142
AA-008-104 32, 143
AA-008-105 31, 32, 142
AB-001-000 25
AB-001-001 161
AB-001-001-7 164
AB-001-001-9 163
AB-001-001-11 163
AB-001-001-13 161
AB-001-002 9
AB-001-003 9
AB-001-004 9, 33
AB-001-005 9
AB-001-006 9, 33
AB-001-007 25
AB-001-010 10
AB-001-011 161
AB-001-011-7 164
AB-001-011-9 163
AB-001-011-11 163
AB-001-012 10
AB-001-013 10
AB-001-014 10
AB-001-015 10
AB-001-016 10
AB-001-022 10
AB-001-023 10
AB-001-024 10
AB-001-025 10
AB-001-026 10
AB-001-040 28
AB-001-041 28
AB-001-042 29
AB-001-043 29
AB-002-001 9
AB-002-002 9
AB-002-003 9
AB-002-004 9
AB-002-005 9
AB-002-010 148
AB-002-011 149
AB-002-012 149
AB-002-013 148
AB-002-014 149
AB-003-001 162
AB-003-002 11
AB-003-003 11, 33
AB-003-004 11, 26
AB-003-005 162
AB-003-006 12
AB-003-007 12, 33
AB-003-008 12
AB-003-009 12
AB-003-010 26
AB-003-011 26
AB-003-012 26
AB-003-013 27
AB-003-014 27
AB-003-015 11
AB-003-020 27
AB-004-001 13, 34, 92, 94
AB-004-002 13
AB-004-003 13, 34, 92, 94
AB-004-004 13
AB-004-005 13
AB-004-006 13
AB-004-007 14, 34
AB-004-008 13, 92, 94
AB-004-009 14
AB-004-009-QL 14
AB-004-010 13
AB-004-010-QL 153
AB-004-011 150

10. INDEX

AB-004-012 150
AB-006-001 115
AB-006-002 115
AB-006-003 115
AB-006-004 115
AB-006-005 115
AB-006-006 115
AB-006-007 115
AB-006-008 115
AB-006-011 118
AB-006-012 118
AB-006-013 118
AB-006-020 118
AB-006-021 118
AB-006-022 118
AB-006-023 118
AB-006-024 118
AB-006-025 118
AB-006-026 118
AB-006-027 118
AB-006-040 119
AB-006-041 119
AB-006-042 119
AB-006-043 119
AB-006-044 119
AB-006-045 119
AB-006-046 119
AB-006-047 119
AB-006-048 119
AB-006-049 119
AB-006-050 119
AB-006-051 119
AB-006-052 118
AB-006-053 118
AB-006-061 117
AB-006-062 117
AB-006-063 117
AB-006-064 117
AB-006-065 117
AB-006-066 117
AB-006-067 117
AB-006-068 117
AB-006-100 112
AB-006-101 112
AB-006-102 112
AB-006-110 112
AB-006-111 112
AB-006-112 112
AB-006-115 112
AB-006-170 109
AB-006-171 109
AB-006-242 113
AB-006-243 113
AB-006-244 113
AB-006-245 113
AB-006-246 113
AB-006-342 113
AB-006-343 113
AB-006-344 113
AB-006-345 113
AB-006-346 113
AB-006-442 113
AB-006-443 113
AB-006-444 113
AB-006-445 113
AB-006-446 113
AB-007-001 87
AB-007-002 87
AB-007-003 87
AB-007-004 88
AB-007-005 88
AB-007-006 88
AB-007-007 88
AB-007-008 87
AC-001-001 15
AC-001-001-R 15
AC-001-002 15, 33
AC-001-002-R 15
AC-001-003 15
AC-001-003-R 15
AC-001-004 15
AC-001-004-R 15
AC-001-005 15
AC-001-005-R 15
AC-001-008 15
AC-001-008-R 15
AC-001-010 130
AC-001-010-R 130
AC-001-011 130
AC-001-011-R 130
AC-001-012 130
AC-001-012-R 130
AC-001-013 130
AC-001-013-R 130
AC-001-014 130
AC-001-014-R 130
AC-001-015 130
AC-001-015-R 130
AC-001-021 165
AC-001-031 144
AC-001-201 15, 130
AC-001-201-OT 19, 135, 137
AC-001-201-T 19, 135, 137
AC-001-202 15, 130
AC-001-202-OT 19, 135, 137
AC-001-202-T 19, 135, 137
AC-001-203 15, 130
AC-001-203-OT 19, 135, 137
AC-001-203-T 19, 135, 137
AC-001-204 15, 130
AC-001-204-OT 19, 135, 137
AC-001-204-T 19, 135, 137
AC-001-205 15, 130
AC-001-206-OT 19, 135, 137
AC-001-206-T 19, 135, 137
AC-001-208 15, 130
AC-001-210 130
AC-001-210-OT 135
AC-001-210-T 135
AC-001-211 130
AC-001-211-OT 135
AC-001-211-T 135
AC-001-212 130
AC-001-212-OT 135
AC-001-212-T 135
AC-001-213 130
AC-001-214 130
AC-001-214-OT 135
AC-001-214-T 135
AC-001-215 130
AC-001-215-OT 135
AC-001-215-T 135
AC-001-221 165
AC-002-001 15
AC-002-001-R 15
AC-002-002 15
AC-002-002-R 15
AC-002-003 15
AC-002-003-R 15
AC-002-004 15
AC-002-004-R 15
AC-002-005 15
AC-002-005-R 15
AC-002-011 130
AC-002-011-R 130
AC-002-012 130
AC-002-012-R 130
AC-002-013 130
AC-002-013-R 130
AC-002-014 130
AC-002-014-R 130
AC-002-015 130
AC-002-015-R 130
AC-002-201 15
AC-002-202 15
AC-002-203 15
AC-002-204 15
AC-002-205 15
AC-002-211 130
AC-002-212 130
AC-002-213 130
AC-002-214 130
AC-002-215 130
AC-003-000 17
AC-003-000-R 17
AC-003-001 17
AC-003-001-R 17
AC-003-002 17
AC-003-002-R 17
AC-003-003 17
AC-003-003-R 17
AC-003-004 17
AC-003-004-R 17
AC-003-005 17
AC-003-005-R 17
AC-003-010 132
AC-003-010-R 132
AC-003-011 132
AC-003-011-R 132
AC-003-012 132
AC-003-012-R 132
AC-003-013 132
AC-003-013-R 132
AC-003-014 132
AC-003-014-R 132
AC-003-015 132

10. INDEX

AC-003-015-R 132
AC-003-021 165
AC-003-200 17
AC-003-200-OT 19
AC-003-200-T 19
AC-003-201 17
AC-003-201-OT 19
AC-003-201-T 19
AC-003-202 17
AC-003-202-OT 19
AC-003-202-T 19
AC-003-203 17
AC-003-204 17
AC-003-204-OT 19
AC-003-204-T 19
AC-003-205 17
AC-003-205-OT 19
AC-003-205-T 19
AC-003-210 132
AC-003-210-OT 134
AC-003-210-T 134
AC-003-211 132
AC-003-211-OT 134
AC-003-211-T 134
AC-003-212 132
AC-003-212-OT 134
AC-003-212-T 134
AC-003-213 132
AC-003-214 132
AC-003-214-OT 134
AC-003-214-T 134
AC-003-215 132
AC-003-215-OT 134
AC-003-215-T 134
AC-003-221 165
AC-003-302 17
AC-003-304 17
AC-003-612-OT 19
AC-003-612-T 19
AC-003-614-OT 19
AC-003-614-T 19
AC-004-004 17
AC-004-004-R 17
AC-004-005 17
AC-004-005-R 17
AC-004-014 132
AC-004-015 132
AC-004-204 17
AC-004-204-OT 19
AC-004-204-T 19
AC-004-205 17
AC-004-205-OT 19
AC-004-205-T 19
AC-004-214 132
AC-004-214-OT 134
AC-004-214-T 134
AC-004-215 132
AC-004-215-OT 134
AC-004-215-T 134
AC-004-622-OT 134
AC-004-622-T 134
AC-004-624-OT 134
AC-004-624-T 134
AC-005-004 17
AC-005-004-R 17
AC-005-014 132
AC-005-014-R 132
AC-005-204 17
AC-005-214 132
AC-006-001 15
AC-006-001-R 15
AC-006-002 15
AC-006-002-R 15
AC-006-003 15
AC-006-003-R 15
AC-006-004 15
AC-006-004-R 15
AC-006-005 15
AC-006-005-R 15
AC-006-011 130
AC-006-011-R 130
AC-006-012 130
AC-006-012-R 130
AC-006-013 130
AC-006-013-R 130
AC-006-014 130
AC-006-014-R 130
AC-006-015 130
AC-006-015-R 130
AC-006-021 165
AC-006-031 144
AC-006-201 15
AC-006-201-OT 19
AC-006-201-T 19
AC-006-202 15
AC-006-202-OT 19
AC-006-202-T 19
AC-006-203 15
AC-006-203-OT 19
AC-006-203-T 19
AC-006-204 15
AC-006-204-OT 19
AC-006-204-T 19
AC-006-205 15
AC-006-205-OT 19
AC-006-205-T 19
AC-006-211 130
AC-006-211-OT 135
AC-006-211-T 135
AC-006-212 130
AC-006-212-OT 135
AC-006-212-T 135
AC-006-213 130
AC-006-213-OT 135
AC-006-213-T 135
AC-006-214 130
AC-006-214-OT 135
AC-006-214-T 135
AC-006-215 130
AC-006-215-OT 135
AC-006-215-T 135
AC-006-221 165
AC-007-001 16
AC-007-001-R 16
AC-007-002 16
AC-007-002-R 16
AC-007-011 131
AC-007-011-R 131
AC-007-012 131
AC-007-012-R 131
AC-007-201 16
AC-007-201-OT 19
AC-007-201-T 19
AC-007-202 16
AC-007-202-OT 19
AC-007-202-T 19
AC-007-211 131
AC-007-211-OT 135
AC-007-211-T 135
AC-007-212 131
AC-007-212-OT 135
AC-007-212-T 135
AC-008-004 16
AC-008-004-R 16
AC-008-005 16
AC-008-005-R 16
AC-008-014 131
AC-008-014-R 131
AC-008-015 131
AC-008-015-R 131
AC-008-204 16
AC-008-205 16
AC-008-214 131
AC-008-215 131
AC-009-001 16
AC-009-001-R 16
AC-009-002 16
AC-009-002-R 16
AC-009-004 16
AC-009-004-R 16
AC-009-011 131
AC-009-011-R 131
AC-009-012 131
AC-009-012-R 131
AC-009-014 131
AC-009-014-R 131
AC-009-201 16
AC-009-201-OT 19
AC-009-201-T 19
AC-009-202 16
AC-009-202-OT 19
AC-009-202-T 19
AC-009-204 16
AC-009-211 131
AC-009-211-OT 135
AC-009-211-T 135
AC-009-212 131
AC-009-212-OT 135
AC-009-212-T 135
AC-009-214 131
AC-012-001 16
AC-012-002 16
AC-012-004 16
AC-012-011 131
AC-012-012 131
AC-012-014 131

10. INDEX

AC-012-015 144
AC-015-002 15
AC-015-002-R 15
AC-015-003-R 15
AC-015-004 15
AC-015-004-R 15
AC-015-012 130
AC-015-012-R 130
AC-015-013-R 130
AC-015-014 130
AC-015-014-R 130
AC-015-021 165
AC-015-202 15
AC-015-203 15
AC-015-204 15
AC-015-212 130
AC-015-213 130
AC-015-214 130
AC-015-221 165
AC-016-001 16
AC-016-001-R 16
AC-016-002 16
AC-016-002-R 16
AC-016-011 131
AC-016-011-R 131
AC-016-012 131
AC-016-012-R 131
AC-016-201 16
AC-016-202 16
AC-016-211 131
AC-016-211-OT 135
AC-016-211-T 135
AC-016-212 131
AC-016-212-OT 135
AC-016-212-T 135
AC-017-001 15
AC-017-001-R 15
AC-017-004 15
AC-017-004-R 15
AC-017-011 130
AC-017-014 130
AC-017-014-R 130
AC-017-201 15
AC-017-204 15
AC-017-211 130
AC-017-214 130
AC-018-001 17
AC-018-001-R 17
AC-018-002 17
AC-018-002-R 17
AC-018-011 132
AC-018-011-R 132
AC-018-012 132
AC-018-012-R 132
AC-018-201 17
AC-018-202 17
AC-018-211 132
AC-018-212 132
AC-019-001 16
AC-019-001-R 16
AC-019-002 16
AC-019-002-R 16
AC-019-011 131
AC-019-011-R 131
AC-019-012 131
AC-019-012-R 131
AC-019-021 165
AC-019-201 16
AC-019-202 16
AC-019-211 131
AC-019-212 131
AC-020-002 16
AC-020-002-R 16
AC-020-004 16
AC-020-004-R 16
AC-020-012 131
AC-020-012-R 131
AC-020-014 131
AC-020-014-R 131
AC-020-202 16
AC-020-204 16
AC-020-212 131
AC-020-214 131
AC-100-001 19, 20, 134, 135, 136, 137
AC-101-201-OT 20
AC-101-201-T 20
AC-101-202-OT 20
AC-101-202-T 20
AC-101-203-OT 20
AC-101-203-T 20
AC-101-204-OT 20
AC-101-204-T 20
AC-101-206-OT 20
AC-101-206-T 20
AC-101-210-OT 137
AC-101-210-T 137
AC-101-211-OT 137
AC-101-211-T 137
AC-101-212-OT 137
AC-101-212-T 137
AC-101-214-OT 137
AC-101-214-T 137
AC-101-215-OT 137
AC-101-215-T 137
AC-103-200-OT 20
AC-103-200-T 20
AC-103-201-OT 20
AC-103-201-T 20
AC-103-202-OT 20
AC-103-202-T 20
AC-103-204-OT 20
AC-103-204-T 20
AC-103-205-OT 20
AC-103-205-T 20
AC-103-210-OT 136
AC-103-210-T 136
AC-103-211-OT 136
AC-103-211-T 136
AC-103-212-OT 136
AC-103-212-T 136
AC-103-214-OT 136
AC-103-214-T 136
AC-103-215-OT 136
AC-103-215-T 136
AC-103-612-OT 20
AC-103-612-T 20
AC-103-614-OT 20
AC-103-614-T 20
AC-104-204-OT 20
AC-104-204-T 20
AC-104-205-OT 20
AC-104-205-T 20
AC-104-214-OT 136
AC-104-214-T 136
AC-104-215-OT 136
AC-104-215-T 136
AC-106-201-OT 20
AC-106-201-T 20
AC-106-202-OT 20
AC-106-202-T 20
AC-106-203-OT 20
AC-106-203-T 20
AC-106-204-OT 20
AC-106-204-T 20
AC-106-205-OT 20
AC-106-205-T 20
AC-106-211-OT 137
AC-106-211-T 137
AC-106-212-OT 137
AC-106-212-T 137
AC-106-213-OT 137
AC-106-213-T 137
AC-106-214-OT 137
AC-106-214-T 137
AC-106-215-OT 137
AC-106-215-T 137
AC-107-201-OT 20
AC-107-201-T 20
AC-107-202-OT 20
AC-107-202-T 20
AC-107-211-OT 137
AC-107-211-T 137
AC-107-212-OT 137
AC-107-212-T 137
AC-109-201-OT 20
AC-109-201-T 20
AC-109-211-OT 137
AC-109-211-T 137
AC-109-212-OT 137
AC-109-212-T 137
AC-114-212-OT 136
AC-114-212-T 136
AC-116-211-OT 137
AC-116-211-T 137
AC-116-212-OT 137
AC-116-212-T 137
AD-001-001 23, 140
AD-001-002 23, 140
AD-001-003 23, 34, 140
AD-001-004 23, 140
AD-001-005 23, 140
AD-001-006 23, 140
AD-001-007 23, 140
AD-001-008 23, 140
AD-001-011 23, 140

10. INDEX

AD-001-012 23, 140
AD-001-013 23, 140
AD-001-014 23, 140
AD-001-015 23, 140
AD-001-016 23, 140
AD-001-017 23, 140
AD-001-018 23, 140
AD-001-021 165
AD-002-001 23, 140
AD-002-002 23, 140
AD-002-003 23, 34, 140
AD-002-005 23, 140
AD-002-011 23, 140
AD-002-012 23, 140
AD-002-013 23, 140
AD-002-015 23, 140
AD-002-021 165
AD-003-003 23, 140
AD-003-005 23, 140
AD-003-007 23, 140
AD-003-008 23, 140
AD-003-013 23, 140
AD-003-015 23, 140
AD-003-017 23, 140
AD-003-018 23, 140
AD-010-001 24, 57, 92, 141, 166, 173
AD-010-002 24, 57, 141, 166, 173
AD-010-003 24, 34, 57, 141, 166, 173
AD-010-006 24, 57, 141, 166, 173
AD-010-007 24, 57, 141, 166, 173
AD-010-008 24, 57, 141, 166, 173
AD-010-011 75
AD-010-012 75
AD-010-013 75, 94
AD-010-015 75
AD-010-024 24, 57, 141, 166, 173
AD-010-025 58, 174
AD-010-026 58, 174
AE-001-001 46, 91
AE-001-002 46
AE-001-003 46
AE-001-004 46, 91
AE-001-005 46
AE-001-006 47
AE-001-007 47
AE-001-008 47
AE-001-009 47
AE-001-010 46
AE-001-011 46
AE-001-012 46
AE-001-013 46
AE-001-014 47
AE-001-015 47
AE-001-016 47, 91
AE-001-017 47
AE-001-030 63
AE-001-031 63, 93
AE-001-032 63, 93
AE-001-033 63
AE-001-034 64, 93
AE-001-035 64
AE-001-036 63
AE-001-037 63
AE-001-038 63
AE-001-039 64
AE-001-040 63
AE-001-041 63
AE-001-042 63
AE-001-043 63
AE-001-044 64
AE-001-045 64
AE-001-046 64
AE-001-047 64
AE-001-048 63
AE-002-001 46
AE-002-002 46
AE-002-003 46
AE-002-004 46
AE-002-005 46
AE-002-006 47
AE-002-007 47
AE-002-010 46
AE-002-011 46
AE-002-012 46
AE-002-013 46
AE-002-014 47
AE-002-015 47
AE-002-016 47
AE-002-031 64
AE-003-001 55
AE-003-003 55
AE-003-004 55
AE-003-010 55
AE-003-014 55
AE-003-016 55
AE-003-031 69
AE-003-032 69
AE-003-033 69
AE-003-034 69
AE-003-036 69
AE-003-040 69
AE-003-041 74
AE-003-043 74
AE-003-044 74
AE-003-045 74
AE-003-046 74
AE-004-002 169
AE-004-004 169
AE-004-016 169
AE-004-018 169
AE-004-019 169
AE-005-002 172
AE-005-004 172
AE-005-016 172
AE-005-018 172
AE-005-019 172
AE-006-002 172
AE-006-004 172
AE-006-016 172
AE-006-018 172
AE-006-019 172
AF-001-001 21, 138
AF-001-002 21, 138
AF-001-003 21, 138
AF-001-004 21, 138
AF-001-005 21, 138
AF-001-006 21, 138
AF-001-007 22, 139
AF-001-008 22, 34, 139
AF-001-009 22, 139
AF-001-010 22, 139
AF-001-011 22, 139
AF-001-012 22, 139
AF-001-013 22, 139
AF-001-014 22, 139
AF-001-015 21, 138
AF-001-016 21, 138
AF-001-017 21, 138
AF-001-018 22, 139
AF-001-019 21, 138
AF-001-021 21, 22, 138, 139
AF-001-022 21, 138
AF-001-023 21, 138
AF-001-024 21, 22, 138, 139
AF-001-025 22, 139
AG-001-001 38, 91, 101, 106
AG-001-001-HF 38, 101
AG-001-002 38, 101
AG-001-003 39
AG-001-003-HF 39
AG-001-004 39
AG-001-006 77
AG-001-007 40
AG-001-008 40
AG-001-009 77
AG-001-010 104
AG-001-011 60, 93
AG-001-012 60
AG-001-017 41
AG-001-018 41
AG-002-001 49
AG-002-002 49
AG-002-003 50
AG-002-004 50
AG-002-006 51
AG-002-007 80
AG-002-008 80
AG-002-011 66
AG-002-012 66
AG-002-013 71
AG-002-015 71
AG-002-016 84
AG-003-001 168, 176
AG-004-001 171, 179
AG-005-001 171
AG-100-006 77
AG-100-007 77
AH-001-001 44
AH-001-002 44
AH-001-011 45, 62, 78, 90, 91, 93, 155
AH-001-012 45, 62
AH-001-013 45, 62, 155
AH-001-015 45, 62
AH-001-016 45, 62
AH-001-017 45, 62

10. INDEX

AH-001-018 45, 62
AH-001-101 44
AH-001-102 44
AH-001-201 42
AH-001-202 42
AH-001-203 42
AH-001-204 42
AH-001-205 42
AH-001-211 42
AH-001-212 42
AH-001-301 43, 61, 78, 90, 91, 93, 154, 155
AH-001-302 43, 61, 78, 91, 154
AH-001-303 43, 61
AH-001-304 43, 61
AH-001-305 43, 61
AH-001-306 43, 61
AH-001-308 43, 78, 155
AH-001-311 43, 61
AH-001-312 43, 61
AH-002-001 54
AH-002-011 56, 68, 82, 158
AH-002-013 56, 68, 82, 158
AH-002-015 73, 86
AH-002-016 73, 86
AH-002-201 52
AH-002-211 72, 85
AH-002-212 52
AH-002-213 72, 85
AH-002-216 72, 85
AH-002-301 53, 67, 81, 157
AH-002-302 53, 67, 81, 157
AH-003-001 167
AH-003-011 167
AH-004-001 170
AH-004-011 170
AH-005-001 170
AH-005-011 170
AI-001-001 98, 106
AI-001-002 99, 106
AI-001-003 98, 106
AI-001-004 100, 106
AI-001-005 100
AI-001-006 99
AI-001-007 102, 106
AI-001-008 102
AI-001-009 102
AI-001-011 102
AI-001-012 100
AI-002-001 97
AI-002-002 97
AI-003-001 175
AI-003-003 168, 176
AI-003-004 175
AI-003-007 177
AI-003-009 177
AI-003-010 177
AI-003-011 177
AI-004-001 178
AI-004-004 178
AI-004-007 180
AI-004-009 180
AI-004-010 180
AI-004-011 180
AI-005-001 103
AI-005-003 103
AI-005-004 103
AI-005-007 105
AI-005-008 105
AI-005-011 105
AJ-001-000 89
AJ-001-001 89
AJ-001-002 89, 92, 94
AJ-001-003 89
AJ-001-010 89
AJ-001-011 89
AJ-001-012 89
AJ-001-013 89
AJ-002-001 89, 92
AJ-002-002 89
AJ-002-003 89
AJ-002-004 89
AJ-002-005 89
AJ-003-001 183
AJ-003-002 183
AJ-003-003 183
AJ-003-004 183
AJ-003-005 183
AJ-003-006 183
AJ-003-007 183
AJ-003-008 183
AJ-003-009 183
AJ-003-010 183
AJ-003-011 183
AJ-003-012 183
AJ-003-013 183
AJ-003-014 183
AJ-003-015 183
AJ-003-016 183
AJ-003-021 183
AJ-003-022 183
AJ-003-023 183
AJ-003-024 183
AJ-003-025 183
AJ-003-026 183
AJ-003-027 183
AJ-003-028 183
AJ-003-029 183
AJ-003-030 183
AJ-003-031 183
AJ-003-032 183
AJ-003-033 183
AJ-003-034 183
AJ-003-035 183
AJ-003-036 183
AJ-003-037 183
AJ-003-038 183
AJ-003-039 183
AJ-003-040 183
AJ-003-041 183
AJ-003-042 183
AJ-003-043 183
AJ-003-044 183
AJ-003-045 183
AJ-003-046 183
AJ-003-050 183
AJ-003-051 183
AJ-003-052 183
AJ-003-053 183
AJ-003-054 183
AJ-003-055 183
AJ-003-056 183
AJ-003-057 183
AJ-003-058 183
AJ-003-059 183
AJ-003-060 183
AJ-003-061 183
AJ-003-062 183
AJ-003-063 183
AJ-003-064 183
AJ-003-065 183
AJ-003-066 183
AJ-003-067 183
AJ-003-068 183
AJ-003-069 183
AJ-003-070 183
AJ-003-071 183
AJ-003-072 183
AJ-003-074 183
AJ-003-081 185
AJ-003-082 185
AJ-003-084 185
AJ-003-086 185
AJ-003-088 185
AJ-003-090 185
AJ-003-092 185
AJ-003-094 185, 190
AJ-003-096 185
AJ-003-098 185
AJ-003-100 185
AJ-003-102 185
AJ-003-104 185
AJ-003-106 185
AJ-003-108 185
AJ-003-110 185
AJ-003-112 185
AJ-003-114 185
AJ-003-116 185
AJ-003-118 185
AJ-003-120 185
AJ-003-121 185
AJ-003-122 185
AJ-003-123 185
AJ-003-124 185
AJ-003-125 185
AJ-003-126 185
AJ-003-127 185
AJ-003-128 185
AJ-003-140 186
AJ-003-141 186
AJ-003-142 186
AJ-003-143 186
AJ-003-151 184
AJ-003-152 184
AJ-003-153 184
AJ-003-154 184

10. INDEX

AJ-003-155 184
AJ-003-156 184
AJ-003-157 184
AJ-003-158 184
AJ-003-170 184
AJ-003-171 184
AJ-003-172 184
AJ-003-173 184
AJ-003-174 184
AJ-003-175 184
AJ-003-176 184
AJ-003-177 184
AJ-003-178 184
AJ-003-179 184
AJ-003-180 184
AJ-003-181 184
AJ-003-182 184
AJ-003-183 184
AJ-003-184 184
AJ-003-185 184
AJ-003-186 184
AJ-003-187 184
AJ-003-188 184
AJ-003-189 184
AJ-003-190 184
AJ-003-191 184
AJ-003-192 184
AJ-003-193 184
AJ-003-194 184
AJ-003-195 184
AJ-003-196 184
AJ-003-197 184
AJ-003-198 184
AJ-003-199 184
AJ-003-200 184
AJ-003-201 184
AJ-003-202 184
AJ-003-203 184
AJ-006-012 186
AJ-010-001 11, 26, 111, 122, 123, 124, 125, 126, 127, 128, 129
AJ-010-003 26, 111, 122, 123, 124, 125, 126, 127, 128, 129
AJ-010-006 47, 55, 64, 69, 74, 90, 92, 94, 102, 105, 106, 169, 172, 177, 180
AJ-010-011 187
AJ-010-012 187
AJ-010-013 187
AJ-010-015 114, 116
AJ-010-018 187
AJ-010-020 90
CI-001-030 58, 174
CI-001-040 58, 174
EJ-030-001 190
EJ-030-002 190
EJ-030-003 190
EJ-030-010 190
EJ-030-031 190
EJ-030-032 190
EK-015-060 24, 57, 141, 166, 173
ES-007-001 189
FA-001-000 189
FA-001-001 189
FA-001-001-G 189
FA-001-001-R 189
FA-002-001 189
FA-002-001-G 189
FA-002-001-R 189
FA-010-001 188
FA-010-002 188
FA-010-003 188
FA-010-004 188
FA-010-005 188
FA-010-006 188

WE COVER THE FOLLOWING AREAS IN SURGERY WITH HIGH QUALITY PRODUCTS

ARTHROSCOPY
CARDIAC MIS
CONTAINER
ENT
HF-SURGERY
GASTROSCOPY
GYNAECOLOGY
LAMINECTOMY
LAPAROSCOPY
PROCTOSCOPY
SPINE SURGERY
UNITS
UROLOGY

STEMA Medizintechnik GmbH | Take-off Gewerbepark 74 | 78579 Neuhausen ob Eck
info@stema-medizintechnik.de | www.stema-medizintechnik.de | +49 7467 79630 0

© 2026 STEMA Medizintechnik GmbH. All rights reserved.
Specifications are subject to change without notice.